

Standardized Extractables testing of polymeric materials enabling effective safety Risk Assessment

Patrick EVRARD, Sr Director SLS Single-use Technologies, Pall Corp.

- Introduction and Regulatory expectations
- BPOG standard proposal
- BPOG vs USP <665>
- Risk Assessment using single-use supplier's Datasets
- Conclusions

REASONS FOR RESTRICTING USE OF SUT




Leachables and extractables are a concern



Breakage of bags and loss of production material is a concern



High cost of disposables (consumables)




Legend




Source: 12th Annual Report and Survey of Biopharmaceutical Manufacturing, April 2015, preliminary data, www.bioplanassociates.com/11th

- Antioxidants
- Stabilizers
- Molding agents
- Extrusion agents
- Polymerization aids
- Pore formers
- Colorants
- Lubricants
- Residual solvents
- Unreacted monomers
- Oligomers
- Degradation products



WHY ARE LEACHABLES A CONCERN?



- US CFR 211.65^a

“Equipment shall be constructed so that surfaces that contact components, in-process materials, or drug products shall not be reactive, additive, or absorptive so as to alter the safety, identity, strength, quality, or purity of the drug product beyond the official or other established requirements.”

- EC Eudralex^b

§ 3.39: “Production equipment should not present any hazard to products. Parts of production equipment that come into contact with the product must not be reactive, additive or absorptive to such an extent that it will affect the quality of the product and thus present any hazard”

a. US Code of Federal Regulations, Part 211, “Current Good Manufacturing Practice for Finished Pharmaceuticals”, Sec. 211.65, “Equipment Construction”, 2015

b. European Commission, EUDRALEX Volume 4, “Good Manufacturing Practices, Medicinal Products for Human and Veterinary Use”, Vol 4, Chapter 3, “Premises and Equipment”, 2014

- *FDA - Guidance for Industry
Container Closure Systems for Packaging Human Drugs and Biologics, U.S.
Department of Health and Human Services, FDA, May 1999*
 - Attachment C – “Extraction studies may be conducted ...”
 - **GENERIC EXTRACTABLES** *Tests on plastics (USP <661>) or elastomers (USP <381>)*
 - **EXTRACTABLES PROFILES** *Biological Reactivity Tests (USP <87> and <88>)*
 - **EXTRACTABLES PROFILES** *Relative extraction profiles of plastics or elastomers*
 - **EXTRACTABLES PROFILES** *Relative extraction profiles of plastics or elastomers*
- *EMA - Guideline on plastic immediate packaging materials,
CPMP/QWP/4359/03, May 2005*

*“...with respect to extractable and leachable data ... it is ultimately your responsibility to assess this data and its applicability to your products and process. CBER recommends a **risk-based approach** be taken in evaluating extractables and leachables where you take multiple aspects into account (e.g., indication, safety issues, product characteristics, dosage, formulation, stability profile, etc.).”*

“If there is **no relevant risk** associated with the (material in question), **vendor data can be cross-referenced** and a detailed justification for the applicability of these data and a justification for additional testing should be submitted.”

Destry M. Sullivan, FDA CBER
IBC Single Use Conference, June, 2010

GENERIC




NVR



TOC


PROFILES



Direct Inj.
GC\MS



Headspace
GC\MS



LC\UV\MS

Hg


Cd



ICP\MS



- Industry group
- Issue guidance documents to facilitate industry's topics of concern




STANDARDIZED EXTRACTABLES DATA

Aimed to bracket >80% Biologics Applications

Standardized Extractables Testing Protocol for Single-Use Systems in Bio

by Weibing Ding, Gary Madsen, Ekta Ma



* Duration for testing storage bags necessary to support 3-year storage at 0 °C ($Q_{10}=2.0$)


IMPLEMENTATION OF BPOG PROTOCOL

Component
Manufacture

Extraction

Analytical
Testing


Reporting



PALL BPOG REPORTS



Technical Report



Descriptive Summary

Autoclaved
Gamma irradiated

STERILE CONNECTOR (autoclave)
STERILE FILTER EKV
STERILE FILTER DFL
STERILE CONNECTOR (gamma)
STERILE FILTER EKV
STERILE FILTER DFL

SLS Global Technical Support

Summary Report

Characterization of Extractables from Gamma Irradiated Pail Filter Connectors (APCI)

Extraction Report

Author: SLS, SLS, SLS, SLS
Reviewer: SLS, SLS, SLS, SLS
Date: 08/15/2016
Version: 1.0

Pall Corporation

Project Number: VAL-AB-003136
Issue Date: June 07, 2016
Page 2 of 21

SUMMARY

Pail Allegro biocompatibles are used under a wide range of process conditions in the commercial production of human and animal biopharmaceutical products such as vaccines and human therapeutics. This report summarizes the results from the extractables studies on the Allegro film performed in accordance with published end-user recommendations to support risk assessments of potential impurities that may migrate from process equipment into a fluid process stream. The extractables study is designed and intended to represent typical "best-case" process conditions, which may exaggerate actual "in-process" or "final product" leachables.

To characterize the extractables profile of the Allegro film material, two lots of Pail Allegro 1-liter biocompatibles (PIL 619-21A) were manufactured under identical manufacturing conditions with the exception that the HDPE port was not incorporated (leaked separately), and the biocompatibles were gamma irradiated at 25 x10⁵ Gy. A small portion of the biocompatible was cut to permit filling with solvent, and then heat sealed to close the biocompatible. The biocompatibles were extracted (dissolved) at approximately 40°C in 50% Ethanol and 1% Polysorbate 80 and 4 incubation conditions (30 minutes at 25°C, 24 hours at 40°C, 21 days at 40°C and 70 days at 40°C). The general material properties of the film and the extraction conditions employed in this study are summarized in Table 1.

Organic extractables analyses were performed on all samples by Direct Injection GC/MS, Headspace GC/MS, and LCP/DAMS (ESI+, APCI+), with a grand summary of these data shown in Table 2. All compounds detected at > 0.1 µg/cm² level are reported. Majority of the extractables detected were < 0.1 µg/cm². Extractables found at > 0.1 µg/cm² are mostly related to extractables degradation. In particular, BQZ-2,6-di-tert-butylphenolacetate, an inactive 180 degradation product, was detected at 0.58 µg/cm² 50% Ethanol and 70%PS80 at no more than 0.15 µg/cm² ICH Q3C residual solvents Hexane, Cyclohexane, Ethyl acetate and 1,1-difluoroethane were detected at < 0.18 µg/cm² whereas Biphenyl A (BPA) was not detected at all in any of the samples. All 7-day samples were assayed for 35 elemental impurities including all ICH Q3D elements by ICP-MS. None of the ICH Q3D elements were detected at > 20 ppb level. A grand summary of the ICP-MS results is presented in Table 3.

The detailed summaries of the extractables data and assay parameters for all analytical methods are presented in Tables 4 - 5 and Tables 6 - 9 (Appendix), respectively.

¹ Ding et al. Standardized Extractables Testing Protocol for Single-Use Systems in Biomanufacturing, 2014. Pharmaceutical Engineering 24 (6), pages 1 - 11.

PALL Corporation www.pall.com

We are ISO 9001:2015 certified. The Pall logo is your assurance of quality.


Extraction Parameters

Test Article	Allegro Biocompatible		
Number of Test Articles	48		
Part Number	619-21A (manufactured without HDPE port)		
Lot Numbers	6010363, 6010455		
Pretreatment	Variables	Units	Values
Gamma irradiation	Dose	kGy	50 ± 5
Autoclave	Time, temperature, number of cycles	minutes, °C, #	N/A
Pre-Rush	Fluid identity, duration, temperature, volume	Name, minutes, °C, L	None
Test Article Extraction Conditions	Variable	Units	Values
	Temperature	°C	25, 40, 40, 40
	Duration	minutes, hours, days	30 minutes, 24 hours, 21 days, 70 days
	Solvent contact surface area	sq cm, cm ²	813
	Solvent volume	mL	130
Surface area to volume ratio	ratio		6.3
	Solvent Loss (pre- and post-weight)	% loss	< 15 %
Supporting Information			
Bags	Film thickness	mm	0.325
	Volume (capacity)	L	1
Tubing	Internal diameter	mm	N/A
	Length	mm	N/A
	Internal diameter	mm	N/A
Tubing connectors and aseptic connectors	Length	mm	N/A
	Internal diameter	mm	N/A
Filters and TFF cassettes	EPA	sq m (µm ²)	N/A
	Internal diameter	mm	N/A
Filling Numbers	Time between film manufacturing and extraction	days	< 60
	Time between gamma irradiation and extraction	days	< 35
All gamma-irradiated components	Type of dose range during normal manufacturing	kGy	27 - 37

Analytical Method Details

Standards	2-Propanol and Methyl ethyl ketone (MEK)		
Limit of Detection	Signal-to-noise ratio ≥ 3 (all standards)		
Limit of Quantification	Signal-to-noise ratio ≥ 10 (all standards)		
Internal Standard	None		
Precision	RSD ≤ 2 % (Retention Time) and ≤ 20% (Area Counts) (n = 6), 1 - 10 ppm (all standards in all solvents)		
Spike Recovery	80 - 120%, 1 - 8 ppm (all standards in all solvents)		
Sample Bracketing	Every 10 sample injections (all standards), % Difference ≤ 20%		
Column	Agilent DB-624 MS, 60m x 0.25 mm, 1.4 µm		
Injection Port Temperature	250 °C		
Headspace Vial Temperature	70 °C		
Injection Volume	1.0 µL		
GC Temperature Program	Rate (°C/min)	Temperature (°C)	Hold Time (min)
	N/A	40.0	4.00
	5.0	50.0	5.00
	5.0	85.0	5.00
	15.0	200.0	5.00
m/z range, amu	35 - 650		

Images of Extraction Process



Summarized Organic Extractables Data (µg/cm²)

Solvent	Method	Mode	RT (min)	Compound	CAS No.	ID Level	ICH Q3C Class	Standard Used for Semi-Quantification	LOD (µg/mL)	LOQ (µg/mL)	Worst-case, µg/cm ²					
											1/2 hr	24 hrs	21 days	70 days		
Water	HS GC/MS	TIC	8.77	2-Methyl-2-propanol	75-65-0	Confirmed	N/A	2-Methyl-2-propanol	0.205	0.685	<LOD	<LOD	<LOD	<LOD		
	HS GC/MS	TIC	27.43	3,4-Dimethyl-3-hexanol	19550-08-4	Tentative	N/A	2-Methyl-2-propanol	0.205	0.685	<LOD	<LOD	<LOD	<LOD		
	LCP/DAMS	ES+, TIC	8.99	Unknown (mass ion 587.6)	N/A	Unknown	N/A	Encrucamide	0.000	0.301	<LOD	<LOD	<LOD	<LOD		
	HS GC/MS	TIC	12.86	Ethyl Acetate	141-78-6	Confirmed	3	Ethyl acetate	0.041	0.138	<LOD	<LOD	<LOD	<LOD		
	HS GC/MS	TIC	19.64	Pentanal	110-62-3	Confirmed	N/A	Pentanal	0.032	0.107	<LOD	<LOD	<LOD	<LOD		
	DI GC/MS	TIC	N/A	None Detected	N/A	N/A	N/A	n-Decane	0.032	0.106	<LOD	<LOD	<LOD	<LOD		
	HS GC/MS	TIC	8.74	2-Methyl-2-propanol	75-65-0	Confirmed	N/A	2-Methyl-2-propanol	0.017	0.057	<LOD	0.0362	0.0378	0.0530		
	HS GC/MS	TIC	12.76	Ethyl acetate	141-78-6	Confirmed	3	Ethyl Acetate	0.026	0.088	<LOD	0.0168	<LOD	<LOD		
	5M NaCl	LCP/DAMS	PDA (210-400 nm), ES+, APCI +/-	N/A	None Detected	N/A	N/A	N/A	BPA (PDA), DEHP (PDA), Encrucamide (ES+), Irganox 1010 (ES-)	0.036, 0.120, 0.211, 0.703, 0.014, 0.047, 0.018, 0.060, Irganox 1010 (APC+APCI-)	0.489, 1.79	1.63, 5.98	<LOD	<LOD	<LOD	<LOD
		DI GC/MS	TIC	N/A	None Detected	N/A	N/A	N/A	n-Decane	0.039, 0.132, 0.025	0.094	<LOD	<LOD	<LOD	<LOD	
HS GC/MS		TIC	8.79	2-Methyl-2-Propanol	75-65-0	Confirmed	N/A	2-Methyl-2-propanol	0.147	0.490	<LOD	<LOD	<LOD	0.0811		
0.1M H ₂ O ₂	DI GC/MS	TIC	6.87	Hexanoic acid	142-62-1	Confident	N/A	2-Octanone	0.001	0.002	<LOD	0.00829	0.0185	0.0191		
	HS GC/MS	TIC	19.63	Pentanal	110-62-3	Confirmed	N/A	Pentanal	0.050	0.166	<LOD	<LOD	<LOD	<LOD		
	LCP/DAMS	PDA (210-400 nm), ES+, APCI +/-	N/A	None Detected	N/A	N/A	N/A	BPA (PDA), DEHP (PDA), Encrucamide (ES+), Irganox 1010 (ES-)	0.161, 0.536, 0.389, 1.30, 0.077, 0.256, 0.062, 0.207, Irganox 1010 (APC+APCI-)	4.68, 8.69	15.6, 28.0	<LOD	<LOD	<LOD	<LOD	





Summarized Inorganic Extractables Data

Element	ICH Q3D Class	Water	Worst case, µg/cm ² (70 days)				
			(M NaCl)	0.1M H ₂ O ₂	0.5M NaCl	50% Ethanol	1% PS80
Ca	3	-	-	-	-	-	
Pb	1	-	-	-	-	-	
Ag	1	-	-	-	-	-	
Cd	1	-	-	-	-	-	
Ni	2A	-	-	-	-	-	
Al	2B	-	-	-	-	-	
As	2B	-	-	-	-	-	
B	2B	-	-	-	-	-	
Br	2B	-	-	-	-	-	
Bi	2B	-	-	-	-	-	
Ce	2B	-	-	-	-	-	
Ag	2B	-	-	-	-	-	
Cr	2B	-	-	-	-	-	
Co	3	-	-	-	-	-	
Fe	3	-	-	-	-	-	
Mn	3	-	-	-	-	-	
Cu	3	-	-	-	-	-	
Sr	3	-	-	-	-	-	
Zn	3	-	-	-	-	-	
Ba	3	-	-	-	-	-	
Na	N/A	0.204	-	0.504	-	-	
Mo	N/A	0.020	-	0.050	-	-	
W	N/A	-	-	-	-	-	
Mg	N/A	-	-	-	-	-	
Al	N/A	-	-	-	-	-	
Cg	N/A	-	-	-	-	-	
Ti	N/A	-	-	-	-	-	
Mn	N/A	-	-	-	-	-	
Pb	N/A	-	-	-	-	-	
Zn	N/A	-	-	-	-	-	
Cr	N/A	-	-	-	-	-	

Notes: "-" indicates < 20 ppb; N/A - not analyzed due to interference with matrix.

Includes all ICH Q3D plus additional elements

IMPACT OF PRE-TREATMENT

Study	Sterilization	Total Compounds
Sterile connector (24 hr)	Autoclaved	
	Gamma irradiated	
PES filter (24 hr)	Autoclaved	
	Gamma irradiated	

↑ ↑

| Common | Unique
compound compound


BPOG RISK ASSESSMENT GUIDANCE




• Risk Assessment model

- Distance along production stream
- Exposure temperature
- Exposure duration
- Process Fluid Interaction
- Dilution ratio

⇒ Risk Level



From "BPOG Best Practices Guide for Evaluating Leachables Risk in Biopharmaceutical Single-Use Systems"



SCALING DATA FOR E/L EVALUATION



Technical Report



Single-use system

$\mu\text{g}/\text{cm}^2$

Scaling to Different Size Biocontainers

B₅₀₀₀ AP
 10⁶ CFU Challenge Level
 0.05 μm

Allegro Film
 Test Type: Extractables
 Revision: 01

Project Number: VAL-AS-003176
 January 27, 2017
 Page 22 of 24

APPENDIX 2 (Supplemental)

In order to facilitate scaling of extractables data to various sized components, extractables data are reported in units of $\mu\text{g}/\text{mL}$ and summarized in units of $\mu\text{g}/\text{cm}^2$ of fluid contact surface area. The tables below show approximate fluid contact surface areas for typical Allegro biocontainers. Example calculations are shown further below.


2D Allegro Biocontainers

Product	Total Volume (L)	Film Contact Area (sq. cm.)	Port Contact Area (sq. cm.)	Part Number
50mL	0.05	180	24	LGR0050ML770
500mL	0.5	523	24	LGR0500ML770
1L	1	894	24	LGR1000ML770
5L	5	2331	44	LGR0005L6600, LGR0005L7700, LGR0005L8800, LGR0005L8800TF
10L	10	3758	44	LGR0010L6600, LGR0010L7700, LGR0010L8800, LGR0010L8800TF
20L	20	5881	44	LGR0020L6600, LGR0020L6600CAP, LGR0020L7700, LGR0020L8800
50L	50	10283	44	LGR0050L6600, LGR0050L6600CAP, LGR0050L7700, LGR0050L8800

3D Allegro Biocontainers

Product	Total Volume (L)	Film Contact Area (sq. cm.)	Port Plate Area (sq. cm.)	Drain Flange Area (sq. cm.)	Part Number
100L	100	14556	117	97	LGR0100L5550, LGR0100L6600, LGR0100L8880, LGR0100L9990
200L	200	21500	117	97	LGR0200L5550, LGR0200L6600, LGR0200L8880, LGR0200L9990
500L	500	38307	117	97	LGR0500L5550, LGR0500L6600, LGR0500L8880, LGR0500L9990
1000L	1000	62966	117	97	LGR1000L5550, LGR1000L6600, LGR1000L8880, LGR1000L9990, LGR1000L8880, LGR1000L9990, LGR1000L9990, LGR1000L9990

Surface/Volume ratio
 Dilution factors



Scaling Examples

Example Calculation Showing Final Concentration of Port Extractables in Filled 1L and 20L Allegro 2D Biocontainers

Compound	Worst-case Amount ($\mu\text{g}/\text{cm}^2$)	Port Contact Area in Biocontainer		Worst-case Amount in Filled Biocontainers ($\mu\text{g}/\text{mL}$) ^a	
		1L	20L	1L	20L
Palmitic acid	3.18	24	44	0.0763	0.0070
1,3-Di-tert-butylbenzene	2.89	24	44	0.0694	0.0064
2-Propanol	1.51	24	44	0.0362	0.0033
Stearic acid	1.44	24	44	0.0346	0.0032
2-(Fornyloxy)-1-phenyl-ethanone	1.23	24	44	0.0295	0.0027
Ba	0.020	24	44	< 0.001	< 0.0001
Mo	0.018	24	44	< 0.001	< 0.0001

^a Worst-case Amount in Filled Biocontainer = [Worst-case Amount ($\mu\text{g}/\text{cm}^2$) x Port Contact Area] / Total Volume

Drug Product

$\mu\text{g}/\text{ml}$

BPOG VS USP <665>



	BPOG
Scope	User requirement specification for single use components
Goal	Intended to standardize and expedite component qualification by bracketing worst case extractables for ~80% of bioprocess applications
Risk-based	No supplier risk assessment. End user responsibility
Time points	3 to 4
Solvents	6
Low pH	0.1 M H ₃ PO ₄ (~pH 1.6)
High pH	0.5 N NaOH (pH 13.5)
Organic	50% EtOH
Aqueous	WFI
Salt	5 M NaCl
Surfactant	1% Polysorbate 80

BPOG VS USP <665>



	BPOG	USP <665> & <1665> (draft)
Scope	User requirement specification for single use components	Standard for extractables profiling of plastics used in pharmaceutical processing equipment
Goal	Intended to standardize and expedite component qualification by bracketing worst case extractables for ~80% of bioprocess applications	Intended to aid supplier selection through meaningful component profiling covering ~80% of applications
Risk-based	No supplier risk assessment. End user responsibility	Yes. Level of data based on risk matrix described in <1665>
Time points	3 to 4	1
Solvents	6	3
Low pH	0.1 M H ₃ PO ₄ (~pH 1.6)	KCl/HCl (pH 3)
High pH	0.5 N NaOH (pH 13.5)	Phosphate (pH 10)
Organic	50% EtOH	50% EtOH
Aqueous	WFI	-
Salt	5 M NaCl	-
Surfactant	1% Polysorbate 80	-

RISK ASSESSMENT USING E&L DATASET



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment




Product Map

- Delivery route
- Dose
- Dose regimen

Product Description

- Organic content
- Surfactants
- pH

Process Map




Process Description

Manufacturing	
Quantity	Yes
Batch	Regen
Material	Locked
Item	Active
Product Ref/Level	Basic
Manufacturing Location	Excipients
	Search
	Errors
Process description	
Please click the starting column and push the scrolling function	
Process Description Light source: <input type="checkbox"/> Temperature control: <input type="checkbox"/> Recipe number: <input type="checkbox"/> Controller: <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> Volume/Weight condition: <input type="checkbox"/> Product ID: <input type="checkbox"/> Product strength: <input type="checkbox"/>	
Product Description Name (in compliance with ICH Q10, and other local, European and US FDA): <input type="text"/> Organical No.: <input type="text"/>	
Title: <input type="text"/>	
Product Data Name: <input type="checkbox"/> Raw Name: <input type="checkbox"/> Name: <input type="checkbox"/> Is it a final product? <input type="checkbox"/> Type: <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
Product Description pH: <input type="text"/> Buffer: <input type="text"/> <input type="text"/> <input type="text"/> Product is a liquid: <input type="checkbox"/> Key Name: <input type="text"/> Excipients: <input type="text"/> Size: <input type="text"/> Excipients: <input type="text"/> Product: <input type="text"/> Excipients: <input type="text"/>	
Successful link refers to corresponding base file name	


RISK ASSESSMENT USING E&L DATASET

1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment

BPOG Component Extractables Data



Compendial Extractables Baseline Data

- USP <87>, <88>, <661>
- EP monographs

RISK ASSESSMENT USING E&L DATASET



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment

Component	Surface Area	Standard Solvents Data Available ($\mu\text{g}/\text{cm}^2$)							Action		Reported Compounds (μg)				
		pH 1.6	WFI	pH 13.5	50% EtOH	5% PS80	5M NaCl	Other	Risk Assess	Test	Compound A	Compound B	Compound C	Compound D	Compound E
PES Filter (KA3EKVP1G)	1530	●	●	●	●	●	●	○	●	○	0.1	0.1	0.1	12	1.4
Allegro Biocontainer	225	●	●	●	●	●	●	○	●	○	0.6	0.6	12	3.1	1.1
Tubing	122	○	●	○	●	○	○	●	●	○	0.1	1	0.4	0.5	0.6
Kleenpak Presto Connector	44	●	●	●	●	●	●	○	○	○	-	0.1	0.1	-	-
Fittings	34	○	●	○	●	○	○	●	○	●	0.2	-	-	0.1	-
Filling Needle	6	○	●	○	●	○	○	○	●	○	-	-	-	-	-
Total Extractable Levels per SUS											0.9	1.8	13	16	3.2

RISK ASSESSMENT USING E&L DATASET



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment

Component	Surface Area	Standard Solvents Data Available ($\mu\text{g}/\text{cm}^2$)							Action		Reported Compounds (μg)				
		pH 1.6	WFI	pH 13.5	50% EtOH	5% PS80	5M NaCl	Other	Risk Assess	Test	Compound A	Compound B	Compound C	Compound D	Compound E
PES Filter (KA3EKVP1G)	1530	●	●	●	●	●	●	○	●	○	0.1	0.1	0.1	12	1.4
Allegro Biocontainer	225	●	●	●	●	●	●	○	●	○	0.6	0.6	12	3.1	1.1
Tubing	122	○	○	○	○	○	○	○	○	○	0.1	1	0.4	0.5	0.6
Kleenpak Presto Connector	44	●	●	●	●	●	●	○	●	○	-	0.1	0.1	-	-
Fittings	34	○	○	○	○	○	○	○	○	○	0.2	-	-	0.1	-
Filling Needle	6	○	○	○	○	○	○	○	○	○	-	-	-	-	-
Total Extractable Levels per SUS											0.9	1.8	13	16	3.2

Application	Total Daily Exposure ($\mu\text{g}/\text{day}$)		AET				
	Route of Administration	Final API Concentration	0.7	1.4	10	13	2.6
	IV injection	1 mg/mL					
	Dosage Concentration	250 $\mu\text{g}/\text{mL}$					
	Maximum Dosage per day	2 mg					
	Duration of Exposure	Lifetime					
	Patient Population	All					

Threshold of Toxicological Concern (TTC) a level of exposure for all chemicals below which there would be no appreciable risk to human health

Analytical Evaluation Threshold (AET) threshold at or above which a leachable should be characterized and reported for toxicological assessment

RISK ASSESSMENT USING E&L DATASET



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment

Component	Surface Area	Standard Solvents Data Available ($\mu\text{g}/\text{cm}^2$)							Action		Reported Compounds (μg)				
		pH 1.6	WFI	pH 13.5	50% EtOH	5% PS80	5M NaCl	Other	Risk Assess	Test	Compound A	Compound B	Compound C	Compound D	Compound E
PES Filter (KA3EKVP1G)	1530	●	●	●	●	●	●	○	●	○	0.1	0.1	0.1	12	1.4
Allegro Biocontainer	225	●	●	●	●	●	●	○	●	○	0.6	0.6	12	3.1	1.1
Tubing	122	○	●	○	○	○	○	○	●	○	0.1	1	0.4	0.5	0.6
Kleenpak Presto Connector	44	●	●	●	●	●	●	○	●	○	-	0.1	0.1	-	-
Fittings	34	○	●	○	●	○	○	○	○	●	0.2	-	-	0.1	-
Filling Needle	6	○	●	○	●	○	○	○	●	○	-	-	-	-	-
Total Extractable Levels per SUS											0.9	1.8	13	16	3.2

Application		0.7	1.4	10	13	2.6
Total Daily Exposure ($\mu\text{g}/\text{day}$)						
Route of Administration	IV injection					
Final API Concentration	1 mg/mL					
Dosage Concentration	250 $\mu\text{g}/\text{mL}$					
Maximum Dosage per day	2 mg					
Duration of Exposure	Lifetime					
Patient Population	All					

Tox		II	I	I	I	I
Toxicology Cramer Class						
Genotoxicity (ICH M7)		N	N	N	N	N
Protein or Cell Interactions		N	N	N	N	N
ICH Q3C (Class)		III	-	-	-	-
Permitted Daily Exposure (PDE)		✓	✓	✓	✓	-

AET



RISK ASSESSMENT USING E&L DATASET



1. Product & Process map
2. Single-Use System
3. Component Quality Information
4. Aggregate and Scale Relevant Data
5. Application Information
6. Toxicological Risk Assessment

Component	Surface Area	Standard Solvents Data Available ($\mu\text{g}/\text{cm}^2$)							Action		Reported Compounds (μg)				
		pH 1.6	WFI	pH 13.5	50% EtOH	5% PS80	5M NaCl	Other	Risk Assess	Test	Compound A	Compound B	Compound C	Compound D	Compound E
PES Filter (KA3EKVP1G)	1530	●	●	●	●	●	●	○	●	○	0.1	0.1	0.1	12	1.4
Allegro Biocontainer	225	●	●	●	●	●	●	○	●	○	0.6	0.6	12	3.1	1.1
Tubing	122	○	○	○	○	○	○	○	○	○	0.1	1	0.4	0.5	0.6
Kleenpak Presto Connector	44	●	●	●	●	●	●	○	●	○	-	0.1	0.1	-	-
Fittings	34	○	○	○	○	○	○	○	○	○	0.2	-	-	0.1	-
Filling Needle	6	○	○	○	○	○	○	○	○	○	-	-	-	-	-
Total Extractable Levels per SUS											0.9	1.8	13	16	3.2


Application		0.7	1.4	10	13	2.6
Total Daily Exposure ($\mu\text{g}/\text{day}$)						
Route of Administration	IV injection					
Final API Concentration	1 mg/mL					
Dosage Concentration	250 $\mu\text{g}/\text{mL}$					
Maximum Dosage per day	2 mg					
Duration of Exposure	Lifetime					
Patient Population	All					

Tox		II	I	I	I	I
Toxicology Cramer Class						
Genotoxicity (ICH M7)		N	N	N	N	N
Protein or Cell Interactions		N	N	N	N	N
ICH Q3C (Class)		III	-	-	-	-
Permitted Daily Exposure (PDE)		✓	✓	✓	✓	-

Component Safety Risk Assessment

AET

SUPPORT LEVELS TO PHARMA END-USER



Standardized extractable data



Consulting services related to core data

Help to build Risk Assessment: identification of compounds to be evaluated, facilitation of Toxicological assessment, ...



Laboratory Services for process specific data
e.g. specific Extractables

Suppliers can provide different levels of support to E&L topics depending on customer's needs

- Standardized extractable data facilitate adoption of SUT
- They enable cGMP extractable and toxicological risk assessment
- Extractable data can also be used as input (alert signal) for product impact risk assessment
- It becomes very important to have high impact components well characterized by standardized extractables protocols



- Special thanks to my colleagues James Hathcock for helping setting up this presentation

