


Robotic Tub Decontamination System - RTDS

Christophe Riedel, Claranor SA



2020 PDA EUROPE

Parenteral Packaging

25-26 FEBRUARY 2020
BASEL | SWITZERLAND

CONNECTING
PEOPLE
SCIENCE AND
REGULATION®

1. Background of the Project
2. Pulsed Light Technology
3. Presentation of the RTDS
4. Validation Procedure
 - a) First validation trial
 - b) Why spray inoculation
 - c) Determination of reference micro-organism
 - d) Development of a BI
5. Features of the RTDS
6. Integration in a filling line

The application of Ready to Use (RTU) container for fill finish technologies in the pharmaceutical industry increased in the recent years.

- Easy to use for smaller batch sizes
- Easy to perform frequent product/container changes
- Increasing availability of RTU container

Steriline, manufacturer of aseptic fillers, was looking to a simple and safe way of TUB decontamination



Challenges of current decontamination processes

Robotic Tub Decontamination
System RTDS

- Suitable for high output 2 to 4 tubs/min
- No decontamination step after un-packing
- Possible contamination through pin-holes not covered

Unpacking devices
with no-touch
transfer

- Suitable for high output 2 to 3 tubs/min
- Batch process
- Residual VHP content in containers possible

Decontamination
with gases : VHP,
NO2

- Suitable for high output up to 6 tubs/min
- Cost intensive
- High weight

E-Beam
Technology

Project triggered and conducted by 2 companies

- **Steriline**, manufacturer of aseptic filling lines for the pharma industry
- **Claranor**, manufacturer of Pulsed Light Decontamination equipments in the food industry

2009 : first ideas about tub decontamination with pulsed light, attempt to start with a global pharma company

2015 : Steriline introduced the idea of manipulating the tub with a robot

2017 : internal validation at 6log on a prototype, and start of a joint project with pharma expert, potential customer, manufacturer of BI's, Steriline, Claranor


2018 : first order of RTDS by a European pharmacy company

2020 : January : FAT with 4log BIs, May: delivery


The pulsed light process consists in producing short duration emissions of intense, broad-spectrum white light from a Xenon lamp to sterilize surfaces.

- The sterilizing effect was discovered first in Japan 1970
- Developed in the US by Pure Pulse 1985-2000
- Approved by FDA for food decontamination in 1996
- Industrially developed:
 - At low power and with UVB for hair removal by many companies
 - At high power and with UVC for packaging decontamination in the food and beverage industry first by Claranor since 2004

2. How to produce Pulsed Light (“flashes”)




2. Pulsed Light : A technology with power



- Very short flash: 300 micro-second
- Intensity: 50,000 x solar intensity on earth
- Spectrum between 200 and 1100 nm

15% UV rich in UVC
50% visible
35% near IR

$$\frac{\text{Lamp energy}}{\text{Flash time}} = \frac{300 \text{ J}}{0,3 \text{ ms}} = 1 \text{ MW}$$




UV emission in the range of 1 kW / cm²

2. Pulsed Light Decontamination Efficiency

Microbiological decontamination of surfaces:

Up to 6 log reduction on the most resistant microorganisms (bacteria, mold, spore)



*Inactivation profiles of A. brasiliensis, B. subtilis, B. pumilus, B. atrophaeus depending on medium Fluence
Results obtained by C. Levy, PhD student, on Claranor 3 lamp pilot equipment, 2012*

2. Pulsed Light Decontamination of Packaging



Cap decontamination: 3-4 log on *Aspergillus brasiliensis* or *Bacillus atrophaeus*, currently up to 82,000 caps/hour

Cup decontamination: 3-5 log in < 0.5s

400 units operating worldwide

Hard operating conditions :


- High intensive industry (>6000h production per year)
- Humidity, vibrations...
- Low educated personnel



Cup decontamination on a low-fat margarine filling line (2013)
3log decontamination (*Aspergillus brasiliensis*) in 2 flashes at 60 strokes/min




3. Presentation of the RTDS



- a) First Validation Trial
- b) Determination of a reference micro-organism
- c) Why spray inoculation ?
- d) Production of a bio-indicator

Qualification Protocol: Count Reduction Test

- Experimental microorganism: **Bacillus pumilus DSM 492**
- **Spray inoculation** of Petri dishes
- Dry overnight under laminar flow
- Average count of positive samples: **$1,6 \cdot 10^8$**
- Treatment of the tubs with Petri dishes stuck on 5 sides through tunnel (transparent UV film protection on petri dishes)
- Recovery of surviving microorganism
- Calculation of the killing rate: **$\log (N / N_0)$**



4.a. First Validation Trial

Attempt A: 8 flashes / side

LOG REDUCTION						
Investigated face	1 st Exp.	2 nd Exp.	3 rd Exp	4 th Exp	5 th Exp	Mean Log-Red. (per side)
Short side #1	6.7	6.2	6.4	6.1	6.7	6.4
Short side #2	6.6	7.5	6.2	6.7	6.5	6.7
Large side #3	7.0	7.2	6.3	7.0	6.5	6.8
Bottom #4	7.2	7.0	6.2	7.4	7.0	7.0
Top # 5	6.8	6.9	6.2	6.5	6.2	6.5
Mean	6.9	7.0	6.3	6.7	6.6	

4.a. First Validation Trial

Further points and questions

- **Results of first validation trial**
 - Test A: 6.7 log with 8 flashes
 - Test B: 7.1 log with 14 flashes
 - No sample < 6log (out of 50 in total)
- **Further steps with perspective of industrial validation**
 - Confirmation of reference micro-organism, use of spray inoculation
 - Development of a BI with the French laboratory Icare for industrial use for end-point test at 6log.
 - FAT

Efficiency assessment

- **Spray inoculation on Petri dish**
 - Monolayer of spores
 - Homogeneous spore distribution
- **Drying overnight**
- **Pulsed Light treatment 1,1J/cm²**
 - Lamps at end of lifetime
 - 1 to X flashes
 - 5 samples for each modality
- **Recovery** of survivors (count reduction test or end-point test)





Deposition of *Bacillus subtilis* spores using an airbrush-spray or spots to study surface decontamination by pulsed light

Caroline Levy^{a,b,c}, Isabelle Bornard^d, Frédéric Carlin^{b,c,*}

^a Clammar S.A. Chemin de la Bollandaie, Agrignon, BP 21 531, F-84916 Avignon cedex 9, France

^b INRA, UMR108 Sécurité et qualité des produits d'origine végétale, F-84000 Avignon, France

^c Université d'Avignon et des Pays de Vaucluse, UMR 408 Sécurité et qualité des produits d'origine végétale, F-84000 Avignon, France

^d INRA, UR407, Unité de pathologie végétale, F-84143 Montfavet, France




Fig. 2. Scanning Electron Microscopy pictures of 10 μL spots of a *B. subtilis* suspension in sterile distilled water containing 10⁶ spores in a 10 μL spot at ×100 (A), and ×8000 (B) magnifications, and 10 μL spots containing 10⁴ spores in a 10 μL spot at ×1000 (C) and ×8000 (D) magnifications.

4.b. Why spray inoculation




Fig. 4. Inactivation of spores of *B. subtilis* DSM 402 inoculated on LB Agar media with (A) spray (●) and with a manual spreading with rakes (○), and (B) inoculated by spot at concentrations of 3×10^6 (●), 3×10^5 (□) and 3×10^4 (▲) spores per 10 μl spot on a polystyrene surface, as a function of fluence (in J cm⁻²). Bars represent standard error.

- Pulsed Light is efficient on spots up to 10^{E4}/10μl.
- Above 10^{E4}, reproducibility is reduced.
- No difficulty to spray 10^{E8} spores on a petri dish





Fig. 1. Dispersion of spores of *B. subtilis* sprayed on polystyrene Petri dishes. Colonies were formed after 24 h incubation at 30 °C from a spore suspension sprayed on the dish, dried for 24 h, and covered with molten LB agar. Volumes of 100 µL of the 10^{-3} (A) and 10^{-4} (B) dilutions of a spore suspension containing 10^8 *B. subtilis* spores mL^{-1} were sprayed. The pictures C and D show the Scanning Electron Microscopy of 100 µL of *B. subtilis* spores at a concentration of 10^8 spores mL^{-1} sprayed on polystyrene surface at $\times 500$ (C), and $\times 20,000$ (D) magnifications.

And what about real life contamination ?

Example from the dairy industry :

- Empty PP cup stored 3 days in our warehouse without protection
- 10 positive samples, with 160 to 880 cfu/cup, **in average 610 cfu/sample** (2,8log)
- Decontamination test with pulsed light and compared to a high output UVC lamp (100W/cm)

	Pulsed Light 1 flash	UVC 2sec	UVC 4sec
Log red	2,1	0,1	0,3



Pulsed Light is efficient against natural contamination

- Pulsed light is efficient against **natural contamination**
- Inoculation in **spots** is not a regulatory requirement
- Spot inoculation does not give reproducible results above 10^E4 to 10^E5 , because pulsed light is not efficient against **clusters** of micro-organisms for high concentrations.
- Spray inoculation enables to inoculate small surfaces with high quantities of spores, and to demonstrate in a reproducible way a local 6log reduction


Conclusion :

Spray inoculation can be considered as acceptable for BI's validating surface decontamination in the pharmacy industry.

Log10 reduction in spores after PL treatment at 3600 V
(similar to industrial configuration)

Strain	<i>Bacillus pumilus</i> ATCC 27142				<i>Geobacillus stearothermophilus</i> ATCC 7953				
	Flash	Exp. 1	Exp. 2	Exp. 3	Mean	Exp. 1	Exp. 2	Exp. 3	Mean
1		3,1	3,9	3,8	3,6	2,4	1,9	2,9	2,4
2		3,7	4,2	4,8	4,2	4,1	3,5	3,8	3,8
3		4,7	>4,7	5,3	5,0	5,1	3,5	4,5	4,4
4		6,1	4,7	>5,4	5,4	> 5,3	> 5,1	4,4	4,9
5		6,4	5,7	>6,0	6,1	6,0	> 5,9	> 5,3	5,7

The study was conducted with 4 microorganisms : *A. brasiliensis*, *B. atrophaeus*, *G. stearothermophilus*, *B. pumilus*
***Bacillus pumilus* ATCC 27142 and *Geobacillus stearothermophilus* ATCC 7953** emerged as the two more resistant



Geobacillus stearothermophilus ATCC 7953 was chosen because having a comparable sensitivity, and easier lab handling.

4.d. BI development


- Definition of shape, packaging, handling: BI in a UV-transparent film
- Stability of micro-organism resistance to pulsed light
- Reproducibility of spray inoculation
- Stability during storage and transport

- Already validated at 4log
- Ongoing validation at 6log




Handling of the BI's

- A dummy Tub was 3D printed with spaces available for BI's.
- This Tub should be used for any qualification



5. Features of the RTDS




- Lamp change every 100,000 tubs
- Energy consumption <3kW
- Traceability on flashes
- Light intensity measurement
- Access to reflectors and lamps from outside, enabling intervention on lamps without breaking sterility of the chamber
- Simple lamp change by user
- Power electronics in changeable racks
- Low maintenance


- Isolator
- AirTight System
- 6-log Reduction with VHP
- Short VHP Decontamination Cycle



- 3-Chamber-System:
 - Inlet Chamber
 - Decontamination Chamber
 - Outlet Chamber
- Uni-directional airflow through HEPA-Filter
- Airflow 0,45 m/s
- Pressure difference of 2 Pa between chambers



- Automatic De-Bagging System
- Pulsed-Light Decontamination
- Automatic De-Liner
- Filling Machine for 360 Syringes/min under Isolator



Questions ?



2020 PDA EUROPE

Parenteral Packaging

CONNECTING
PEOPLE
SCIENCE AND
REGULATION®