



# Laser Based Headspace Analysis for CCIT

*Josine Wilmer, Study Manager at LIGHTHOUSE*



# Overview

## Part 1: Laser based headspace analysis

- Technology overview
- Application to CCI testing of sterile products

## Part 2: Theoretical background

- How does gas ingress work?
- How can theory be applied?

## Part 3: Case studies and examples

- Process optimization
- Typical use

Laser light matches frequency of target molecule.

Amount of absorbed laser light is dependent on concentration of target molecule in headspace.




**Nondestructive**

## What gases can be measured?

- Headspace oxygen
- Headspace carbon dioxide
- Headspace moisture (water vapor)
- Headspace total pressure levels



# Application to CCI testing-Products packaged under a modified atmosphere



## Product

Stoppered under a modified atmosphere

Lyophilized product:

Partial vacuum with nitrogen or full vacuum

Oxygen sensitive liquids:

Purged with nitrogen

## Storage Condition

Ambient air

## Leak Indicating Measurement


Increase in oxygen level and/or increase in pressure.

## Application

Laboratory instruments or fully automated inspection

# Application to CCI testing-Products packaged under a non-modified atmosphere

Chamber that can be purged or evacuated



## Product

Stoppered under a non-modified atmosphere.  
Headspace is air at 1 atm  
Typical of many liquids

## Storage Conditions


Must put container in a chamber to drive a headspace change (similar to Blue Dye testing, but HSA is more sensitive and doesn't involve dye). After a pre-determined conditioning time then remove and test in a headspace analyzer.

## Leak Indicating Measurement

Decrease in oxygen level/increase in CO<sub>2</sub>/  
decrease in pressure level.

## Application

Laboratory instruments



## What type of product-packages?

- Clear or amber glass
- Transparent plastics
- Liquid, lyophilized or powder filled
- Vials, syringes, ampoules, cartridges
- Nominal volume ranging from 0.2mL to 250mL



Laboratory and At-line  
Instruments and accessories



Automated Inspection Machines




Strategic partnership with Bosch for CCI machines with Lighthouse laser measurement technology inside.

Instrument and machine qualification using NIST traceable standards.



N=100	Headspace Oxygen (% atm)				
	Standard Label	Known Value	Meas. Mean	Error	St. Dev.
	0.0	0.000	0.01	0.01	0.02
	1.0	1.005	0.96	-0.04	0.03
	2.0	2.004	1.98	-0.03	0.03
	4.0	3.998	4.02	0.02	0.04
	8.0	7.999	8.13	0.13	0.03
	20.0	20.00	19.93	-0.06	0.04

↑ Accuracy      ↑ Precision







## For every delivered system:

- Full validation package:
  - Functional Requirements (FR)
  - Design Specification (DS)
  - Traceability Matrix (TM)
  - Installation Qualification (IQ)
  - Operational Qualification (OQ)
- Users and data managed in a database solution for 21-CFR-11 compliance and full audit trail
- Certificates of NIST traceable calibration standards
- Optional yearly re-certification of standards





## Part 2

### Gas diffusion theory

In general, there are two different ways by which gas can flow through a defect in and out of a pharmaceutical container:

- **Effusion:** gas flow that is generated by a total pressure difference across the container defect
- **Diffusion:** gas flow of a particular gas that is generated by a partial pressure difference of that gas across the container defect

**Understanding this gas flow enables the development of CCI test methods based on the measurement of gas ingress**

$$\vec{J} = -D\vec{\nabla}n \quad \text{Fick's 1<sup>st</sup> Law}$$

$$\frac{\partial P_i(t)}{\partial t} = \frac{-D \cdot A_0}{V} \frac{\partial P_i(z, t)}{\partial z}$$

$$\%Oxygen(t) = 20.9\% \cdot \left[ 1 - \exp\left(-\frac{\alpha_{Diff}}{V}t\right) \right]$$

New USP <1207> states:


“Mathematical models appropriate to leak flow dynamics may be used to predict the time required for detecting leaks of various sizes or rates.”

**The change in oxygen concentration will be exponential with respect to time**

**Diffusion Parameter**

$$\alpha_{Diff} \left[ \frac{cm^3}{s} \right] = \frac{D \cdot A_0}{L}$$

**The Diffusion Parameter is a function of the Diffusion Coefficient,  $D$ , the defect cross-sectional Area,  $A_0$ , and Depth,  $L$ .**



With fixed values for:

$$D = 0.22 \text{ cm}^2/\text{s}$$

$$A_0 = 20 \mu\text{m}^2 \text{ (} 5 \mu\text{m } \varnothing \text{)}$$


$$V = 18 \text{cc (15R)}$$

Obtain an empirical  
depth parameter value:

$$L = 6 \mu\text{m}$$

Model matches the data  $\pm 0.3$  %-atm oxygen at every point

## Predicted oxygen concentration versus time for ideal defects





## Part 3

# Case studies and examples



## Overview Case studies

1. Filling line CCI qualification
2. Method development
3. Package development syringe CCI and permeation
4. Cold Storage CCI Study
5. 100% inspection of lyophilized product






# Case Study 1: Filling line CCI qualification

- **Objective: Generate data demonstrating that the filling process produces good CCI for a specific vial-stopper combination**
- **CCI Study:**
  - Produce empty stoppered & crimped vials with the process. Initial headspace is 1 atm of air (20.9% oxygen).
  - Use headspace gas ingress model to design a sample chamber evacuation, nitrogen backfill, and sample hold cycle.
  - Measure samples for changes in headspace composition. Include positive controls having a 10  $\mu\text{m}$  micro-capillary through the stopper.

	Headspace Oxygen Level after defined evacuation, backfill, and storage cycle [% atm]		
Sample Vial	1 $\mu\text{m}$ ideal defect	0.6 $\mu\text{m}$ ideal defect	0.5 $\mu\text{m}$ ideal defect
2R	1	4.1	7.1
6R	3.4	12.1	15.0
20R	11.1	18.8	19.5

# Case Study 1: Filling line CCI qualification

## Results headspace CCI test




**Example: Gas ingress measurements straightforwardly produced robust CCI data**

## Objective

Detection of 5 micron leak within 30 minutes

## Sample set

- 6R DIN clear tubing vial – 1.5mL product
- Positive controls: 2µm, 5µm, 10µm and 15µm laser drilled defects
  - Glass defects
  - Metal plate defects



**Nominal hole size 5 µm**




Image provided by Lenox Laser



## Case Study 2: Method development

### **Study 1: Manufacturing conditions**

- Determine purge quality

### **Study 2: API reactivity**


- Oxidation rate

### **Study 3: CCI method development**

- Diffusion tests with vials with known defects
- Effusion test with vials with known defects
- Method protocol

## Study 1: Manufacturing conditions


- 50 product, water-filled and empty samples



## Study 2: API reactivity

- 50 product samples opened to air and followed over time


**Mean measured headspace oxygen level monitored over time**



	Oxygen (% atm)
Start	19.59
2 hours	19.50
26 hours	19.18
36 hours	17.63
144 hours	15.76


## Study 3: CCI method development

- Diffusion tests with vials with known defects



## Study 3: CCI method development

- Effusion tests with vials with known defects





# Case Study 3: Package Development Syringe CCI and Permeation


## Sample set

- Two syringe types: Staked needle and Luer cone
- Positive controls: 2 $\mu$ m, 5 $\mu$ m, 10 $\mu$ m laser drilled glass defects

Step 1: Establishing the inherent leak rate of a package.


Step 2: Determine leak rate for positive controls

Step 3: Model MALL (0.2  $\mu$ m for sterility assurance)




# Case Study 4: CCI testing for vials stored on dry ice (CO<sub>2</sub>)

### Headspace oxygen



### Headspace CO<sub>2</sub>



## Case


- Air headspace vials stored on dry ice (CO<sub>2</sub>)
- Storage on dry ice increases risk of container closure integrity loss,
- Conventional rubber stoppers lose elasticity at -80°C risking CO<sub>2</sub> ingress

## Result

- 3 containers revealed decreased oxygen levels
- Same vials revealed increased CO<sub>2</sub> levels




# Case Study 5: 100% Inspection of lyo product



Total batch size: 29048  
Number rejected: 16  
Reject rate: 0.06%




# Case Study 5: 100% Inspection of lyo product



Total batch size: 29156  
Number rejected: 568  
Reject rate: 1.95%


# Case Study 5: 100% Inspection of lyo product



Results of 6  
chronological  
batches

**Not a robust  
process**

## Case Study 5: 100% Inspection of lyo product



### Case 100% inspection

4 years of manufacturing data:

- 156 lots
- Total 1.6 million vials

### Results

44-lots (28%) with zero rejects

3-lots had > 2% reject rate

Average reject rate was 0.27%

**It is difficult to manufacture a perfect batch**

## Pro's

- Non-destructive
- Rapid
- Quantitative results
- Deterministic method
- Operator independent
- Applicable over whole leak range
- Permanent & temporary leaks detectable

## Con's

- Not all fill levels
- Sample needs to be transparent to laser
- Inline production inspection needs modified headspace





Thank you!