

Ensuring the Product Robustness, the Science and the Technologies for **Proven Single-Use Container Closure Integrity**

Case Study: Helium Integrity Test



*Marc Hogreve
Senior Engineer Integrity Testing
Sartorius Stedim Biotech*



Single use process integrity is a key industry & regulatory challenge

Validated SU Container Closure Integrity (CCI) required to:

- Strengthen regulatory compliance
- Improve patient safety

SU container integrity testing (SU-CCIT) is meaningful when it correlates to liquid leaks & microbial ingress

- Lack of understanding of defect size causing liquid leaks and microbial ingress under real life conditions
- Current CCIT not correlated to leaks & microbial ingress

||| Bag failures cost
~\$100K to
\$1M per bag

R. Wong, Bayer

||| Leaks have been
reportedly
responsible for up to
\$20M worth of
products per year at
some larger
facilities.

Bioplan

Users cannot show real proof of process integrity to regulatory organizations

Why are concerns about Integrity a key industry challenge for SUS?



SUS Used in Drug Substance & Drug Product Process Steps
Require Higher Assurance of Container Closure Integrity (CCI)



“We’re very much concerned about the integrity of SUS used in final filling applications.”

FDA Inspector 02|2018

FDA-ASTM Workshop - Oct 2016 ⁽¹⁾ Examples of SU validation activities and Integrity testing

Validation for SU fill finish sterile product manufacturing

- **Physical integrity test correlated to microbial ingress**
- Packaging integrity at the supplier, post shipping & post-installation



Microbiological challenge testing

- Identification of defect size that would allow ingress of bacteria
under process conditions

⁽¹⁾ SUS: A Microbiology Product Quality Perspective - Patricia F. Hughes, Ph.D., Branch Chief (Acting), CDER/OPQ/OPF/DMA

FDA-ASTM Workshop - Apr 2018 ⁽¹⁾ SUS fill finish assemblies for sterile products

Must meet requirements to ensure flow-path sterility and integrity

- Supporting sterilization validation summary data and information on the gamma irradiation process
- **Integrity tests are also reviewed during the review of the BLA and on inspection**



An example of an information request sent to an applicant is as follows:

“You have submitted supplier’s leak testing data and method qualification information for the single use assembly. Implement a leak integrity test for the single use assembly prior to use for filling at the facility.”

⁽¹⁾ Sterile drug product, New developments - Patricia F. Hughes, Ph.D., Branch Chief (Acting), CDER/OPQ/OPF/DMA



Industry Challenges



- **USP<1207> 2016** *package integrity evaluation – sterile products*



- **ASTM E55 WK64337 2018** *Standard Practice for Integrity Assurance & Testing of Single-use Systems*
- **ASTM E55 WK64975 2018** *Test Method for Microbial Ingress Testing on Single-use Systems*



- **New PDA TR** *on Pharmaceutical package integrity (to replace existing TR27 from 1998)*



- **BPSA 2017** *Design, Control, and Monitoring of SUS for Integrity Assurance*

The challenge is to ensure that no product leakage or microbial ingress can occur during the process

- QbD & validation to reach consistent robustness
- Process control, quality control & integrity testing



SU Product
Robustness

- Understand liquid leakage & bacteria ingress mechanisms
- Correlate detection limit to liquid leaks & microbial ingress



SU Container
Integrity Science


- Implement 100% physical supplier integrity testing
- Provide point of use testing



Integrity Testing
Technology

Combined with the robustness of our self-deploying bags, our integrity tests ensure that no product leakage or microbial ingress can occur during the process

QbD, process control and quality control to reach consistent robustness



~ **3000 robustness tests**

Tensile Strength: 930


Elongation at break: 930

Puncture tests: 930

Flex durability: 186

Water burst tests : 26

Most stringent standards:
ASTM D4169



Passed
**ASTM Shipping
Test**

Self-deploying bag
after installation



Avoid mishandling
&
reduce bag failures

Process & Quality Controls Ensure CCI along our **Entire Production Cycle**


Process control reduces the risks of defects on film, seals, welds & connections

- Extrusion temperature & speed,
- Weld & seal temperature & time

“ SU-CCI is the Result of Process Qualification and Quality Controls Along the Entire Product Life Cycle

Quality controls confirm the absence of leaks


- Film & seal strength integrity,
- Microbial Ingress Testing – ISO15747 - Immersion in B. Diminuta
- 100% Bag Chamber Leak Test (BC-LT)



Intermediate process control, quality control & final product integrity testing



Bag chamber
leak test
40 – 90 μm



Discard 400ppm defects from
manufacturing operations
Integrity Assurance Level (IAL) 10^{-4}



Assembly
Helium SIT
2 μm



Discard 500ppm additional defects
from manufacturing operations
Integrity Assurance Level (IAL) 10^{-5}

But, why 2 μm ?



PoU-IT
10 μm
FlexAct® BT



Discard 20ppm additional defects
from transportation & handling
Integrity Assurance Level* (IAL) 10^{-6}

Existing data in peer-reviewed publications on **critical leak size for sterile package using micro tubes**

Liquid leak & bacterial ingress depend on


- Process & applications conditions
- Liquid properties (mainly surface tension)
- Size of defect

Correlation between liquid leak & bacterial Ingress

- Microbial ingress cannot occur without liquid flow (Kirsch & Morton, 1989)

Gibney, 2000
$$P_o > P_{atm} + \left[\left(\frac{4\sigma}{D_h} - \rho g L \right) \times 0.390 \right]$$

Treshold Leak Diameter for DI Water vs Microbial Ingress



2µm is the maximum allowable leakage limit that never produce liquid leak or microbial ingress in any process conditions using micro tubes

Understand liquid leakage and microbial ingress mechanisms on film materials used in single use biomanufacturing

- Demonstrate the relation between microbial ingress and liquid leak
- Establish a predictive model to determine the Maximum Allowable Leakage Limit (MALL*) under any process conditions
- Establish a correlation between liquid leak/microbial ingress and physical integrity testing
- Develop and validate the physical test methods with detection limits that guarantee the absence of liquid leak/microbial ingress in SUS

* MALL: is the greatest leak size tolerable that poses no risk to product safety, USP<1207>

Process Conditions Impact Liquid Leaks & Microbial Sizes


Hydrostatic pressure in storage



Flexboy[®] 20 L

Flexsafe[®] 500 L

Shipping can generate up to 20g acceleration



Application

Pressure range

Storage

~10 – 100 mbar

Shipping

~100 – 250 mbar

Interim conclusion from results obtained on **liquid leak & microbial ingress studies** on SUS


- Results obtained on film materials fall into the existing models for sterile package using micro tubes
- Both studies, liquid leak and microbial ingress ones, tend to confirm that leak sizes for liquid flow are not significantly different from leak sizes for sterility loss
- The conditions of use (pressure) does significantly impact the Maximum Allowable Leakage Limit (MALL)
- 2 μm is the MALL under any conditions and both, for liquid leak and microbial ingress

Main Requirements

- Test of **finished product** (bag chamber incl. lines & connectors)
- **100% non-destructive** testing (deterministic approach)
- **Sensitivity** that can be correlated to **bacterial ingress**
- Capable to be implemented **in-line in routine production**


Test Concept

- Porous spacers avoid masking effect
- Restraining plates reduce stress on bag
- Test can be used as “Integrity Test” covering the whole SUS
- Gas leak rate can be correlated to a leak size of 2 μ m




small inflation volume + high test pressure = best sensitivity

Feasibility Tests



Feasibility Tests

- Feasibility tests done with three different test gases
- Results based on 20 samples per leak size (50l bags)
- One decade differentiation in the time frame between 4-10 s
- Helium gas tracer was the best test method to achieve a level of detection **down to 2 μm**



Challenges for Industrialization

- High permeability of helium through the thin bag film
- Outgassing of volatile components
- Vacuum level vs. tact time
- Helium background concentration & absorption

Solutions for Industrialization

- Machine concept with fast helium filling to measure before the start of permeation
=> dynamic leak rate evaluation
- High vacuum to get lowest background level
- High vacuum pump capacity to reach vacuum in a reasonable time frame
- 2 chamber design for alternating operation to reach target tact time
- Flexible chamber design to cover different product families

Finished Product SIT to Test **100%** of 2D Bag Assemblies Used in Critical DS & DP Applications

- **Helium Gas Tracer in Vacuum Chamber**
- Restraining plates
 - Reduce stress on bag
 - Provide small inflation volume
 - Allow high test pressure:
300 mbar
- Porous spacer avoid masking effect of potential leaks
- With low background noise and dynamic leak rate measurement helium gas tracer is the best method to achieve reliably highest sensitivity
- Detection limit: **2 µm correlated** to microbial ingress and liquid leaks at 300mbar pressure (shipping pressure conditions)



Validation Approach

- A **reliable** test method using samples from multiple production lots
- A **reproducible, accurate** and **robust** test method using a 6 sigma confidence interval for the differentiation between defective and non-defective samples
- A **sensitive** test method using defective samples, calibrated on a pinhole equivalent
- A **One-for-all** validation procedure, whatever the test method is



Definition of Defect

- Aerosol BCT & liquid leak studies performed on different bag film materials
- Identification of MALL for microbial ingress & liquid leak under **different process conditions**
- 2µm estimated as MALL for worst-case (shipping) conditions -> target for the helium test validation
- Deliberate defects manufactured and confirmed with flow calibration

Definition of Defect


- Different types of defects tested: patches and capillaries
- Defects placed at worst-case positions



Parameter Study

- Evaluate impact of different product designs
 - Volume
 - Components
 - Material
- Determine optimum test parameter set(s)
- Tests on a limited number of test samples


Parameter Study Results on 50L Bags



Validation Study


- Verification of selected test parameter set(s)
- Determine pass/fail criteria (time & leak rate limit) => **dynamic leak rate**
- Tests on a statistically significant number of test samples

Validation Study Results on 50L Bags



Validation Study Result Overview

- Measurement of defective & non-defective samples
- Use of different types of artificial defect, like film patches & capillaries
- Defects deliberately calibrated
- Reliability proven using samples from multiple batches
- **2µm detection validated applying a 6 sigma confidence interval**



Validation study with >380 tests to establish test specifications

Proven Integrity of your process improves patient safety, regulatory compliance, production costs and drug availability

- Enhance patient and operator safety
- Prevents any risk of high value product loss at commercial phase
- Meets cGMP regulatory expectations for SU container closure integrity
- Speeds up drug product manufacturing capacity and market availability with SU proven process integrity





- Questions?