



# **CCIT - Container Closure Integrity Testing following USP 1207**

**Tony Stauffer/ Tatiana Simental CCIT s.a./Pti**

**St Prex, Switzerland**

**2019-02**

# CCIT - New USP 1207 guidance document

Documents specifies and differentiates acceptable test methods:

- **Deterministic Leak Test Methods**
  - **Vacuum (ASTM F2338-09)**
  - **Electrical Conductivity and Capacitance tests (HVLD)**
  - **Head Space analysis**
- **Seal Quality test methods**
  - **Airborne Ultrasound (ASTM F3004-13)**
- **Probabilistic leak test methods (not considered good practice)**
  - **Water bath**
  - **Blue Dye Ingress**
  - **Microbial ingress test**
  - **Burst test**




- Measures and verifies container closure integrity
- Detects vapor or gas release
- Test sensitivity down to 0.01 cc/min. (1 – 1.5 micron)
- Case studies prove more reliable than dye ingress
- ASTM F2338-09



# VeriPac Differential Vacuum Test Method

(simplified scPti hematics)


ASTM F2338-09



*US Patents 6,513,366 and Pend.*



## Leak Rate Distribution










# Electrical Capacitance and High Voltage Detection also known as HVLD

**Pre-filled Syringes**  
**Vials & Ampoules**  
**BFS**  
**Bottles and Cups**  
**Pouches & IV Bags**



# E-Scan based on HVLD<sup>mc</sup> for liquid filled parenteral products



**Off-line laboratory system to  
Inspect liquid filled  
Vials  
Syringes  
Ampules**

**Lower HV application**


**mc: micro current**

**Improved SNR**


**Good for low conductivity liquids  
Incl distilled water**

**Product not exposed to HV**

**Negligible Ozone production**




## Good Sample



2 capacitors

## Leak



1 capacitor

$V$  – High Voltage Source


$R$  – Electric Resistance of the product

$C_1$  – Capacitor 1: Glass between the inspection electrode and product


$C_2$  – Capacitor 2: Glass between the detection electrode and product

$I_1$  – current produced when product container is sealed

$I_2$  – current produced when product container is defective



- Signal Detection 20ml Vial – Integrator – Signal has amplitude changes
- High Voltage:16kV
- **Red - 3.6um leakage**
- **Black sample without defect**
- **350 RPM, 5mm/s, filled with distilled water**




# Test with E-Scan<sup>®</sup>


## (positive controls)

Sample id.	Nominal leak	Positive	
		1 ml syringe	
		Water	Alb.
		Volt	Volt
44	5	6.75	
41	10	8.38	
38	20	8.19	
39	?	10.00	
40	?	7.76	
42	?	8.94	
43	?	9.89	
45	?	8.52	
46	?	8.60	
23	5		5.79
35	10		8.48
7	20		9.06
8	?		8.29
9	?		9.82
10	?		9.45
11	?		5.65
31	?		9.89
32	?		7.64
33	?		7.83
34	?		6.28
19	?		7.47
20	?		10.00
21	?		7.17
22	?		6.29

Sample id.	Nominal leak	Positive	
		2.25 ml syringe	
		Water	Alb.
		Volt	Volt
M2	5	8.19	
M3	10	9.46	
M	20	10.00	
N2	?	9.22	
O2	?	10.00	
N3	?	9.22	
O3	?	7.95	
N	?	10.00	
O	?	10.00	
G2	5		5.38
G3	10		6.08
G	20		9.71
H2	?		7.09
I2	?		4.70
J2	?		5.45
K2	?		3.38
H3	?		7.06
I3	?		4.87
J3	?		5.83
K3	?		6.35
H	?		7.77
I	?		9.03
J	?		8.74
K	?		8.60


# Voltage for Good and Defective samples









# High Conductivity Liquid 1'000uS





# High Voltage - Technology Comparison

( distilled Water w/ 1ml Syringe )



	Conventional HVLD	E-Scan 655 ( HVLD <sup>mc</sup> )
Product Exposure	7,000V with 18.5kV input	300V with 18.5kV input
Ozone production	0.305ppm ozone in 1 minute	0.004ppm ozone in 10 minutes

# Ozone Creation

**HVLD - 0.305ppm Ozone  
after 1 min**



**HVLD<sup>mc</sup> – 0.001ppm Ozone  
after 6 minutes**



# Feasibility Study – PFS Samples tested

- 1 ml and 2.25 ml syringes
- Prefilled with H<sub>2</sub>O and Albumin (175 mg/ml, Cp  $\approx$  2.5)



# Summary of results (negative controls)

Negative controls		# samples	VeriPac VP-455 ( vacuum decay)		E-Scan 655 ( HVLD)	
			found negative	found negative	found negative	found negative
1 ml	Water	15	15	100 %	15	100 %
	Albumin	24	24	100 %	24	100 %
2.25 ml	Water	15	15	100 %	15	100 %
	Albumin	24	24	100 %	24	100 %

- All negative samples are identified as such with both VeriPac® and E-Scan® instruments
- No false positives

# Summary of results (positive controls)

Positive Controls			# samples	VeriPac VP-455 (vacuum decay)		E-Scan 655 (HVLD)	
				Found positive	Found positive	Found positive	Found positive
5 $\mu$ m	1 ml	Water	3	0	0 %	3	100%
		Albumin	5	0	0 %	5	100%
	2.25 ml	Water	3	0	0 %	3	100 %
		Albumin	5	0	0 %	4	80 %
10 $\mu$ m	1 ml	Water	3	3	100 %	3	100%
		Albumin	5	0	0 %	5	100%
	2.25 ml	Water	3	0	0 %	3	100 %
		Albumin	5	0	0 %	5	100 %
20 $\mu$ m	1 ml	Water	3	3	100 %	3	100%
		Albumin	5	0	0 %	5	100%
	2.25 ml	Water	3	3	100 %	3	100 %
		Albumin	5	0	0 %	5	100 %

- No albumin prefilled positive sample could be detected with Vacuum Decay
- E-Scan® allows to identify all positive samples except one

**Non-Destructive, Non-Invasive, No Sample preparation**

**Follows USP 1207 Guidelines**

**ASTM test method > ISO/FDA recognized**

**Repeatable as well as Reproducible**

**Deterministic Quantitative (Informative)**

**Accurate and Reliable**

**Simple and Robust**

**Cost Effective**

**Zero Waste**



# Case study 3 (Extra)


Scratch



Heat



Water Cool






- 16kV
- Sensitivity 50%
- Speed 15mm/sec
- 250RPM
- Exposure voltage ~300V
- Negligible Ozone creation



**0-10V  
Range**



Controlled crack produced.  
 Certified flow measurement with Helium mass spec.  
 Peak signal response is recorded.  
 Voltage signal response with each rotation.



**Sample 1**



He Certified Leak Rate	
95.0	10 <sup>-5</sup> mbar·l/s
0.05700	sccm
2.67	~μm



Sample	leak rate			Visual	Size um
	mbar·l/s	10 <sup>-5</sup> mbar·l/sec	sccm		
1	0.00095	95	0.05700	Large crack	2.67
2	0.0000021	0.21	0.00013	Medium crack	0.13
3	0.000014	1.4	0.00084	Medium crack	0.32
4	0	0	0.00000	Small Scratch	0.00
6	0.00067	67	0.04020	Large crack	2.24
7	0.015	1500	0.90000	Large crack	10.61
9	0.00000029	0.029	0.00002	Small crack	0.05
10	0.00029	29.0	0.01740	Large crack	1.47
11	0.074	7400	4.44000	Large crack	23.56
12	0.055	5500	3.30000	Large crack	20.31
16	0.0014	140	0.08400	Small crack	3.24
18	0	0	0.00000	Small Scratch	0.00
19	0	0	0.00000	Small Scratch	0.00
20	0.016	1600	0.96000	Large crack	10.95

Sample 2



Sample 6



Sample 3



Sample 16




**Sample 6**

**He Certified Leak Rate**

67.0	$10^{-5}$ mbar·l/s
0.04020	sccm
2.24	~ $\mu$ m



**10.0V**




**Sample 3**


**He Certified Leak Rate**

29.0	$10^{-5}$ mbar·l/s
0.01740	sccm
1.47	~ $\mu$ m



**10.0V**






**Case Study 3**

**Sample 2**

He Certified Leak Rate	
0.21	$10^{-5}$ mbar·l/s
0.00013	sccm
0.13	~ $\mu$ m




**10.0V**




**Sample 3**

He Certified Leak Rate	
1.40	$10^{-5}$ mbar·l/s
0.00084	sccm
0.32	~ $\mu$ m



**7.6V**



CCIT methodology requires consideration of multiple factors.

Wide range of positive control methods; use them responsibly.

HVLD is not time critical (time zero vs. day 66).

MicroCurrent HVLD generally effective for wide range of product conductivities.

Naturally occurring defects below the 1 $\mu$ m level can be readily detected using HVLD.

If a technology can accurately and reliably test for defects below the MALL, the debate over the MALL becomes negligible.



# Global Quality Solutions



**Thank You!**

**Tony Stauffer/Tatiana Ponce-Simental**  
**contact@ccit.com**

**USA: +1 914 337 2005**

**Switzerland: +41 21 805 0020**

**[www.ptiusa.com](http://www.ptiusa.com)**

**[www.ccit.com](http://www.ccit.com)**



## **CCIT Laboratory – Pti-Europe**

### **Address:**

**CCIT s.a.  
Chemin du Glapin 4  
1162 St-Prex  
Switzerland**



# Acknowledgements to staff members of Pti and CCIT