

Case Study: Test Method Development and Validation for 2 mL Glass Vials Containing 0.70 mL of Drug Product using HVLD

Turan Gül, Ph.D; Eurofins Proxy Laboratories; turangul@eurofins.com

With significant contribution from Jennifer Roark, Eurofins Medical Device Testing

Venice, March 21-22, 2019



- Method Development Strategy
- Aim of the study
- Methods and Materials
- Experimental Setup
- Results
- Conclusion

➤ Container and Product Details

- Container Type: Vial, Prefilled Syringe, Ampoules
- Container Material: Glass, Plastic, other etc.
- Closure Type: Rubber stopper, Aluminum Cap
- Drug product detail: Formulation, Liquid small or large molecule (Protein),
- Size of the container,
- The liquid fill volume of the container
- Materials
 - Negative controls (non-defective containers)
 - Positive controls: Laser drilled holes at 2, 5, 10, 15 and 20 μm



Instrumental (HVLD) Details

- Equipment Qualification (IQ, OQ and PQ)
- Container specific equipment:
 - sample holder
 - resistor block
- Instrumental parameters:
 - voltage
 - gain sensitivity
 - speed, rotation
 - motion profile

- **Validation Parameters <USP 1207.1>**
 - Specificity
 - LOD
 - Accuracy and Range
 - Precision (Repeatability, Intermediate Precision)

To develop and validate a test method for an optimal CCIT for 2 mL vials containing approximately 1.0 mL fill of drug product using HVLD





2 mL vial filled with 0.70 mL DP



PTI E-scan 655 HVLD

- Drug Product (DP)/ Formulation buffer
- 2 mL empty vials (45 units)
- 2 mL vial cap (45 units)
- 2 mL vial stopper (45 units)
- Negative controls: 2 mL Vials (30 units)
- Positive Controls: 2 mL Vials with laser drilled holes at 2, 5 and 10 μm (10 units each)


- Negative Control: 2 mL Vials (30 units)
- Positive Control: 2 mL Vials with holes at 2, 5 and 10 μm (10 units each)
- In total 90 Positive and 90 Negative Control measurements
- 3 experimental replicates, 3 different days, 2 analysts



Defect Location at the Shoulder of the Vial

Parameter	Description
High voltage	13 kV
Sensitivity	55%
Rotation	300 RPM
Speed	8 mm/s

Motion Profile				
Step	X	Y	V	W
1	43.5	50.0	32.0	50.0
2	43.5	5.1	32.0	8.4
3	43.5	5.1	31.5	8.4
4	43.5	5.1	32.0	8.4
5	41.8	5.1	32.0	8.4
6	41.8	8.0	32.0	8.4
7	32.3	8.0	32.0	8.4
8	21.2	8.0	16.2	8.4
9	12.7	8.0	24.8	8.4
10	12.7	5.3	24.8	8.4
11	12.7	5.3	25.3	8.4
12	12.7	5.3	24.8	8.4






LOD: the smallest breach in the vials allowing a “consistent differentiation” from voltage observed in non-defective vials

- LOD was established as 2 μm
- Confidence 95%

Specificity: the ability to differentiate between leaking and non-leaking vials, despite interfering factors that may cause false detection.



Example of HVLD-scanning profile result

n=3, 90 measurements*	
Parameters	Negative controls
Average	2.66 VDC
StdDev	0.26
Mean + StdDev	2.92 VDC
Mean + 6xStdDev	4.19 VDC
Mean + 10xStdDev	5.22 VDC
Confidence 95%	0.05 VDC
Confidence Low	2.61 VDC
Confidence High	2.71 VDC
Minimum	2.36 VDC
Maximum	3.97 VDC

Rejection Threshold: 4.19 VDC

LOD Criterion: Smallest leakage rate (or leak size) that a test method can reliably detect.

*All negative controls passed

Precision: the ability of the method to yield reliable and repeatable data (n=3, 2 analysts, in total 90 measurements for positive and negative vials)


Accuracy: the ability to differentiate samples with defects greater than or equal to the established LOD from samples with defects less than the established LOD, using voltage measurements.

Range: the range over which the method has shown to have acceptable precision accuracy from the established LOD up to the greatest defect size tested

Acceptance criteria:

- The voltage results for every positive control and gross leak vial tested is \geq Rejection threshold (4.19 VDC) (for Accuracy and Precision)
- The voltage result for every negative control vial tested is less than 4.19 VDC (for Accuracy and Precision)

Summary of Development Data



n=3, total 90 measurements*			
Parameters	2 μ Defect (n=30)	5 μ Defect (n=30)	10 μ Defect (n=30)
Average	8.22 VDC	8.49 VDC	7.24 VDC
StdDev	1.38	1.27	1.20
Mean - StdDev	6.84 VDC	7.22 VDC	6.04 VDC
Mean - 2xStdDev	5.46 VDC	5.94 VDC	4.84 VDC
Confidence 95%	0.49 VDC	0.46 VDC	0.43 VDC
Confidence Low	7.73 VDC	8.03 VDC	6.81 VDC
Confidence High	8.71 VDC	8.94 VDC	7.67 VDC
Minimum	5.63 VDC	5.00 VDC	4.96 VDC
Maximum	10.00 VDC	10.00 VDC	10.00 VDC

*All tested positive controls passed

- The method was specific to differentiate non-defective (negative controls) vials from defective (positive) ones
- The method LOD was determined to be 2 μm
- The rejection threshold was determined to be 4.19 VDC
- The data showed acceptable accuracy and precision
- For gross leak testing, stopper defects (fiber leaks) are not recommended, instead cracks on the body of the vials are recommended to be used as gross leaks to evaluate the range of the method.

<1207> Package integrity evaluation –Sterile products

<1207.1> Package integrity in the product life cycle –Test method selection and validation

<1207.2> Package integrity leak test technologies

- Jennifer Roark, B.S., Manager Chemistry & Container Testing, Eurofins Lancaster Laboratories
- Lei Li, Ph.D., Associate Engineer Advisor Delivery and Device R&D, Eli Lilly
- Bert Al, Ph.D., Senior Scientist, Eurofins Proxy Laboratories

THANK YOU!