



# Elastomeric Closures for Lyophilization Applications

**Francis Carroll**  
**Specialist TCS, EU**

---

Date 27<sup>th</sup> Oct 2022

PDA Freeze Drying in Practice, Osterode, Germany





- Considerations for Elastomeric Stoppers for Lyophilization
- Challenges & Solutions: Stopper pop-up
- Challenges & Solutions: Stoppers sticking to the Lyo Shelf
- Challenges & Solutions: Moisture Uptake and Transmission
- Vial System Considerations
- New Stopper for Lyophilization Applications

# Considerations for Lyophilization Closures: Elastomer Formulation



- Low absorption and adsorption characteristics
- Low extractable volatiles
- Low oxygen transmission
- Low moisture absorption
- Low moisture permeation
- Easy to dry
- Low particulates
  
- Low coring/fragmentation
- Good resealing
- Optimum container closure integrity
- Good machineability



# Considerations for Lyophilization Closures: Design



## Legged Design

- More flexible during insertion
- Symmetric: keeps horizontal position during freeze drying
- Twining effect of stoppers possible



## Igloo Design

- More stable in the vented position due to more contact with the glass vial and a less flexible plug part
- Asymmetric balance point: can get out of vertical axis during stoppering
- Igloo design prevents twinning

# Challenges & Solutions: Stopper Pop-Up/Pop-Off



## Check list:

- ✓ Blowback features
- ✓ Interference fit of the stopper plug and vial neck
- ✓ Siliconization of the stopper and vial

# Challenges & Solutions: Stopper Pop-Up/Pop-Off

## Interference Fit of stopper plug and vial neck



**Interference Fit** : Analysis which shows the dimensional interference between stopper plug and vial inner neck




✓	Stopper plug outer diameter [mm]	13.10
✓	Vial neck inner diameter* [mm]	12.60
✓	Nominal Interference Fit [%]	<b>3.8</b>

**Industry Standard is 2 - 10%**

Consider min and max values!

Stopper plug outer diameter

Vial neck inner diameter




$$Interference\ fit\ [\%] = 100 \cdot \frac{[(stopper\ plug\ outer\ diameter) - (vial\ inner\ neck\ diameter)]}{(stopper\ plug\ outer\ diameter)}$$

\*\* According to ISO 8362-1:2009 Injection containers and accessories — Part 1: Injection vials made of glass tubing  
 Holistic Considerations in Optimizing a Sterile Product Package to Ensure Container Closure Integrity. Fran L. DeGrazio. PDA Journal of Pharmaceutical Science and Technology Jan 2018, 72 (1) 15-34; DOI: 10.5731/pdajpst.2017.007658


# Challenges & Solutions: Stopper Pop-Up/Pop-Off Blowback features



**Non-blowback stopper**



**Blowback stopper**



protruding nubs or rings  
for stopper positioning



Blowback feature


# Challenges & Solutions: Stoppers Sticking to the Lyo Shelf



The closure surface can stick to the lyo shelf and lift the whole vial when closing the vials after lyophilization

## Solution: FluroTec<sup>®</sup> film to prevent sticking

- Eliminate sticking to the lyo chamber shelf
- Reduce clumping/sticking and twinning issues
- Minimize adsorption and absorption
- Minimize extractable volatiles from the elastomer
- Increase lubricity and thus ensuring a smoother transport in the filling line, leading to higher throughput
- Base elastomer formulation remains the same - regulatory burden for change is smaller





# Challenges & Solutions: Moisture Uptake and Transmission



- Moisture Vapor Transmission Rate (MVTR) can be measured
- For most drugs, the target residual moisture is < 1%  
lyophilized product
- Contributions to moisture include
  - Washing procedure
  - Sterilization cycle
  - Drying parameters dependent
  - Migration through the stopper




# Challenges & Solutions: Moisture Uptake and Transmission



→ Container closure integrity (CCI):  
Moisture ingress via the stopper-vial interface  
when lacking intact seal

→ Moisture Vapor Transmission Rate (MVTR):  
Water vapor permeates from the  
environment through the stopper during  
long-term storage

→ Residual moisture in the closure:  
Desorbs or releases from the stopper after  
processing (washing and steam sterilization)




# Challenges & Solutions: Moisture Uptake and Transmission



Reference: ISO 15106-3:2005, ASTM F1927

- **Critical Parameters:** temperature, humidity, thickness, flow rate, gas concentration
- Test item: Web
- Equilibrium needed before analysis
- MVTR is a material constant but is dependent on stopper design
- Typical data range:  $0.06 - 6 \text{ g/m}^2 \cdot \text{day}$




# Challenges & Solutions: Moisture Uptake and Transmission



## Case Study: MVTR and moisture uptake by the lyophilization cake

- Autoclaving: 60 min, 121°C  
Drying: at 105°C for 8 hours
- Filling: 2 mL of 5% lactose in water solution
- Storage of lyophilized product vials for 3 years at 25°C/60% RH for 36 months




**Moisture uptake by the lyo cake correlates with the MVTR of the rubber formulation**

# Vial System Considerations: Container Closure Integrity



- For an optimal CCI - use interference fit and stack up assessment
- The land seal is the main sealing surface
- An integral container closure system supports a low MVTR
- Methods to measure CCI, e.g., in USP <1207> CCI Testing



# Vial System Considerations: Stack-Up Analysis of Vial, Stopper and Crimp Seal




**Stack-Up Analysis:** Shows the stack-up of vial, stopper and seal



✓	Seal skirt length [mm]	7.5
✓	Stopper flange thickness [mm]	3.43
✓	Vial crown height * [mm]	3.60
✓	Nominal excess skirt length ** [mm]	<b>1.04</b>

**Industry Standard is 0.76 – 1.3 mm**




$$\text{Excess skirt length} = (\text{seal skirt length}) - (\text{aluminum thickness}) - (\text{vial crown height}) - [(\text{stopper flange thickness}) \cdot (1 - \% \text{ stopper compression})]$$

\* According to ISO 8362-1:2009 Injection containers and accessories — Part 1: Injection vials made of glass tubing  
 \*\* at 20% compression, aluminum thickness 0.2 mm

# Vial System Considerations: Crimping Process



- Capper parameter variables for sealing rail capper
- Head height
- Pressure block
- Vial rest position
- Pre-compression force
- Sealing rail vertical position
- Sealing rail lateral position
- Sealing rail angles and angle contour
- Compression zone
- Applied top spring force




# Vial System Considerations: Residual Seal Force



Residual Seal Force (RSF) is the stress a compressed elastomeric closure flange continues to exert on a vial land seal after crimping

- The stopper acts like a “compressed spring”
- Quantifying the RSF is a test method for the indirect estimation of elastomeric closure compression
- Sufficient compression is essential to seal integrity

Helium Leak Rate versus Residual Seal Force (20 mm Stoppers)



RSF can be correlated to CCI, but needs to be evaluated for each specific container/closure/seal combination



# Vial System Considerations: Volatile Leachables



Reconstituted solutions might show haze formation

- Volatiles from the rubber composition can adsorb to the solid cake
- Antioxidants, oils, waxes, oligomers, low molecular weight PDMS, ...
- This is more likely with legacy formulations than with modern formulations
- FluroTec<sup>®</sup> barrier film can help prevent leachables issues
- A “low extracting” rubber formulation might help as well ...



Blue area: FluroTec<sup>®</sup> film coverage  
Blue is for illustration purposes only

# New Stopper for Lyophilization Application: 4040/40



## Formulation:

Chlorobutyl 4040/40

## Design:

V-50-I (13mm), S-87-I (20mm), Igloo LyoTec®

Stopper

FluroTec® on top

B2 Coating on plug part

Westar® washed

Ready-to-Sterilize and Ready-to-Use

Steam and Gamma compatible


Envision™ verification




# Quality by Design Approach: 4040/40 Development



## Historical Approach



## 4040/40 Portfolio Approach




**West's deep elastomer expertise & Quality by Design approach to development of the 4040/40 portfolio has focused on mitigating risk from raw material selection of the formulation to the final product**


# 4040/40 Leachables Over Time with Comparators



## Volatile Leachables Lyo Cake Gas Chromatography–Mass Spectrometry (GC-MS) Headspace Reporting Threshold $\geq 0.1 \mu\text{g/g}$




# 4040/40 Optimized Properties: Low Moisture Content



*Sterilization of Stoppers*

- › Ready-to-sterilize product
- › Steam sterilized @ 121°C for 1 hr
- › Oven-dried @105°C
- › Improve steam sterilization throughput with optimized drying properties
- › Total Residual Moisture  $\leq 0.2\%$  H<sub>2</sub>O for Lyophilization Stoppers



# 4040/40 Optimized Properties: Compatible Across ISO Glass Vial Blowback Types



TESTING PROGRAM



ISO standard tested	13mm (V-50-I)	20mm (S-87-I)
European blowback (EBB)	Compatible	Compatible
No blowback (NBB)	Compatible	Compatible
American blowback (ABB)	Compatible	Compatible



**Simplify global supply chain management with stopper designs compatible with various glass blowback geometries**



# Thank you!

This presentation is for informational purposes only. West's products and solutions are sold on the basis that it is the customer's responsibility to evaluate and test the West product or solution to determine its compatibility with other materials and fitness for any end use. WEST MAKES NO WARRANTIES, WHETHER EXPRESS, IMPLIED OR STATUTORY, INCLUDING, WITHOUT LIMITATION, THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE, RELATING TO THE INFORMATION IN THIS DOCUMENT/PRESENTATION.

West & the diamond logo, Westar, STERILIZABLEBAG, Envision, FluroTec<sup>®</sup>, B2 and LyoTec are trademarks and registered trademarks of West Pharmaceutical Services, Inc. in the U.S. & other jurisdictions. FluroTec technology is licensed from Daikyo Seiko, Ltd.

© 2022 by West Pharmaceutical Services, Inc., Exton, PA. All rights reserved.

This material is protected by copyright. No part of it may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying or otherwise, without written permission of West Pharmaceutical Services, Inc.