

Helium Leak Testing: Pharmaceutical CCIT Methods and Applications

Presenter:

Jeff Morrow-Lucas


Product Manager – Helium Solutions
Packaging Technologies & Inspection



Helium Leak: What is it and How does it work?

- Gas molecules flow into an ionization chamber/analyzer cell.
- Filament electron beam ionizes molecules.
- Ions travel according to an acceleration voltage, and a magnetic field causes deflection relative to mass/charge.
- The analyzer cell is specifically tuned for the mass/charge of helium.
- Stream of helium ions hitting the target is reported as leak rate.
- Mass Spectroscopy is a four-step process: ionize, accelerate, deflection, & detect.
- The helium flow rate must be normalized to 100% helium.
- Leak rate is expressed as a flow rate, i.e., mbar-L/sec

He Leak Testing: Gas Flow Animation



Package Formats & Applications Most Suitable for He Leak Testing

- Glass Vials/Syringes/Cartridges
- Foil Blistercards & Foil Pouches/Sachets
- Pharma Foil Tubes – Centrifuge Tube Containers
- Plastic Packaging – Ophthalmic Bottles, Vials, Syringes
- IV Bags – Sniffer Mode & Hard Vacuum Methods
- Medical Devices – Pacemakers, Injectors
- Industrial Applications – Weld Joints, Tube Connections, Valves, Pumps

He Leak Testing: The 3 Principle Methods

1. *Hard Vacuum Method-Vacuum Chamber*

- Samples filled/flushed with helium, placed in vacuum chamber, chamber evacuated, and global helium leak rate obtained.
- Examples: Vials, Blistercards, pouches
- Need Helium concentration to normalize Helium Leak rate obtained from the vacuum chamber test.

2. *100% Helium Flow*

- Area of leak interest is placed under vacuum using custom designed fixtures.
- 100% helium gas flow is maintained on atmosphere pressure side of test fixture. Any increase over background leak rate is a breach or leak pathway
- Ex: Any Plastic container: vials, bottles, syringes, cartridges. Glass containers with pre-drilled hole for helium gas flow.
- No need for leak rate correction since 100% helium concentration is used during test.
- Fast, simple method. Most cycles < 30 seconds

3. *Sniffer Mode Method*

- Site specific probe used to scan for helium leakage in packages/objects containing helium gas.
- Ex: IV bag seals & port fittings, large vessels & tubing
- Not quantitative – pass/fail only – helium in atmosphere limits detection

He Leak Testing: Hard Vacuum Apparatus



Rectangular Chambers
Blistercards/Pouches



Cylindrical Chambers
Vials

He Leak Testing: 100% He Flow Apparatus



Plastic Vial Test Fixture



Plastic Vial Fixture
(test ready)



Foil Tube - Closure
Leak Test Fixture

He Leak Testing: 100% He Flow Apparatus



Induction Seal
Vacuum Fixture



Syringe Barrel
Needle Fuse Leak Test Fixture

He Leak Testing: 100% He Flow Apparatus



Glass Jar – Closure Fixture



Helium Flow Meter

- Device allows a measured flow of helium into a container for package for continuous 100% helium flow method testing.
- Flowmeter calibrated specific to Helium gas.

He Leak Testing: Sniffer Mode



Helium Sniffer Probe



Foil Bag Sniff Leak Test

Cold Temp Helium Leak Testing



- Enables low temperature vial & syringe leak testing at temperatures from 0°C to ~ -150 °C
- Provides real-time sample temperature monitoring via TC connection; vacuum chamber attaches to inlet port of Helium Detector

He Leak Testing: Validation Hardware



External Helium Leak Standards

- NIST traceable Helium Leak standards are available in a wide range of helium leak rates
- Validates the detector response of the SIMS Helium Leak Testing systems up to the test port – point of sample testing.
- Direct connection to Helium Leak Detector inlet port.



For additional information, please contact me:

j.morrowlucas@pti-ccit.com

Or visit our site:

<https://www.pti-ccit.com/technologies/helium-leak-detection>

