

Introduction to current sterilization methods

Hydrogen peroxide decontamination

Lucia Maita, Technical Customer Support, Fedegari Group

Agenda

Introduction

Decontamination, Gas sterilization & Vapour phase sterilization

Hydrogen peroxide

Chemical properties & Sporicidal Mechanism

Decontamination technologies

Fedegari Isolation technology

Development & Validation of a decontamination process

Normative references & practical case studies

Decontamination

The removal of microorganisms by **disinfection** or **sterilization**.



Disinfectant—A chemical or physical agent that destroys or removes vegetative forms of harmful microorganisms when applied to a surface.

Sterilant—An agent that destroys all forms of microbial life including fungi, viruses, and all forms of bacteria and their spores. Sterilants are liquid or vapor-phase agents.

USP43 <1072> Disinfectant and Antiseptics


Decontamination

...“decontaminated” refers to an **item** or **surface** that has been subjected to a **process that eliminates viable bioburden**.

USP43 (1208) Sterility Testing – Validation of Isolator System

Decontamination

The overall process of **removal or reduction of any contaminants** (chemical, waste, residue or microorganisms) from an area, object, or person. The method of decontamination used (**e.g. cleaning, disinfection, sterilisation**) should be chosen and validated to achieve a level of cleanliness appropriate to the intended use of the item decontaminated. See also Bio-decontamination.



Bio - Decontamination: A process that eliminates viable bioburden via use of sporicidal chemical agents.

EudraLex Vol. 4 - Annex 1 – Glossary

Gas sterilization (vapour phase sterilization)

Gas sterilization of surfaces may be used for the sterilization of *primary packaging materials, equipment* and some *pharmaceuticals*.

European Pharmacopoeia 10.0 – 5.1.1 Methods of preparation of sterile products- Gas sterilisation (vapour phase sterilization)



Sterilizing agent: physical or chemical entity or combination of entities, having sufficient microsporidical activity to achieve sterility under specified conditions.

ISO 11139:2018, 3.288

Sporicidal Vapor Phase Decontamination

The destruction or inactivation of microbial spores using a **vapor or gaseous agent**.

PDA TR. No. 51, Glossary of Terms

Sterilization

A process used to render a product, environment or pieces of equipment free from viable microorganisms with a specific probability.

PDA TR. No. 51, Glossary of Terms

Technical Report No. 51
Biological Indicators for
Gas and Vapor-Phase
Decontamination Processes:
Specification, Manufacture,
Control and Use




2010


Vapor Phase Sterilization

Sterilization can be accomplished using *sporicidal agents suspended in air* (i.e., vapor). Sterilizing agents that operate in this fashion include *hydrogen peroxide* (H_2O_2), *peracetic acid* (CH_3CO_3CH), *formaldehyde* (CH_2O), and *glutaraldehyde* [$CH_2(CH_2CHO)_2$] in aqueous solution.


USP43 (1229 .11) Vapor phase sterilization




Hydrogen Peroxide
 H_2O_2



Chlorine dioxide
 ClO_2



Nitrogen dioxide
 NO_2



Ozone
 O_3

General Classification of Antiseptics, Disinfectants, and Sporicidal Agents

| Chemical Entity | Classification | Example |
|----------------------------------|--|---|
| Aldehydes | Sporicidal agent | 2% Glutaraldehyde |
| Alcohols | General purpose disinfectant, antiseptic, antiviral agent | 70% Isopropyl alcohol, 70% alcohol |
| Chlorine and sodium hypochlorite | Sporicidal agent | 0.5% Sodium hypochlorite |
| Phenolics | General purpose disinfectant | 500 µg per g Chlorocresol, 500 µg per g chloroxylenol |
| Ozone | Sporicidal agent | 8% Gas by weight |
| Hydrogen peroxide | Vapor phase sterilant, liquid sporicidal agent, antiseptic | 4 µg per g H ₂ O ₂ vapor, 10%–25% solution, 3% solution |
| Substituted diguanides | Antiseptic agent | 0.5% Chlorhexidine gluconate |
| Peracetic acid | Liquid sterilant, vapor phase sterilant | 0.2% Peracetic acid, 1 µg per g peracetic acid |
| Ethylene oxide | Vapor-phase sterilant | 600 µg per g Ethylene oxide |
| Quaternary ammonium compounds | General purpose disinfectant, antiseptic | Concentration dependent on application, Benzalkonium chloride |
| β-Propiolactone | Sporicidal agent | 100 µg per g β-Propiolactone |

USP43 (1072) Disinfectants and antiseptics

The Resistance of Some Clinically Important Microorganisms to Chemical Disinfectants

(Listed in Order of Decreasing Resistance)

| Type of Microorganisms | Examples |
|--|---|
| Bacterial spores | <i>Bacillus subtilis</i> and <i>Clostridium sporogenes</i> |
| Mycobacteria | <i>Mycobacterium tuberculosis</i> |
| Nonlipid-coated viruses | Poliovirus and rhinovirus |
| Fungal spores and vegetative molds and yeast | <i>Trichophyton</i> , <i>Cryptococcus</i> , and <i>Candida</i> spp. |
| Vegetative bacteria | <i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , and <i>Salmonella</i> spp. |
| Lipid-coated viruses | Herpes simplex virus, hepatitis B virus, and human immunodeficiency virus |

USP43 (1072) Disinfectant and Antiseptics

Sporicidal process

A gaseous, vapour or liquid treatment applied to **surfaces**, using an agent that is recognised as capable of killing bacterial and fungal spores.


..... The process is applied to *internal surfaces of the isolator* and **external surfaces of materials inside the isolator**, when conventional sterilization methods are not required.



The application of a sporicidal process to isolators is not considered to be a sterilization process in the same way as, for example, a sealed container subjected to a validated dry heat, moist heat or irradiation process.

PIC/S - Recommendation used for aseptic processing and sterility testing, Glossary

Hydrogen Peroxide – H₂O₂



PubChem



Oxidizer

Corrosive

Irritant

Hydrogen Peroxide is a peroxide and **oxidizing agent** with **disinfectant**, **antiviral** and **anti-bacterial** activities. Upon rinsing and gargling or topical application, hydrogen peroxide exerts its **oxidizing activity** and produces *free radicals which leads to oxidative damage to proteins and membrane lipids*. This may inactivate and destroy pathogens and may prevent spreading of infection.

H₂O₂ Decontamination: when?

Heat-sensitive materials (including electronic devices) that should be transferred between classified areas (class C,D → class A, B) in order to minimize the risk of contamination

H₂O₂ Decontamination: where?

Fedegari Isolator technology




FCDV

*Fedegari Cabinet
Decontamination Vapor*



FCTS

*Fedegari Sterility Test
Isolator*



Customizable solutions


*Flexible Filling Line
Isolator*

An effective H₂O₂ decontamination: how?

The biocide should reach **all the surfaces** that should be decontaminated at a **specific concentration** for a pre-determined **amount of time**.

This means that the biocide must have:

- Good and complete **distribution**
- Sufficient contact **time** at a specified **concentration**




Items exposed to the process should have their surfaces exposed to the greatest extent possible


... There are several effective approaches to hydrogen peroxide (H_2O_2) injection, including **continuous**, **intermittent**, or **all at once**. Some of the systems utilize an **evacuation or drying step** prior to introduction of the hydrogen peroxide (H_2O_2) to allow for increased concentration without excess condensation. Alternatively, hydrogen peroxide (H_2O_2) can be introduced as a **liquid, followed by target heating**. Following the exposure period, the chamber or target is **aerated** to an acceptable level for further processing of materials and/or personnel exposure (whichever is lowest) prior to opening and removing the sterilized article.

USP43 (1229.11) – Vapor Phase sterilization


Typical process conditions



700 –1300 ppm H₂O₂ vapor during decontamination



≤ 5 – 90 % Humidity




≤ 25°C Temperature during the process



Total cycle time: based on the product type

Decontamination Cycle



An effective H₂O₂ decontamination: how?



A combination of **physical** and **biological methods** is used to determine the **optimum decontamination conditions**.



Biocide concentration
Process temperature
Relative humidity
Load configuration
Exposure time




Chemical Indicator
Biological Indicator

H₂O₂ distribution

... **chemical indicators** can provide a limited indication of sterilant distribution.


USP 43 (1229.11) Vapor phase sterilization



Chemical Indicator: Test system that reveals change in one or more pre-defined process variables based on a chemical or physical change resulting from exposure to a process

UNI EN ISO 17665, Terms and Definition 3.6

H₂O₂ distribution: Empty chamber mapping



Exposed



Not Exposed



H₂O₂ distribution: BIs & H₂O₂ Process Validation

Biological indicators are required to **demonstrate the effectiveness of the process** they play a crucial role during the cycle development and validation.

USP 43 (1229.5) and PIC/S Recommendation clause 9.4.13



Biological Indicator: It is a well-characterized preparation of a specific microorganism that has known resistance to a specific sterilization process.



USP 43 (1229.5) Biological Indicators For Sterilization

The use of **multiple BIs** at each test location is recommended to more adequately support the process lethality.

USP 43 (1229.11) Vapor phase sterilization

H₂O₂ distribution: BIs & H₂O₂ Process Validation

The BIs should be evenly distributed throughout the chamber and the load



Key


X locations for BIs, PCDs or temperature sensors

NOTE The diagram shows examples for locations in typical chamber usable space. Different chamber sizes can require more or fewer locations however a similar distribution pattern can be used.

ISO 22442 – Annex K (informative) Recommended validation test procedures


BI placement & handling

Proper BI Placement During VHP Decontamination Cycles - Kurt McCauley Spore News MesaLabs Volume 9, No. 5



BIs & CIs Layout

- T1 Left corner of the sealing front side
- T2 Left corner of the sealing backside
- T3 Left corner of the bottom front side
- T4 Left corner of the bottom backside
- T5 Middle of the left wall
- T6 Right corner of the sealing front side
- T7 Right corner of the sealing backside
- T8 Right corner of the bottom front side
- T9 Right corner of the bottom backside
- T10 Middle of the right wall
- T11 Between the bottle and Petri package - top shelf
- T12 Between the canisters - top shelf
- T13 Middle of the middle shelf
- T14 Middle of the bottom shelf
- T15 Between the bottle of the bottom shelf, the left side
- T16 Middle of the bench
- T17 Middle of the ceiling
- T18 Middle of the front wall (transfer unit door)
- T19 Middle of the back wall




BIs & CIs Mapping

- T1 Left corner of the sealing front side
- T2 Left corner of the sealing backside
- T3 Left corner of the bottom front side
- T4 Left corner of the bottom backside
- T5 Middle of the left wall
- T6 Right corner of the sealing front side
- T7 Right corner of the sealing backside
- T8 Right corner of the bottom front side
- T9 Right corner of the bottom backside
- T10 Middle of the right wall
- T11 Between the bottle and Petri package - top shelf
- T12 Between the canisters - top shelf
- T13 Middle of the middle shelf
- T14 Middle of the bottom shelf
- T15 Between the bottle of the bottom shelf, left side
- T16 Middle of the bench
- T17 Middle of the ceiling
- T18 Middle of the front wall (transfer unit door)
- T19 Middle of the back wall



BIs & CIs Mapping

3 BIs and 1 CI at each chosen position





PRACTICAL CASE STUDY


Simulation of a

Passbox reduced load

qualification



1. **Chamber loading** according to the layout
2. **CIs placement**
3. **BIs placement**



KEEP CALM
and wait for ...
...1 ppm!

Thank you for your attention

