

# Headspace Moisture and Water Activity Applications for Lyophilized Product

*Rapid non-destructive moisture determination using headspace analysis*



Connecting People, Science and Regulation®




Connecting People, Science and Regulation®

*Dr. Derek Duncan  
Director Product Line  
Lighthouse Instruments  
dduncan@lighthouseinstruments.com*

CONNECTING  
PEOPLE  
SCIENCE  
AND  
REGULATION®

- Thirty years ago the gravimetric or loss-on-drying method was the standard moisture determination method for freeze-dried pharmaceutical product.
  - Karl-Fischer (KF) titration used more frequently as a new method starting in the early 90's.
  - Total moisture content measured by both methods was different – so people started with correlation studies.

[Ref.] *Measurement of Final Container Residual Moisture in Freeze-Dried Biological Products*, International Symposium on Biological Product Freeze-Drying and Formulation, Bethesda, USA, 1990; *Developm. Biol. Standard.*, Vol. 74, pp. 153-164 (Karger, Basel, 1991)



- Disadvantages of the standard moisture determination methods:
  - Risk of moisture contamination is high - methods dependent on operator skill
  - Methods are time & resource intensive
  - Methods are destructive meaning high value samples are lost



Very difficult (impossible) to analyze statistical sample sets of finished product:


- **Insight into the drying process is limited**
- **Accuracy of moisture stability specifications also limited** due to destructive analysis and assumptions about identical stability samples

- Headspace moisture analysis - Introduction
  - Comparison to standard moisture methods
  - Relation to water activity
  
- Industry case study examples
  - Formulation development
  - Process optimization and scale up
  - Manufacturing & QC
  
- Q & A

## Where is the water?

*Headspace measures free water. Other methods do not distinguish between different forms of water.*

**Total cake moisture**  
*(Karl-Fischer titration)*




**Water activity ( $a_w$ ) is the partial vapor pressure of water at a fixed temperature**

**Headspace moisture**  
*(laser-based headspace analysis)*

**Partial cake moisture**  
*(loss on drying, ...NIR)*

# Headspace Moisture – correlation with total moisture content




## Water activity definition – thermodynamic property

- **Water activity**, or  $a_w$  is the **partial vapor pressure of water** in a substance divided by the standard state partial vapor pressure of water.
  - Indicates how much ‘free’ water is available for chemical reaction
  - A water activity of 0.80 means the vapor pressure is 80% of that of pure water

Laser light matches absorption frequency of water vapor molecule.

Amount of absorbed laser light is dependent on partial water vapor pressure in headspace.



**Non-destructive** (product not lost, sample can be re-measured)

**Rapid** (~1 sec)

- Measurement of the water vapor pressure in a sealed container at a defined temperature is a water activity determination.
- Measurement concept can be applied to general containers containing dry material.
- **There will be a trend to replace some total moisture measurements with a water activity determination – USP <922>.**






The same moisture method can be scaled up from manual benchtop platforms, to semi-automated platforms, to automated production scale inspection machines.


# Example: Robust efficient lyo cycle optimization

- Headspace water vapor correlates to product dryness
- 100% headspace moisture determination to characterize and optimize freeze drying cycle

Headspace moisture as a function of tray position Initial lyo cycle



Headspace moisture as a function of tray position Optimized lyo cycle



**Total analysis time (both batches): ~ 1 day!**



- Model: 5% sucrose solution, 6R vials, 2.5 ml fill
- Full headspace moisture analysis of shelf:
  - Water vapor range of 0.65 to 2.04 mbar correlated to 0.8 to 1.4 % (KF)
  - Headspace results clearly show an inhomogeneous batch
- Edge effects front and side result in lower values at front and side positions (acrylic door, covered with stainless steel plate)

FD02      5% sucrose

Sh2

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0,958	0,811	0,812	0,848	0,852	0,859	0,844	0,876	0,849	0,868	0,831	0,897	0,995	1,032	1,106	
0,843	0,866	0,764	0,911	0,807	0,825	0,832	0,853	0,91	0,96	0,831	0,88	0,918	0,868	1,01	
0,731	0,895	0,865	0,802	1,126	0,837	0,823	1,062	1,061	0,914	1,221	0,83	0,981	0,967	0,945	
0,783	0,825	0,818	1,065	0,808	0,868	1,09	0,93	0,901	1,446	0,837	0,952	0,843	0,929	0,898	
0,719	0,847	1,06	1,132	0,821	1,116	0,954	1,141	0,976	1,137	1,027	0,888	1,087	0,835	0,807	
0,986	0,791	0,845	0,774	1,6	1,43	0,941	1,236	0,966	0,888	1,05	0,899	0,997	0,837	0,771	
0,718	0,766	0,839	0,794	0,817	1,135	0,921	0,767	1,504	0,892	2,035	1,109	0,921	0,77	0,689	
0,71	0,72	0,917	0,845	1,461	0,988	1,238	1,323	1,073	1,373	1,107	0,949	0,75	0,801	0,683	
0,704	0,922	0,946	0,724	1,647	1,015	1,308	1,466	1,213	0,787	0,881	0,797	0,855	0,687	0,752	
0,681	0,709	0,747	0,84	0,705	1,097	1,165	0,982	1,868	1,265	1,213	1,099	1,083	0,742	0,677	
0,76	0,879	0,734	1,399	1,297	0,808	1,351	1,31	1,007	0,93	1,093	1,275	1,014	0,762	0,743	
0,706	0,724	0,823	0,79	0,767	1,087	1,431	1,277	1,081	0,789	1,44	0,965	0,76	0,717	0,657	
0,69	0,725	1,015	1,073	1,471	0,735	1,331	1,126	0,999	1,016	1,155	0,918	0,662	0,785	0,721	
0,649	0,724	0,719	0,755	1	1,226	1,686	0,723	0,991	0,702	1,087	1,017	0,877	0,809	0,859	
0,677	1,165	1,158	0,835	0,832	1,394	0,846	1,432	1,314	0,903	0,774	0,987	0,792	0,747	0,684	
0,657	0,808	0,769	0,878	0,806	1,108	1,27	1,313	0,829	0,714	1,261	0,832	1	0,709	0,703	
0,692	0,902	0,756	0,921	1,156	0,846	0,933	0,928	0,766	0,749	0,835	1,265	0,748	0,685	0,726	
0,654	0,699	0,721	0,886	0,884	0,723	0,832	0,785	1,173	0,775	0,837	0,745	0,745	0,734	0,764	
0,695	0,672	0,783	0,754	0,799	0,876	0,765	0,81	0,986	0,728	0,69	0,771	0,678	0,843	0,813	
0,656	0,711	0,665	0,737	0,853	0,729	0,772	1,037	0,827	0,733	0,764	0,781	0,725	0,758	0,769	
0,707	0,695	0,659	0,698	0,657	0,754	0,746	0,813	0,797	0,67	0,736	0,794	0,791	0,787	0,69	
0,746	0,713	0,844	0,758	0,773	0,7	0,699	0,793	0,814	0,733	0,718	0,728	0,802	0,693	0,755	

- Model: 5% sucrose solution
- Full headspace moisture analysis of shelf:
  - 0.18 to 3.79 mbar, i.e. 0.6% to 2.2% KF
  - More inhomogeneous batch
- Cake appearance acceptable, minor shrinkage in some vials
- Edge effects front, back and side result in lower values (stainless steel door)

FD04 5% sucrose

SH2

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0.379	0.358	0.434	0.389	0.361	0.378	0.366	0.357	0.184	0.465	0.37	0.366	0.38	0.388	0.399	
0.362	0.342	0.481	0.361	0.724	0.401	0.483	0.433	0.369	0.647	0.369	0.433	0.483	0.469	0.339	
0.321	0.514	0.411	1.765	0.691	0.556	0.942	1.684	1.448	0.867	1.779	0.76	1.103	0.398	0.369	
0.414	1.166	2.279	0.701	0.936	1.52	1.245	1.912	0.644	0.599	0.944	1.194	0.734	0.627	0.351	
0.314	0.434	1.781	2.45	1.222	0.978	1.355	0.609	2.214	1.23	2.338	1.149	1.935	0.746	0.448	
0.435	0.851	1.363	0.657	1.437	0.727	1.364	2.411	1.397	2.607	2.705	1.178	1.808	0.552	0.32	
0.39	0.439	0.515	2.332	1.569	2.863	3.314	0.796	0.815	0.671	2.009	2.124	0.806	0.435	0.356	
0.419	1.044	1.344	0.868	2.009	2.887	2.36	1.791	2.045	2.154	1.073	0.881	0.738	0.548	0.343	
0.324	0.51	2.236	1.741	0.861	1.579	3.786	2.423	0.558	2.298	0.857	3.308	1.881	0.669	0.408	
0.408	0.421	1.537	2.209	2.348	0.604	3.19	1.626	1.831	2.212	0.795	2.845	1.278	0.493	0.325	
0.305	1.504	1.257	1.083	1.698	1.799	2.104	1.63	1.624	2.199	2.284	0.896	2.069	1.529	0.395	
0.488	0.804	0.504	3.424	0.602	3.228	1.459	2.468	2.233	1.008	1.876	0.701	2.491	0.747	0.483	
0.336	0.694	0.935	1.43	1.388	2.623	0.913	1.828	2.009	0.793	1.93	1.001	2.921	1.918	0.443	
0.4	0.515	2.229	1.095	1.452	2.353	2.917	2.318	2.467	2.45	0.804	2.731	1.434	0.53	0.357	
0.437	0.4	0.481	1.381	1.302	2.294	2.495	3.326	3.278	0.658	2.368	1.241	0.785	0.504	0.39	
0.366	0.906	3.334	3.688	1.201	2.144	1.051	1.202	2.837	1.399	2.307	2.965	1.256	0.384	0.363	
0.301	0.369	2.257	1.301	1.004	2.114	1.919	2.571	2.187	2.045	0.63	2.032	0.707	0.751	0.355	
0.335	0.45	0.519	2.507	1.388	1.109	1.045	1.202	0.615	2.912	1.841	2.991	1.494	0.455	0.359	
0.305	0.359	0.782	1.541	1.199	0.578	1.669	0.49	1.303	2.892	0.815	2.355	1.55	0.416	0.424	
0.34	0.378	0.675	0.483	0.584	0.657	0.667	0.705	0.97	0.792	1.29	0.398	0.591	0.417	0.431	
0.338	0.399	0.452	0.34	0.388	0.395	0.36	0.353	0.516	0.473	0.38	0.394	0.378	0.374	0.486	
0.342	0.393	0.356	0.379	0.376	0.351	0.38	0.367	0.361	0.392	0.363	0.401	0.352	0.364	0.358	

## 5% sucrose, lyophilized with **non-optimized** cycle on two different freeze-dryers

### Freeze-dryer Nr 1

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0.9580.8110.8120.8480.8520.8590.8440.8760.8490.8680.8310.8970.9951.0321.106															
0.8430.8660.7640.9110.8070.8250.8320.8530.910.980.8310.880.9180.8681.01															
0.7310.8950.8650.8021.1260.8370.8231.0621.0610.9141.2210.830.9810.9670.945															
0.7830.8250.8181.0650.8080.8681.090.930.9011.4460.8370.9520.8430.9290.898															
0.7190.8471.061.1320.8211.1160.9541.1410.9761.1371.0270.8881.0870.8350.807															
0.9860.7910.8450.7741.61.430.9411.2360.9660.8881.050.8990.9970.8370.771															
0.7180.7660.8390.7940.8171.1350.9210.7671.5040.8922.0351.1090.9210.770.689															
0.710.720.9170.8451.4610.9881.2381.3231.0731.3731.1070.9490.750.8010.683															
0.7040.9220.9460.7241.6471.0151.3081.4661.2130.7870.8810.7970.8550.6870.752															
0.6810.7090.7470.840.7051.0971.1650.9821.8681.2651.2131.0991.0830.7420.677															
0.760.8790.7341.3991.2970.8081.3511.311.0070.931.0931.2751.0140.7620.743															
0.7060.7240.8230.790.7671.0871.4311.2771.0810.7891.440.9650.760.7170.657															
0.690.7251.0151.0731.4710.7351.3311.1260.9991.0161.1550.9180.6620.7850.721															
0.6490.7240.7190.7551.2261.6860.7230.9910.7021.0871.0170.8770.8090.859															
0.6771.1651.1580.8350.8321.3940.8461.4321.3140.9030.7740.9870.7920.7470.684															
0.6570.8080.7690.8780.8061.1081.271.3130.8290.7141.2610.8321.07090.703															
0.6920.9020.7560.9211.1560.8460.9330.9280.7660.7490.8351.2650.7480.6850.726															
0.6540.6990.7210.8860.8840.7230.8320.7851.1730.7750.8370.7450.7450.7340.764															
0.6950.6720.7830.7540.7990.8760.7650.810.9860.7280.690.7710.6780.8430.813															
0.6560.7110.6650.7370.8530.7290.7721.0370.8270.7330.7640.7810.7250.7580.769															
0.7070.6950.6590.6980.6570.7540.7460.8130.7970.670.7360.7940.7910.7870.689															
0.7460.7130.8440.7580.7730.70.6990.7930.8140.7330.7180.7280.8020.6930.755															

### Freeze-dryer Nr 2

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0.3790.3580.4340.3890.3610.3780.3660.3570.1840.4650.370.3660.380.3880.399															
0.3620.3420.4810.3610.7240.4010.4830.4330.3690.6470.3690.4330.4830.4690.339															
0.3210.5140.4111.7650.6910.5560.9421.6841.4480.8671.7790.761.1030.3980.369															
0.4141.1662.2790.7010.9361.521.2451.9120.6440.5990.9441.1940.7340.6270.351															
0.3140.4341.7812.451.2220.9781.3550.6092.2141.232.3381.1491.9350.7460.448															
0.4350.8511.3630.6571.4370.7271.3642.4111.3972.6072.7051.1781.8080.5520.32															
0.390.4390.5152.3321.5692.8633.3140.7960.8150.6712.0092.1240.8060.4350.356															
0.4191.0441.3440.8682.0092.8872.361.7912.0452.1541.0730.8810.7380.5480.343															
0.3240.512.2361.7410.8611.5793.7862.4230.5582.2980.8573.3091.8810.6690.408															
0.4080.4211.5372.2092.3480.6043.191.6261.8312.2120.7952.8451.2780.4930.325															
0.3051.5041.2571.0831.6981.7992.1041.631.6242.1992.2840.8962.0691.5290.395															
0.4880.8040.5043.4240.6023.2281.4592.4682.2331.0081.8760.7012.4910.7470.483															
0.3360.6940.9351.431.3882.6230.9131.8282.0090.7931.931.0012.9211.9180.443															
0.40.5152.2291.0951.4522.3532.9172.3182.4672.450.8042.7311.4340.530.357															
0.4370.40.4811.3811.3022.2942.4953.3263.2780.6582.3681.2410.7850.5040.39															
0.3660.9063.3343.6881.2012.1441.0511.2022.8371.3992.3072.9651.2560.3840.363															
0.3010.3692.2571.3011.0042.1141.9192.5712.1872.0450.632.0320.7070.7510.355															
0.3350.450.5192.5071.3881.1091.0451.2020.6152.9121.8412.9911.4940.4550.359															
0.3050.3590.7821.5411.1990.5781.6690.491.3032.8920.8152.3551.550.4160.424															
0.340.3780.6750.4830.5840.6570.6670.7050.970.7921.290.3980.5910.4170.431															
0.3380.3990.4520.340.3880.3950.360.3530.5160.4730.380.3940.3780.3740.486															
0.3420.3930.3560.3790.3760.3510.380.3670.3610.3920.3630.4010.3520.3640.358															

5% sucrose, **2.5% BSA**, lyophilized with an **optimized cycle** on the two different freeze-dryers

## Freeze-dryer Nr 1

FD02 5% sucr 25 mg/ml BSA

Sh3

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0.32	0.32	0.34	0.35	0.34	0.37	0.35	0.33	0.33	0.40	0.44	0.38	0.39	0.59	0.54	
0.36	0.33	0.33	0.34	0.39	0.35	0.35	0.41	0.39	0.46	0.35	0.39	0.42	0.40	0.31	
0.35	0.40	0.35	0.49	0.37	0.35	0.44	0.37	0.51	0.45	0.36	0.35	0.41	0.43	0.43	
0.41	0.33	0.33	0.43	0.36	0.55	0.40	0.48	0.48	0.41	0.44	0.40	0.40	0.35	0.35	
0.31	0.46	0.35	0.39	0.39	0.43	0.39	0.46	0.53	0.47	0.54	0.40	0.38	0.38	0.34	
0.46	0.34	0.34	0.36	0.35	0.52	0.41	0.39	0.51	0.40	0.47	0.38	0.37	0.37	0.29	
0.35	0.39	0.35	0.41	0.40	0.45	0.37	0.43	0.42	0.50	0.41	0.42	0.38	0.37	0.36	
0.41	0.37	0.38	0.53	0.35	0.38	0.45	0.45	0.57	0.41	0.38	0.37	0.37	0.35	0.33	
0.35	0.38	0.35	0.40	0.51	0.46	0.38	0.37	0.43	0.54	0.38	0.39	0.47	0.31	0.28	
0.31	0.34	0.42	0.35	0.40	0.60	0.38	0.37	0.42	0.40	0.38	0.44	0.33	0.30	0.31	
0.30	0.32	0.33	0.34	0.37	0.39	0.65	0.61	0.42	0.66	0.47	0.45	0.34	0.31	0.27	
0.31	0.46	0.35	0.34	0.38	0.40	0.47	0.37	0.45	0.59	0.42	0.32	0.32	0.26	0.29	
0.30	0.32	0.43	0.51	0.55	0.54	0.35	0.44	0.38	0.39	0.36	0.33	0.42	0.29	0.29	
0.39	0.45	0.31	0.51	0.42	0.37	0.44	0.40	0.62	0.44	0.47	0.33	0.29	0.30	0.31	
0.37	0.33	0.31	0.47	0.36	0.46	0.65	0.37	0.44	0.50	0.42	0.33	0.30	0.30	0.30	
0.37	0.41	0.41	0.34	0.41	0.60	0.43	0.49	0.45	0.55	0.34	0.44	0.34	0.30	0.26	
0.39	0.34	0.34	0.46	0.38	0.48	0.43	0.34	0.44	0.34	0.38	0.33	0.30	0.29	0.30	
0.35	0.45	0.43	0.57	0.38	0.52	0.59	0.38	0.48	0.46	0.49	0.35	0.30	0.29	0.26	
0.34	0.28	0.35	0.34	0.39	0.34	0.41	0.51	0.50	0.43	0.40	0.39	0.34	0.30	0.31	
0.33	0.31	0.41	0.35	0.39	0.38	0.39	0.31	0.36	0.32	0.32	0.38	0.32	0.31	0.25	
0.28	0.33	0.34	0.28	0.31	0.26	0.30	0.32	0.33	0.31	0.29	0.30	0.27	0.28	0.26	
0.27	0.27	0.27	0.30	0.35	0.30	0.28	0.27	0.31	0.29	0.31	0.26	0.36	0.35	0.25	

## Freeze-dryer Nr 2

FD04 5% sucrose 25 mg/ml BSA

Sh3

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
0.201	0.23	0.198	0.285	0.221	0.212	0.246	0.225	0.246	0.242	0.189	0.273	0.233	0.205	0.195	
0.223	0.299	0.238	0.317	0.267	0.335	0.254	0.243	0.354	0.191	0.312	0.237	0.211	0.185	0.204	
0.234	0.244	0.372	0.297	0.402	0.311	0.435	0.512	0.385	0.297	0.255	0.359	0.371	0.243	0.218	
0.323	0.272	0.439	0.5	0.373	0.361	0.398	0.407	0.405	0.376	0.458	0.325	0.377	0.295	0.261	
0.245	0.259	0.384	0.446	0.379	0.51	0.298	0.323	0.387	0.272	0.21	0.382	0.274	0.185	0.265	
0.297	0.271	0.456	0.324	0.331	0.401	0.458	0.469	0.464	0.482	0.323	0.43	0.302	0.279	0.211	
0.212	0.313	0.44	0.367	0.428	0.368	0.248	0.266	0.367	0.397	0.377	0.499	0.355	0.274	0.293	
0.284	0.252	0.368	0.246	0.408	0.353	0.434	0.347	0.433	0.393	0.388	0.285	0.416	0.291	0.266	
0.217	0.316	0.306	0.434	0.363	0.499	0.361	0.32	0.406	0.461	0.303	0.409	0.312	0.429	0.325	
0.27	0.383	0.397	0.355	0.422	0.377	0.466	0.551	0.294	0.429	0.32	0.282	0.38	0.278	0.26	
0.19	0.32	0.239	0.495	0.531	0.359	0.412	0.108	0.478	0.328	0.358	0.306	0.32	0.273	0.336	
0.387	0.305	0.375	0.335	0.455	0.63	0.486	0.381	0.398	0.424	0.418	0.214	0.406	0.369	0.208	
0.346	0.435	0.513	0.371	0.382	0.487	0.484	0.475	0.319	0.337	0.299	0.361	0.284	0.315		
0.209	0.258	0.318	0.583	0.503	0.372	0.495	0.341	0.448	0.451	0.387	0.349	0.287	0.384	0.205	
0.195	0.328	0.463	0.549	0.377	0.581	0.379	0.415	0.563	0.432	0.449	0.334	0.279	0.314	0.252	
0.289	0.279	0.32	0.482	0.476	0.646	0.442	0.317	0.527	0.398	0.404	0.408	0.359	0.301	0.349	
0.23	0.336	0.36	0.372	0.361	0.372	0.369	0.482	0.374	0.357	0.344	0.432	0.205	0.285	0.306	
0.365	0.306	0.403	0.515	0.36	0.438	0.492	0.396	0.484	0.453	0.355	0.377	0.492	0.385	0.21	
0.279	0.339	0.372	0.452	0.319	0.42	0.349	0.763	0.434	0.363	0.421	0.377	0.335	0.295	0.235	
0.32	0.26	0.313	0.514	0.296	0.371	0.481	0.443	0.474	0.389	0.329	0.353	0.329	0.3	0.231	
0.165	0.206	0.231	0.262	0.324	0.287	0.231	0.233	0.308	0.222	0.319	0.238	0.304	0.22	0.225	
0.166	0.18	0.203	0.172	0.17	0.168	0.188	0.172	0.176	0.213	0.225	0.186	0.167	0.199		

Headspace moisture analysis of statistical sample sets enabled efficient cycle optimization and clear demonstration of freeze dryer equivalence

# USP <922> Headspace moisture as a formal method



Connecting People, Science and Regulation®

- The previous case studies demonstrated that headspace moisture analysis gives very accurate and efficient insight into process and product quality during process development and scale-up.
- USP <922> describes headspace moisture (water activity) as a formal moisture determination method that in many cases is more useful than a total moisture content determination.

- **What do you mean by ‘water content’?**

- adsorbed water?
- chemically bound water?
- free water?
- hydration shells (e.g., of proteins)?
- water of crystallization?

**Not just pragmatically, but fundamentally the moisture content paradigm being used by the pharmaceutical industry (total water content) is not the correct one.**



## USP<922> Water Activity implemented in Q1 2021.

- “Total water content is an important quality attribute...However, water may be allocated in more than one compartment within these materials. Some of it may be tightly bound and not available to participate in chemical, biochemical, or physicochemical reactions (e.g., as hydrate salts), whereas some of the water may be more loosely bound and more freely available to participate in reactions... **It is important to establish what fraction of the total water is contained in the latter category, and the determination of water activity ( $a_w$ ) provides this information.**”
- “The determination of  $a_w$  aids in decisions during ingredient and product processes design, ingredient selection, packaging selection, and product storage. These include:
  - ...
  - Reducing the degradation of active ingredients within product formulations, especially those susceptible to chemical hydrolysis
  - ...
  - Providing a complementary method to the Karl Fischer titration for monitoring changes in water content
  - Controlling and monitoring physical, chemical, and microbial product stability
  - ...”

Headspace moisture


Total cake moisture  
(Karl-Fischer)



Formulation A



Formulation B



Formulation C

INCREASING WATER ACTIVITY

Total water content in the cake is equal in all product vials but the water activity is very different (why? - ratio of bound to unbound water is different).


- Assess correlation between headspace water vapor pressure (water activity) and degradation of a hydrolytic small molecule API (mannitol excipient)
- Measure initial headspace moisture of several lyo sample groups
- Since measurement is non-destructive, the same samples can be put on stability
- Degradation of the API measured on stability with HPLC


Moisture	0.50%	1.00%	1.50%	2.00%	3.00%
25°C	X	X	X	X	X
40°C	X	X	X	X	X
60°C	X	X	X	X	X

<b>Targeted moisture content groups:</b>	<b>0.50 %</b>	<b>1.00 %</b>	<b>1.50 %</b>	<b>2.00 %</b>	<b>3.00%</b>
<b>Corresponding water vapor levels:</b>	<b>4 torr</b>	<b>6 torr</b>	<b>8 torr</b>	<b>10 torr</b>	<b>14 torr</b>



## Product degradation results





$$k = \frac{-(\text{Ln}(\frac{[At]}{[Ao]}))}{t}$$

- Correlation of product stability to water activity is exceptional
- **Potential to replace KF with rapid non-destructive technique** for lyo moisture determination



ELSEVIER

Contents lists available at ScienceDirect

## Journal of Pharmaceutical Sciences

journal homepage: [www.jpharmsci.org](http://www.jpharmsci.org)

Pharmaceutics, Drug Delivery and Pharmaceutical Technology

### Noninvasive Moisture Detection in Lyophilized Drug Product Using NIR Spectrometer and Headspace Moisture Analyzer

Yu-Elaine Pu\*, Ying Li, Dong Xiang

*Drug Product Development, Bristol-Myers Squibb, Summit, NJ, United States*

#### ARTICLE INFO

**Article history:**

Received 2 August 2022

Revised 8 November 2022

Accepted 8 November 2022

Available online xxx

**Keywords:**

NIR

Laser headspace


Moisture detection

Lyophilization

#### ABSTRACT

Utilization of laser headspace and near-infrared (NIR) methods provides rapid and non-destructive approaches for moisture detection of lyophilized drug products to facilitate lyophilization formulation characterization and process development. In the present study, the NIR method was developed based on a partial least square regression (PLSR) model calibrated and validated with Karl Fisher (KF) data, whereas the laser headspace method was developed with aid of dynamic vapor sorption (DVS) method so that the water vapor pressure measured from the headspace of a lyophilized drug product vial can be converted directly to water content value through the water vapor sorption isotherm of the lyophilized drug product bypassing KF calibration. The water contents of lyophilized samples obtained from both methods agreed well with KF data, with a root mean squared error of prediction (RMSEP) of less than 0.15%. The pros and cons of NIR and laser headspace method were evaluated. The results suggest that traditional off-line KF method can be potentially replaced by at-line laser headspace method combined with water sorption isotherm data from DVS. Further studies may be needed to evaluate the quantitation limit and generality of this method to a variety of lyophilized formulations.

© 2022 American Pharmacists Association. Published by Elsevier Inc. All rights reserved.




DVS system used to generate a moisture adsorption isotherm for a lyophilized biopharmaceutical formulation.

**Isotherm correlates water weight gain (Wt % Water) and relative humidity conditions (%ERH water activity) measured in the headspace**

After generation of a moisture isotherm, a rapid non-destructive water activity measurement gives both the water activity and water content of a formulation.

**KF titration not needed anymore.**

Pu, Yu-Elaine, Ying Li, and Dong Xiang. "Noninvasive Moisture Detection in Lyophilized Drug Product Using NIR Spectrometer and Headspace Moisture Analyzer." *Journal of Pharmaceutical Sciences* (2022). <https://doi.org/10.1016/j.xphs.2022.11.009>



Book | © 2023

## Principles and Practices of Lyophilization in Product Development and Manufacturing

[Home](#) > [Book](#)**Editors:** [Feroz Jameel](#)

In depth coverage of lyophilization processes with case studies

Presents engineering antibodies as therapeutic agents and new FDA requirements

Includes cryopreservation of biologicals and CFD modeling

**Part of the book series:** [AAPS Advances in the Pharmaceutical Sciences Series](#) (AAPS, volume 59)[Home](#) > [Principles and Practices of Lyophilization in Product Development and Manufacturing](#) > [Chapter](#)

## Laser-Based Headspace Moisture Analysis for Rapid Nondestructive Moisture Determination of Lyophilized Products

[Derek Duncan](#) , [James R. Veale](#), [Ken Victor](#) & [Adriaan H. de Goeij](#)Chapter | [First Online: 25 April 2023](#)**215** AccessesPart of the [AAPS Advances in the Pharmaceutical Sciences Series](#) book series (AAPS, volume 59)

### Abstract


Residual moisture content is a critical quality parameter for freeze-dried pharmaceutical product and exists in different forms. Only the relatively “free” – in contrast to the more “bound” – water is capable of inducing degradation events within the product. However, traditional moisture analysis techniques employing gravimetric or Karl Fischer (KF) titration methods cannot distinguish between the different water states. Furthermore, these methods



35 years of using and studying freeze drying: what did we learn?


Gerhard Winter, LMU München

08.05.2024, Torino ISL-FD




## Benefits and interpretation of water activity

- $a_w$  is a comparable value
- $a_w$  enables comparison of environments for the antibody
- Close to  $a_w$  of native antibody in liquid formulation
- Degradation reactions dependent on  $a_w$ <sup>[1]</sup>



$0.10 < a_w < 0.25$

Various excipient (combinations) and residual moistures



$a_w$

Q How is  $a_w$  modeled?

Thermodynamic aspects in stabilisation of biopharmaceuticals in lyophilized formulations



Christoph Brandenbusch, Gabriele Sadowski

[1] Sandulachi, 2012



- Current moisture determination methods for freeze-dried pharmaceutical product (KF titration, LOD) give limited insight into process and product quality.
- Headspace moisture analysis (water activity) is an analytical tool useful for formulation development, process optimization, freeze dryer qualification, and as a formal water activity test method to determine product stability.
- USP <922> is recommending changing the moisture paradigm from total water content to water activity.

**Goal: A robust science-based understanding of process and product quality.**



Connecting People, Science and Regulation®

**Thank you  
for your attention!**

