PDA Week 2025

Palm Springs, CA 2025

Material and Packaging Selection and Compatibility with Chlorine Dioxide Gas Sterilization

Chlorine Dioxide Gas Properties

Chlorine dioxide gas (CD) is a unique single electron-transfer-oxidizing agent and has been termed a "stable free radical." CD has the microbiocidal properties of chlorine and is highly soluble in water, but unlike chlorine, its reaction chemistry does not lead to the formation of chlorinated organic products (trihalomethane and chloramine). CD is a respiratory irritant with an 8-hour timeweighted average exposure limit of 0.1 ppm; the 15-minute short-term exposure limit is 0.3 ppm. At use concentrations, CD is not flammable nor is it carcinogenic or ozone depleting (Long, Battisti, Elyath 2000). In addition, chlorine dioxide gas breaks down into chlorite, chlorate, and chloride which have relatively low toxicities in humans. CD gas and its derivatives, chlorite and chlorate, on the other hand, have relatively low toxicities in humans (11-13, 18, 24). They are not mutagenic (16, 17).

Chlorine Dioxide Gas Generation

Chlorine Dioxide (CD) gas by ClorDiSys is generated with the CD Cartridge (EPA label # 80802-1) and is registered with the US-Environmental Protection Association as a sterilant gas. The US-EPA label specifies chlorine dioxide can be used to sterilize prefilled syringes, cold chain products, medical devices, as well as other applications. Prefilled syringes and medical devices are typically sterilized in an ambient temperature vacuum sterilizer.

Chlorine dioxide is not sufficiently stable to be stored, it must be generated onsite at the point of use. There are many ways to generate CD gas and one method is using, a low-level chlorine gas (2%) which is passed over solid sodium chlorite CD Generating Cartridges which convert the chlorine to pure chlorine dioxide (>99% yield).

$Cl_2(g) + 2NaClO_{2(S)}$ to $2ClO_{2(g)} + 2NaCl_{(S)}$

The 2% chlorine gas is mixed with nitrogen. This stoichiometrically limits the chlorine dioxide gas concentration to 4% which is significantly below the 10% level which is potentially explosive. Chlorine purity levels should be greater than 99%. The CD Generating Cartridges have a 1-year shelf life and generate CD gas for up to 300 minutes before they require replacement. The shelf life of reagent gas tanks are manufacturer specific and must be followed. The sterilizer control systems automatically prevent the CD Cartridges from being used for more than the allowable capacity. This process produces a pure chlorine dioxide gas and is injected into the target chamber for a fixed time or until a dosage setpoint has been met. The chlorine dioxide flow rate is 20 liters per minute +/- 4 liters per minute. After the exposure is completed, the gas is typically passed through a carbon scrubber to fully eliminate any chlorine dioxide gas. Alternatively, chlorine dioxide may be exhausted to the outside environment via house exhaust systems. Federal state and local regulations must be verified prior to exhausting CD gas.

Chlorine Dioxide Gas Properties and Cycle Monitoring

Sterilization chambers for chlorine dioxide include ambient pressure chambers or vacuum pressure chambers. The target chamber choice depends upon the product requiring sterilization. If the product is simple in its geometry and has no small openings, then an ambient pressure chamber can be used. If the product is complex in its geometry and has many lumens or tubing, or small openings, then a vacuum process may be required. A vacuum sterilizer will remove most of the air particles from both the chamber as well as internal portions of the device itself, then replace the air with humidity and CD gas. This allows the moisture, as well as the CD, to penetrate tight areas of a device or packaging.

The product must be wrapped in suitable packaging to maintain sterility. Tyvek packaging is often utilized to serve as the sterile barrier primary packaging since it allows CD gas to penetrate and does not allow organisms to contaminate the device inside. Primary packaging may also include blister packaging or other packaging alongside a gas permeable material, which includes Tyvek or medical grade paper. Product to be sterilized may also be placed into secondary packaging which can include unit cartons, display boxes, pouches, trays, etc. Additionally, items may be placed within tertiary packaging or shippers. Chlorine dioxide gas is compatible with the sterilization of cellulosic materials which permits the incorporation of paper labels on devices as well as fiberboard or corrugated packaging. This allows for the simplification of the sterilization process and flexibility on the packaging of a device.

Chlorine dioxide gas is considered a true ambient temperature process since the temperature is not manipulated by the process. Chlorine dioxide gas a yellow-green color and this allows the concentration to be monitored by a UV-Vis spectrophotometer. Throughout the cycle, samples of gas are pulled from the chamber and passed through the UV-Vis spectrophotometer to give a precise concentration reading. its concentration can be monitored and thus it has the ability to utilize product parametric release. Chlorine dioxide gas also has a boiling point between -20° C and -40° C at use concentrations, thereby making it a true gas at room temperatures (15-25° C). Being a true gas, it will not condense at use temperatures like a vapor would. Instead, a gaseous process is a more reliable and consistent process with better penetration. Chlorine dioxide has a molecular weight of 67.5 but stratification is not an issue with chlorine dioxide. Testing at ambient pressure was performed demonstrating that with minimal circulation, gas concentration is still evenly distributed throughout a chamber.

Sterilization Chambers

Sterilization chambers for chlorine dioxide include ambient pressure chambers or vacuum pressure chambers. The target chamber choice depends upon the product requiring sterilization. If the product is simple in its geometry and has no small openings, then an ambient pressure chamber can be used. If the product is complex in its geometry and has many lumens or tubing, or small openings, then a vacuum process may be required. A vacuum sterilizer will remove most of the air particles from both the chamber as well as internal portions of the device itself, then replace the air with humidity and CD gas. This allows the moisture, as well as the CD, to penetrate tight areas of a device or packaging.



For chlorine dioxide gas the typical pressure/temperature chamber sterilization cycle consists of the following steps

Step 1- Pull an initial vacuum to remove much of the air out of the chamber (typically 5-10 kPa).

Step 2- Perform a leak test to ensure that there is no significant air leakage into the chamber.

Step 3- Precondition. This step raises the relative humidity to the setpoint (generally 65-75%).

Step 4- Condition. This step holds the relative humidity at the setpoint for typically 30-90

Step 5- Charge. This step raises the chlorine dioxide gas concentration to the setpoint (typically 1-30 mg/L).

Step 6- Exposure. This step holds the chlorine dioxide gas concentration at the desired CD concentration, expressed in mg/L, until either the dosage or the exposure time is attained. If the concentration drops for any reason the Steridox-100 will inject more gas into the system until concentration returns to its desired setpoint.

Step 7- Aeration. This step removes the chlorine dioxide gas from the chamber by pulling vacuum and breaking with filtered air. Typically, 12 vacuum/break cycles are required to remove the gas to safe levels.

Step 8- Introduction of filtered air into the chamber to bring the chamber back up to ambient pressure (HEPA filtered air is an option).

Sample Material Compatibility with Chlorine Dioxide Gas Exposure

Testing from Syensqo involved the following materials that are commonly used in healthcare and pharmaceutical applications. Chlorine dioxide exposure was highly robust and represented an atypically high dose cycle to mimic a challenging product's requirements. All testing indicates


- IXEF GS-1022 GY02
- Udel P1700 NT15
- Veradel 3300 WH001
- Radel R-5100 NT15
- Avaspire AV-651 BG15
- Avaspire AV-651 BG20 GF30

ł	nigh material co
Syensqo Material	Change in MW (%) from co hour expos
Ixef© PARA	< 1.5%
Udel© PSU	<1.3%
Veradel© PPSU	<1.1%
Radel© PSU	<0.7%
Avaspire© PAEK	<1.2%
Avaspire© GF PAEK	Re-analysis re
ensure there was no chemical interaction bet vzed each sample after exposure for changes	in molecular weight. A signific

Testing:

Four Delrin grades – SC655, SC631, PC652, SC698


- **Properties Evaluated:**
- Tensile properties
- Notched Charpy impact
- Color change
- Weight change
- Delrin grade Grade description SC655 NC010 Medium viscosity acetal homopolyme Medium-high viscosity, nucleated, acetal SC631 NC010 Medium viscosity, lubricated acetal homopolyme Low viscosity, lubricated acetal homopolymer SC698 NC010 (with improved processing and surface finish)





(Above) Steridox-100 Chlorine Dioxide Gas Sterilizer

SYENSQO



*Nominal strain at break reported

MD&M West • SPE MPD MiniTec 2025


Test outline:

ClorDiSys sterilized Oliver Packaging's proprietary films used in pouch configuration to determine compatibility of the packaging after exposure to chlorine dioxide gas sterilization and vacuum levels as deep as 5 kPa

Cycle Parameters:

3500 ppm-hour chlorine dioxide gas setpoint 1 mg/L chlorine dioxide gas concentration setpoint 70% relative humidity for 60 minutes

Biaxially Orier



Elongation %

MD	Non-sterile	CIO2
	Initial	Initial
	%	%
1	125	129
2	123	120
3	124	116
4	118	111
5	142	113
6	123	105
Avg	126	116
CD	Non-sterile	CIO2
CD	Non-sterile Initial	CIO2 Initial
CD		
CD 1	Initial	Initial
	Initial %	Initial %
1	Initial % 84	Initial % 114
1 2	Initial % 84 87	Initial % 114 106
1 2 3	Initial % 84 87 87	Initial % 114 106 116
1 2 3 4	Initial % 84 87 87 87 87	Initial % 114 106 116 117

Puncture Resistance

	Non- sterile	CIO2
Thru PET	Initial	Initial
1	6.27	6.52
2	6.43	6.55
3	6.31	6.50
4	6.3	6.35
5	5.74	6.46
AVG	6.21	6.38

4. RESULTS	
Acceptance Criteria	The test procedure was consi At the 48-hour timepoint, the cytotoxic grades of 0. At the 48-hour timepoint, the 4. All replicates of a sample rece
The device packaging was remove	
was added to submerge the device. The test article was extract	
Pre-Extraction Observations	The test article appeared nor The test extract appeared nor
Post-Extraction	The test article appeared norr
Observations	The test extract appeared nor
After extraction, the cells were dosed with 1.0 ± 0.1 ml of study a	
and ≥ 70% humidity for 24 hours. The triplicate wells were grad	
	Test Artic
Replicate	1
24-hour timepoint	0
48-hour timepoint	0
Highest Reactivity	None
Highest Classification	Non-cytotoxic
Deviations or Exclusions:	No deviations or exclusions
The final results relate only to the te	est items as received by the lab
	11 C

(Above) Cytotoxicity results for pre-filled syringe blister tray and Tyvek® lid, concluding that it is non-cytotoxic after sterilization

By Emily Lorcheim, PMP VP, Sterilization Technologies ClorDiSys

Packaging Testing after Chlorine Dioxide Exposure Soliver Type/Format Material Description

UT-73/ LF-1250AV	Pouch	Uncoated 1073B Tyvek sealed to peelable PET/PE film (adhesive lam)
ST-7382C / UF-1250AV	Pouch	Coated 1073B Tyvek sealed to PET/PE film (adhesive lam)
UT-73 / LF-2550AV	Pouch	Uncoated 1073B Tyvek sealed to peelable Nylon/PE film (adhesive lam)
UT-73 / LF-2550	Pouch	Uncoated 1073B Tyvek sealed to peelable Nylon/PE film (extrusion lam)
ST-7382C / PETG	Coated Tyvek sealed to PETG strip	Water-based adhesive coated 1073B Tyvek lidding sealed to PETG strip
XT-73NP34 / PETG	Coated Tyvek sealed to PETG strip	Hot melt adhesive coated 1073B Tyvek lidding sealed to PETG strip
HDPE CleanCut Cards	HDPE Card	28mil HDPE

Sample Test Report from Study

F-1250AV Peelable PET Film Lamination	
12μm (48ga) BOPET Adhesive 50μm (2mil) Peelable Sealant	



Left: Unsterilized Control Sample

Right: Sterilized Sample

Tensile Strength		
MD	Non-sterile	CIO2
	Initial	Initial
	lb f	lb f
1	16.42	17.0
2	16.38	16.5
3	16.52	15.9
4	16.38	16.8
5	17.18	16.1
Avg	16.6	16.4
	PSI	PSI
	6568	6792
	6552	6596
	6608	6352
	6552	6700
	6872	6448
Avg	6630	6578
CD	Non-sterile	CIO2
	Initial	Initial
	lb f	lb f
1	17.41	18.9
2	17.82	19.1
	17.70	19.5
3		17.0
3 4	18.10	19.6
	18.10 18.30	
4		19.6
4 5	18.30	19.6 18.8
4 5 6	18.30 18.09	19.6 18.8 19.4
4 5 6	18.30 18.09 17.90	19.6 18.8 19.4 19.2
4 5 6	18.30 18.09 17.90 PSI	19.6 18.8 19.4 19.2 PSI
4 5 6	18.30 18.09 17.90 PSI 6964	19.6 18.8 19.4 19.2 PSI 7564
4 5 6	18.30 18.09 17.90 PSI 6964 7128	19.6 18.8 19.4 19.2 PSI 7564 7564
4 5 6	18.30 18.09 17.90 PSI 6964 7128 7080	19.6 18.8 19.4 19.2 PSI 7564 7564 7564
4 5 6	18.30 18.09 17.90 PSI 6964 7128 7080 7240	19.6 18.8 19.4 19.2 PSI 7564 7564 7640 7792


Peel Strength Non-sterile CIO2 Initial Initial os 3 1.99 Pos 3 2.23 1.98 Pos 4 2.89 2.30 Pos 4 2.08 2.58 Avg 2.63 2.61

Property	Non-Sterile	CIO ₂
Elongation	MD: 126% CD: 90%	MD: 116% CD: 111%
Puncture Resistance	6.211bf	6.38 lbf
Tensile Strength	MD: 6,630 psi CD: 7,161 psi	MD: 6,578 psi CD: 7,657 psi
Peel Strength	2.63 lb/in	2.61 lb/in

esult: no loss in physical performance properties

Cytotoxicity Testing

A study included variety of products and packaging have been tested for cytotoxicity. Packaging is to be tested for cytotoxicity to comply with ASTM F2475-20. All devices and packaging were determined to by non-cytotoxic after sterilization.



4. RESULTS he test procedure was considered valid At the 48-hour timepoint, the triplicate wells of the vehicle and negative controls received Acceptance Criteria At the 48-hour timepoint, the triplicate wells of the positive control received cytotoxic grades of e interior components of the pre-filled syringe primary packaging were removed and excluded. 35.0 ml of extraction media wa lded to the inside of the packaging. The test article was extracted at 37 °C for 24 hours with agitation set at 60 RPM. e test article appeared norma bservations The test extract appeared normal in triplicate wells and incubated in at 37 °C, 5% Co ter extraction, the cells were dosed with 1.0 \pm 0.1 ml of study article extract nd ≥ 70% humidity for 24 hours. The triplicate wells were graded at 24- and 48-hour timepoints. **Test Article Results** Replicate 24-hour timepoint 48-hour timepoint Highest Reactivity Highest Classification Non-cytotoxic Deviations or Exclusions: No deviations or exclusions e final results relate only to the test items as received by the laboratory and only to the items specific to this projection

(Above) Cytotoxicity results for pre-filled syringe concluding that it is non-cytotoxic after sterilization