

## Chloé Stoll, PhD

Clinical & Human Factors Senior Specialist – Medical & Scientific Affairs – BD Medical - Pharmaceutical Systems

### Background

User interface design is a critical component of generic drug-device combination product development that must be considered from an early stage. The FDA guidance on **Abbreviated New Drug Applications (ANDA)** advises pharmaceutical companies to carefully consider and seek to minimize differences from the **Reference Listed Drug (RLD)**. To assist pharmaceutical companies in identifying user interface differences between an RLD and the Generic Drug-Device Combination Product (GDDCP), the FDA guidance on ANDA detailed how to perform a **Threshold Analysis**:

#### Threshold Analysis

##### Labelling Comparisons

Side-by-side, line-by-line comparison of the relevant sections of the prescribing information, instructions for use, and descriptions of the delivery device constituent parts of the GDDCP and its RLD

##### Comparative Task Analysis

Systematically analyze and compare the sequential activities required for the end-users to use the device and administer the drug product

##### Physical Comparison of Delivery device Constituent Part

Visual, auditory, tactile examination of the physical features (size, shape, feedback) of the RLD, compared to those of the delivery device constituent part of the proposed generic combination product

In the threshold analysis, points of comparisons are classified in **3 different outcomes**

	Outcome	FDA expectation
✓	No design difference	No additional data needed
✓	Minor design difference	Rational
✗	Other design difference	Additional information/data needed.

A **CUHF study** assesses whether the differences identified introduce a risk that might impact the clinical effect or safety profile of the generic combination product

### Delivery Device Comparison

A **Drug Delivery Device Developer (DDDD)** can support pharmaceutical companies in their **ANDA submission** of final GDDCP by authoring a portion of the Threshold Analysis. A **comparative analysis** is a partial threshold analysis where a DDDD compares the delivery devices on the aspects that are not drug-related.

#### Comparative Analysis

RLD & GDDCP (as component) Device comparisons

**Comparative Task Analysis**

- ✓ Overview of design differences & potential associated risks
- ✓ Support to pharmaceutical companies to build HF strategy

\*Labelling comparison is not in scope of DDDD

- ✓ If no 'other design differences' are identified, pharmaceutical companies can leverage DDDD data to complete their Threshold Analysis with labelling comparisons.
- ✗ If 'other design differences' are observed, pharmaceutical companies will have to provide additional data to FDA by performing, for example, a final CUHF study with final drug device combination product.

#### Abbreviations


ANDA : Abbreviated New Drug Application  
 RLD: Referenced Listed drug  
 GDDCP: Generic Drug Device Combination Product  
 CUHF : Comparative Use Human Factors  
 DDDD : Drug Delivery Device Developer

#### Reference

FDA Draft Guidance: CUHF Studies for a Drug-Device Combination Product Submitted in an ANDA, January 2017.

### Formative CUHF Study

A DDDD can conduct a **formative CUHF** study (with GDDCP filled with water) with a low number of users to obtain insight on the user errors rate of a final CUHF against the RLD.




#### A Formative CUHF Study:

- Does not replace the final CUHF study
- Provides insight on the use error rate and potential related risks
- Is a basis to design the final CUHF study with final combination product as per the FDA guidance (e.g., sample size, non inferiority margin)

### Case Study: A 3-Step vs 4-Step Autoinjector Comparison

BD Physioject™ Disposable Autoinjector (GDDCP) compared to Tirzepatide Autoinjector (RLD)

#### Comparative Analysis:

Similar features	Potential consequence
<ul style="list-style-type: none"> <li>✓ Activation button</li> <li>✓ 360° viewing window</li> </ul>	
	
Other differences identified	Potential consequence
<ul style="list-style-type: none"> <li>x Different unlocking mechanism RLD: Turn lock ring to unlock GDDCP: press needle cover to unlock</li> <li>x Different end of dose indicators RLD: visual and audio GDDCP: visual</li> </ul>	<ul style="list-style-type: none"> <li>Cannot unlock the autoinjector Potential risk: delayed treatment</li> <li>Premature removal of the device Potential risk: underdosing</li> </ul>

#### Formative CUHF Study:


Two tasks were impacted by an "other design difference" and assessed in the study:

- Unlock the device
- Wait until the injection is complete

#### Results:

- ✓ Same use error was observed: premature removal of the device
- ✓ No use error observed with the unlocking mechanism
- ✓ No new risk has been identified in this formative study

	RLD	DDCP
Use error rate	0.13	0.2
Within subject correlation	0.8	



### Conclusion

- In this formative study, the BD Physioject™ Disposable Autoinjector appears to be a suitable candidate to be a GDDCP of the compared RLD.
- DDDDs can support pharmaceutical companies to de-risk their ANDA submission by proactively developing supportive documentation and conducting formative CUHF studies to support final CUHF study design and risk evaluation.