


VHP as a Utility



john_klostermyer@steris.com

VHP as a Utility

Outline

- VHP 101
 - Advantages
 - Process
 - Vapor vs. Mist
 - Kill Curve
 - Emissions
- Portable or Modular
 - Why Modular
- Installation
 - Modular Schematic
 - Options for Integration
 - Single Pass
 - Recirculating
- Examples
 - Large Room
 - Pass-Throughs
 - HEPA's
 - BSL Labs
 - Automated Sequences

Why Use VHP?

✓ Consistency & Distribution

- Wet surfaces / minimal contact times -not an issue
- Passes through HEPA filters
- Decontaminates biosafety cabinets and HEPAs during room decon
- Kills airborne and surface microbes

✓ Labor

- Minimal labor required
- Easy to validate

✓ Environmental

- Excellent material compatibility
- Low toxicity
- No residues
- EPA approved



Reduce



Boiling Points:

H₂O 100 C


H₂O₂ 150°C



...if you can see it,
it's not a vapor


VHP[®] Kill Matrix

1 mgH₂O₂/liter air = 720ppm



G. stearothermophilus spores inoculated on
Stainless Steel Coupons at 30°C

VIRUS III



**300m3 (10600ft3)
Application Example**

Which VHP System?

portable



- ✓ Spaces not yet defined
- ✓ Uses in different buildings
- ✓ Typically less than 10,000ft³
- ✓ Cycle time not a constraint
- ✓ Use of fans not an issue
- ✓ Less frequent use

modular



- ✓ Large and small spaces up to ~80,000 ft³
- ✓ Same enclosures repeatedly
- ✓ Frequent use (chamber)
- ✓ Short cycle times
- ✓ Automated sequenced decontamination of multiple rooms

Why Modular?

Keep Equipment Outside Space

Save space within room / Pass Through

Avoid cross contamination

Keep maintenance activities outside

No Set Up

Decon at the Press of a Button

Run Sequential Decons via BMS*

Reduced Handling of Peroxide

Excellent Distribution

Cost

Less expensive than multiple portables

Save on labor

The easier to use – the more frequent
the use – the cleaner the space

*BMS = Building Management System



Pharmacy Compounding

Syringe Decon







Chambers

VHP Modular Integration Schematic


Single Pass





Utilities & Communications

Single Pass System



Large Room Biopharmaceutical Fermentation Suite

Volume: 32,000ft³ (900m³)

Ceiling height: 28ft (8,5m)

Single pass, No fans

6-log reduction

Cycle time 6 hours

New construction

Cycle Phase	Time min.	Airflow	Injection g/min
Dehumidification	30	6 A.E./ hour	-
Condition	30	120 cfm	96
Decontamination	90	120 cfm	60
Aeration	210	40 A.E./ hour	-



A.E. = Air Exchange

Pass-Through Chambers



Shown above

Enclosure v Volume ft3	Enclosure surface material	Injection rate Condition g/min.	Injection rate Decon g/min.	Decon time Min 6 log	Decon airflow ft3/min	PPM	Aeration airflow ft3/min	Total Cycle Time min.
460 (6x8x9.5'L)	Stainless	32	23	12	120	1000	765	45
175 (4x6x7'L)	Epoxy paint	12	9	8	40	950	1750	30

Automated Sequential Zone BSL3 Lab Decontamination Single Pass



EMERGENCY BUTTON FOR INITIATION OF DECONTAMINATION SYSTEM, SIMILAR TO FLEE MODEL, ABOVE, PROVIDED WITH 5021-57 WALL PLATE, PROVIDE ENCLOSURE WITH FLIP LID SIMILAR TO FLEE MODEL, 501-1000, PROVIDE COVER WITH CUSTOM LABEL TO READ "EMERGENCY DECONTAMINATION LIFT LID AND PUSH BUTTON". MOUNT BUTTONS AT 42" ABOVE FLOOR TO CENTER OF BUTTON. (TYPICAL OF 8 LOCATIONS)



Zone	
	1
	2
	3
	4
	5
	6
	7

BSL Lab & BSC Decon Single Pass

Simultaneous Decon. of Primary Containment

A2 type biosafety cabinets can be decontaminated together with the room

- exhaust dampers above cabinets are closed
- cabinet blowers left on




HEPA Filter Decontamination



RABS – (Restricted Access Barrier Systems)

Modular VHP systems can rapidly decontaminate both Active & Passive RABS and the rooms they are housed in




U.R.S.

- VHP is treated as a building utility
- Decon can be performed after each
 - product change
 - product batch
 - during an incident
 - after maintenance
- No sealing of perimeter required, surrounding areas can continue to operate
- Zone 1 (AHUs 2-4, 2-5, and 2-6) will be decontaminated together
- Zone 2 (AHUs 2-9, 2-10 and 2-11) will be decontaminated together
- Multiple zones will be decontaminated sequentially
- Decontamination time per zone – about 2 hours *
- Aeration time per zone about 2 hours * (at 10 AE/hour)

* Typical timing only, actual timing depends on room load and configuration

Scope of work

- Conceptual design and specifications
 - Material, size, maximum lengths of VHP cPVC piping, connection to AHU ducting, etc
 - Type and sizing of dampers and valves
 - Recommended air exchanges
 - Utility requirements
- Equipment skid (1 x T4, 1 x Munters, handshake dry contacts, etc)
- SS braided tubing from T4 to bulk Vaprox (up to 10m)
- Sensors for H₂O₂ safety monitoring – 1 per zone
- Installation supervision



Process Concept Distribution

Automatización

VIRUS-I

HR 41 HR 42 Prog. Horario **Horarios**

Desinf.

- Deshumidificación
- En cola
- Maniobra 1
- Fumigación
- Maniobra 2
- Ventilación
- Maniobra 3

Consignas

VIRUS-II

HR 14 HR 15 Prog. Horario **Horarios**

Desinf.

- Deshumidificación
- En cola
- Maniobra 1
- Fumigación
- Maniobra 2
- Ventilación
- Maniobra 3

Consignas


VIRUS-III

HR 10 Prog. Horario **Horarios**

Desinf.

- Deshumidificación
- En cola
- Maniobra 1
- Fumigación
- Maniobra 2
- Ventilación
- Maniobra 3

Consignas



CÉLULAS 1

HR 47 HR 48 Prog. Horario **Horarios**

Desinf.

- Deshumidificación
- En cola
- Maniobra 1
- Fumigación
- Maniobra 2
- Ventilación
- Maniobra 3

Consignas


CONCENTRACIÓN

Prog. Horario **Horarios**

Desinf.

- Deshumidificación
- En cola
- Maniobra 1
- Fumigación
- Maniobra 2
- Ventilación
- Maniobra 3

Consignas



MandoGeneral Desinfecciones

Consignas de desinfección

Humedad

Tiempo aireación





Direct injection setup

Technical Area




VHP M1000-T4

Proposed area of installation
on AHU mechanical floor




Zone 1

Zone 2



Recirculation concept




BI's incubated at 55° C for 7 days



Application Cycle Developpement





F & B

**Horizontal Filling Lines
Flash Decon Theory**



F&B application



Selected Installations

Pharma /

Animal Health

- GSK
- Intervet
- Fresenius
- Sanofi
- Merial
- Alcon
- Pfizer
- Boehringer
- Hisun

Public Health Labs

- Indiana
- New Jersey
- West Virginia

Others

- Tripler- US Army
- Univ. Nebraska
- Lawrence Livermore
National Labs
- INRS
- NCI
- WP-AFB



B & V TESTING and STERIS Advanced Biodecontamination Solutions (ABS)



Advanced Biodecontamination
Solutions

STERIS®



- **ABS: a Flexible service offering matching the field service expertise of B & V TESTING with STERIS VHP technology and EPA-registered consumables**

Contract VHP biodecontamination services



- National leader in biodecontamination and contamination control technologies testing, certification and maintenance services:
 - 30 years experience performing gaseous Biodecontamination Services (formaldehyde, VHP, CD)
 - Testing, Certification and Maintenance of Cleanrooms, Biological Safety Cabinets, and HEPA-filtered systems

STERIS®



- STERIS is a global leader in infection prevention, contamination control, surgical and critical care technologies, and more. Manufacturers of EPA registered Vaprox® 35% and 59% hydrogen peroxide (EPA reg. no. 58779-4), Spor-Klenz (EPA reg. no. 52252-4-1043 and VHP® technology