

Laser Induced Fluorescence

LIF Fundamentals and Applications


TSI Inc. Model 9510-BD
BioTrak® Real-time viable particle counter



UNDERSTANDING, ACCELERATED

What is fluorescence?

- + Molecule absorbs energy from light source
- + Some energy is lost to molecular vibrations
- + Remaining energy released as light at a higher wavelength






What is Autofluorescence?

- + Cells contain many molecules that fluoresce when exposed to UV light
- + Metabolites of microorganisms
- + Look for autofluorescence to detect viable particles

Molecule	Approximate Peak Fluorescence (nm)
NAD(P)H	450
Retinol	500
Riboflavin	550
Folic acid	450
Pyridoxine	400
Tyrosine	305
Tryptophan	325
Flavin	540




Detecting Autofluorescence



Environmental Monitoring using LIF

- + Biofluorescent Particle Counters (BFPC)
- + Sampling method identical to particle counting

- Fully automated
- Real-time results





- + Viable Particle Count or Autofluorescing Unit (AFU)



Environmental Monitoring using Traditional Microbiology

- + Collect an air sample onto growth media passively or actively



- + Incubate plates to allow for the formation of visible colonies (≥ 2 days)
- + Count the colony forming units (CFU)



Method Comparison

No method detects everything

- + Each provides an estimate of airborne viable contamination based on detection capabilities
- + Traditional methods estimate based on growth under the chosen incubation conditions
- + LIF methods estimate based on fluorescence consistent with what is expected for viable microorganisms



Method Comparison

Timeliness of results can be the difference between preventing an adverse event or being able to react to one

- + Traditional methods take days to obtain results – can only try to determine the possible impact to product of the viable contamination detected
- + LIF methods obtain results in real-time – can be proactive and immediately work to mitigate the issue to prevent an impact to your product




Method Comparison

Manual methods can have data integrity concerns

- + Traditional methods rely on manual plate counts
 - susceptible to human error
- + LIF methods are fully automated



LIF Applications




LIF Applications

Root Cause Analysis



Real-Time Monitoring



Continuous Monitoring



Scanning



LIF Applications

Monitoring in Grade C/D Areas:

- + Continuous monitoring supplements current EM
- + Reduce risk by being proactive instead of reactive
- + Reduce sampling



LIF Applications

Monitoring Aseptic Processes:

- + Continuous monitoring system
- + Only requires an isokinetic probe in Grade A area
- + Immediate actions can be taken in the event of an excursion



In-process Environmental Monitoring

1. No interventions • No: product waste, nozzle clogging, downtime, extra sterility testing...
2. Increase process knowledge • real-time, continuous data
3. Real-time corrective action
4. Readily integrated & automated • Samples like total particle counter
5. Automated operation, alarms, analysis, etc... (FMS)

Benefits:

- **Increase Throughput**
- **Reduce Risk**



Implementation

Challenges

- + Understanding the results
- + Regulatory scrutiny
 - Will depend on application
 - Agencies encourage the use of new technologies
 - How will inspectors react?
 - BioPhorum Operations Group (BPOG)
 - Process and Environmental Monitoring Method (PEMM)



Implementation

BPOG Article in the PDA Journal April 2019:
Continuous Microbiological Environmental
Monitoring for Process Understanding and
Reduced Interventions in Aseptic Manufacturing
Jeffrey Weber, James Hauschild, Pieta Ijzerman-
Boon, et al.

PDA J Pharm Sci and Tech 2019, 73 121-134

Access the most recent version at
[doi:10.5731/pdajpst.2018.008722](https://doi.org/10.5731/pdajpst.2018.008722)



Implementation

- + Process and Environmental Monitoring Methods (PEMM) group shares information about the technology and applications for real-time biofluorescent particle detection methods and how they help to address major industry concerns. To find the recent article go to: <https://www.tsi.com/products/cleanroom-particle-counters/real-time-viable-particle-counter/biotrak-real-time-viable-particle-counter/>



Questions?

