




Theory 4:

- Mechanistic principles of (Parvo-) Virus retention
- Challenges of implementing virus filtration into continuous manufacturing
- Virus filters as bioprocess subject – current hot topics (ATMPS, facility segregation)


Mechanistic Principles of Parvovirus Retention

How does it work ?



Mechanistic Principles of Parvovirus Retention

Size exclusion enabled by pore structures across the membrane that trap viruses and allow protein molecules to pass through




Mechanistic Principles of Parvovirus Retention



Planova™ Filter



Hollow fiber




Membrane Cross section



Hollow fibers have a three-dimensional network structure consisting of voids connected by capillaries.


Mechanistic Principles of Parvovirus Retention

Parvovirus






~18 – 26 nm

Retrovirus



~80 – 110 nm

Mechanistic Principles of Parvovirus Retention



Cross sectional view

Porcine Parvovirus (PPV): 18 - 22 nm



Load: 12.32 log₁₀ (TCID₅₀/m²)

Immuno fluorescent Staining for PPV
FITC conjugated anti porcine parvovirus

5 g/L IgG, 100 mM NaCl, pH 4.5

Mechanistic Principles of Parvovirus Retention

Planova 20N





Virus load:12.32 log (TCID₅₀/m²)

Immuno fluorescent Staining for PPV
FITC conjugated anti porcine parvovirus

5 g/L IgG, 100 mM NaCl, pH 4.5

Mechanistic Principles of Parvovirus Retention

Planova BioEX




Virus load: 12.32 log (TCID₅₀/m²)

Immuno fluorescent Staining for PPV
FITC conjugated anti porcine parvovirus


5 g/L IgG, 100 mM NaCl, pH 4.5

Mechanistic Principles of Parvovirus Retention



All nanofilters on the market have a similar size exclusion mechanism. Different filter brand, pore structure & thickness: different permeability, protein loading capacity, parvovirus removal.

Mechanistic Principles of Parvovirus Retention




u ; Flow velocity is controlled by filtration pressure
 U_{cr} ; Critical velocity to overcome the Brownian motion

Mechanism 1: Size exclusion
Mechanism 2: Hydrodynamic force

*A Yamamoto, THongo-Hirasaki, YUchi, HHayashida and FNagoya.
 Effect of hydrodynamic forces on virus removal capability of Planova™ filters,
 AIChE Journal, 2014, 60(6): 2286–2297*

Mechanistic Principles of Parvovirus Retention

Virus behavior simulation under different level of flow velocity




The degree of virus mobility was controlled by flow velocity.
The movable distance of virus increased with pressure decrease.

*A Yamamoto, THongo-Hirasaki, YUchi, H Hayashida and FNagoya.
Effect of hydrodynamic forces on virus removal capability of Planova™ filters,
AIChE Journal, 2014, 60(6): 2286–2297*

Mechanistic Principles of Parvovirus Retention

Mechanism 3: physicochemical property & effect of solution condition (hypothesis)



Tomoko Hongo, Asahi, PDA 2013

Mechanistic Principles of Parvovirus Retention


Mechanism 3: physicochemical property & effect of solution condition (hypothesis)

	pl	pH	
		4	7
PPV	5-5.5	(+)	(-)
IgG (poly)	6.8-10	(++)	(+)
	PPV-protein interaction	PPV-IgG; repulsive	PPV-IgG; attractive (complex)

Tomoko Hongo, Asahi, PDA 2013

Mechanistic Principles of Parvovirus Retention

Mechanism 3: physicochemical property & effect of solution condition (hypothesis)



0 or 10 mg/ml IgG/ 100 mM NaCl, pH 7,
0.5 vol% serum-free PPV spiking, 230~250 L/m²

Pressure (bar)	PPV LRV(pool)		Relative Distance of Peak Position (%)	
	IgG(+)	IgG(-)	IgG(+)	IgG(-)
0.2	≥ 5.6	3.9	42	55
0.8	≥ 5.6	≥ 5.3	42	49


Pressure (bar)	IgG (-)	PPV LRV		
		pH 4	pH 5	pH 6
0.2		3.7	3.7	3.8
0.8		≥ 4.5	≥ 4.7	≥ 5.0

- Peak position without IgG shifted to outer side under low pressure.
- PPV behavior with/ without IgG was different at pH 7.
- At pH7, PPV with IgG may behave larger size.

Tomoko Hongo, Asahi, PDA 2013

Mechanistic Principles of Parvovirus Retention

Virus filters have a finite virus capture capacity



10 mg/mL human IgG
Serum-free PPV spiking


Tomoko Hongo-Hirasaki, Asahi, 2014 Planova Workshop (adapted)

Mechanistic Principles of Parvovirus Retention


VIRUS RETENTIVE FILTRATION

Virus Break Through

Dependency of Residual infectivity and total viral load
All filter // all viruses



Dependency of Residual infectivity and total viral load
All filter // all viruses < 30 nm



Questions ?



AsahiKASEI

Challenges of Implementing Virus Filtration into Continuous Manufacturing


Daniel Strauss, PhD
Principal Scientist
Asahi Kasei Bioprocess America

Agenda

- 1) Virus filtration integration
- 2) Fluctuating solution conditions
- 3) Viral clearance validation studies

What are Continuous Processes


Fully Integrated Continuous Process



(Konstantinov and Cooney, 2015)


Hybrid Processes

Upstream



(Konstantinov and Cooney, 2015)

Downstream



(Shamashkin, et al., 2015)

Amenable to Continuous Processing

- Flow-through process
- Constant flow-rate operation
- Viral clearance is robust
- Can be run for long times

There is very little data published for continuous virus filtration processes!

Reviews:

REVIEW

BIOTECHNOLOGY
and
BIOENGINEERING

Adapting Viral Safety Assurance Strategies to Continuous Processing of Biological Products

Sarah A. Johnson, Matthew R. Brown, Scott C. Lute, Kurt A. Brorson
DBRRII, Office of Biotechnology Products, Office of Pharmaceutical Quality, Center for Drug Evaluation and Research, Food and Drug Administration, Silver Spring, Maryland 20993; telephone: 240-402-5730; e-mail: sarah.johnson1@fda.hhs.gov

Evolving Needs For Viral Safety Strategies in Continuous Monoclonal Antibody Bioproduction

Andrew Clutterbuck,¹ Michael A. Cunningham,² Cedric Geyer,¹ Paul Genest,² Mathilde Bourguignat,¹ and Helge Berg¹


¹Technology Management, Millipore SAS, 39 Route Industrielle de la Hardt, 67124 Molsheim, France

²Technology Management, EMD Millipore Corporation, 290 Concord Road, Billerica, MA 01821, USA


Viral Filtration Processes

Example Continuous Virus Filtration Process


- Constant flow – set by process
- Initial transmembrane pressure based on vendor recommendation
- Switch to new filter:
 - Validated throughput, OR
 - Maximum TMP



Feedstock Variation




Output from a Continuous CEX Unit Operation




Solution Variations

- Product Concentration
- Salt Concentration
- pH
- Impurities

Feedstock Variation



Output from a Flow Through Q Membrane



Solution Variations

- Product Concentration
- Salt Concentration
- pH
- Impurities

Viral Clearance Robustness

Protein Concentration

IgG Conc. (g/L)	PPV LRV	
	Planova 20N	Planova BioEX
1	≥ 5.67	≥ 5.42
5	≥ 5.37	≥ 5.78
10	≥ 6.00	≥ 5.35
30	≥ 5.58	≥ 5.28
50	≥ 5.67	≥ 5.10

Salt Concentration

NaCl Conc. (mM)	PPV LRV	
	Planova 20N	Planova BioEX
1	≥ 5.84	N/A
50	N/A	≥ 5.48
100	≥ 6.00	≥ 5.28
200-250	≥ 5.67	≥ 5.28
500	≥ 6.00	≥ 5.92


(Hongo-Hirasaki, PDA Virus and TSE Safety Forum, 2011)

- Virus filters provide excellent viral clearance over wide ranges of conditions
- But what about the effects of fluctuating conditions during the filtration?

Collaboration between Asahi Kasei and FDA:

Do fluctuating solution conditions impact virus removal by VF?

Simulated peaks running virus filters using an AKTA:



Buffer A: baseline condition

Buffer B: Same as A except
one variable

Buffer A (Baseline condition):


- 1 g/L Human Gamma Globulin (HGG, SeraCare)
- 20 mM Acetate, 50 mM NaCl, pH 6.0
- Spiking with PP7 bacteriophage at ~ 7 log PFU/mL

Buffer B:

- 10 g/L Human Gamma Globulin
- 20 mM Acetate, 50 mM NaCl, pH 6.0
- Spiking with PP7 bacteriophage at ~ 7 log PFU/mL

Protein Peak

Planova 20N with High Protein Peak



- 0.001 m² filters
- Planova 20N at 0.5 mL/min = 30 LMH
- Planova BioEX at 1.0 mL/min = 60 LMH

Effect of Protein Peak

Planova 20N

Sample	PP7 Titer (log PFU/mL)	LRV _{Instantaneous}
Average Load	6.0	N/A
Pre-peak	≤ 0.0	≥ 6.0
Peak F1	≤ 1.0	≥ 5.0
Peak F2	≤ 1.0	≥ 5.0
Peak F3	≤ 1.0	≥ 5.0
Peak F4	≤ 1.0	≥ 5.0
Peak F5	≤ 1.0	≥ 5.0
Post-Peak F1	≤ 0.0	≥ 6.0
Post-Peak F2	≤ 0.0	≥ 6.0

Protein peaks have little or no impact on virus removal

Effect of Protein Peak

Planova BioEX

Sample	PP7 Titer (log PFU/mL)	LRV _{Instantaneous}
Average Load	6.7	N/A
Pre-peak	≤ 0.0	≥ 6.7
Peak F1	≤ 1.0	≥ 5.7
Peak F2	≤ 1.0	≥ 5.7
Peak F3	≤ 1.0	≥ 5.7
Peak F4	≤ 1.0	≥ 5.7
Peak F5	≤ 1.0	≥ 5.7
Post-Peak F1	≤ 0.0	≥ 6.7
Post-Peak F2	≤ 0.0	≥ 6.7

Protein peaks have little or no impact on virus removal


Effect of Protein Peak

Planova 20N – Run 2 (Higher than recommended load titer)

Sample	PP7 Titer (log PFU/mL)	LRV _{Instantaneous}
Average Load	8.1	N/A
Pre-peak	1.5	6.6
Peak F1	2.4	5.7
Peak F2	2.5	5.6
Peak F3	3.4	4.8
Peak F4	2.2	5.9
Peak F5	2.0	6.2
Post-Peak F1	2.6	5.6
Post-Peak F2	3.3	4.8

Effect of Protein Peak

Planova 20N – Run 2 (Higher than recommended load titer)



Protein peaks may impact virus removal under challenging conditions


Effect of Conductivity Peak

Buffer A:

- 1 g/L HGG in 20 mM Acetate, pH 6.0, ~7 log pfu/mL PP7
- 10 mM NaCl

Buffer B:

- Same with 500 mM NaCl




Effect of Conductivity Peak

	Planova 20N		Planova BioEX	
Sample	PP7 Titer (log PFU/mL)	LRV _{Inst.}	PP7 Titer (log PFU/mL)	LRV _{Inst.}
Average Load	7.8	N/A	7.6	N/A
Pre-peak	≤ 0.0	≥ 7.8	≤ 0.0	≥ 7.6
Peak F1	≤ 1.0	≥ 6.8	≤ 1.0	≥ 6.6
Peak F2	≤ 1.0	≥ 6.8	≤ 1.0	≥ 6.6
Peak F3	≤ 1.0	≥ 6.8	≤ 1.0	≥ 6.6
Peak F4	≤ 1.0	≥ 6.8	≤ 1.0	≥ 6.6
Peak F5	≤ 1.0	≥ 6.8	≤ 1.0	≥ 6.6
Post-Peak F1	≤ 0.0	≥ 7.8	≤ 0.0	≥ 7.6
Post-Peak F2	≤ 0.0	≥ 7.8	≤ 0.0	≥ 7.6


Conductivity peaks have little or no impact on virus removal

Impact of Feedstock Variation


B&E Chrom.



Virus Filter



Mixing Tank



Mixing Vessel


- Sufficient volume to mix multiple elution peaks
- Advantages:
 - Better utilization of filter area – no pressure spikes
 - Avoid fluctuating backpressure on upstream steps
 - Potential impact on validation strategy

Viral Clearance Validation Challenges:

- In-line spiking
- Complex equipment setup
- Long filtration runs (multiple days)
 - Laborious
 - Virus and product stability issues
- Startup and stop of process may differ from steady-state
- Product profile through full process may change




In-Line Spiking



Product Conc.

Virus Titer
(inline spike)

Virus Titer
(actual)



Throughput

In-Line spiking captures protein fluctuations, but not virus fluctuations


Effect of Virus Peak

Planova 20N

- Load A: 6.7 log PFU/mL
- Load B: 8.8 log PFU/mL

Planova BioEX


- Load A: none
- Load B: 7.8 log PFU/mL



All pools and fractions had complete clearance!


Effect of Protein, Salt, and Virus Peak

Load A	1 g/L HGG	50 mM NaCl	7.0 log PFU/mL
Load B	10 g/L HGG	500 mM NaCl	7.9 log PFU/mL



All pools and fractions had complete clearance!


Viral Clearance Validation



Product Conc.


Virus Titer
(spike)

Virus Titer
(actual)



With sufficient mixing, batch spiking may be representative

Planova 20N



Conditions:

0.001m² Planova 20N
0.15 g/L HGG
50 mM Acetate, pH 6.0
20 mM NaCl
0.5 mL/min = 30 LMH

- Load titers decreased over the course of the run, limiting potential LRVs

Summary

- **Continuous processes are coming!!!**

- Virus filtration implementation into a fully integrated continuous process has significant challenges associated
 - Integration strategy
 - Viral clearance validation
 - Integrity testing
 - Automation

- Virus filtration itself is highly robust...
 - But we still need to get better at demonstrating its capabilities

Virus Filters as Bioprocess Subject - Current Hot Topics

Asahi**KASEI**
BIOPROCESS

RAW MATERIAL SAFETY CONCEPTS – CASE STUDIES FOR MEDIA TREATMENT

Sebastian B. Teitz, PhD
Product Manager & Scientific Coordinator
Asahi Kasei Bioprocess Europe

CAACB Meeting on Raw Material Treatment,
April 2018, Cambridge, USA



www.ak-bio.com

© 2018 Asahi Kasei Medical Co., Ltd.

Introduction

- Raw Material Safety

Case Studies

- High Volume Media filtration
- Porcine Cirovirus
- Mycoplasma Treatment by Nanofiltration
- Nanofiltration of Microbial Fermentation Media Components

Conclusion

- Considerations



Introduction

Contamination by blood-borne pathogens
 → 10.000s of patients affected

1989	1992		1993	1994	1995/1996/1997
HIV	B19	HAV	HCV	HBV	HAV
PPSB	F-VIII	F-VIII	Ivlg	PPSB	F-VIII

the Why? - 2) Recombinant Proteins

- Impactful events!
- Shortage in drug supply to patient.
- Competitors product fast-tracked.

	Virus /Host Cell	Events
1985-1989	Orbivirus /CHO EHDV /CHO	2
1990-1994	MMV /CHO MMV /CHO	2
1995-1999	Reovirus /Hu 1° Kidney Vesivirus /CHO CVV /CHO	3
2000-2004	CVV /Unknown CVV /Unknown	2
2005-2009	Vesivirus /CHO Vesivirus /CHO MMV /CHO CVV /CHO MMV /CHO Vesivirus /CHO	6
2010+	MMV /CHO PCV-1 /Vero	2
Unknown	MMV /BHK-21 Human Adenovirus /HEK293 Reovirus /Unknown	3
Total		20

Mike Wiebe, CAACB, IBC Viral Safety, Huntington, Feb 25 2013

2) Recombinant proteins



Associated / suspected contamination, e.g. by rodents at some point in the supply chain.



A single mouse feces can contain
>10.000 i.u. of parvovirus.
(Besselsen et al. Comp Med (2008) 58: 140)

the Why? - 3) Viral Vaccines contamination events by Mycoplasma*

- Initially detected by elevated total DNA content in final doses
- Investigation confirmed *M.arginini* (+ some *M.fermentans* & *M.hyorhinis*) as contaminants
- Suspected sources: bovine (BSA, *M.arginini*), porcine (Trypsin, *M.hyorhinis*), human (commensal, *M.ferementans*)



*Eric Sarcey, Sanofi Pasteur, CAACB Workshop on Contamination with Difficult to Detect Bacteria, Boston, April 2016


the Why? - 4) Viral Vaccines contamination event with Porcine Circovirus*

- PCV1 identified as contaminant of the paediatric Rotarix vaccine through MPS
- 100.000s of children were exposed to a live virus.
- All Vero-cell banks back to MCB (1983) found positive for PCV1
- Suspected entry point: porcine Trypsin, used during MCB generation.




*Delwart et al., Viral nucleic acids in live-attenuated vaccines: detection of minority variants and an adventitious virus, J Virol. 2010 Jun;84(12):6033-40

Bioreactor protection by barrier filter



Barrier & Beyond | Damon Asher | 2014 PDA Europe Mycoplasma | Sep 2014


Bioreactor Virus Safety Assurance Model



Barrier & Beyond | Damon Asher | 2014 PDA Europe Mycoplasma | Sep 2014

- Specific nanofilter for USP ? → NO. Same virus as in DSP to be removed!
- “Low cost” nanofilter? → NO. High quality nanofilter required
- Much higher flux ? → NO → 1) Nano-pore limitation
2) or with sacrificing LRV
- High volumetric loads?
feasible → YES! Unlike in DSP, higher L/m² loads
↓ ↓ ↓
Longer filtration times required


1) CD-CHO Medium Filtration



- ✓ No impact of the virus spike on Filtration Volume
- ✓ **Consistent** performance
- ✓ **20N:**
2 000 L/m² in 1 day
5 000 L/m² in 3 days
- ✓ **BioEX:**
same as 20N
+ 10 000 L/m² in 7 days

Konstantin Agolli, Asahi Kasei, BioInnovation 2016, Berlin, February 10th, 2016

1) CD-CHO Medium Filtration




- ✓ **No virus detected (↑)**
- ✓ Difference in PPV LRV is due to differences in assay sensitivity

Konstantin Agolli, Asahi Kasei, BioInnovation 2016, Berlin, February 10th, 2016

1) CD-CHO Medium Filtration

Filter cost (€/L medium) vs. Filtration time (days)



- ✓ Price of CD-CHO:
< 30 €/L
- ✓ **BioEX:**
After 1.5 days, NF cost
< 10 % medium cost
- ✓ **20N:**
After 1 day, NF cost < 10
% medium cost
- ✓ **The longer the filtration
time, the more cost
effective !**

Konstantin Agolli, Asahi Kasei, BioInnovation 2016, Berlin, February 10th, 2016

2) Porcine Circovirus Nanofiltration

Planova 12.5 nm & 10 nm qPCR - removal data in DMEM

Planova 12.5 nm Filtration

Sample Code	PCV Loads- Run 1 (log ₁₀)			PCV Loads- Run 2 (log ₁₀)		
	-DNase	+DNase	RF*	-DNase	+DNase	RF*
SSM	7.17	6.85	-	7.17	6.86	-
PreF	7.00	6.83	-	7.00	6.83	-
NF2	6.64	4.59	2.24	6.69	4.54	2.29

Planova 10 nm Filtration

Sample Code	PCV Loads- Run 1 (log ₁₀)			PCV Loads- Run 2 (log ₁₀)		
	-DNase	+DNase	RF*	-DNase	+DNase	RF*
SSM	7.13	6.71	-	7.13	6.71	-
PreF	6.86	6.76	-	6.86	6.76	-
NF2 with conc.	NA	≤2.70	≥4.06	NA	≤2.70	≥4.06

* RF calculated relative to PreF

2) Porcine Circovirus Nanofiltration

Planova 15N qPCR - removal data in DMEM

Single 15N Filtration

Sample ID	PCV Loads- Run 1 (log ₁₀)			PCV Loads- Run 2 (log ₁₀)		
	-DNase	+DNase	RF	-DNase	+DNase	RF
SSM	7.07	6.66	-	7.07	6.66	-
PreF	6.86	6.71	0.00	6.86	6.71	0.00
NF2 (1x 15N)	6.57	5.15	1.56	6.34	5.24	1.47

Serial 2x 15N Filtration

Sample ID	PCV Loads- Run 1 (log ₁₀)			PCV Loads- Run 2 (log ₁₀)		
	-DNase	+DNase	RF	-DNase	+DNase	RF
SSM	6.99	6.85	-	6.99	6.85	-
PreF	6.86	7.04	0.00	6.86	7.04	0.00
NF2 (2x 15N)	6.50	4.25	2.79	6.65	4.25	2.79

2) Porcine Circovirus Nanofiltration


Planova qPCR - removal data in a mAb preparation

Sample ID	Planova 20N		Planova 15N		BioEX	
	+DNase	RF*	+DNase	RF*	+DNase	RF*
SSM	6,62	-	5.33	-	6.42	-
Hold control	6.20	-	5.05	-	6.39	-
PreF	6.59	-	5.58	-	6.21	-
NF1	≤4.18	≥2.41	≤3.58	≥2.00	≤ 4.18	≥2.44
NF2	≤3.98	≥2.61	≤3.38	≥2.20	≤ 3.98	≥2.64
NF3	≤3.78	≥2.81	≤3.18	≥2.40	≤ 3.78	≥2.84

Andy Bailey, ViruSure, 14th Planova Workshop, Cologne, Nov 09th, 2011

3) Mycoplasma Nanofiltration

Filtration Profile of diluted *Acholeplasma laidlawii* preparation (incl. beef heart broth, yeast extract, horse serum...)



3) Mycoplasma Nanofiltration

Removal of *Acholeplasma laidlawii* preparation by Planova filters

Filter	Titer (Log ₁₀ cfu/mL)		LRV
	Load	Filtrate	
P20N	7.24	≤1.65	≥5.59
PBioEX	7.24	≤1.65	≥5.59

Masayasu Takahara, Asahi Kasei, 19th Planova Workshop, Philadelphia, Sept 22nd, 2016

4) Nanofiltration of Microbial Fermentation Media Components

Summary of Nanofiltration experiments

Nr.	Media type	Concentration [g/L]	Volume per 4000 L scale [L]	PN20			BioEX		
				Flow	Average flux [L/h/m ²]	Area for 4000L scale [m ²]	Flow	Average flux [L/h/m ²]	Area for 4000L scale [m ²]
1	Glucose Feed	>100	>200	decrease	<10	>50	constant	10-100	>10
2	Vitamin solution	<50	<20	constant	10-100	<0.1	constant	>100	<0.1
3	Salt solution	>100	20-200	constant	10-100	0.1-0.5	constant	10-100	<0.1
4	Amino acid stock	<50	20-200	constant	10-100	<0.1	constant	10-100	<0.1
5	Tetracycline- alcohol	<50	<20	decrease	10-100	<0.1	blocked	n.a.	n.a.
6	Tetracycline- water	<50	<20	constant	10-100	<0.1	decrease	>100	<0.1
7	IAA solution	<50	20-200	constant	10-100	0.1-0.5	blocked	n.a.	n.a.
8	Trace elements solution	>100	<20	constant	10-100	<0.1	decrease	>100	<0.1
9	Kanamycin Solution	50-100	<20	constant	10-100	<0.1	constant	>100	<0.1
10	Fe-sulfate stock	50-100	<20	constant	10-100	<0.1	constant	>100	<0.1
11	Inducer	50-100	<20	constant	10-100	0.1-0.5	constant	10-100	<0.1
12	Media solution	<50	20-200	constant	10-100	0.1-0.5	constant	>100	<0.1
13	Sterile addition	>100	>200	decrease	10-100	>0.5	constant	>100	>0.5
14	Fe-chloride stock	>100	<20	constant	10-100	<0.1	constant	>100	<0.1

- Flow rates were higher with BioEX filter (pressure was also ~ 3 fold higher)
- No filter blocking observed for PN20 filter; BioEX: filter blocked for 2 media
- Not feasible for “glucose feed” (1) and “sterile addition” (13)

- Flow rates were constant and comparable for most media
- Application of PN20 seems feasible for most media
- Limitations due to organic solvents, viscosity and large volumes

Considerations

Considerations - Commercial Aspects Media Treatment

Case Examples – Results and Discussions

Low Flow Rate – 5 L/min

- Assumptions
 - Minimal capital costs
 - No flux decay
 - 100 LMH/bar, 3 bar, \$10000/m²
- Costs per year
 - Total = \$2.8M / year
 - \$23000 / batch

High Flow Rate – 100 L/min

- Assumptions
 - Automated VFC
 - No flux decay
 - 100 LMH/bar, 3 bar, \$10000/m²
- Costs per year
 - Total = \$55.3M / year
 - \$460000 / batch

Comments and Discussions

- Feasible for small scale batch and perfusion production
- Large scale batch production options: 1) Filter only a portion of total volume, 2) Use longer process (e.g. 24 hrs), 3) Wait for better filters

Considerations - Media Treatment - Way Forward?

Large scale recombinants:

HTST (or UVC)



Insensitive Media Components



Virus filtration



Sensitive Media Components



ATMPs – Gene Therapy

- Virus-mediated delivery: more (mixed) history ...
 - **NLE, small**
Adeno-associated virus
traditional down-stream methods applicable (NF, SD..)
 - **NLE, large**
Adenovirus
traditional down-stream methods applicability limited (no NF)
 - **LE, large**
Retrovirus / Lentivirus, Herpes simplex, Vaccinia
traditional down-stream methods applicability limited (no NF, SD)


→ differences that matter !



Virus Filtration as Upstream Barrier

Advanced Therapy Medicinal Products, ATMPs

- Cell-based therapies
 - Gene therapy vectors
- upstream barrier as the only option (?)



Acknowledgements

Virusure: Andy Bailey, Natascha Hodosi

Boehringer Ingelheim: Simon Haidinger

Shire: Andreas Wieser, Thomas Kreil

BioReliance: Kate Smith, Isobel Radigan, Doli Patel, Hidenori Coy Meiseki,
Yuki Shimazoe, Steven McDade

Asahi Kasei: Masayasu Takahara, Bixente Martirene, Konstantin Agolli,
Tomoko Hongo Hirasaki, Daniel Strauss



Virus Filters as Bioprocess Subject - Current Hot Topics

Facility Segregation


Why?

→ Regulatory pressure - companies are asked to segregate pre- & post-virus processes by agencies

Although...


No clear regulatory guidance! *“...the nature of the product as well as the equipment used will determine the level of segregation needed to avoid cross--contamination.”* European GMP Guidance, Annex 2

Facility Segregation




Facility Segregation

VF is often THE point of segregation.



Ballroom facilities:
 everything in one room
 → need for closed systems and
 aseptic assembly options



Drivers:

Process economics


Ease of use

Change over time and safety for multi-product facilities
(fewer multiple-use equipment)

Flexibility

Facility Segregation → Single-Use

Single-use system anatomy



(PDA Technical Report No. 66 - Application of Single-Use Systems in Pharmaceutical Manufacturing)

Virus Filters as Bioprocess Subject - Current Hot Topics

Industry Approaches to Facility Segregation for Viral Safety

Paul W. Barone
MIT Center for Biomedical Innovation

2017 CAACB Spring Workshop | April 11 - 12 | Boston, MA

Industry Approaches to Facility Segregation for Viral Safety

Paul W. Barone
MIT Center for Biomedical Innovation

2017 CAACB Spring Workshop | April 11 - 12 | Boston, MA

Questions ?

