

#### Debates and Challenges Concerning Disinfectant Validation and Methods to Circumvent Them



PDA/ISPE Australia September 20, 2019 Jim Polarine Jr. MA. Sr. Technical Service Manager

# Agenda



Disinfectant Validation Practices

## **End-User Disinfectant Validation Components**



- In vitro testing
  - Suspension testing (also called Time Kill Study)
  - Carrier Testing (also called Coupon Testing)
- In situ testing
- Environmental monitoring
  - Data trending (6-12 months, reviewed monthly)
  - Identification of organisms (mold, yeast, and bacteria)



## Disinfectant Qualification Procedure Recommendations

- USP 42 <1072> Disinfectants and Antiseptics
  - Use-dilution tests
  - Surface Challenge tests
- ASTM E2614-15 Guide for evaluation of Cleanroom Disinfectants
- ISO 14698 (parts1-3)
  - Surface evaluation, focus on cleaning
- PDA TR No. 70 on Cleaning and Disinfection (October, 2015)

#### **In Vitro Options for Testing**

- AOAC
  - Use-dilution Test Methods (955.14, 955.15, 964.02)
  - Sporicidal Activity of Disinfectants (966.04)
  - Germicidal Spray Products as Disinfectants
- ASTM
  - Time Kill Method
  - Spray Slide
  - Sanitizer method (E1153)
  - Wipe method
  - Quantitative Carrier Method (E2111-11 & E2197-11)
  - Biofilm Method (E1427)
  - Viral Testing (Suspension E1052-11)
  - Viral Testing (Carrier E1053-11)
  - Standard Guide for Evaluation of Cleanroom Disinfectants (E2614-15)
- · Variations of all of the above

Copyright © 2018 STERIS Corporation. All Rights Reserved.

STERIS STERIS

#### **More In Vitro Options**



- EN

  - 1276 (bacterial suspension test)
    1040 (bacterial suspension test)
    1650 (fungal suspension test)
    13704 (sporicidal suspension test)
    13697 (Carrier test)-Revised 2015
  - 14476 (Viral Testing)
  - 14348 (TB Testing)
  - 14885:2015
  - 16777 (Viral Hard Surface test)
- AFNOR (France)
  - NFT 72-150 Suspension
    NFT 72-190 Carrier Test
- VAH (DGHM) (Germany, Carrier & Suspension Tests)
- TGA (Australia)

## EN 13697: 2015 Revisions



- Some of the key revisions include:
  - Obligatory soiling conditions for *P. aeruginosa* changed from bovine albumin to skimmed milk
  - An evaluation for 75% mature (spiny)
    - A. brasiliensis spores prior to testing
  - Nearly all method verification acceptance criteria modified (now including a neutralizer toxicity evaluation)
  - Calculation of weighted mean for counts

#### **Spiny Spores**



Life Sciences



## EN 13697



- Being a prescriptive test method allows for consistency across European facilities
- Video





#### In Vitro Carrier Comparison

#### EN 13697

#### Inoculum



Copyright © 2018 STERIS Corporation. All Rights Reserved.

#### Test Product





#### In Vitro Carrier Comparison

#### **ASTM E 2197**

#### Inoculum



#### **Test Product**





#### USP 42 <1072> 2"x2" Coupons?



- USP 42 <1072> does not provide specific guidance on recovery methods
- Established reference methods that specify recovery methods, utilize smaller coupons
- Using larger coupons can negatively impact some recovery methods
- The volume of inoculum and test product used in prescriptive reference methods obviates the need for larger coupons
- USP <1072> was really meant to be a "wipe method" Tony Cundell

#### USP 42 <1072> 2"x2"

• Necessary?



#### **Coupon Size Debate**



- USP 42 <1072> Calls for 2" x 2" (5.08 cm x 5.08 cm) coupons-no other operatic details specified
- PDA TR # 70 Calls for 3.8 cm X 3.8 cm
- ASTM E2197-11 Calls for 1 cm disc
- EN 13697 (2015) Calls for 2 cm disc
- Some End Users 28 X 28cm and 5 X 5cm
- Larger coupons can limit possible recovery methods
- Having scientifically sound method, more important than arbitrary size

## Key Considerations for In Vitro Testing



- Use-dilution
- Temperature (hot WFI drops, use in cold room?)
- Technique
  - Suspension vs. carrier
  - Substrates
  - Neutralization/dilution
  - Subculture techniques
- Microorganisms
- Efficacy requirements

#### **Substrates for Carrier Testing**

- Traditional methods (AOAC and ASTM)
  - Stainless steel disks, penicylinders or coupons
  - · Watch glasses or glass slides
  - · Porcelain penicylinders and silk suture loops
- Cleanroom disinfectant qualifications representative materials
  - Stainless steel (416, 316, 316L, 306, 304)
  - Various plastics and elastomers
  - · Lexan curtains
  - Kydex (thermoplastic alloy used for ceilings and walls)
  - Bodycote aluminum wall
  - Epoxy-coated flooring
  - Polymeric flooring
  - MMA Flooring
  - Vinyl Flooring
  - Terrazo Flooring
  - Acyrlic and Grout
  - Saniflex
  - Paints (Epoxy and Water Based) & Sealants
  - Gaskets (EPDM, Teflon)
  - Rubber or Nitrile gloves





#### **Suspension Testing**

- Often called "Time Kill" study
- Estimates the in vitro activity of the biocide
- Often used for preliminary evaluation of several different biocides
- Not required, but useful screening tool





Copyright © 2018 STERIS Corporation. All Rights Reserved.



#### **Carrier Testing**

- Simulates practical conditions of disinfectant use and application
- Test organisms are dried on coupons made of varied substrates
- End-user required to perform carrier tests to qualify disinfectants



Copyright © 2018 STERIS Corporation. All Rights Reserved.





#### **Neutralization Methods**

- Elimination of inhibitory residual disinfectant activity
  - Chemical neutralization of the active
  - Dilution generally not effective alone (alcohols)
  - Filtration + Rinsing separating the active from the organism
- Issues
  - Antimicrobial activity of neutralizer (toxicity)
    - Thioglycollate, thiosulfate, and sodium sulfite can be toxic
  - If ineffective, contact time is inaccurate
- Validation of neutralization is required



## **Microorganism Selection**



- Environmental isolates must be considered
  - Broad spectrum
  - Most frequently occurring
  - High levels in the environment
  - Demonstrated decontamination difficulty at the facility
  - "Worst Case"
- USP (ATCC or USDA) challenge organisms may also be considered but environmental isolates are the most critical

#### **Microorganism Selection**

	Microorganism	Examples			
More Resistant	Prions	Scrapie, Creutzfeld-Jacob disease, Chronic wasting disease			
	Bacterial Spores	Bacillus, Geobacillus, Clostridium			
	Protozoal Oocysts	Cryptosporidium			
	Helminth Eggs	Ascaris, Enterobius			
	Mycobacteria	Mycobacterium tuberculosis, M. terrae, M. chelonae			
	Small, Non-Enveloped Viruses	Poliovirus, Parvoviruses, Papilloma viruses			
	Protozoal Cysts	Giardia, Acanthamoeba			
	Fungal Spores	Aspergillus, Penicillium			
	Gram negative bacteria	Pseudomonas, Providencia, Escherichia			
	Vegetative Fungi and Algae	Aspergillus, Trichophyton, Candida, Chlamydomonas			
	Vegetative Helminths and Protozoa	Ascaris, Cryptosporidium, Giardia			
	Large, non-enveloped viruses	Adenoviruses, Rotaviruses			
	Gram positive bacteria	Staphylococcus, Streptococcus, Enterococcus			
Less Resistant	Enveloped viruses	HIV, Hepatitis B virus, Herpes Simplex virus			

Bacillus cereus / sphaericus Bacillus subtilis / G. stearothermop hilus Clostridium spp.

From McDonnell, "Antisepsis, Disinfection, and Sterilization: Types, Action, and Resistance" 2007, ASM Press



#### **Debate Regarding Coupon Testing**



Pros for not testing

- Reduce testing and resources costs significantly
- Have one centralized coupon study as a reference
- BPOG and PQRI

#### Cons for not testing

 There are in fact more resistant strains of bacterial spores such as Bacillus cereus that do not conform

## **General Efficacy Recommendations** STERIS



Life Sciences

- Suspension acceptance criteria
  - 4-5 log reduction
- Carrier acceptance criteria USP 42 <1072>
  - 2 log reduction bacterial spores
  - 3 log reduction vegetative bacteria
  - PDA TR #70
    - 1-5min disinfectant and sporicide >1 log reduction
    - 90sec sanitizer >1 log reduction

#### PDA TR # 70: Table 5.2.2-1



Antimicrobial Chemical Agent	Organism Type	Suggested Contact Time	Suggested Minimum Reduction
Sanitizer	Non-spore formers	max. 90 sec	>1 Log
Disinfectant/Sporicide	Non-spore formers	1 - 5 min	>1 Log
Disinfectant/Sporicide	Mycoplasma	1 - 5 min	>1 Log
Sporicide	Mold Spores	1 - 5 min	>1 Log
Sporicide	Bacterial Spores	1 - 5 min	>1 Log



#### STERIS

Life Sciences



## In Vitro Testing Considerations Contributors to Test Failures



- Recovery issues post-drying (P. aeruginosa)
- Inoculum prep (e.g. fungal spores)
- Coupon prep (autoclaving peeling Saniflex)
- Improper dilution of Concentrate
- Inappropriate biocide for spores
- Insufficient contact time should match SOP
- US vs. EU requirements

## Testing Against Fungal Spores



- Trichophyton mentagrophytes is US EPA standard (easily killed)
- Cleanroom users test Aspergillus brasiliensis (typically the most difficult to kill mold)



#### **Case Study on Substrates**



Life Sciences

Efficacy (log reduction) of Low pH phenolic: (1:256) against test microorganisms on representative surfaces

Surface	Staphylococcu s epidermidis	Pseudomonas aeruginosa	Corynebacteri um glutamicum	Candida albicans	Aspergillus brasiliensis	Penicillium chrysogenum
Stainless Steel	6.62	>6.10 <sup>b</sup>	4.18	>4.31 <sup>b</sup>	<3.00 <sup>c</sup>	4.95
Glass	6.85	6.42	5.26	>5.80 <sup>b</sup>	2.98	5.11
Aluminum	6.35	5.69	5.14	>3.93 <sup>b</sup>	<3.00 <sup>c</sup>	3.48
Ероху	4.36	4.45	4.48	3.19	<3.00 <sup>c</sup>	<3.00 <sup>c</sup>
Enamel	>6.05 <sup>b</sup>	>5.72 <sup>b</sup>	5.45	>3.92 <sup>b</sup>	<3.00 <sup>c</sup>	2.83
Acrylic	4.53	6.06	4.49	2.92	<3.00 <sup>c</sup>	<3.0 c
Mipolam	4.36	3.87	4.29	4.37	<3.00 <sup>c</sup>	3.25
Vinyl	4.08	3.68	3.93	2.61	<3.00 <sup>c</sup>	2.1
Hardwood	5.18	>4.54 <sup>b</sup>	5.26	3.2	<3.00 <sup>c</sup>	2.59
Melamine Covered Wood	>5.38 <sup>b</sup>	>5.64 <sup>b</sup>	>5.09 <sup>b</sup>	>5.12 <sup>b</sup>	3.65	3.95
Plastic	>5.73 <sup>b</sup>	>5.32 <sup>b</sup>	>5.05 <sup>b</sup>	>4.04 <sup>b</sup>	<3.00 <sup>c</sup>	2.44
Plexiglas	>5.90 <sup>b</sup>	5.62	4.83	>4.40 <sup>b</sup>	<3.00 <sup>c</sup>	3.85
Chromium	6.55	5.95	6.63	4.08	<3.00 <sup>c</sup>	2.61

<sup>a</sup> Disinfectant Efficacy = (Log MSP<sub>(positive control)</sub> - Log MSP<sub>(test coupons</sub>)), where MSP<sub>(Positive Control)</sub> = Mean surviving population on positive control coupons; MSP<sub>(test coupon</sub>) = Mean surviving population on test coupons after disinfectant treatment; <sup>b</sup> Each of triplicate coupons showed no growth after disinfectant treatment; <sup>c</sup> Each of triplicate coupons showed TNTC growth

#### **Hard Surface Test Results**



Life Sciences



#### **Environmental Isolate Testing**





# Most Common Causes for Failures in Efficacy Testing



Life Sciences

General	<ul> <li>Testing biocide against inappropriate microbes</li> <li>Using inappropriate methods</li> <li>Inadequate planning</li> <li>Insufficient contact time</li> </ul>
Neutralization	<ul> <li>Inadequate neutralization</li> <li>Neutralizer toxicity</li> </ul>
Inoculum	<ul> <li>Poor viability of inoculum suspensions</li> <li>Fungal and bacterial spore suspensions prepared incorrectly</li> </ul>
Surfaces	<ul> <li>Porous surfaces</li> <li>Coupons not amenable to steam sterilization</li> <li>Uneven inoculation or product coverage due to curvature or surface tension</li> </ul>
Recovery	<ul> <li>Lethality after drying (e.g. <i>P. aeruginosa)</i></li> <li>Setting artificially high log reduction targets</li> <li>Final plates are not countable</li> <li>Recovery method not validated</li> </ul>

Copyright © 2018 STERIS Corporation. All Rights Reserved.



31/

#### **Neutralizers**



PDA TR # 70: Table 5.2.1-1

Antimicrobial Chemical Agent	Neutralizing Agent			
Alcohols	Dilution or Polysorbate 80			
Sodium Hypochlorite	Sodium Thiosulfate			
Quaternary Ammonium Compounds	Polysorbate 80 and Lecithin			
Phenolic Compounds	Dilution or Polysorbate 80 and Lecithin			
Hydrogen Peroxide/Peracetic Acid and Hydrogen Peroxide	Catalase			

#### **Common Chemical Neutralizers**



Life Sciences

Neutralizer	Biocide Class			
Bisulfate	Gluteraldehyde			
Catalase	Hydrogen Peroxide			
Glycine	Aldehydes			
Lecithin	Quats, Phenolics, Bis-biguanides			
Letheen	Quats			
Mg+2 or Ca+2 ions	EDTA			
Polysorbate (Tween)	Quats, Phenolics, Iodine			
Sodium Thiosulfate	Sodium Hypochlorite, Iodine			

#### **Neutralizing Broths**



Life Sciences

Ingredient	AOAC	DEB	LET	NIH	TAT	TPL
Beef extract	5.0		5.0			
Casitone				15.0		
Cystine				0.5		
Dextrose		10.0		5.5		2.5
Lecithin		7.0	0.7		5.0	0.7
Peptamin	10.0		10.0			
Polysorbate 20					43.2	
Polysorbate 80		5.0	5.0			15.0
Sodium bisulfite		2.5				
Sodium chloride	5.0		5.0	2.5		
Sodium thioglycollate		1.0		0.5		
Sodium thiosulfate		6.0				
Soytone						3.0
Tryptone		5.0			20.0	17.0
Yeast extract		2.5		5.0		

Sutton, SW et al. 2002. Validation of Microbial Recovery From Disinfectants. PDA J Pharma. Sci. Technol. 56(5):255-266.

## **Viability of Inoculum**



- Making sure the bacteria at the right phase of growth
- Making sure to isolate the fungal spores with a glass gauze fritted filter or glass wool (testing spores and not mycelia or mycelial mat)
- Checking the viability of the culture and making sure no cross contamination is present



#### Inoculum Preparation—Fungal Spores



Incubate cultures for a sufficient length of time before harvesting spores





Copyright © 2018 STERIS Corporation. All Rights Reserved.

Courtesy Dan Klein

#### **Aspergillus Spores**



Life Sciences



Copyright © 2018 STERIS Corporation. All Rights Reserved. Copyright © 2017 STERIS Corporations. All Rights Reserved.

#### **Cleanroom Fungi**



Life Sciences



Courtesy Dan Klein Copyright © 2018 STERIS Corporation. All Rights Reserved.

#### **Aspergillus Spores**



Aspergillus Brasiliensis Spores 4.49 µm 4.49 µm 4.49 µm 101 100 µm 100 µm 101 100 µm 100 µ

Copyright © 2018 STERIS Corporation. All Rights Reserved.

**Courtesy Bruce Ritts** 

#### Aspergillus brasiliensis







#### **Cladosporium** Spores



Courtesy Bruce Ritts





## Leptosphaerulina Spores





1194 1.50kV 6.1mm x3.00k SE

10.0μm

# SEM: Pseudomonas 5,000X magnification





#### **Bacillus Subtilis**





Copyright © 2018 STERIS Corporation. All Rights Reserved. Copyright © 2017 STERIS Corporations. All Rights Reserved.

#### **Surface/Coupon Issues**



- Surface type and condition can have a huge impact on efficacy
- Preparation of surfaces prior to testing
  - Autoclaving may not be acceptable for some surfaces
  - Residues must be removed
  - No rusting or pitting of surfaces
- Some surfaces pose a challenge during qualification studies:
  - Peeling after sterilization
  - Surface tension (issue on Epoxy, Vinyl, and Terrazzo)
  - Paints and Glove Materials

## **Surface Type and Condition**



- Visually smooth surfaces can be irregular
- Older or damaged surfaces can be more challenging
- Glass and stainless steel typically the least challenging



## **Residue Analysis: SEM**





Copyright © 2018 STERIS Corporation. All Rights Reserved.

**Courtesy Bruce Ritts** 

#### Surface Conditions Effect Performance





Courtesy Bruce Ritts Copyright © 2018 STERIS Corporation. All Rights Reserved.



#### **Surface Preparation**



Autoclaving may not be acceptable for some surfaces, gypsum board with paint (Saniflex)



#### **Surface Tension Issue**



Life Sciences





## **Surface Creation Issue**



Coupon creation led to unrepresentative texture



Courtesy of Erin Kruesi, STERIS Laboratories



## **Surface Degradation Issue**





Courtesy of Erin Kruesi, STERIS Laboratories



## **Recovery Method Issues**



- Contact plates (rarely used)
- Swabs
- Direct inoculation of coupons into neutralizing media
  - Requires sterile coupons
  - May include manual or automated dislodging
- Stomacher bags (Food Industry)
- Recovery method must be validated/verified
- Sonication, vortexing, and glass beads.
- Final plates must be countable to calculate log reduction



# **Testing Bacterial Spore Formers**



- Bacterial spores are significantly more challenging to inactivate than vegetative cells (worst case for a strain)
- Bacterial spores generally require sporicidal agent to inactivate
- Spore form is most likely form of *Bacillus*, *Paenibacillus*, etc. to be encountered in classified area



# **Testing Bacterial Spore Formers**

Can you consistently test against "vegetative" *Bacillus*?





# **Testing Bacterial Spore Formers in DET**

- •Bacterial spore formation is a stochastic event
- Sporulation occurs within sub-populations (i.e. the entire suspension does not sporulate at the same time)
- Sporulation kinetics is complicated, depending on a multitude of environmental factors-NOT just nutrient scarcity

#### **Paenibacillus 24 Hour Culture**



Life Sciences



Courtesy of Jamie Knutson, STERIS Laboratories

#### **Paenibacillus 18 Hour Culture**





Life Sciences

Courtesy of Jamie Knutson, STERIS Laboratories

61/

#### 18 Hour Paenibacillus Culture





Life Sciences

Courtesy of Jamie Knutson, STERIS Laboratories



# **Disinfectant Qualification Study Tips**

- AOAC methods are inappropriate for this testing (but some procedures such as inoculum prep, etc. can be of value)
- EN-13697 (2015) and ASTM E2197-11 offer valuable insight into quantitative surface testing
- PDA TR #70 (2015) is useful in determining log reductions
- Up-front proactive planning is extremely important
- Combining physical removal and chemical kill in one study is not recommended
- Consistency is crucial to a positive outcome
- Reading the product labels to understand product claims and limitations is necessary
- Incorporate expiry dating specified in internal SOPs into the study
- Using a contract lab to perform testing sounds easy but still requires time, effort, and vigilance
- Auditing the contract lab is very useful Copyright © 2018 STERIS Corporation. All Rights Reserved.

# Keys to a Successful Qualification



- Effective Antimicrobial agents
- Effective and repeatable testing protocol
- Effective sanitization procedures
- Effective change control procedures

# Requalification



- Review annually to assess risk/ whether changes have occurred
- If new bioburden appears at high levels or inherently resistant organisms
- Re-evaluate every five years to determine if any repeat testing is needed due to testing deficiencies

# **Summary Slide**



Current Industry Best Practice in Disinfectant Validation

# Acknowledgements



- Special Thanks to Bruce Ritts and Stacey Gish for SEM work.
- Special Thanks to Dan Klein and David Shields.