



# Testing of Pre-filled Syringes

*Erik Berndt*

*Zwick, Ulm*

## Erik Berndt

### Zwck Industry Manager Medical & Pharmaceutical



#### Background/ Experiences:

- Mechanical Engineer / Textile Engineer
- 1997-2003: R&D ITA, RWTH Aachen University  
(Biomaterials/Implants)
- 2003-2008: R&D Paul Hartmann AG, Heidenheim  
(Compression/Immobilization products)
- Materials, production, testing, regulatory, launching, patents

**ISO 11040-4:** Prefilled syringes — Part 4: Glass barrels for injectables and sterilized subassembled

**ISO 80369-7:** Small-bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications

**ISO 80369-20:** Small-bore connectors for liquids and gases in healthcare applications – Part 20: Common test methods

**(ISO 594-2:** Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment – Part 2: Lock Fittings)



Pre-filled syringes



- C1** Flange breaking resistance
- C2** Luer cone breaking resistance
- E** Glide force test method to evaluate syringe lubrication
- F** Needle penetration test
- G1** Needle pull-out force
- G2** Closure system liquid leakage test
- G3** Luer lock adaptor collar pull-off force
- G4** Luer lock adaptor collar torque resistance
- G5** Luer lock rigid tip cap unscrewing torque
- G6** Pull-off force of the tip cap or the needle shield



## C1 - Flange breaking resistance








## C2 - Luer cone breaking resistance







# F - Needle penetration test








# G1 - Needle pull-out force



## G3 - Luer lock adaptor collar pull-off force




## G4 - Luer lock adaptor collar torque resistance







## G5 Luer lock rigid tip cap unscrewing torque







# G6 Method 1 Pull-off force of the tip cap or the needle shield





## G6 Method 2 Pull-off force of the tip cap or the needle shield








Testing device following DIN EN ISO 7886-1



Universal testing rig for syringes and cartridges

- Recording of pictures (setups, samples before/after testing) is integrated in testXpert III.
- Option Videocapture to synchronize testing results with recorded movies of testing sequences (variety of camera systems usable)

Movie 






Synchronized measuring of the thickness of the plunger while performing a continuous/single ejection test.



Cartridge with plunger



Plunger compression measurement



## Advantages



- Compact and robust design, space saving and relatively low-priced
- Single and superimposed tensile/compression/torsion testing, synchronization of both test axes
- Test methods of the two axes can be used in any combination
- High-resolution rotation angle and travel measurement
- Operation with standard PC or laptop (no additional interface card required) and testXpert III testing software
- Optional CE-compliant safety device for use with hazardous specimens and / or higher rotational speeds



Questions?