

Haselmeier: Injection Pens for Clinical Trials

Terence O'Hagan General Manager Haselmeier, Inc



Headline



"Working together to inspire your patients"

Fast Facts

- Established in 1920
- Self Injection Systems
- Standard Platform devices
- **Customized Solutions**
- Proprietary Technology
- Pharma Packaging
- Connected Devices







Global Design & Manufacturing

- Manual & Automated Assembly
- Final Assembly & Drug Packaging
- IP Development & Ownership
- Factories
 - Germany
 - Czech Republic
 - India









A personalized injection device with...



- ... one or more preset doses in one device?
- ... flexible dose volume setting not available with fixed-dose injectors?
- ... a mechanism to prevent the setting of unintended doses?
- ... an easy, rapid and cost-effective customization of the preset doses?

Haselmeier has a new solution for you:

It takes Just one. The new





Clinical Trial Challenges

> Complex

Manufacturing

Ethical/Regulatory

Trial Protocols

Supplier Management

Scheduling

Reporting

Contracts & agreement





D-Flex Self Injection Device



- Flexibility
- Improved Patient Safety
- Clinic to Commercial Simplicity








Flexibility



- In D-Flex the **Dose Selector** defines the intended doses
- The dose regimen can be selected upon assembly by choosing the appropriate Dose Selector!
- The **same** assembly **equipment** can be used for different variants
- 3 Versions: **Max Dose** of 0.1 ml, 0.3 ml, 0.6 ml


A truly **flexible** platform!





Patient Safety

Patients cannot set the dose in between the preset fixed doses











Clinic to Commercial Simplicity





Clinic to Commercial Simplicity

The D-Flex simplifies process from clinical trial to commercial launch



for DVT

Other drugs




Benefits of D-Flex For Clinical Trials

- Simplify Dosing Studies
 - Parallel Dose Comparison
 - **Crossover Study**
 - **Dose Titration**
 - Dose escalation
- Improve Patient Compliance
- Facilitate Low Cost Phase Transfers










A flexible and powerful platform.

D-Flex adds value to



Patients



A single device instead of several devices.

Preset doses instead of variable dose device

Patient cannot set other than the preset doses.

Pharma companies



Improved time-tomarket and CAPEX

One platform for an easy customization of preset doses.

The same platform to cover clinical doseranging studies up to commercial product.₁₃



Thank you!

Terence O'Hagan

Haselmeier Inc. 126 John Street, Suite 11 Lowell, MA 01852

+1 469.386.6950

t.ohagan@haselmeier.com

Haselmeier.com

Visit us at Booth X35

