

Desorption

Georg Frinke - *Process Engineer*

Volunteer for PDA


Overview

- Theory of Desorption
- Practical aspects of Desorption
- PAT

Theory of Desorption

Residual moisture

- When leaded into process, the water content of the product was defined by water mass percentage to full mass
- The residual moisture is calculated by the water mass to dry mass => moisture levels above 100% are possible
- For long-term-storage the moisture content is of interest => moisture levels of 2...5% are common, sometimes a minimum moisture level is also required e.g. viruses or proteins



Universität Dortmund, 1976, S. 10.3)


Sorption Isotherm

- Residual moisture level depends on Temperature and surrounding (partial) vapor pressure
- Moisture equilibrium is shown in sorption isotherm
- Conditions above/below isotherm lead to desorption/adsorption
- Water molecules are directly bound to dry product

Sorption Isotherm of dried potatoe pieces

Sorption Isotherm

- Residual moisture level depends on Temperature and surrounding (partial) vapor pressure
- Convergence against the sorption isotherm



Adsorption / Desorption

- Residual moisture (detected by Karl Fischer-Method) consists of desorbable and non-desorbable water
- The integration of the desorption rate versus infinity allows – at least – theoretically the prediction of the residual moisture
- Overdried products would adsorb water at same conditions (“remoistering”) e.g. freeze dried coffee

Arrhenius is everywhere

- W_B molar Bonding Energy (fractions...several eV)
- K Sorption Constant
- R general Gas Constant
- T absolute Temperature

- $K < 0 \Rightarrow$ Desorption
- $K > 0 \Rightarrow$ Adsorption

$$\text{Sorption Rate} = K \cdot e^{\frac{-W_B}{R \cdot T}}$$

The Oetjen Idea

$$dW_{(t)} = \int_t^{\infty} \left(K \cdot e^{\frac{-W_B}{R \cdot T}} dt \right)$$

- Prediction of Residual Moisture is possible when Equilibrium moisture of pressure and temperature is known

Desorption

- The Desorption Rate is linear depending on the Residual Moisture
- An empiric model (RM vs. DR) can be defined during cycle development for more effective cycle control with reduced development effort

$$DR \sim \dot{R}M_{KF}(t)$$

Practical Aspects of Desorption

Desorption Phase

- Process Temperature profile equilibrates during desorption phase
- Shelf Temperature at Desorption should be near the max. storage temperature, not above
- Process vacuum is of minor interest, higher pressure level allows good moisture equilization, but inhibits desorption

Optimization of Desorption Phase

- A general hint to estimate duration can not be given
- The use of PAT Tools allows the safe control of the right step time
- Frequently accomplished Desorption Rate Measurements enable safe detection of end point (historical Freeze Dryers can use the right calculated pressure rise)

$$DR = \frac{\Delta p / [\mu\text{bar}]}{\Delta t / [\text{s}]} \cdot \frac{V_c / [\text{m}^3]}{v_s / [\text{K}]} \cdot \frac{1}{m_{Tr} / [\text{kg}]} \cdot 0,7803 = \left[\frac{1}{\text{h}} \right]$$

Further Optimization of desorption phase

- As shown before the desorption rate is the first derivation of the desorbable Water
- Desorption time depends on concentration and crystal structure during freezing, impact of layer height is marginal
- Remark that the pharmaceutical standard for moisture measurement detects all Water content, including non-desorbable Water



PAT


PAT @ Desorption:

Process Control for Desorption

- Shelf Temperature
- Chamber Pressure

Process Feedback for Desorption

- Desorption Rate
 - TDLAS
 - Analysis of Pressure Rise Measurement



PAT @ Desorption:

Further Process Feedback for Desorption

- Vapor Concentration
 - Cold Plasma
 - NIR Gas Analysis
 - Dew Point Detection (minor accuracy)
 - Inert Gas Flow Monitoring (minor accuracy)
 - Comparative Pressure Measurement (minor accuracy)

The material shown in this presentation does not represent the official opinion of Bayer Pharma or any linked subsidiary company

Bayer Pharma or any linked subsidiary company are not responsible for any content of this presentation

Some Questions of mine:

- ✓ Relation between residual Moisture and Desorption Rate?
- ✓ What is most important, when achieving a minimum and maximum moisture level
- ✓ PAT Tools for desorption?

Thank you for your attention!

Now its correctly dried

Questions?