

Martin Christ GmbH
USB flash drive
Batch data 21112303.00



Recording file: 21112303.00
Recording started by:
Begin recording: 2021-11-23 14:59:40
End recording: 2021-11-26 8:22:40
Recording time (hhh:mm:ss): 65:22:58
Unit type : Epsilon 2-6D (LSCplus)

- Channel01: Vac.Pirani
- Channel02: SetShelf
- Channel03: Ice cond.1
- Channel04: Vac.cap.
- Channel05: WTMplus #1
- Channel06: WTMplus #2
- Channel07: WTMplus #3
- Channel08: WTMplus #4
- Channel09: WTMplus #5
- Channel10: WTMplus #6
- Channel11: WTMplus #7
- Channel12: WTMplus #8

Product no.: divers
Product name:
Batch no.:
Batch name:
Batch size: 20R
Batch volume per unit: 5-10
Note:
PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17

LPCplus

Document issued by:

Signature:


Page: 1 / 209

Martin Christ GmbH
 USB flash drive
 Process graphics 21112303.00



Vac.Pirani SetShelf Ice cond.1 Vac.cap. WTMplus #1 WTMplus #2 WTMplus #3 WTMplus #4 WTMplus #5 WTMplus #6 WTMplus #7 WTMplus #8

Intermediate Valve 1



Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	2 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



Date (yyyymmdd)	Time (hh:mm:ss)	Duration (hhh:mm:ss)	Event												Note
2021-11-23	14:59:41	000:00:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	14:59:41	000:00:01	mBar	°C	°C	mBar	°C	°C	°C	°C	°C	°C	°C	°C	
	14:59:41	000:00:01	1000	20,0	-67,4	1,18	---	---	-0,6	---	---	1,2	---	---	
	15:00:41	000:01:01	1000	20,0	-59,8	1,18	---	---	-0,8	---	---	1,0	---	---	
	15:01:42	000:02:02	1000	15,0	-53,4	1,18	---	---	-0,9	---	---	0,8	---	---	
	15:02:43	000:03:03	1000	14,0	-47,9	1,18	---	---	-1,0	---	---	0,6	---	---	
	15:03:41	000:04:01	1000	13,1	-45,8	1,18	---	---	-1,0	---	---	0,4	---	---	
	15:04:42	000:05:02	1000	12,0	-82,0	1,18	---	---	-1,0	---	---	0,4	---	---	
	15:05:43	000:06:03	1000	11,0	-82,4	1,18	---	---	-1,0	---	---	0,4	---	---	
	15:06:43	000:07:03	1000	10,0	-82,8	1,18	---	---	-0,7	---	---	0,5	---	---	
	15:07:40	000:08:00	1000	9,1	-81,1	1,18	---	---	-0,2	---	---	0,9	---	---	
	15:08:40	000:09:00	1000	8,1	-82,9	1,18	---	---	0,4	2,6	---	1,5	---	---	
	15:09:42	000:10:02	1000	7,0	-84,3	1,18	---	---	1,0	2,7	---	2,2	---	---	
	15:10:42	000:11:02	1000	6,0	-85,1	1,18	---	---	1,6	2,9	---	2,8	---	---	
	15:11:42	000:12:02	1000	5,0	-85,7	1,18	---	---	2,1	3,1	---	3,5	---	---	
	15:12:40	000:13:00	1000	4,1	-86,0	1,18	---	---	2,5	3,3	---	4,0	---	---	
	15:13:41	000:14:01	1000	3,0	-86,3	1,18	---	---	2,8	3,4	---	4,4	---	---	
	15:14:42	000:15:02	1000	2,0	-86,6	1,18	---	---	3,1	3,6	---	4,6	---	---	
	15:15:44	000:16:04	1000	1,0	-86,7	1,18	---	---	3,2	3,6	---	4,7	---	---	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	3 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



15:16:40	000:17:00	1000	0,1	-86,8	1,18	---	3,5	3,3	3,6	---	4,7	---	---
15:17:41	000:18:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:17:41	000:18:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:17:41	000:18:01	1000	-1,0	-87,0	1,18	---	3,4	3,3	3,5	---	4,6	---	---
15:18:43	000:19:03	1000	-2,0	-87,2	1,18	---	3,4	3,1	3,3	---	4,5	---	---
15:19:40	000:20:00	1000	-2,9	-87,3	1,18	---	3,2	3,0	3,2	---	4,3	---	---
15:20:42	000:21:02	1000	-4,0	-87,5	1,18	---	3,0	2,7	2,9	---	4,0	---	---
15:21:44	000:22:04	1000	-5,0	-87,6	1,18	---	2,8	2,4	2,7	---	3,7	---	---
15:22:40	000:23:00	1000	-5,0	-87,6	1,18	---	2,6	2,0	2,5	---	3,4	---	---
15:23:41	000:24:01	1000	-5,0	-87,6	1,18	---	2,3	1,6	2,4	---	3,0	---	---
15:24:42	000:25:02	1000	-5,0	-87,7	1,18	---	2,0	1,2	---	---	2,7	---	---
15:25:40	000:26:00	1000	-5,0	-87,8	1,18	---	1,9	0,8	---	---	2,3	---	---
15:26:42	000:27:02	1000	-5,0	-87,8	1,18	---	1,1	0,4	---	---	1,9	---	---
15:27:40	000:28:00	1000	-5,0	-87,8	1,18	---	1,0	0,1	---	---	1,5	---	---
15:28:40	000:29:00	1000	-5,0	-88,0	1,18	---	0,9	-0,3	---	---	1,1	---	---
15:29:42	000:30:02	1000	-5,0	-88,1	1,18	---	0,7	-0,6	---	---	0,8	---	---
15:30:44	000:31:04	1000	-5,0	-88,1	1,18	---	---	-0,9	---	---	0,4	---	---
15:31:41	000:32:01	1000	-5,0	-88,2	1,18	---	-0,2	-1,2	---	---	0,1	---	---
15:32:44	000:33:04	1000	-5,0	-88,3	1,18	---	-0,2	-1,4	---	---	-0,1	---	---
15:33:41	000:34:01	1000	-5,0	-88,5	1,18	---	-0,2	-1,7	---	---	-0,4	---	---
15:34:43	000:35:03	1000	-5,0	-88,6	1,18	---	-0,4	-1,9	---	---	-0,6	---	---

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 4 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



15:35:43	000:36:03	1000	-5,0	-88,8	1,18	---	-0,5	-2,0	-0,8	---	-0,8	---	---
15:36:42	000:37:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:36:42	000:37:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:36:42	000:37:02	1000	-5,0	-88,9	1,18	---	---	-2,2	-0,9	---	-1,0	---	---
15:37:42	000:38:02	1000	-5,0	-88,9	1,18	---	---	-2,4	-1,0	---	-1,3	---	---
15:38:44	000:39:04	97,3	-5,0	-89,6	1,18	---	-1,2	-2,6	-1,1	---	-1,7	---	---
15:39:41	000:40:01	13,1	-5,0	-90,3	1,18	---	-1,2	-2,7	-1,2	---	-2,0	-1,0	---
15:40:42	000:41:02	1000	-5,1	-90,6	1,18	---	-1,4	-2,9	-1,4	---	-2,0	---	---
15:41:40	000:42:00	1000	-5,6	-89,3	1,18	-0,4	-1,4	-1,9	-1,5	---	-1,7	---	---
15:42:40	000:43:00	1000	-6,1	-89,5	1,18	-0,5	---	-0,8	-1,7	---	-1,4	---	---
15:43:40	000:44:00	1000	-6,6	-89,4	1,18	-0,5	---	-0,4	-1,8	---	-1,5	---	---
15:44:40	000:45:00	1000	-7,1	-89,3	1,18	-0,7	---	-0,4	-1,8	---	-1,6	---	---
15:45:42	000:46:02	1000	-7,6	-89,2	1,18	-0,8	---	-0,4	-1,9	---	-1,8	---	---
15:46:43	000:47:03	1000	-8,1	-89,1	1,18	-0,9	---	-0,4	---	---	-1,9	---	---
15:47:41	000:48:01	1000	-8,6	-88,7	1,18	-1,1	---	-0,4	---	---	-2,1	---	---
15:48:42	000:49:02	1000	-9,1	-89,2	1,18	-1,2	---	-0,4	---	---	-2,3	---	---
15:49:42	000:50:02	1000	-9,6	-88,1	1,18	-1,4	---	-0,4	---	---	-2,5	---	---
15:50:43	000:51:03	1000	-10,1	-88,6	1,18	-1,6	---	-0,4	-0,4	---	-2,7	---	---
15:51:44	000:52:04	1000	-10,6	-88,5	1,18	-1,7	---	-0,4	-0,4	---	-3,0	---	---
15:52:41	000:53:01	1000	-11,1	-88,5	1,18	-1,8	---	-0,4	-0,4	---	-3,2	---	---
15:53:41	000:54:01	1000	-11,6	-88,2	1,18	-2,0	---	-0,5	-0,4	---	-3,5	---	---

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 5 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



15:54:44	000:55:04	1000	-12,1	-88,8	1,18	-2,2	---	-0,5	-0,4	---	-3,8	---	---
15:55:44	000:56:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:55:44	000:56:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:55:44	000:56:04	1000	-12,6	-88,9	1,18	-2,4	---	-0,5	-0,4	---	-4,1	---	---
15:56:44	000:57:04	1000	-13,2	-88,6	1,18	-2,5	---	-0,5	-0,4	---	-4,4	---	---
15:57:40	000:58:00	1000	-13,6	-88,5	1,18	-2,7	---	-0,5	-0,5	---	-4,7	---	---
15:58:43	000:59:03	1000	-14,1	-89,3	1,18	-3,0	---	-0,6	-0,5	---	-5,0	---	---
15:59:43	001:00:03	1000	-14,6	-88,9	1,18	-3,2	---	-0,6	-0,5	---	-5,4	---	---
16:00:44	001:01:04	1000	-15,2	-88,5	1,18	-3,4	---	-0,6	-0,6	---	-5,7	---	---
16:01:41	001:02:01	1000	-15,6	-89,3	1,18	-3,6	---	-0,7	-0,6	---	-6,1	---	---
16:02:43	001:03:03	1000	-16,1	-89,3	1,18	-3,7	---	-0,7	-0,6	---	-6,4	---	---
16:03:43	001:04:03	1000	-16,6	-89,2	1,18	---	---	-0,8	-0,7	---	-6,8	---	---
16:04:41	001:05:01	1000	-17,1	-89,1	1,18	-4,3	---	-0,8	-0,7	---	-7,2	---	---
16:05:41	001:06:01	1000	-17,6	-89,5	1,18	-4,4	---	-0,9	-0,8	---	-7,5	---	---
16:06:42	001:07:02	1000	-18,1	-89,5	1,18	-4,5	---	-0,9	-0,8	---	-7,9	---	---
16:07:43	001:08:03	1000	-18,6	-89,8	1,18	-4,7	---	-1,0	-0,9	---	-8,3	---	---
16:08:44	001:09:04	1000	-19,2	-89,7	1,18	-4,9	---	-1,1	-1,0	---	-8,7	---	---
16:09:41	001:10:01	1000	-19,6	-90,0	1,18	-5,0	---	-1,2	-1,1	---	-9,1	---	---
16:10:43	001:11:03	1000	-20,1	-89,9	1,18	-5,2	---	-1,3	-1,2	---	-9,5	---	---
16:11:44	001:12:04	1000	-20,7	-90,0	1,18	-5,5	---	-1,4	-1,2	---	-9,8	---	---
16:12:44	001:13:04	1000	-21,2	-89,9	1,18	-5,7	---	-1,5	-1,3	---	-10,2	---	---

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	6 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	16:13:40	001:14:00	1000	-21,6	-89,9	1,18	-5,9	---	-1,6	-1,5	---	-10,6	---	---	
	16:14:40	001:15:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	16:14:40	001:15:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	16:14:40	001:15:00	1000	-22,1	-90,0	1,18	-6,1	---	-1,7	-1,6	---	-11,0	---	---	
	16:15:40	001:16:00	1000	-22,6	-89,7	1,18	-6,3	---	-1,9	-1,7	---	-11,4	---	---	
	16:16:41	001:17:01	1000	-23,1	-90,0	1,18	-6,6	---	-2,0	-1,8	---	-11,8	---	---	
	16:17:41	001:18:01	1000	-23,6	-90,0	1,18	-6,8	---	-2,2	-2,0	---	-12,2	---	---	
	16:18:41	001:19:01	1000	-24,1	-90,3	1,18	-7,1	---	-2,4	-2,1	---	-12,5	---	---	
	16:19:42	001:20:02	1000	-24,6	-90,3	1,18	-7,3	---	-2,6	-2,3	---	-12,9	---	---	
	16:20:43	001:21:03	1000	-25,1	-89,7	1,18	---	---	-2,9	-2,4	---	-13,3	---	---	
	16:21:43	001:22:03	1000	-25,6	-85,6	1,18	---	---	-3,2	-2,6	---	-13,6	---	---	
	16:22:43	001:23:03	1000	-26,1	-80,0	1,18	---	---	-3,7	-2,7	---	-14,0	---	---	
	16:23:44	001:24:04	1000	-26,7	-83,3	1,18	---	---	-4,4	-2,9	---	-14,4	---	---	
	16:24:41	001:25:01	1000	-27,1	-84,0	1,18	---	---	-5,4	-3,1	---	-14,7	---	---	
	16:25:41	001:26:01	1000	-27,6	-80,0	1,18	---	---	-6,6	-3,3	---	-15,1	---	---	
	16:26:41	001:27:01	1000	-28,1	-83,9	1,18	---	---	-7,9	-3,4	---	-15,5	---	---	
	16:27:42	001:28:02	1000	-28,6	-84,9	1,18	---	---	-9,3	-3,6	---	-15,8	---	---	
	16:28:42	001:29:02	1000	-29,1	-80,0	1,18	---	---	-10,7	-3,8	---	-16,2	---	---	
	16:29:42	001:30:02	1000	-29,6	-83,6	1,18	---	---	-12,1	-4,1	---	-16,6	---	---	
	16:30:40	001:31:00	1000	-30,1	-81,5	1,18	---	---	-13,3	-4,3	---	-17,0	---	---	
	16:31:43	001:32:03	1000	-30,6	-78,6	1,18	---	---	-14,5	-4,5	---	-17,4	---	---	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 7 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



16:32:44	001:33:04	1000	-31,2	-83,7	1,18	---	---	-15,6	-4,7	---	-17,8	---	---
16:33:40	001:34:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
16:33:40	001:34:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
16:33:40	001:34:00	1000	-31,6	-81,2	1,18	---	---	-16,7	-5,0	---	-18,1	---	---
16:34:42	001:35:02	1000	-32,1	-81,2	1,18	---	---	-17,6	-5,2	---	-18,5	---	---
16:35:43	001:36:03	1000	-32,6	-77,3	1,18	---	---	-18,6	-5,5	---	-18,9	-8,4	---
16:36:40	001:37:00	1000	-33,1	-82,9	1,18	---	---	-19,4	-5,7	---	-19,3	-8,6	---
16:37:40	001:38:00	1000	-33,6	-80,7	1,18	---	---	-20,2	-6,0	---	-19,7	-8,9	---
16:38:42	001:39:02	1000	-34,1	-78,4	1,18	---	---	-20,9	-6,3	---	-20,1	-9,4	---
16:39:42	001:40:02	1000	-34,6	-78,4	1,18	---	---	-21,6	-6,5	---	-20,5	-10,0	---
16:40:42	001:41:02	1000	-35,1	-75,3	1,18	---	---	-22,2	-6,8	---	-20,9	-10,6	---
16:41:43	001:42:03	1000	-35,6	-75,9	1,18	---	---	-22,9	-7,1	---	-21,3	-11,3	---
16:42:40	001:43:00	1000	-36,1	-78,8	1,18	---	---	-23,5	-7,4	---	-21,7	-12,0	---
16:43:42	001:44:02	1000	-36,6	-79,3	1,18	---	-8,2	-24,1	-7,7	---	-22,1	-12,7	---
16:44:42	001:45:02	1000	-37,1	-78,8	1,18	---	-8,2	-24,6	-8,0	---	-22,5	-13,5	---
16:45:43	001:46:03	1000	-37,6	-76,1	1,18	---	-8,4	-25,1	-8,2	---	-22,8	-14,4	---
16:46:41	001:47:01	1000	-38,1	-74,4	1,18	---	---	-25,6	-8,5	---	-23,3	-15,3	---
16:47:42	001:48:02	1000	-38,6	-71,1	1,18	---	---	-26,0	-8,8	---	-23,6	-16,4	---
16:48:42	001:49:02	1000	-39,1	-69,3	1,18	---	---	-26,5	-9,0	---	-24,0	-17,5	---
16:49:43	001:50:03	1000	-39,6	-63,2	1,18	---	---	-27,0	---	---	-24,4	-18,7	---
16:50:44	001:51:04	1000	-40,0	-56,5	1,18	---	---	-27,5	---	---	-24,8	-19,9	-21,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature: Page: 8 / 209
Document issued by:

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



16:51:40	001:52:00	1000	-40,0	-71,2	1,18	---	---	-28,0	---	---	-25,2	-21,1	-21,7
16:52:40	001:53:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
16:52:40	001:53:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
16:52:40	001:53:00	1000	-40,0	-72,1	1,18	---	---	-28,6	---	---	-25,7	-22,3	-22,1
16:53:40	001:54:00	1000	-40,0	-74,8	1,18	---	---	-29,1	---	---	-26,1	-23,4	-23,0
16:54:41	001:55:01	1000	-40,0	-74,5	1,18	---	---	-29,6	---	---	-26,6	-24,4	-23,7
16:55:41	001:56:01	1000	-40,0	-73,6	1,18	---	---	-30,0	---	---	-27,0	-25,2	-24,4
16:56:41	001:57:01	1000	-40,0	-71,3	1,18	---	---	-30,4	---	---	-27,3	-26,1	-25,1
16:57:42	001:58:02	1000	-40,0	-74,9	1,18	---	-21,8	-30,8	---	-22,9	-27,7	-26,8	-25,8
16:58:42	001:59:02	1000	-40,0	-79,0	1,18	---	-22,3	-31,2	---	-23,2	-28,0	-27,5	-26,5
16:59:44	002:00:04	1000	-40,0	-78,1	1,18	---	-22,7	-31,6	---	-23,3	-28,4	-28,3	-27,1
17:00:46	002:01:06	1000	-40,0	-78,4	1,18	---	-23,6	-32,0	---	-26,0	-28,7	-28,9	-27,8
17:01:42	002:02:02	1000	-40,0	-79,3	1,18	---	-24,7	-32,4	---	-26,4	-29,0	-29,5	-28,4
17:02:43	002:03:03	1000	-40,0	-81,8	1,18	-9,7	-26,0	-32,7	---	-26,8	-29,3	-30,0	-28,9
17:03:45	002:04:05	1000	-40,0	-82,8	1,18	-10,1	-26,9	-33,1	---	-27,8	-29,6	-30,5	-29,4
17:04:40	002:05:00	1000	-40,0	-79,6	1,18	-10,7	-27,7	-33,4	---	-28,6	-29,9	-31,0	-29,9
17:05:42	002:06:02	1000	-40,0	-85,4	1,18	-11,5	-28,5	-33,7	---	-29,4	-30,1	-31,4	-30,3
17:06:42	002:07:02	1000	-40,0	-83,8	1,18	-12,4	-29,1	-34,0	---	-30,1	-30,4	-31,8	-30,8
17:07:44	002:08:04	1000	-40,0	-79,3	1,18	-13,3	-29,7	-34,3	---	-30,8	-30,6	-32,2	-31,3
17:08:44	002:09:04	1000	-40,0	-85,6	1,18	-14,3	-30,3	-34,5	---	-31,3	-30,8	-32,5	-31,7
17:09:44	002:10:04	1000	-40,0	-83,5	1,18	-15,4	-30,9	-34,8	---	-31,9	-31,0	-32,8	-32,0

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 9 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



17:10:44	002:11:04	1000	-40,0	-81,9	1,18	-16,5	-31,4	-35,0	-31,3	-32,3	-31,2	-33,1	-32,4	
17:11:42	002:12:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
17:11:42	002:12:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
17:11:42	002:12:02	1000	-40,0	-85,0	1,18	-17,2	-31,8	-35,3	-31,3	-32,8	-31,4	-33,4	-32,7	
17:12:43	002:13:03	1000	-40,0	-84,8	1,18	---	-32,2	-35,5	-31,4	-33,0	-31,6	-33,6	-32,9	
17:13:41	002:14:01	1000	-40,0	-83,1	1,18	---	-32,6	-35,7	-31,9	-34,0	-31,8	-33,9	-33,2	
17:14:42	002:15:02	1000	-40,0	-86,3	1,18	---	-33,0	-35,8	-32,5	-34,2	-31,9	-34,1	-33,4	
17:15:43	002:16:03	1000	-40,0	-86,8	1,18	---	-33,3	-36,0	-32,8	---	-32,1	-34,3	-33,7	
17:16:44	002:17:04	1000	-40,0	-82,1	1,18	---	-33,6	-36,2	-33,1	-35,0	-32,2	-34,5	-33,9	
17:17:40	002:18:00	1000	-40,0	-83,2	1,18	---	-33,9	-36,3	-33,4	-35,1	-32,4	-34,7	-34,1	
17:18:43	002:19:03	1000	-40,0	-87,5	1,18	---	-34,2	-36,5	-33,6	-35,2	-32,5	-34,8	-34,3	
17:19:40	002:20:00	1000	-40,0	-86,3	1,18	---	-34,4	-36,6	-33,9	---	-32,6	-35,0	-34,4	
17:20:41	002:21:01	1000	-40,0	-81,4	1,18	---	-34,7	-36,7	-34,2	-35,9	-32,7	-35,2	-34,6	
17:21:41	002:22:01	1000	-40,0	-87,5	1,18	---	-34,9	-36,8	-34,5	-36,0	-32,9	-35,3	-34,7	
17:22:42	002:23:02	1000	-40,0	-85,9	1,18	---	-35,1	-37,0	-34,7	-36,1	-33,0	-35,4	-34,9	
17:23:43	002:24:03	1000	-40,0	-80,3	1,18	---	-35,3	-37,1	-34,9	---	-33,1	-35,5	-35,0	
17:24:45	002:25:05	1000	-40,0	-86,8	1,18	---	-35,4	-37,1	-35,1	---	-33,2	-35,6	-35,2	
17:25:40	002:26:00	1000	-40,0	-88,4	1,18	---	-35,6	-37,2	-35,3	---	-33,3	-35,7	-35,3	
17:26:42	002:27:02	1000	-40,0	-81,8	1,18	---	-35,8	-37,3	-35,5	-37,0	-33,4	-35,8	-35,4	
17:27:42	002:28:02	1000	-40,0	-85,2	1,18	---	-35,9	-37,4	-35,8	-37,0	-33,4	-35,9	-35,5	
17:28:43	002:29:03	1000	-40,0	-88,1	1,18	---	-36,0	-37,5	-36,0	-37,1	-33,5	-36,0	-35,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 10 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



17:29:44	002:30:04	1000	-40,0	-85,6	1,18	---	-36,2	-37,6	-36,1	-37,2	-33,6	-36,1	-35,7
17:30:40	002:31:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
17:30:40	002:31:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
17:30:40	002:31:00	1000	-40,0	-84,7	1,18	---	-36,3	-37,6	-36,3	-37,3	-33,7	-36,2	-35,8
17:31:40	002:32:00	1000	-40,0	-86,6	1,18	---	-36,4	-37,7	-36,5	-37,4	-33,7	-36,2	-35,8
17:32:42	002:33:02	1000	-40,0	-87,0	1,18	---	-36,5	-37,7	---	-37,4	-33,8	-36,3	-35,9
17:33:43	002:34:03	1000	-40,0	-88,5	1,18	---	-36,6	-37,8	-36,9	-37,5	-33,9	-36,3	-36,0
17:34:45	002:35:05	1000	-40,0	-82,3	1,18	-33,1	-36,7	-37,9	-36,9	-37,7	-33,9	-36,4	-36,1
17:35:40	002:36:00	1000	-40,0	-84,6	1,18	-33,2	-36,8	-37,9	-37,0	-37,8	-34,0	-36,4	-36,1
17:36:41	002:37:01	1000	-40,0	-89,1	1,18	-33,3	-36,8	-37,9	-37,1	-37,8	-34,0	-36,5	-36,2
17:37:41	002:38:01	1000	-40,0	-88,4	1,18	-33,5	-36,9	-38,0	-37,2	-37,8	-34,1	-36,5	-36,2
17:38:42	002:39:02	1000	-40,0	-85,0	1,18	-33,6	-37,0	-38,0	-37,3	-37,9	-34,1	-36,6	-36,3
17:39:44	002:40:04	1000	-40,0	-85,2	1,18	-33,8	-37,1	-38,0	-37,3	---	-34,2	-36,6	-36,3
17:40:44	002:41:04	1000	-40,0	-88,5	1,18	-33,9	-37,1	-38,1	-37,4	-38,0	-34,2	-36,6	-36,4
17:41:40	002:42:00	1000	-40,0	-87,6	1,18	-34,0	-37,2	-38,1	-37,5	-38,0	-34,3	-36,7	-36,4
17:42:43	002:43:03	1000	-40,0	-83,2	1,18	-34,1	-37,2	-38,1	-37,6	---	-34,3	-36,7	-36,5
17:43:43	002:44:03	1000	-40,0	-86,4	1,18	-34,2	-37,3	-38,2	-37,6	-38,1	-34,3	-36,7	-36,5
17:44:44	002:45:04	1000	-40,0	-87,2	1,18	-34,3	-37,3	-38,2	-37,6	-38,2	-34,4	-36,8	-36,5
17:45:44	002:46:04	1000	-40,0	-88,8	1,18	-34,4	-37,4	-38,2	-37,7	-38,2	-34,4	-36,8	-36,6
17:46:42	002:47:02	1000	-40,0	-82,7	1,18	-34,5	-37,4	-38,3	-37,8	-38,2	-34,4	-36,8	-36,6
17:47:43	002:48:03	1000	-40,0	-85,6	1,18	-34,6	-37,4	-38,3	-37,8	-38,3	-34,5	-36,9	-36,6

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 11 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	17:48:44	002:49:04	1000	-40,0	-88,4	1,18	-34,7	-37,5	-38,3	-37,8	-38,3	-34,5	-36,9	-36,6	
	17:49:41	002:50:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	17:49:41	002:50:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	17:49:41	002:50:01	1000	-40,0	-88,2	1,18	-34,7	-37,5	-38,3	-37,9	---	-34,6	-36,9	-36,7	
	17:50:43	002:51:03	1000	-40,0	-84,6	1,18	-34,8	-37,5	-38,3	-37,9	---	-34,6	-36,9	-36,7	
	17:51:45	002:52:05	1000	-40,0	-86,1	1,18	-34,9	-37,5	-38,3	-38,0	-38,4	-34,6	-36,9	-36,8	
	17:52:41	002:53:01	1000	-40,0	-88,8	1,18	-34,9	-37,6	-38,4	-38,0	-38,4	-34,6	-37,0	-36,8	
	17:53:42	002:54:02	1000	-40,0	-88,9	1,18	-35,0	-37,6	-38,4	-38,0	-38,4	-34,7	-37,0	-36,8	
	17:54:42	002:55:02	1000	-40,0	-87,4	1,18	-35,0	-37,6	-38,4	-38,1	-38,4	-34,7	-37,0	-36,8	
	17:55:43	002:56:03	1000	-40,0	-83,0	1,18	-35,1	-37,7	-38,4	-38,1	-38,4	-34,7	-37,0	-36,8	
	17:56:45	002:57:05	1000	-40,0	-86,0	1,18	-35,1	-37,7	-38,4	-38,1	-38,4	-34,7	-37,0	-36,8	
	17:57:40	002:58:00	1000	-40,0	-88,5	1,18	-35,1	-37,7	-38,4	-38,1	-38,4	-34,8	-37,0	-36,8	
	17:58:41	002:59:01	1000	-40,0	-87,7	1,18	-35,2	-37,7	-38,4	-38,1	-38,5	-34,8	-37,0	-36,9	
	17:59:43	003:00:03	1000	-40,0	-83,2	1,18	-35,2	-37,7	-38,4	-38,1	-38,5	-34,8	-37,0	-36,9	
	18:00:40	003:01:00	1000	-40,0	-88,4	1,18	-35,2	-37,7	-38,4	-38,1	-38,5	-34,8	-37,0	-36,9	
	18:01:41	003:02:01	1000	-40,0	-90,0	1,18	-35,3	-37,7	-38,4	-38,1	-38,5	-34,8	-37,1	-36,9	
	18:02:42	003:03:02	1000	-40,0	-87,1	1,18	-35,3	-37,8	-38,5	-38,2	-38,5	-34,8	-37,1	-36,9	
	18:03:40	003:04:00	1000	-40,0	-81,1	1,18	-35,3	-37,8	-38,5	-38,2	-38,5	-34,9	-37,1	-36,9	
	18:04:42	003:05:02	1000	-40,0	-86,5	1,18	-35,4	-37,8	-38,5	-38,2	-38,5	-34,9	-37,1	-36,9	
	18:05:44	003:06:04	1000	-40,0	-87,0	1,18	-35,4	-37,8	-38,5	-38,2	-38,5	-34,9	-37,1	-37,0	
	18:06:40	003:07:00	1000	-40,0	-90,6	1,18	-35,4	-37,8	-38,5	-38,2	-38,5	-34,9	-37,1	-37,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 12 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



18:07:40	003:08:00	1000	-40,0	-82,8	1,18	-35,4	-37,8	-38,5	-38,2	-38,5	-34,9	-37,1	-37,0
18:08:40	003:09:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:08:40	003:09:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:08:40	003:09:00	1000	-40,0	-86,9	1,18	-35,5	-37,8	-38,5	-38,2	-38,5	-35,0	-37,1	-37,0
18:09:42	003:10:02	1000	-40,0	-89,2	1,18	-35,5	-37,8	-38,5	-38,2	-38,5	-35,0	-37,1	-37,0
18:10:43	003:11:03	1000	-40,0	-87,1	1,18	-35,5	-37,8	-38,5	-38,2	-38,6	-35,0	-37,1	-37,0
18:11:45	003:12:05	1000	-40,0	-88,7	1,18	-35,5	-37,8	-38,5	-38,2	-38,6	-35,0	-37,2	-37,0
18:12:41	003:13:01	1000	-40,0	-82,2	1,18	-35,5	-37,8	-38,5	-38,3	-38,6	-35,0	-37,2	-37,0
18:13:41	003:14:01	1000	-40,0	-88,2	1,18	-35,5	-37,8	-38,5	-38,3	-38,6	-35,0	-37,1	-37,0
18:14:43	003:15:03	1000	-40,0	-89,4	1,18	-35,6	-37,9	-38,5	-38,3	-38,6	-35,0	-37,2	-37,0
18:15:43	003:16:03	1000	-40,0	-86,5	1,18	-35,6	-37,9	-38,5	-38,3	-38,6	-35,1	-37,2	-37,1
18:16:44	003:17:04	1000	-40,0	-83,4	1,18	-35,6	-37,9	-38,6	-38,3	-38,6	-35,1	-37,2	-37,1
18:17:40	003:18:00	1000	-40,0	-84,9	1,18	-35,6	-37,9	-38,5	-38,3	-38,6	-35,1	-37,2	-37,1
18:18:41	003:19:01	1000	-40,0	-88,9	1,18	-35,6	-37,9	-38,6	-38,3	-38,6	-35,1	-37,2	-37,1
18:19:42	003:20:02	1000	-40,0	-88,8	1,18	-35,7	-37,9	-38,6	-38,3	-38,6	-35,1	-37,2	-37,1
18:20:42	003:21:02	1000	-40,0	-81,8	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,1	-37,2	-37,1
18:21:42	003:22:02	1000	-40,0	-88,5	1,18	-35,7	-37,9	-38,6	-38,3	---	-35,1	-37,2	-37,1
18:22:40	003:23:00	1000	-40,0	-90,1	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,1	-37,2	-37,1
18:23:40	003:24:00	1000	-40,0	-87,6	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,2	-37,3	-37,1
18:24:41	003:25:01	1000	-40,0	-89,3	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,2	-37,2	-37,1
18:25:40	003:26:00	1000	-40,0	-84,3	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,2	-37,3	-37,1

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 13 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



18:26:42	003:27:02	1000	-40,0	-88,9	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,2	-37,2	-37,1
18:27:40	003:28:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:27:40	003:28:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:27:40	003:28:00	1000	-40,0	-90,6	1,18	-35,7	-37,9	-38,6	-38,3	-38,7	-35,2	-37,2	-37,1
18:28:48	003:29:08	1000	-40,0	-82,6	1,18	-35,8	-37,9	-38,6	-38,3	-38,7	-35,2	-37,2	-37,1
18:29:44	003:30:04	1000	-40,0	-87,6	1,18	-35,8	-37,9	-38,6	-38,3	-38,7	-35,2	-37,2	-37,1
18:30:44	003:31:04	1000	-40,0	-89,2	1,18	-35,8	-37,9	-38,5	-38,3	-38,7	-35,2	-37,3	-37,1
18:31:40	003:32:00	1000	-40,0	-87,2	1,18	-35,8	-38,0	-38,6	-38,3	-38,7	-35,2	-37,3	-37,2
18:32:42	003:33:02	1000	-40,0	-89,1	1,18	-35,8	-38,0	-38,6	-38,3	-38,6	-35,2	-37,3	-37,2
18:33:43	003:34:03	1000	-40,0	-81,8	1,18	-35,8	-38,0	-38,6	-38,3	-38,6	-35,2	-37,3	-37,2
18:34:41	003:35:01	1000	-40,0	-85,9	1,18	-35,8	-38,0	-38,6	-38,3	-38,7	-35,3	-37,3	-37,2
18:35:41	003:36:01	1000	-40,0	-90,0	1,18	-35,8	-38,0	-38,6	-38,3	-38,7	-35,3	-37,3	-37,2
18:36:43	003:37:03	1000	-40,0	-88,9	1,18	-35,8	-38,0	-38,6	-38,3	-38,7	-35,3	-37,3	-37,2
18:37:40	003:38:00	1000	-40,0	-88,5	1,18	-35,8	-38,0	-38,6	-38,3	---	-35,3	-37,3	-37,2
18:38:40	003:39:00	1000	-40,0	-83,2	1,18	-35,8	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:39:42	003:40:02	1000	-40,0	-85,2	1,18	-35,9	-38,0	-38,6	-38,3	-38,7	-35,3	-37,3	-37,2
18:40:44	003:41:04	1000	-40,0	-89,6	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:41:44	003:42:04	1000	-40,0	-89,3	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:42:40	003:43:00	1000	-40,0	-88,9	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:43:41	003:44:01	1000	-40,0	-81,1	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:44:44	003:45:04	1000	-40,0	-86,5	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 14 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



18:45:40	003:46:00	1000	-40,0	-89,4	1,18	---	-38,0	-38,6	-38,4	-38,6	-35,3	-37,3	-37,2
18:46:42	003:47:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:46:42	003:47:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:46:42	003:47:02	1000	-40,0	-90,3	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,3	-37,3	-37,2
18:47:44	003:48:04	1000	-40,0	-85,1	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:48:41	003:49:01	1000	-40,0	-82,0	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:49:43	003:50:03	1000	-40,0	-87,5	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:50:43	003:51:03	1000	-40,0	-89,3	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:51:43	003:52:03	1000	-40,0	-89,8	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:52:43	003:53:03	1000	-40,0	-86,6	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:53:40	003:54:00	1000	-40,0	-82,5	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:54:43	003:55:03	1000	-40,0	-87,0	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:55:41	003:56:01	1000	-40,0	-89,3	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:56:41	003:57:01	1000	-40,0	-89,4	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
18:57:41	003:58:01	1000	-40,0	-88,2	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,3	-37,3
18:58:42	003:59:02	1000	-40,0	-89,4	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,3
18:59:43	004:00:03	1000	-40,0	-83,9	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,3	-37,3
19:00:40	004:01:00	1000	-40,0	-87,9	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,3	-37,2
19:01:42	004:02:02	1000	-40,0	-90,9	1,18	---	-38,0	-38,6	-38,4	-38,7	-35,4	-37,4	-37,2
19:02:43	004:03:03	1000	-40,0	-89,8	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,2
19:03:40	004:04:00	1000	-40,0	-89,5	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,3

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 15 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



19:04:42	004:05:02	1000	-40,0	-81,6	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,3
19:05:40	004:06:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
19:05:40	004:06:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
19:05:40	004:06:00	1000	-40,0	-87,1	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,2
19:06:42	004:07:02	1000	-40,0	-88,9	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,4	-37,4	-37,3
19:07:43	004:08:03	1000	-40,0	-88,7	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:08:44	004:09:04	1000	-40,0	-89,6	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:09:40	004:10:00	1000	-40,0	-89,3	1,18	---	-38,1	-38,6	-38,4	-38,8	-35,5	-37,4	-37,3
19:10:41	004:11:01	1000	-40,0	-82,3	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:11:42	004:12:02	1000	-40,0	-87,5	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:12:40	004:13:00	1000	-40,0	-90,4	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,3	-37,3
19:13:42	004:14:02	1000	-40,0	-89,8	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:14:42	004:15:02	1000	-40,0	-89,2	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:15:44	004:16:04	1000	-40,0	-88,7	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:16:44	004:17:04	1000	-40,0	-85,8	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:17:44	004:18:04	1000	-40,0	-90,1	1,18	-36,0	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:18:40	004:19:00	1000	-40,0	-89,2	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:19:42	004:20:02	1000	-40,0	-87,2	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:20:42	004:21:02	1000	-40,0	-83,7	1,18	-36,0	-38,1	-38,6	-38,4	---	-35,5	-37,4	-37,3
19:21:43	004:22:03	1000	-40,0	-87,9	1,18	-36,0	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:22:40	004:23:00	1000	-40,0	-89,6	1,18	-36,0	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 16 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



19:23:40	004:24:00	1000	-40,0	-90,0	1,18	-36,0	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:24:41	004:25:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
19:24:41	004:25:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
19:24:41	004:25:01	1000	-40,0	-88,7	1,18	-36,0	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:25:43	004:26:03	1000	-40,0	-90,9	1,18	-36,1	-38,1	-38,6	-38,4	---	-35,5	-37,4	-37,3
19:26:41	004:27:01	1000	-40,0	-82,1	1,18	-36,1	-38,1	-38,6	-38,4	---	-35,5	-37,4	-37,3
19:27:42	004:28:02	1000	-40,0	-87,1	1,18	-36,1	-38,1	-38,6	-38,4	---	-35,5	-37,4	-37,3
19:28:44	004:29:04	1000	-40,0	-90,1	1,18	-36,1	-38,1	-38,6	-38,4	---	-35,5	-37,4	-37,3
19:29:40	004:30:00	1000	-40,0	-89,7	1,18	-36,1	-38,1	-38,6	-38,4	-38,6	-35,5	-37,4	-37,3
19:30:41	004:31:01	1000	-40,0	-89,5	1,18	-36,1	-38,1	-38,6	-38,5	-38,6	-35,6	-37,4	-37,3
19:31:42	004:32:02	1000	-40,0	-91,0	1,18	-36,1	-38,1	-38,6	-38,5	-38,7	-35,5	-37,4	-37,3
19:32:42	004:33:02	1000	-40,0	-87,9	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:33:42	004:34:02	1000	-40,0	-86,1	1,18	-36,1	-38,1	-38,6	-38,4	---	-35,6	-37,4	-37,3
19:34:42	004:35:02	1000	-40,0	-89,0	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:35:44	004:36:04	1000	-40,0	-88,7	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,5	-37,4	-37,3
19:36:44	004:37:04	1000	-40,0	-88,8	1,18	-36,1	-38,1	-38,6	-38,4	-38,8	-35,5	-37,4	-37,3
19:37:42	004:38:02	1000	-40,0	-87,4	1,18	-36,1	-38,1	-38,6	-38,5	-38,8	-35,5	-37,4	-37,3
19:38:42	004:39:02	1000	-40,0	-84,4	1,18	-36,1	-38,1	-38,6	-38,5	-38,8	-35,6	-37,4	-37,3
19:39:43	004:40:03	1000	-40,0	-86,3	1,18	-36,1	-38,1	-38,6	-38,4	-38,6	-35,6	-37,4	-37,3
19:40:41	004:41:01	1000	-40,0	-89,1	1,18	-36,1	-38,1	-38,6	-38,4	-38,6	-35,6	-37,4	-37,3
19:41:42	004:42:02	1000	-40,0	-90,5	1,18	-36,1	-38,1	-38,6	-38,5	---	-35,6	-37,4	-37,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 17 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



19:42:43	004:43:03	1000	-40,0	-89,9	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,3
19:43:43	004:44:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
19:43:43	004:44:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
19:43:43	004:44:03	1000	-40,0	-87,6	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:44:44	004:45:04	1000	-40,0	-87,5	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:45:44	004:46:04	1000	-40,0	-89,9	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,3
19:46:40	004:47:00	1000	-40,0	-88,2	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,3
19:47:41	004:48:01	1000	-40,0	-89,2	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,4
19:48:41	004:49:01	1000	-40,0	-88,2	1,18	---	-38,1	-38,7	-38,4	-38,7	-35,6	-37,4	-37,3
19:49:42	004:50:02	1000	-40,0	-80,9	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:50:44	004:51:04	1000	-40,0	-88,7	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:51:42	004:52:02	1000	-40,0	-88,7	1,18	-36,1	-38,1	-38,7	-38,5	---	-35,6	-37,4	-37,4
19:52:44	004:53:04	1000	-40,0	-90,9	1,18	-36,1	-38,1	-38,7	-38,5	-38,7	-35,6	-37,4	-37,4
19:53:40	004:54:00	1000	-40,0	-88,3	1,18	-36,1	-38,1	-38,7	-38,5	-38,7	-35,6	-37,4	-37,4
19:54:41	004:55:01	1000	-40,0	-89,7	1,18	-36,1	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,3
19:55:41	004:56:01	1000	-40,0	-79,9	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,4
19:56:42	004:57:02	1000	-40,0	-88,0	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,3
19:57:43	004:58:03	1000	-40,0	-89,4	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,4
19:58:41	004:59:01	1000	-40,0	-89,7	1,18	-36,1	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,4
19:59:42	005:00:02	1000	-40,0	-89,7	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,6	-37,4	-37,4
20:00:43	005:01:03	1000	-40,0	-82,2	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 18 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	20:01:40	005:02:00	1000	-40,0	-88,2	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,4	
	20:02:40	005:03:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:02:40	005:03:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:02:40	005:03:00	1000	-40,0	-90,0	1,18	---	-38,1	-38,6	-38,4	-38,7	-35,6	-37,4	-37,4	
	20:03:42	005:04:02	1000	-40,0	-87,8	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,6	-37,4	-37,4	
	20:04:43	005:05:03	1000	-40,0	-87,0	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:05:43	005:06:03	1000	-40,0	-88,0	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:06:40	005:07:00	1000	-40,0	-81,9	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:07:42	005:08:02	1000	-40,0	-85,8	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:08:43	005:09:03	1000	-40,0	-89,8	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:09:44	005:10:04	1000	-40,0	-90,3	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:10:41	005:11:01	1000	-40,0	-90,5	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:11:42	005:12:02	1000	-40,0	-82,3	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:12:42	005:13:02	1000	-40,0	-87,4	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:13:43	005:14:03	1000	-40,0	-89,8	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,7	-37,4	-37,4	
	20:14:40	005:15:00	1000	-40,0	-87,7	1,18	---	-38,1	-38,7	-38,5	-38,8	-35,7	-37,4	-37,4	
	20:15:40	005:16:00	1000	-40,0	-91,1	1,18	---	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4	
	20:16:42	005:17:02	1000	-40,0	-89,7	1,18	---	-38,2	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:17:42	005:18:02	1000	-40,0	-85,7	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:18:43	005:19:03	1000	-40,0	-83,5	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:19:44	005:20:04	1000	-40,0	-88,8	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 19 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



20:20:44	005:21:04	1000	-40,0	-88,9	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:21:41	005:22:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
20:21:41	005:22:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
20:21:41	005:22:01	1000	-40,0	-89,7	1,18	-36,2	-38,1	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:22:42	005:23:02	1000	-40,0	-89,4	1,18	-36,2	-38,1	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:23:41	005:24:01	1000	-40,0	-89,2	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:24:42	005:25:02	1000	-40,0	-82,3	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4
20:25:43	005:26:03	1000	-40,0	-86,1	1,18	---	-38,1	-38,6	-38,5	-38,7	-35,7	-37,4	-37,4
20:26:44	005:27:04	1000	-40,0	-90,0	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,4	-37,4
20:27:41	005:28:01	1000	-40,0	-89,7	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:28:42	005:29:02	1000	-40,0	-91,3	1,18	---	-38,2	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:29:43	005:30:03	1000	-40,0	-83,5	1,18	---	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:30:44	005:31:04	1000	-40,0	-80,2	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4
20:31:44	005:32:04	1000	-40,0	-87,8	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:32:41	005:33:01	1000	-40,0	-90,0	1,18	-36,2	-38,2	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4
20:33:42	005:34:02	1000	-40,0	-90,4	1,18	-36,2	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:34:43	005:35:03	1000	-40,0	-89,1	1,18	-36,2	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:35:43	005:36:03	1000	-40,0	-87,8	1,18	-36,2	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:36:43	005:37:03	1000	-40,0	-89,7	1,18	-36,2	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:37:42	005:38:02	1000	-40,0	-82,9	1,18	-36,2	-38,1	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4
20:38:44	005:39:04	1000	-40,0	-89,2	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,7	-37,5	-37,4

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus			
Document issued by:				Signature:	Page: 20 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:39:40	005:40:00	1000	-40,0	-88,8	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:40:40	005:41:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:40:40	005:41:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:40:40	005:41:00	1000	-40,0	-90,3	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:41:41	005:42:01	1000	-40,0	-90,9	1,18	-36,2	-38,1	-38,6	-38,5	-38,7	-35,7	-37,5	-37,4	
	20:42:43	005:43:03	1000	-40,0	-89,1	1,18	-36,2	-38,1	-38,6	-38,4	-38,7	-35,7	-37,5	-37,4	
	20:43:43	005:44:03	1000	-40,0	-80,5	1,18	-36,2	-38,1	-38,6	-38,4	-38,7	-35,7	-37,5	-37,4	
	20:44:44	005:45:04	1000	-40,0	-86,2	1,18	-36,2	-38,1	-38,6	-38,4	-38,7	-35,7	-37,5	-37,4	
	20:45:41	005:46:01	1000	-40,0	-87,3	1,18	-36,2	-38,1	-38,7	-38,4	-38,7	-35,7	-37,5	-37,4	
	20:46:42	005:47:02	1000	-40,0	-87,0	1,18	-36,2	-38,2	-38,7	-38,5	-38,8	-35,7	-37,5	-37,4	
	20:47:43	005:48:03	1000	-40,0	-91,9	1,18	-36,2	-38,2	-38,7	-38,5	---	-35,7	-37,5	-37,4	
	20:48:41	005:49:01	1000	-40,0	-89,3	1,18	-36,2	-38,2	-38,7	-38,5	-38,7	-35,8	-37,5	-37,4	
	20:49:40	005:50:00	1000	-40,0	-84,2	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,8	-37,5	-37,4	
	20:50:42	005:51:02	1000	-40,0	-83,9	1,18	-36,2	-38,1	-38,7	-38,5	-38,7	-35,8	-37,5	-37,4	
	20:51:44	005:52:04	280	-40,0	-90,9	1,18	-36,2	-38,1	-38,7	-38,4	-38,7	-36,0	-37,5	-37,4	
	20:52:41	005:53:01	41,3	-40,0	-89,8	1,18	-36,3	-38,2	-38,7	-38,4	-38,7	-36,2	-37,5	-37,4	
	20:53:43	005:54:03	8,16	-40,0	-83,0	1,18	-36,3	-38,2	-38,7	-38,5	-38,8	-36,3	-37,6	-37,5	
	20:54:40	005:55:00	2,36	-40,0	-89,7	1,18	-36,4	-38,2	-38,7	-38,5	-38,8	-36,3	-37,6	-37,5	
	20:55:42	005:56:02	0,853	-40,0	-83,2	0,731	-36,4	-38,2	-38,7	-38,5	-38,8	-36,2	-37,6	-37,4	
	20:56:42	005:57:02	0,470	-40,0	-75,6	0,379	-36,4	-38,2	-38,7	-38,5	-38,7	-36,2	-37,6	-37,3	
	20:57:43	005:58:03	0,317	-40,0	-89,6	0,237	-36,4	-38,1	-38,6	-38,4	-38,7	-36,2	-37,5	-37,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 21 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	20:58:44	005:59:04	0,245	-40,0	-81,0	0,165	-36,3	-38,0	-38,6	-38,4	-38,7	-36,1	-37,4	-37,1	
	20:59:44	006:00:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:59:44	006:00:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:59:44	006:00:04	0,212	-39,8	-90,7	0,137	-36,3	-38,0	-38,6	-38,4	-38,7	-36,1	-37,4	-37,1	
	21:00:41	006:01:01	0,204	-39,5	-88,5	0,136	-36,5	-38,2	-38,7	-38,4	-38,7	-36,1	-37,6	-37,2	
	21:01:42	006:02:02	0,203	-39,2	-94,0	0,131	-36,7	-38,3	-38,8	-38,6	-38,8	-36,0	-37,8	-37,3	
	21:02:42	006:03:02	0,206	-38,9	-91,0	0,134	-37,0	-38,5	-38,9	-38,7	-38,9	-36,0	-38,0	-37,4	
	21:03:43	006:04:03	0,210	-38,5	-84,0	0,137	-37,2	-38,6	-39,0	-38,8	-39,0	-36,0	-38,1	-37,4	
	21:04:41	006:05:01	0,201	-38,3	-91,7	0,130	-37,3	-38,7	-39,1	-38,9	-39,1	-36,0	-38,2	-37,5	
	21:05:41	006:06:01	0,200	-37,9	-89,4	0,129	-37,5	-38,7	-39,1	-39,0	-39,1	-35,9	-38,2	-37,5	
	21:06:40	006:07:00	0,200	-37,7	-81,7	0,131	-37,6	-38,8	-39,2	-39,1	-39,2	-35,9	-38,3	-37,5	
	21:07:41	006:08:01	0,202	-37,3	-89,9	0,132	-37,7	-38,9	-39,2	-39,1	-39,2	-35,9	-38,3	-37,5	
	21:08:43	006:09:03	0,203	-37,0	-87,4	0,134	-37,8	-38,9	-39,2	-39,2	-39,3	-35,8	-38,3	-37,5	
	21:09:44	006:10:04	0,206	-36,7	-91,2	0,136	-37,9	-39,0	-39,2	-39,2	-39,3	-35,8	-38,3	-37,5	
	21:10:40	006:11:00	0,206	-36,4	-91,6	0,137	-37,9	-39,0	-39,2	-39,3	-39,3	-35,7	-38,3	-37,5	
	21:11:41	006:12:01	0,199	-36,1	-84,5	0,130	-37,9	-39,0	-39,2	-39,3	-39,3	-35,7	-38,3	-37,4	
	21:12:41	006:13:01	0,199	-35,8	-94,2	0,132	-38,0	-39,0	-39,2	-39,3	-39,3	-35,6	-38,3	-37,4	
	21:13:42	006:14:02	0,201	-35,5	-91,3	0,133	-38,0	-39,0	-39,2	-39,3	-39,3	-35,5	-38,3	-37,3	
	21:14:42	006:15:02	0,201	-35,2	-83,0	0,134	-38,0	-39,0	-39,2	-39,3	-39,3	-35,5	-38,2	-37,3	
	21:15:43	006:16:03	0,203	-34,9	-78,2	0,136	-38,0	-39,0	-39,1	-39,3	-39,3	-35,4	-38,2	-37,2	
	21:16:40	006:17:00	0,203	-34,6	-90,8	0,137	-38,0	-38,9	-39,1	-39,2	-39,3	-35,3	-38,1	-37,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	22 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



21:17:44	006:18:04	0,203	-34,3	-82,9	0,137	-38,0	-38,9	-39,0	-39,2	-39,2	-35,2	-38,1	-37,0
21:18:42	006:19:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
21:18:42	006:19:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
21:18:42	006:19:02	0,196	-34,0	-79,2	0,131	-38,0	-38,9	-39,0	-39,2	-39,2	-35,0	-38,0	-37,0
21:19:44	006:20:04	0,198	-33,7	-91,0	0,133	-38,0	-38,8	-38,9	-39,1	-39,2	-34,9	-37,9	-36,9
21:20:44	006:21:04	0,200	-33,4	-82,8	0,135	-37,9	-38,8	-38,9	-39,1	-39,1	-34,8	-37,9	-36,8
21:21:40	006:22:00	0,202	-33,1	-80,3	0,137	-37,9	-38,8	-38,9	-39,0	-39,1	-34,7	-37,8	-36,7
21:22:43	006:23:03	0,197	-32,8	-91,3	0,131	-37,9	-38,7	-38,8	-39,0	-39,0	-34,5	-37,7	-36,6
21:23:40	006:24:00	0,199	-32,5	-84,3	0,133	-37,8	-38,7	-38,7	-39,0	-39,0	-34,4	-37,6	-36,5
21:24:41	006:25:01	0,203	-32,2	-91,3	0,136	-37,8	-38,6	-38,7	-38,9	-38,9	-34,3	-37,6	-36,4
21:25:43	006:26:03	0,205	-31,9	-88,7	0,138	-37,8	-38,5	-38,6	-38,9	-38,8	-34,1	-37,5	-36,3
21:26:44	006:27:04	0,199	-31,6	-79,8	0,132	-37,7	-38,5	-38,6	-38,8	-38,8	-34,0	-37,4	-36,2
21:27:41	006:28:01	0,201	-31,3	-91,1	0,135	-37,6	-38,4	-38,5	-38,7	-38,7	-33,8	-37,3	-36,1
21:28:42	006:29:02	0,204	-31,0	-87,9	0,137	-37,5	-38,4	-38,4	-38,7	-38,6	-33,6	-37,2	-36,0
21:29:43	006:30:03	0,199	-30,7	-79,4	0,131	-37,5	-38,3	-38,3	-38,6	-38,5	-33,5	-37,1	-35,9
21:30:45	006:31:05	0,201	-30,4	-91,0	0,134	-37,4	-38,2	-38,3	-38,5	-38,4	-33,3	-37,0	-35,8
21:31:41	006:32:01	0,205	-30,1	-87,4	0,135	-37,4	-38,1	-38,2	-38,4	-38,4	-33,1	-36,9	-35,7
21:32:43	006:33:03	0,202	-29,8	-79,1	0,134	-37,3	-38,1	-38,1	-38,4	-38,3	-33,0	-36,8	-35,5
21:33:43	006:34:03	0,201	-29,5	-90,4	0,132	-37,2	-38,0	-38,0	-38,3	-38,2	-32,8	-36,7	-35,4
21:34:40	006:35:00	0,204	-29,2	-88,2	0,134	-37,1	-37,9	-37,9	-38,2	-38,1	-32,6	-36,6	-35,3
21:35:41	006:36:01	0,207	-28,9	-79,5	0,137	-37,0	-37,8	-37,9	-38,1	-38,1	-32,4	-36,5	-35,2

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 23 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:36:44	006:37:04	0,202	-28,5	-89,8	0,132	-37,0	-37,7	-37,8	-38,1	-38,0	-32,3	-36,3	-35,1	
	21:37:42	006:38:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:37:42	006:38:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:37:42	006:38:02	0,203	-28,3	-87,1	0,134	-36,9	-37,6	-37,7	-38,0	-37,9	-32,1	-36,2	-35,0	
	21:38:43	006:39:03	0,207	-27,9	-79,4	0,136	-36,8	-37,6	-37,6	-37,9	-37,8	-31,9	-36,1	-34,9	
	21:39:41	006:40:01	0,200	-27,7	-90,2	0,131	-36,7	-37,5	-37,6	-37,8	-37,7	-31,7	-36,0	-34,8	
	21:40:41	006:41:01	0,203	-27,3	-87,6	0,133	-36,6	-37,4	-37,5	-37,7	-37,6	-31,6	-35,9	-34,6	
	21:41:41	006:42:01	0,206	-27,0	-80,5	0,135	-36,5	-37,3	-37,4	-37,7	-37,5	-31,4	-35,8	-34,5	
	21:42:44	006:43:04	0,204	-26,7	-89,5	0,132	-36,5	-37,2	-37,3	-37,6	-37,4	-31,2	-35,7	-34,4	
	21:43:40	006:44:00	0,202	-26,4	-86,7	0,132	-36,4	-37,1	-37,2	-37,5	-37,4	-31,0	-35,6	-34,3	
	21:44:44	006:45:04	0,206	-26,1	-78,1	0,135	-36,3	-37,0	-37,1	-37,4	-37,3	-30,8	-35,5	-34,2	
	21:45:40	006:46:00	0,209	-25,8	-88,7	0,137	-36,2	-36,9	-37,1	-37,4	-37,2	-30,6	-35,4	-34,1	
	21:46:42	006:47:02	0,203	-25,5	-85,1	0,131	-36,1	-36,9	-37,0	-37,3	-37,1	-30,4	-35,2	-33,9	
	21:47:43	006:48:03	0,206	-25,2	-84,1	0,134	-36,0	-36,8	-36,9	-37,2	-37,0	-30,2	-35,2	-33,8	
	21:48:44	006:49:04	0,210	-24,9	-90,2	0,137	-35,9	-36,7	-36,9	-37,1	-36,9	-30,0	-35,1	-33,7	
	21:49:40	006:50:00	0,202	-24,6	-84,8	0,131	-35,8	-36,6	-36,8	-37,0	-36,8	-29,8	-34,9	-33,6	
	21:50:40	006:51:00	0,205	-24,3	-93,9	0,133	-35,8	-36,5	-36,7	-37,0	-36,7	-29,6	-34,8	-33,5	
	21:51:42	006:52:02	0,209	-24,0	-90,2	0,136	-35,7	-36,4	-36,6	-36,9	-36,7	-29,4	-34,7	-33,3	
	21:52:45	006:53:05	0,202	-23,7	-83,3	0,130	-35,6	-36,3	-36,6	-36,8	-36,6	-29,2	-34,6	-33,2	
	21:53:42	006:54:02	0,205	-23,4	-80,2	0,133	-35,5	-36,3	-36,5	-36,8	-36,5	-29,0	-34,5	-33,1	
	21:54:41	006:55:01	0,208	-23,1	-90,1	0,135	-35,4	-36,2	-36,4	-36,7	-36,4	-28,8	-34,4	-33,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 24 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:55:41	006:56:01	0,200	-22,8	-82,9	0,130	-35,3	-36,1	-36,4	-36,6	-36,3	-28,6	-34,4	-32,9
	21:56:43	006:57:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	21:56:43	006:57:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	21:56:43	006:57:03	0,204	-22,5	-90,9	0,132	-35,2	-36,0	-36,3	-36,5	-36,2	-28,4	-34,3	-32,8
	21:57:40	006:58:00	0,207	-22,2	-89,9	0,135	-35,2	-35,9	-36,2	-36,5	-36,1	-28,2	-34,2	-32,6
	21:58:41	006:59:01	0,211	-21,9	-82,5	0,138	-35,1	-35,9	-36,2	-36,4	-36,1	-28,0	-34,0	-32,5
	21:59:44	007:00:04	0,204	-21,6	-90,6	0,132	---	-35,8	-36,1	-36,3	-36,0	-27,7	-33,9	-32,4
	22:00:44	007:01:04	0,207	-21,3	-88,5	0,135	---	-35,7	-36,0	-36,3	-35,9	-27,5	-33,9	-32,3
	22:01:41	007:02:01	0,211	-21,0	-81,7	0,138	---	-35,6	-36,0	-36,2	-35,8	-27,3	-33,8	-32,2
	22:02:42	007:03:02	0,203	-20,7	-89,6	0,132	---	-35,5	-35,9	-36,1	-35,8	-27,1	-33,7	-32,1
	22:03:40	007:04:00	0,207	-20,4	-87,9	0,134	---	-35,5	-35,8	-36,0	-35,7	-26,9	-33,6	-32,0
	22:04:43	007:05:03	0,211	-20,1	-81,2	0,138	---	-35,4	-35,7	-36,0	-35,6	-26,6	-33,5	-31,8
	22:05:44	007:06:04	0,203	-19,8	-88,4	0,132	---	-35,3	-35,7	-35,9	-35,5	-26,4	-33,4	-31,7
	22:06:40	007:07:00	0,207	-19,5	-85,8	0,134	---	-35,2	-35,6	-35,8	-35,4	-26,2	-33,3	-31,6
	22:07:44	007:08:04	0,211	-19,1	-83,2	0,137	---	-35,1	-35,5	-35,8	-35,3	-26,0	-33,2	-31,5
	22:08:41	007:09:01	0,203	-18,9	-89,8	0,131	-34,1	-35,1	-35,5	-35,7	-35,3	-25,8	-33,1	-31,4
	22:09:44	007:10:04	0,207	-18,5	-83,5	0,134	-34,1	-35,0	-35,4	-35,6	-35,2	-25,6	-33,0	-31,3
	22:10:45	007:11:05	0,212	-18,2	-91,3	0,138	-34,0	-34,9	-35,3	-35,5	-35,1	-25,3	-32,9	-31,2
	22:11:40	007:12:00	0,203	-18,0	-90,0	0,131	-33,9	-34,8	-35,2	-35,5	-35,0	-25,1	-32,9	-31,0
	22:12:41	007:13:01	0,207	-17,7	-82,7	0,135	-33,8	-34,7	-35,2	-35,4	-34,9	-24,9	-32,8	-30,9
	22:13:43	007:14:03	0,212	-17,3	-89,2	0,138	-33,7	-34,6	-35,1	-35,3	-34,8	-24,7	-32,7	-30,8

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 25 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



22:14:40	007:15:00	0,203	-17,1	-88,1	0,131	-33,7	-34,6	-35,0	-35,3	-34,8	-24,5	-32,6	-30,7
22:15:41	007:16:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
22:15:41	007:16:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
22:15:41	007:16:01	0,208	-16,7	-83,9	0,135	---	-34,5	-34,9	-35,2	-34,7	-24,3	-32,5	-30,6
22:16:42	007:17:02	0,205	-16,4	-88,7	0,132	---	-34,4	-34,9	-35,1	-34,6	-24,1	-32,4	-30,4
22:17:43	007:18:03	0,204	-16,1	-85,2	0,132	---	-34,3	-34,8	-35,0	-34,5	-23,8	-32,3	-30,3
22:18:44	007:19:04	0,210	-15,8	-92,7	0,136	---	-34,2	-34,7	-35,0	-34,4	-23,6	-32,3	-30,2
22:19:40	007:20:00	0,200	-15,5	-90,0	0,130	-33,2	-34,1	-34,7	-34,9	-34,3	-23,4	-32,2	-30,1
22:20:41	007:21:01	0,205	-15,2	-83,4	0,133	-33,2	-34,0	-34,6	-34,8	-34,3	-23,1	-32,1	-30,0
22:21:43	007:22:03	0,210	-14,9	-90,6	0,136	---	-33,9	-34,5	-34,7	-34,2	-22,9	-32,0	-29,9
22:22:44	007:23:04	0,201	-14,6	-88,1	0,130	-33,0	-33,9	-34,5	-34,7	-34,1	-22,7	-31,9	-29,8
22:23:44	007:24:04	0,206	-14,3	-81,3	0,134	-33,0	-33,8	-34,4	-34,6	-34,0	-22,5	-31,8	-29,7
22:24:40	007:25:00	0,211	-14,0	-88,6	0,136	---	-33,7	-34,3	-34,5	-33,9	-22,2	-31,8	-29,5
22:25:41	007:26:01	0,202	-13,7	-86,3	0,131	-32,6	-33,6	-34,2	-34,4	-33,8	-22,0	-31,7	-29,4
22:26:43	007:27:03	0,208	-13,4	-92,1	0,134	-32,6	-33,5	-34,2	-34,3	-33,8	-21,8	-31,6	-29,3
22:27:40	007:28:00	0,210	-13,1	-90,4	0,137	-32,5	-33,4	-34,1	-34,3	-33,7	-21,5	-31,5	-29,2
22:28:42	007:29:02	0,203	-12,8	-84,0	0,131	---	-33,3	-34,0	-34,2	-33,6	-21,3	-31,5	-29,1
22:29:43	007:30:03	0,208	-12,5	-88,0	0,135	---	-33,3	-34,0	-34,1	-33,5	-21,1	-31,4	-29,0
22:30:44	007:31:04	0,212	-12,2	-86,1	0,138	---	-33,2	-33,9	-34,0	-33,4	-20,8	-31,3	-28,9
22:31:40	007:32:00	0,204	-11,9	-93,2	0,132	-32,0	-33,1	-33,8	-34,0	-33,4	-20,6	-31,3	-28,8
22:32:41	007:33:01	0,208	-11,6	-89,9	0,135	-32,0	-33,0	-33,7	-33,9	-33,3	-20,4	-31,2	-28,7

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 26 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



22:33:42	007:34:02	0,201	-11,3	-83,4	0,129	---	-32,9	-33,7	-33,8	-33,2	-20,2	-31,1	-28,6
22:34:43	007:35:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
22:34:43	007:35:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
22:34:43	007:35:03	0,205	-11,0	-88,2	0,133	-31,8	-32,8	-33,6	-33,7	-33,1	-20,0	-31,1	-28,4
22:35:40	007:36:00	0,209	-10,7	-87,1	0,136	-31,8	-32,8	-33,5	-33,6	-33,0	-19,8	-31,0	-28,3
22:36:40	007:37:00	0,202	-10,4	-91,7	0,130	---	-32,7	-33,5	-33,6	-33,0	-19,5	-30,9	-28,2
22:37:43	007:38:03	0,206	-10,1	-89,6	0,133	---	-32,6	-33,4	-33,5	-32,9	-19,3	-30,9	-28,1
22:38:43	007:39:03	0,211	-9,8	-82,8	0,137	---	-32,5	-33,3	-33,4	-32,9	-19,1	-30,8	-28,0
22:39:40	007:40:00	0,203	-9,5	-88,9	0,131	---	-32,5	-33,3	-33,3	-32,8	-18,9	-30,7	-27,9
22:40:42	007:41:02	0,207	-9,2	-87,5	0,134	---	-32,4	-33,2	-33,3	-32,7	-18,6	-30,7	-27,8
22:41:42	007:42:02	0,200	-8,9	-92,3	0,129	---	-32,3	-33,2	-33,2	-32,6	-18,4	-30,6	-27,7
22:42:40	007:43:00	0,204	-8,6	-89,3	0,132	-31,2	-32,2	-33,1	-33,1	-32,6	-18,2	-30,5	-27,6
22:43:44	007:44:04	0,209	-8,2	-83,2	0,136	-31,2	-32,2	-33,0	-33,1	-32,5	-17,9	-30,5	-27,5
22:44:40	007:45:00	0,201	-8,0	-89,3	0,130	-31,2	-32,1	-33,0	-33,0	-32,4	-17,7	-30,4	-27,4
22:45:42	007:46:02	0,205	-7,6	-87,5	0,133	---	-32,0	-32,9	-32,9	-32,3	-17,5	-30,4	-27,3
22:46:42	007:47:02	0,211	-7,3	-80,1	0,137	-30,8	-32,0	-32,9	-32,8	-32,3	-17,3	-30,3	-27,2
22:47:41	007:48:01	0,203	-7,0	-89,1	0,131	-30,8	-31,9	-32,8	-32,7	-32,2	-17,0	-30,2	-27,1
22:48:44	007:49:04	0,209	-6,7	-83,7	0,135	---	-31,8	-32,7	-32,7	-32,1	-16,8	-30,2	-27,0
22:49:44	007:50:04	0,201	-6,4	-88,1	0,130	-30,7	-31,8	-32,7	-32,6	-32,1	-16,6	-30,1	-26,9
22:50:40	007:51:00	0,205	-6,1	-87,8	0,133	-30,7	-31,7	-32,6	-32,5	-32,0	-16,3	-30,1	-26,8
22:51:41	007:52:01	0,211	-5,8	-90,5	0,137	---	-31,6	-32,5	-32,4	-31,9	-16,1	-30,0	-26,7

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 27 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	22:52:42	007:53:02	0,202	-5,5	-89,2	0,131	---	-31,5	-32,5	-32,3	-31,9	-15,9	-29,9	-26,7	
	22:53:42	007:54:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:53:42	007:54:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:53:42	007:54:02	0,208	-5,2	-83,5	0,135	---	-31,5	-32,4	-32,3	-31,8	-15,7	-29,9	-26,6	
	22:54:44	007:55:04	0,201	-4,9	-87,9	0,130	---	-31,4	-32,4	-32,2	-31,8	-15,4	-29,8	-26,5	
	22:55:41	007:56:01	0,205	-4,6	-87,4	0,133	---	-31,3	-32,3	-32,1	-31,7	-15,2	-29,8	-26,4	
	22:56:42	007:57:02	0,211	-4,3	-91,6	0,137	-30,1	-31,3	-32,2	-32,0	-31,6	-15,0	-29,7	-26,3	
	22:57:43	007:58:03	0,202	-4,0	-89,4	0,131	-30,1	-31,2	-32,2	-31,9	-31,5	-14,7	-29,7	-26,2	
	22:58:44	007:59:04	0,208	-3,7	-90,4	0,135	---	-31,1	-32,2	-31,8	-31,5	-14,5	-29,6	-26,1	
	22:59:40	008:00:00	0,201	-3,4	-88,0	0,129	---	-31,1	-32,1	-31,8	-31,4	-14,3	-29,6	-26,0	
	23:00:40	008:01:00	0,206	-3,1	-85,3	0,133	-29,8	-31,0	-32,0	-31,7	-31,4	-14,0	-29,5	-25,9	
	23:01:41	008:02:01	0,211	-2,8	-90,4	0,137	-29,8	-30,9	-32,0	-31,6	-31,3	-13,8	-29,5	-25,8	
	23:02:41	008:03:01	0,204	-2,5	-88,5	0,131	---	-30,9	-32,0	-31,5	-31,3	-13,6	-29,4	-25,8	
	23:03:43	008:04:03	0,209	-2,2	-92,1	0,136	-29,6	-30,8	-31,9	-31,4	-31,2	-13,3	-29,4	-25,7	
	23:04:44	008:05:04	0,201	-1,9	-89,1	0,130	-29,6	-30,8	-31,8	-31,3	-31,1	-13,1	-29,3	-25,6	
	23:05:42	008:06:02	0,206	-1,6	-83,9	0,134	---	-30,7	-31,8	-31,3	-31,1	-12,9	-29,3	-25,5	
	23:06:42	008:07:02	0,211	-1,3	-87,9	0,137	---	-30,6	-31,7	-31,2	-31,0	-12,6	-29,2	-25,4	
	23:07:43	008:08:03	0,204	-1,0	-86,5	0,132	---	-30,6	-31,7	-31,1	-31,0	-12,4	-29,2	-25,4	
	23:08:44	008:09:04	0,209	-0,7	-90,1	0,136	-29,3	-30,5	-31,7	-31,1	-30,9	-12,1	-29,1	-25,3	
	23:09:40	008:10:00	0,201	-0,4	-89,8	0,130	-29,2	-30,4	-31,6	-31,0	-30,9	-11,9	-29,1	-25,2	
	23:10:44	008:11:04	0,209	-0,1	-91,2	0,134	-29,2	-30,4	-31,6	-30,9	-30,8	-11,6	-29,0	-25,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 28 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:11:45	008:12:05	0,212	0,2	-89,5	0,138	-29,1	-30,3	-31,5	-30,8	-30,8	-11,4	-29,0	-25,0	
	23:12:43	008:13:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:12:43	008:13:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:12:43	008:13:03	0,206	0,5	-80,6	0,133	---	-30,3	-31,5	-30,7	-30,7	-11,2	-28,9	-25,0	
	23:13:44	008:14:04	0,210	0,9	-88,7	0,137	---	-30,2	-31,5	-30,6	-30,6	-11,0	-28,9	-24,9	
	23:14:44	008:15:04	0,203	1,2	-84,1	0,132	---	-30,1	-31,4	-30,6	-30,6	-10,7	-28,9	-24,8	
	23:15:41	008:16:01	0,207	1,4	-87,5	0,135	---	-30,1	-31,3	-30,5	-30,5	-10,5	-28,8	-24,7	
	23:16:44	008:17:04	0,201	1,8	-86,6	0,130	---	-30,0	-31,3	-30,4	-30,5	-10,3	-28,8	-24,7	
	23:17:40	008:18:00	0,206	2,0	-89,2	0,134	-28,7	-30,0	-31,3	-30,3	-30,4	-10,0	-28,7	-24,6	
	23:18:42	008:19:02	0,211	2,4	-88,3	0,137	-28,7	-30,0	-31,2	-30,2	-30,4	-9,8	-28,7	-24,5	
	23:19:43	008:20:03	0,205	2,7	-90,7	0,133	-28,6	-29,9	-31,2	-30,2	-30,3	-9,5	-28,6	-24,4	
	23:20:44	008:21:04	0,209	3,0	-89,5	0,136	-28,6	-29,8	-31,1	-30,1	-30,3	-9,3	-28,6	-24,4	
	23:21:44	008:22:04	0,203	3,3	-90,2	0,132	---	-29,8	-31,1	-30,0	-30,2	-9,1	-28,6	-24,3	
	23:22:42	008:23:02	0,208	3,6	-89,3	0,135	---	-29,8	-31,1	-29,9	-30,2	-8,9	-28,5	-24,2	
	23:23:40	008:24:00	0,201	3,9	-88,0	0,130	---	-29,7	-31,0	-29,8	-30,2	-8,6	-28,5	-24,1	
	23:24:44	008:25:04	0,207	4,2	-89,2	0,134	---	-29,7	-31,0	-29,8	-30,1	-8,4	-28,5	-24,1	
	23:25:40	008:26:00	0,211	4,5	-93,4	0,138	---	-29,7	-31,0	-29,7	-30,1	-8,2	-28,4	-24,0	
	23:26:42	008:27:02	0,204	4,8	-89,8	0,132	---	-29,7	-31,0	-29,6	-30,0	-7,9	-28,4	-23,9	
	23:27:44	008:28:04	0,210	5,1	-87,3	0,137	---	-29,6	-30,9	-29,5	-30,0	-7,7	-28,4	-23,8	
	23:28:41	008:29:01	0,202	5,4	-89,2	0,131	---	-29,6	-30,9	-29,4	-29,9	-7,4	-28,3	-23,8	
	23:29:42	008:30:02	0,208	5,7	-84,7	0,136	-27,9	-29,5	-30,9	-29,3	-29,9	-7,2	-28,3	-23,7	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 29 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	23:30:43	008:31:03	0,202	6,0	-87,0	0,131	-27,9	-29,5	-30,8	-29,3	-29,9	-6,9	-28,3	-23,6	
	23:31:40	008:32:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:31:40	008:32:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:31:40	008:32:00	0,207	6,3	-87,0	0,134	---	-29,5	-30,8	-29,2	-29,8	-6,7	-28,3	-23,6	
	23:32:41	008:33:01	0,201	6,6	-88,6	0,129	---	-29,4	-30,8	-29,1	-29,8	-6,5	-28,2	-23,5	
	23:33:43	008:34:03	0,206	6,9	-90,6	0,133	---	-29,4	-30,7	-29,0	-29,7	-6,2	-28,2	-23,4	
	23:34:44	008:35:04	0,212	7,2	-88,1	0,137	---	-29,3	-30,7	-28,9	-29,7	-6,0	-28,2	-23,3	
	23:35:41	008:36:01	0,205	7,5	-88,5	0,132	---	-29,3	-30,7	-28,8	-29,6	-5,7	-28,1	-23,3	
	23:36:42	008:37:02	0,210	7,8	-89,6	0,137	---	-29,2	-30,6	-28,7	-29,6	-5,5	-28,1	-23,2	
	23:37:44	008:38:04	0,204	8,1	-88,2	0,132	---	-29,2	-30,6	-28,6	-29,5	-5,3	-28,1	-23,1	
	23:38:40	008:39:00	0,209	8,4	-90,8	0,136	---	-29,1	-30,5	-28,5	-29,5	-5,0	-28,0	-23,1	
	23:39:44	008:40:04	0,203	8,7	-88,9	0,131	---	-29,1	-30,5	-28,5	-29,4	-4,8	-28,0	-23,0	
	23:40:40	008:41:00	0,208	9,0	-90,5	0,135	---	-29,0	-30,5	-28,4	-29,4	-4,5	-28,0	-23,0	
	23:41:41	008:42:01	0,202	9,3	-87,5	0,130	---	-29,0	-30,5	-28,3	-29,3	-4,3	-27,9	-22,9	
	23:42:44	008:43:04	0,207	9,6	-86,8	0,135	---	-28,9	-30,4	-28,2	-29,3	-4,1	-27,9	-22,8	
	23:43:40	008:44:00	0,201	9,9	-87,7	0,130	---	-28,9	-30,4	-28,1	-29,2	-3,8	-27,9	-22,8	
	23:44:42	008:45:02	0,207	10,0	-86,7	0,134	---	-28,8	-30,4	-28,0	-29,2	-3,6	-27,9	-22,7	
	23:45:43	008:46:03	0,210	10,0	-87,6	0,135	---	-28,8	-30,3	-27,9	-29,2	-3,4	-27,9	-22,7	
	23:46:40	008:47:00	0,205	10,0	-87,4	0,133	---	-28,7	-30,3	-27,8	-29,1	-3,1	-27,8	-22,6	
	23:47:40	008:48:00	0,211	10,0	-89,4	0,138	---	-28,7	-30,3	-27,7	-29,1	-2,9	-27,8	-22,5	
	23:48:43	008:49:03	0,204	10,0	-90,1	0,132	---	-28,7	-30,3	-27,7	-29,1	-2,7	-27,8	-22,5	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 30 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



			0,210	10,0	-88,6	0,137	---	-28,6	-30,2	-27,6	-29,0	-2,5	-27,8	-22,5	
			Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
			mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:49:44	008:50:04	0,210	10,0	-88,6	0,137	---	-28,6	-30,2	-27,6	-29,0	-2,5	-27,8	-22,5	
	23:50:44	008:51:04	0,203	10,0	-91,1	0,132	---	-28,6	-30,2	-27,5	-29,0	-2,2	-27,8	-22,4	
	23:51:41	008:52:01	0,209	10,0	-88,4	0,136	---	-28,6	-30,2	-27,4	-29,0	-2,0	-27,8	-22,4	
	23:52:43	008:53:03	0,202	10,0	-91,7	0,131	---	-28,6	-30,2	-27,4	-29,0	-1,8	-27,8	-22,4	
	23:53:44	008:54:04	0,208	10,0	-87,5	0,135	---	-28,6	-30,2	-27,3	-29,0	-1,6	-27,8	-22,3	
	23:54:40	008:55:00	0,201	10,0	-93,0	0,130	---	-28,5	-30,2	-27,3	-29,0	-1,5	-27,8	-22,3	
	23:55:41	008:56:01	0,206	10,0	-88,6	0,134	---	-28,5	-30,2	-27,2	-29,0	-1,3	-27,8	-22,3	
	23:56:44	008:57:04	0,204	10,0	-91,1	0,131	---	-28,5	-30,2	-27,1	-29,0	-1,1	-27,9	-22,3	
	23:57:41	008:58:01	0,204	10,0	-90,2	0,133	---	-28,5	-30,2	-27,1	-29,0	-1,0	-27,9	-22,3	
	23:58:42	008:59:02	0,211	10,0	-91,1	0,138	---	-28,5	-30,3	-27,1	-29,0	-0,9	-27,9	-22,3	
	23:59:43	009:00:03	0,203	10,0	-89,8	0,132	---	-28,5	-30,3	-27,0	-29,0	-0,7	-27,9	-22,3	
2021-11-24	00:00:44	009:01:04	0,210	10,0	-89,8	0,137	---	-28,5	-30,3	-27,0	-29,0	-0,6	-28,0	-22,3	
	00:01:44	009:02:04	0,203	10,0	-91,1	0,131	---	-28,5	-30,3	-26,9	-29,0	-0,5	-28,0	-22,3	
	00:02:44	009:03:04	0,208	10,0	-88,3	0,135	---	-28,5	-30,3	-26,9	-29,0	-0,4	-28,0	-22,4	
	00:03:40	009:04:00	0,201	10,0	-92,8	0,130	---	-28,5	-30,3	-26,8	-29,0	-0,3	-28,0	-22,4	
	00:04:41	009:05:01	0,206	10,0	-88,6	0,134	---	-28,5	-30,3	-26,8	-29,0	-0,2	-28,0	-22,4	
	00:05:44	009:06:04	0,200	10,0	-88,5	0,129	---	-28,5	-30,3	-26,8	-29,0	-0,1	-28,1	-22,4	
	00:06:40	009:07:00	0,207	10,0	-87,3	0,133	---	-28,5	-30,3	-26,7	-29,0	0,0	-28,1	-22,4	
	00:07:40	009:08:00	0,210	10,0	-87,5	0,137	---	-28,5	-30,3	-26,7	-29,0	0,1	-28,1	-22,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus			
Document issued by:			Signature:		Page: 31 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	00:08:42	009:09:02	0,204	10,0	-87,6	0,132	---	-28,6	-30,3	-26,6	-29,0	0,2	-28,1	-22,4	
	00:09:42	009:10:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:09:42	009:10:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:09:42	009:10:02	0,208	10,0	-88,7	0,136	---	-28,6	-30,3	-26,6	-29,0	0,2	-28,2	-22,5	
	00:10:43	009:11:03	0,202	10,0	-88,0	0,131	---	-28,6	-30,3	-26,6	-29,1	0,3	-28,2	-22,5	
	00:11:44	009:12:04	0,208	10,0	-92,8	0,135	---	-28,6	-30,3	-26,5	-29,1	0,4	-28,2	-22,5	
	00:12:40	009:13:00	0,200	10,0	-87,6	0,129	---	-28,6	-30,3	-26,5	-29,1	0,5	-28,3	-22,5	
	00:13:40	009:14:00	0,205	10,0	-92,3	0,133	---	-28,6	-30,3	-26,5	-29,1	0,5	-28,3	-22,5	
	00:14:41	009:15:01	0,210	10,0	-89,6	0,137	---	-28,6	-30,3	-26,4	-29,1	0,6	-28,3	-22,5	
	00:15:42	009:16:02	0,204	10,0	-90,8	0,132	---	-28,6	-30,3	-26,4	-29,1	0,6	-28,4	-22,6	
	00:16:44	009:17:04	0,209	10,0	-89,9	0,136	---	-28,6	-30,3	-26,4	-29,1	0,7	-28,4	-22,6	
	00:17:44	009:18:04	0,204	10,0	-90,0	0,132	-26,4	-28,6	-30,3	-26,3	-29,1	0,8	-28,4	-22,6	
	00:18:41	009:19:01	0,207	10,0	-88,7	0,135	-26,4	-28,6	-30,3	-26,3	-29,1	0,8	-28,5	-22,6	
	00:19:41	009:20:01	0,202	10,0	-92,3	0,130	-26,4	-28,6	-30,3	-26,2	-29,1	0,9	-28,5	-22,6	
	00:20:42	009:21:02	0,206	10,0	-88,8	0,134	-26,5	-28,6	-30,3	-26,2	-29,1	0,9	-28,5	-22,6	
	00:21:43	009:22:03	0,207	10,0	-91,2	0,134	-26,5	-28,6	-30,3	-26,2	-29,2	1,0	-28,5	-22,6	
	00:22:44	009:23:04	0,204	10,0	-89,8	0,133	-26,5	-28,6	-30,3	-26,1	-29,2	1,0	-28,6	-22,7	
	00:23:41	009:24:01	0,210	10,0	-90,2	0,137	-26,5	-28,6	-30,4	-26,1	-29,2	1,1	-28,6	-22,7	
	00:24:42	009:25:02	0,203	10,0	-91,1	0,131	-26,5	-28,6	-30,4	-26,1	-29,2	1,1	-28,6	-22,8	
	00:25:44	009:26:04	0,208	10,0	-88,2	0,135	-26,5	-28,6	-30,4	-26,1	-29,2	1,1	-28,6	-22,8	
	00:26:42	009:27:02	0,202	10,0	-87,5	0,131	-26,6	-28,6	-30,4	-26,1	-29,2	1,2	-28,7	-22,8	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 32 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



00:27:42	009:28:02	0,207	10,0	-87,1	0,135	-26,6	-28,6	-30,4	-26,1	-29,2	1,2	-28,7	-22,8
00:28:44	009:29:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
00:28:44	009:29:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
00:28:44	009:29:04	0,201	10,0	-93,2	0,130	-26,5	-28,6	-30,4	-26,0	-29,2	1,3	-28,7	-22,9
00:29:44	009:30:04	0,205	10,0	-89,3	0,133	-26,5	-28,6	-30,4	-26,0	-29,3	1,3	-28,7	-22,9
00:30:41	009:31:01	0,210	10,0	-91,3	0,137	-26,5	-28,6	-30,4	-26,0	-29,3	1,3	-28,8	-22,9
00:31:42	009:32:02	0,204	10,0	-89,3	0,132	-26,5	-28,6	-30,4	-25,9	-29,3	1,4	-28,8	-22,9
00:32:42	009:33:02	0,209	10,0	-89,5	0,136	-26,6	-28,7	-30,4	-25,9	-29,3	1,4	-28,9	-22,9
00:33:42	009:34:02	0,202	10,0	-88,8	0,131	-26,6	-28,7	-30,4	-25,9	-29,3	1,4	-28,9	-23,0
00:34:42	009:35:02	0,208	10,0	-90,9	0,136	-26,6	-28,7	-30,5	-25,9	-29,3	1,5	-28,9	-23,0
00:35:44	009:36:04	0,202	10,0	-88,4	0,130	-26,6	-28,7	-30,5	-25,8	-29,3	1,5	-28,9	-23,0
00:36:44	009:37:04	0,207	10,0	-93,1	0,134	-26,6	-28,7	-30,4	-25,8	-29,3	1,5	-29,0	-23,0
00:37:41	009:38:01	0,211	10,0	-88,4	0,138	-26,6	-28,7	-30,4	-25,8	-29,3	1,5	-29,0	-23,0
00:38:41	009:39:01	0,205	10,0	-89,7	0,133	-26,6	-28,7	-30,4	-25,8	-29,3	1,6	-29,0	-23,1
00:39:41	009:40:01	0,209	10,0	-87,9	0,136	-26,7	-28,7	-30,5	-25,7	-29,3	1,6	-29,1	-23,1
00:40:42	009:41:02	0,204	10,0	-93,8	0,132	-26,6	-28,7	-30,5	-25,7	-29,4	1,6	-29,1	-23,1
00:41:43	009:42:03	0,208	10,0	-87,5	0,135	-26,7	-28,7	-30,4	-25,7	-29,4	1,7	-29,1	-23,1
00:42:43	009:43:03	0,202	10,0	-92,0	0,131	-26,7	-28,7	-30,4	-25,6	-29,4	1,7	-29,1	-23,2
00:43:44	009:44:04	0,207	10,0	-87,7	0,134	-26,7	-28,7	-30,4	-25,6	-29,4	1,7	-29,2	-23,2
00:44:41	009:45:01	0,211	10,0	-86,4	0,138	-26,7	-28,7	-30,4	-25,6	-29,4	1,7	-29,2	-23,2
00:45:43	009:46:03	0,204	10,0	-86,8	0,133	---	-28,7	-30,4	-25,6	-29,4	1,7	-29,2	-23,2

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 33 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	00:46:40	009:47:00	0,210	10,0	-93,0	0,136	---	-28,7	-30,4	-25,6	-29,4	1,8	-29,3	-23,3	
	00:47:40	009:48:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:47:40	009:48:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:47:40	009:48:00	0,202	10,0	-89,0	0,131	---	-28,7	-30,4	-25,5	-29,4	1,8	-29,3	-23,3	
	00:48:40	009:49:00	0,206	10,0	-89,9	0,134	---	-28,7	-30,4	-25,5	-29,4	1,8	-29,3	-23,3	
	00:49:41	009:50:01	0,211	10,0	-91,5	0,138	---	-28,7	-30,4	-25,5	-29,4	1,9	-29,3	-23,4	
	00:50:43	009:51:03	0,205	10,0	-88,1	0,133	---	-28,7	-30,4	-25,5	-29,4	1,9	-29,4	-23,4	
	00:51:43	009:52:03	0,210	10,0	-91,5	0,137	---	-28,8	-30,4	-25,4	-29,5	1,9	-29,4	-23,4	
	00:52:42	009:53:02	0,203	10,0	-89,2	0,131	---	-28,8	-30,4	-25,4	-29,5	1,9	-29,4	-23,4	
	00:53:43	009:54:03	0,209	10,0	-89,3	0,136	-26,8	-28,8	-30,4	-25,4	-29,5	1,9	-29,4	-23,5	
	00:54:43	009:55:03	0,203	10,0	-87,4	0,131	-26,8	-28,8	-30,4	-25,4	-29,5	2,0	-29,5	-23,5	
	00:55:43	009:56:03	0,208	10,0	-93,2	0,135	-26,8	-28,8	-30,5	-25,4	-29,5	2,0	-29,5	-23,5	
	00:56:45	009:57:05	0,200	10,0	-88,6	0,129	-26,8	-28,8	-30,5	-25,3	-29,5	2,0	-29,5	-23,5	
	00:57:41	009:58:01	0,206	10,0	-88,4	0,133	-26,8	-28,8	-30,4	-25,3	-29,5	2,0	-29,6	-23,6	
	00:58:43	009:59:03	0,210	10,0	-88,7	0,137	-26,8	-28,8	-30,4	-25,3	-29,5	2,1	-29,6	-23,6	
	00:59:44	010:00:04	0,205	10,0	-90,8	0,132	-26,9	-28,8	-30,5	-25,3	-29,5	2,1	-29,6	-23,6	
	01:00:40	010:01:00	0,208	10,0	-89,7	0,135	-26,9	-28,8	-30,5	-25,3	-29,5	2,1	-29,6	-23,6	
	01:01:41	010:02:01	0,201	10,0	-89,1	0,130	-26,9	-28,8	-30,5	-25,2	-29,5	2,1	-29,6	-23,7	
	01:02:40	010:03:00	0,206	10,0	-88,7	0,133	-26,9	-28,8	-30,4	-25,2	-29,6	2,1	-29,7	-23,7	
	01:03:40	010:04:00	0,210	10,0	-87,8	0,137	-26,9	-28,8	-30,4	-25,2	-29,6	2,1	-29,7	-23,7	
	01:04:43	010:05:03	0,204	10,0	-93,3	0,132	-26,9	-28,8	-30,4	-25,2	-29,6	2,1	-29,7	-23,8	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 34 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	01:05:44	010:06:04	0,208	10,0	-88,6	0,136	-26,9	-28,9	-30,4	-25,2	-29,6	2,1	-29,7	-23,8
	01:06:41	010:07:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	01:06:41	010:07:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	01:06:41	010:07:01	0,202	10,0	-87,3	0,131	-26,9	-28,9	-30,5	-25,1	-29,6	2,2	-29,7	-23,8
	01:07:42	010:08:02	0,207	10,0	-86,5	0,135	-26,9	-28,9	-30,4	-25,1	-29,6	2,2	-29,8	-23,8
	01:08:42	010:09:02	0,204	10,0	-93,4	0,131	-27,0	-28,9	-30,4	-25,1	-29,6	2,2	-29,8	-23,9
	01:09:44	010:10:04	0,204	10,0	-88,0	0,132	-27,0	-28,9	-30,4	-25,1	-29,6	2,2	-29,8	-23,9
	01:10:40	010:11:00	0,208	10,0	-93,0	0,136	-27,0	-28,9	-30,4	-25,1	-29,6	2,2	-29,9	-23,9
	01:11:40	010:12:00	0,202	10,0	-88,6	0,131	-27,0	-28,9	-30,5	-25,0	-29,6	2,3	-29,9	-23,9
	01:12:41	010:13:01	0,207	10,0	-88,9	0,135	-27,0	-28,9	-30,5	-25,0	-29,7	2,3	-29,9	-24,0
	01:13:43	010:14:03	0,205	10,0	-88,6	0,132	-27,0	-28,9	-30,5	-25,0	-29,7	2,3	-29,9	-24,0
	01:14:40	010:15:00	0,206	10,0	-93,5	0,133	-27,1	-28,9	-30,5	-25,0	-29,7	2,3	-29,9	-24,0
	01:15:41	010:16:01	0,210	10,0	-87,8	0,137	-27,1	-28,9	-30,4	-25,0	-29,7	2,3	-29,9	-24,0
	01:16:42	010:17:02	0,203	10,0	-91,5	0,132	-27,1	-28,9	-30,4	-24,9	-29,7	2,3	-29,9	-24,1
	01:17:44	010:18:04	0,207	10,0	-88,8	0,135	-27,1	-28,9	-30,4	-24,9	-29,7	2,3	-30,0	-24,1
	01:18:40	010:19:00	0,201	10,0	-91,4	0,130	-27,1	-28,9	-30,4	-24,9	-29,7	2,4	-30,0	-24,1
	01:19:42	010:20:02	0,204	10,0	-89,3	0,133	-27,1	-28,9	-30,4	-24,9	-29,7	2,4	-30,0	-24,1
	01:20:44	010:21:04	0,210	10,0	-90,4	0,137	-27,1	-29,0	-30,5	-24,9	-29,7	2,4	-30,0	-24,2
	01:21:41	010:22:01	0,203	10,0	-90,1	0,132	-27,2	-29,0	-30,5	-24,9	-29,7	2,4	-30,1	-24,2
	01:22:43	010:23:03	0,208	10,0	-89,5	0,136	-27,2	-29,0	-30,5	-24,8	-29,8	2,4	-30,1	-24,2
	01:23:44	010:24:04	0,200	10,0	-88,3	0,130	-27,2	-29,0	-30,5	-24,8	-29,8	2,4	-30,1	-24,2

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature: Page: 35 / 209
Document issued by:

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:24:40	010:25:00	0,207	10,0	-90,2	0,135	-27,2	-29,0	-30,4	-24,8	-29,8	2,4	-30,1	-24,2
	01:25:41	010:26:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	01:25:41	010:26:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	01:25:41	010:26:01	0,210	10,0	-89,8	0,137	-27,2	-29,0	-30,4	-24,8	-29,8	2,5	-30,1	-24,3
	01:26:43	010:27:03	0,203	10,0	-87,7	0,132	-27,2	-29,0	-30,4	-24,8	-29,8	2,5	-30,1	-24,3
	01:27:44	010:28:04	0,208	10,0	-86,9	0,135	-27,2	-29,0	-30,5	-24,8	-29,8	2,5	-30,2	-24,3
	01:28:40	010:29:00	0,202	10,0	-87,5	0,130	-27,3	-29,0	-30,5	-24,7	-29,8	2,5	-30,2	-24,3
	01:29:41	010:30:01	0,207	10,0	-91,3	0,134	-27,3	-29,0	-30,5	-24,7	-29,8	2,5	-30,2	-24,4
	01:30:40	010:31:00	0,210	10,0	-87,6	0,137	-27,3	-29,0	-30,5	-24,7	-29,8	2,5	-30,2	-24,4
	01:31:42	010:32:02	0,204	10,0	-93,0	0,132	-27,3	-29,0	-30,5	-24,7	-29,9	2,5	-30,2	-24,4
	01:32:43	010:33:03	0,207	10,0	-89,3	0,135	-27,3	-29,1	-30,5	-24,7	-29,9	2,5	-30,3	-24,4
	01:33:40	010:34:00	0,202	10,0	-91,2	0,130	-27,3	-29,1	-30,5	-24,7	-29,9	2,6	-30,3	-24,5
	01:34:41	010:35:01	0,206	10,0	-90,2	0,133	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,3	-24,5
	01:35:40	010:36:00	0,211	10,0	-90,1	0,138	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,3	-24,5
	01:36:41	010:37:01	0,204	10,0	-90,5	0,132	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,3	-24,5
	01:37:41	010:38:01	0,209	10,0	-89,1	0,136	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,4	-24,6
	01:38:40	010:39:00	0,201	10,0	-88,5	0,130	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,4	-24,6
	01:39:41	010:40:01	0,206	10,0	-89,4	0,134	-27,4	-29,1	-30,5	-24,6	-29,9	2,6	-30,4	-24,6
	01:40:40	010:41:00	0,210	10,0	-89,4	0,137	-27,5	-29,1	-30,5	-24,6	-29,9	2,7	-30,4	-24,6
	01:41:40	010:42:00	0,205	10,0	-89,4	0,133	-27,5	-29,1	-30,5	-24,5	-29,9	2,7	-30,4	-24,6
	01:42:42	010:43:02	0,208	10,0	-89,8	0,135	-27,5	-29,1	-30,5	-24,5	-30,0	2,7	-30,4	-24,6

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 36 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	01:43:42	010:44:02	0,202	10,0	-88,6	0,131	-27,5	-29,1	-30,5	-24,5	-30,0	2,7	-30,5	-24,7	
	01:44:44	010:45:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:44:44	010:45:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:44:44	010:45:04	0,205	10,0	-87,4	0,134	-27,5	-29,2	-30,5	-24,5	-30,0	2,7	-30,5	-24,7	
	01:45:43	010:46:03	0,211	10,0	-88,2	0,138	-27,5	-29,2	-30,5	-24,4	-30,0	2,7	-30,5	-24,7	
	01:46:41	010:47:01	0,203	10,0	-90,3	0,131	-27,5	-29,2	-30,5	-24,4	-30,0	2,7	-30,5	-24,7	
	01:47:42	010:48:02	0,208	10,0	-87,4	0,135	-27,6	-29,2	-30,5	-24,4	-30,0	2,7	-30,6	-24,8	
	01:48:42	010:49:02	0,200	10,0	-91,8	0,129	-27,6	-29,2	-30,5	-24,4	-30,0	2,7	-30,6	-24,8	
	01:49:44	010:50:04	0,204	10,0	-89,9	0,133	-27,6	-29,2	-30,5	-24,4	-30,0	2,8	-30,6	-24,8	
	01:50:40	010:51:00	0,210	10,0	-88,9	0,137	-27,6	-29,2	-30,5	-24,3	-30,0	2,8	-30,6	-24,8	
	01:51:42	010:52:02	0,203	10,0	-89,3	0,131	-27,6	-29,2	-30,5	-24,3	-30,1	2,8	-30,6	-24,9	
	01:52:44	010:53:04	0,208	10,0	-89,2	0,135	-27,6	-29,2	-30,5	-24,3	-30,1	2,8	-30,6	-24,9	
	01:53:41	010:54:01	0,212	10,0	-89,9	0,138	-27,6	-29,2	-30,5	-24,3	-30,1	2,8	-30,6	-24,9	
	01:54:42	010:55:02	0,206	10,0	-89,6	0,133	-27,7	-29,2	-30,5	-24,3	-30,1	2,8	-30,7	-24,9	
	01:55:42	010:56:02	0,209	10,0	-89,3	0,136	-27,7	-29,2	-30,5	-24,3	-30,1	2,8	-30,7	-24,9	
	01:56:40	010:57:00	0,203	10,0	-89,8	0,131	-27,7	-29,2	-30,4	-24,3	-30,1	2,9	-30,7	-25,0	
	01:57:45	010:58:05	0,206	10,0	-89,2	0,134	-27,7	-29,2	-30,5	-24,3	-30,1	2,9	-30,7	-25,0	
	01:58:42	010:59:02	0,201	10,0	-88,5	0,129	-27,7	-29,3	-30,5	-24,2	-30,1	2,9	-30,7	-25,0	
	01:59:48	011:00:08	0,205	10,0	-92,5	0,133	-27,7	-29,3	-30,5	-24,2	-30,1	2,9	-30,7	-25,0	
	02:00:44	011:01:04	0,208	10,0	-87,3	0,136	-27,7	-29,3	-30,5	-24,2	-30,1	2,9	-30,8	-25,0	
	02:01:42	011:02:02	0,203	10,0	-91,9	0,130	-27,8	-29,3	-30,5	-24,2	-30,1	2,9	-30,8	-25,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 37 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	02:02:44	011:03:04	0,206	10,0	-89,8	0,134	-27,8	-29,3	-30,5	-24,2	-30,2	2,9	-30,8	-25,1	
	02:03:40	011:04:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:03:40	011:04:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:03:40	011:04:00	0,210	10,0	-90,3	0,137	-27,8	-29,3	-30,5	-24,2	-30,2	2,9	-30,8	-25,1	
	02:04:45	011:05:05	0,203	10,0	-88,6	0,131	-27,8	-29,3	-30,5	-24,2	-30,2	2,9	-30,8	-25,1	
	02:05:45	011:06:05	0,208	10,0	-90,6	0,136	-27,8	-29,3	-30,5	-24,1	-30,2	2,9	-30,9	-25,2	
	02:06:40	011:07:00	0,200	10,0	-88,7	0,130	-27,8	-29,3	-30,5	-24,1	-30,2	2,9	-30,9	-25,2	
	02:07:40	011:08:00	0,206	10,0	-91,6	0,134	-27,8	-29,3	-30,5	-24,1	-30,2	2,9	-30,9	-25,2	
	02:08:42	011:09:02	0,210	10,0	-89,0	0,137	-27,9	-29,3	-30,5	-24,1	-30,2	3,0	-30,9	-25,2	
	02:09:40	011:10:00	0,203	10,0	-90,6	0,132	-27,9	-29,3	-30,5	-24,1	-30,2	3,0	-30,9	-25,2	
	02:10:43	011:11:03	0,208	10,0	-89,7	0,134	-27,9	-29,4	-30,5	-24,1	-30,2	3,0	-31,0	-25,2	
	02:11:45	011:12:05	0,201	10,0	-87,8	0,130	-27,9	-29,4	-30,5	-24,1	-30,2	3,0	-31,0	-25,2	
	02:12:41	011:13:01	0,204	10,0	-90,3	0,133	-27,9	-29,4	-30,5	-24,1	-30,3	3,0	-31,0	-25,3	
	02:13:42	011:14:02	0,211	10,0	-87,1	0,136	-27,9	-29,4	-30,5	-24,1	-30,3	3,0	-31,1	-25,3	
	02:14:42	011:15:02	0,202	10,0	-92,1	0,131	-28,0	-29,4	-30,5	-24,1	-30,3	3,0	-31,2	-25,3	
	02:15:42	011:16:02	0,207	10,0	-88,8	0,134	-28,0	-29,4	-30,5	-24,1	-30,3	3,0	-31,2	-25,3	
	02:16:44	011:17:04	0,204	10,0	-87,6	0,131	-28,0	-29,4	-30,5	-24,0	-30,3	3,0	-31,2	-25,4	
	02:17:42	011:18:02	0,204	10,0	-87,9	0,132	-28,0	-29,4	-30,5	-24,0	-30,3	3,0	-31,3	-25,4	
	02:18:43	011:19:03	0,209	10,0	-93,1	0,136	-28,0	-29,4	-30,5	-24,0	-30,3	3,1	-31,3	-25,4	
	02:19:41	011:20:01	0,202	10,0	-88,1	0,130	-28,0	-29,4	-30,5	-24,0	-30,3	3,1	-31,4	-25,4	
	02:20:42	011:21:02	0,205	10,0	-90,3	0,133	-28,0	-29,4	-30,5	-24,0	-30,3	3,1	-31,4	-25,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 38 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	02:21:43	011:22:03	0,208	10,0	-90,5	0,136	-28,1	-29,5	-30,5	-24,0	-30,3	3,1	-31,4	-25,5	
	02:22:44	011:23:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:22:44	011:23:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:22:44	011:23:04	0,203	10,0	-88,7	0,131	-28,1	-29,5	-30,5	-24,0	-30,3	3,1	-31,5	-25,5	
	02:23:40	011:24:00	0,207	10,0	-88,1	0,134	-28,1	-29,5	-30,5	-23,9	-30,3	3,1	-31,5	-25,5	
	02:24:40	011:25:00	0,200	10,0	-89,6	0,129	-28,1	-29,5	-30,5	-23,9	-30,4	3,1	-31,5	-25,5	
	02:25:42	011:26:02	0,204	10,0	-88,7	0,132	-28,1	-29,5	-30,5	-23,9	-30,4	3,1	-31,6	-25,5	
	02:26:43	011:27:03	0,208	10,0	-89,3	0,136	-28,1	-29,5	-30,5	-23,9	-30,4	3,2	-31,6	-25,6	
	02:27:44	011:28:04	0,201	10,0	-88,6	0,130	-28,2	-29,5	-30,5	-23,9	-30,4	3,2	-31,6	-25,6	
	02:28:42	011:29:02	0,206	10,0	-89,9	0,134	-28,2	-29,5	-30,5	-23,9	-30,4	3,2	-31,7	-25,6	
	02:29:45	011:30:05	0,210	10,0	-89,0	0,137	-28,2	-29,5	-30,5	-23,9	-30,4	3,2	-31,7	-25,6	
	02:30:44	011:31:04	0,203	10,0	-89,6	0,132	-28,2	-29,5	-30,5	-23,8	-30,4	3,2	-31,7	-25,7	
	02:31:44	011:32:04	0,207	10,0	-90,4	0,134	-28,2	-29,5	-30,5	-23,8	-30,4	3,2	-31,8	-25,7	
	02:32:43	011:33:03	0,211	10,0	-88,0	0,138	-28,2	-29,5	-30,5	-23,9	-30,4	3,2	-31,8	-25,7	
	02:33:44	011:34:04	0,203	10,0	-90,3	0,132	-28,2	-29,5	-30,5	-23,9	-30,4	3,2	-31,8	-25,7	
	02:34:44	011:35:04	0,208	10,0	-88,1	0,135	-28,2	-29,6	-30,5	-23,8	-30,4	3,2	-31,9	-25,7	
	02:35:40	011:36:00	0,200	10,0	-91,7	0,130	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-31,9	-25,7	
	02:36:40	011:37:00	0,204	10,0	-89,4	0,132	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-31,9	-25,7	
	02:37:42	011:38:02	0,210	10,0	-87,7	0,137	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-32,0	-25,8	
	02:38:44	011:39:04	0,202	10,0	-88,6	0,131	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-32,0	-25,8	
	02:39:44	011:40:04	0,208	10,0	-91,8	0,135	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-32,0	-25,8	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 39 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	02:40:40	011:41:00	0,210	10,0	-88,8	0,137	-28,3	-29,6	-30,5	-23,8	-30,5	3,3	-32,1	-25,8
	02:41:40	011:42:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	02:41:40	011:42:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	02:41:40	011:42:00	0,204	10,0	-90,6	0,132	-28,3	-29,6	-30,5	-23,7	-30,5	3,3	-32,1	-25,8
	02:42:42	011:43:02	0,207	10,0	-89,3	0,135	-28,4	-29,6	-30,5	-23,7	-30,5	3,3	-32,1	-25,9
	02:43:42	011:44:02	0,201	10,0	-88,9	0,130	-28,4	-29,6	-30,5	-23,8	-30,5	3,3	-32,2	-25,9
	02:44:40	011:45:00	0,204	10,0	-90,4	0,132	-28,4	-29,6	-30,5	-23,8	-30,5	3,3	-32,2	-25,9
	02:45:41	011:46:01	0,209	10,0	-87,2	0,136	-28,4	-29,6	-30,5	-23,7	-30,5	3,4	-32,2	-25,9
	02:46:42	011:47:02	0,202	10,0	-89,6	0,131	-28,4	-29,6	-30,5	-23,7	-30,5	3,4	-32,3	-25,9
	02:47:43	011:48:03	0,207	10,0	-88,2	0,134	-28,4	-29,6	-30,5	-23,7	-30,5	3,4	-32,3	-26,0
	02:48:42	011:49:02	0,210	10,0	-91,8	0,137	-28,4	-29,6	-30,5	-23,7	-30,5	3,4	-32,3	-26,0
	02:49:42	011:50:02	0,203	10,0	-88,4	0,131	-28,5	-29,7	-30,5	-23,7	-30,5	3,4	-32,4	-26,0
	02:50:44	011:51:04	0,207	10,0	-90,0	0,135	-28,5	-29,7	-30,5	-23,7	-30,5	3,4	-32,4	-26,0
	02:51:44	011:52:04	0,209	10,0	-88,3	0,135	-28,5	-29,7	-30,5	-23,7	-30,6	3,4	-32,4	-26,1
	02:52:45	011:53:05	0,204	10,0	-92,8	0,133	-28,5	-29,7	-30,5	-23,7	-30,6	3,4	-32,4	-26,1
	02:53:40	011:54:00	0,208	10,0	-87,9	0,135	-28,5	-29,7	-30,5	-23,7	-30,6	3,4	-32,5	-26,1
	02:54:41	011:55:01	0,201	10,0	-92,0	0,130	-28,6	-29,7	-30,5	-23,7	-30,6	3,4	-32,5	-26,1
	02:55:42	011:56:02	0,204	10,0	-89,6	0,133	-28,6	-29,7	-30,5	-23,6	-30,6	3,4	-32,5	-26,1
	02:56:40	011:57:00	0,210	10,0	-90,3	0,136	-28,6	-29,7	-30,5	-23,6	-30,6	3,4	-32,6	-26,2
	02:57:40	011:58:00	0,202	10,0	-90,4	0,130	-28,6	-29,7	-30,5	-23,6	-30,6	3,4	-32,6	-26,2
	02:58:41	011:59:01	0,207	10,0	-88,3	0,134	-28,6	-29,7	-30,4	-23,6	-30,6	3,4	-32,6	-26,2

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 40 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:59:44	012:00:04	0,210	10,0	-90,3	0,137	-28,7	-29,7	-30,4	-23,6	-30,6	3,5	-32,7	-26,2	
	03:00:40	012:01:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:00:40	012:01:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:00:40	012:01:00	0,204	10,0	-88,4	0,132	-28,7	-29,7	-30,4	-23,6	-30,6	3,5	-32,7	-26,2	
	03:01:44	012:02:04	0,208	10,0	-91,9	0,135	-28,7	-29,7	-30,5	-23,6	-30,6	3,5	-32,7	-26,2	
	03:02:43	012:03:03	0,200	10,0	-87,1	0,130	-28,7	-29,8	-30,5	-23,6	-30,6	3,5	-32,7	-26,3	
	03:03:41	012:04:01	0,204	10,0	-89,3	0,132	-28,7	-29,8	-30,5	-23,6	-30,6	3,5	-32,8	-26,3	
	03:04:44	012:05:04	0,208	10,0	-87,1	0,136	-28,8	-29,8	-30,5	-23,6	-30,7	3,5	-32,8	-26,3	
	03:05:41	012:06:01	0,202	10,0	-92,0	0,130	-28,8	-29,8	-30,5	-23,6	-30,7	3,5	-32,8	-26,3	
	03:06:45	012:07:05	0,206	10,0	-88,8	0,133	-28,8	-29,8	-30,4	-23,6	-30,7	3,5	-32,9	-26,3	
	03:07:42	012:08:02	0,210	10,0	-91,1	0,137	-28,8	-29,8	-30,5	-23,6	-30,7	3,5	-32,9	-26,3	
	03:08:44	012:09:04	0,202	10,0	-89,4	0,131	-28,8	-29,8	-30,5	-23,5	-30,7	3,5	-32,9	-26,4	
	03:09:40	012:10:00	0,207	10,0	-91,1	0,134	-28,8	-29,8	-30,5	-23,5	-30,7	3,5	-32,9	-26,4	
	03:10:41	012:11:01	0,212	10,0	-89,4	0,137	-28,9	-29,8	-30,5	-23,5	-30,7	3,5	-33,0	-26,4	
	03:11:44	012:12:04	0,204	10,0	-89,6	0,132	-28,9	-29,8	-30,4	-23,5	-30,7	3,6	-33,0	-26,4	
	03:12:41	012:13:01	0,207	10,0	-90,6	0,135	-28,9	-29,8	-30,4	-23,5	-30,7	3,6	-33,0	-26,5	
	03:13:44	012:14:04	0,200	10,0	-88,2	0,130	-28,9	-29,8	-30,4	-23,5	-30,7	3,6	-33,0	-26,5	
	03:14:40	012:15:00	0,204	10,0	-92,0	0,132	-28,9	-29,9	-30,4	-23,5	-30,7	3,6	-33,0	-26,5	
	03:15:42	012:16:02	0,208	10,0	-89,9	0,135	-28,9	-29,9	-30,5	-23,5	-30,7	3,6	-33,0	-26,5	
	03:16:43	012:17:03	0,202	10,0	-90,1	0,130	-29,0	-29,9	-30,4	-23,5	-30,8	3,6	-33,0	-26,5	
	03:17:43	012:18:03	0,205	10,0	-91,3	0,133	-29,0	-29,9	-30,4	-23,5	-30,7	3,6	-33,1	-26,5	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 41 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:18:43	012:19:03	0,210	10,0	-87,7	0,136	-29,0	-29,9	-30,4	-23,5	-30,7	3,6	-33,1	-26,6	
	03:19:41	012:20:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:19:41	012:20:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:19:41	012:20:01	0,202	10,0	-90,1	0,131	-29,0	-29,9	-30,4	-23,5	-30,7	3,6	-33,1	-26,6	
	03:20:41	012:21:01	0,208	10,0	-87,5	0,134	-29,0	-29,9	-30,4	-23,5	-30,8	3,6	-33,1	-26,6	
	03:21:42	012:22:02	0,210	10,0	-87,9	0,138	-29,0	-29,9	-30,4	-23,5	-30,8	3,6	-33,1	-26,6	
	03:22:42	012:23:02	0,204	10,0	-86,6	0,132	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,1	-26,6	
	03:23:44	012:24:04	0,208	10,0	-90,8	0,135	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,1	-26,7	
	03:24:40	012:25:00	0,200	10,0	-86,9	0,129	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,7	
	03:25:42	012:26:02	0,204	10,0	-90,2	0,132	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,7	
	03:26:43	012:27:03	0,208	10,0	-89,2	0,135	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,7	
	03:27:40	012:28:00	0,200	10,0	-89,9	0,130	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,8	
	03:28:42	012:29:02	0,204	10,0	-90,7	0,132	-29,1	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,8	
	03:29:43	012:30:03	0,208	10,0	-87,6	0,136	-29,2	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,8	
	03:30:43	012:31:03	0,201	10,0	-88,3	0,130	-29,2	-29,9	-30,4	-23,4	-30,8	3,7	-33,2	-26,8	
	03:31:44	012:32:04	0,207	10,0	-87,7	0,134	-29,2	-30,0	-30,4	-23,4	-30,8	3,7	-33,2	-26,8	
	03:32:44	012:33:04	0,210	10,0	-86,5	0,137	-29,2	-30,0	-30,3	-23,4	-30,8	3,7	-33,2	-26,8	
	03:33:45	012:34:05	0,203	10,0	-87,4	0,131	-29,2	-30,0	-30,3	-23,4	-30,8	3,8	-33,2	-26,8	
	03:34:42	012:35:02	0,206	10,0	-92,7	0,134	-29,2	-30,0	-30,4	-23,4	-30,8	3,8	-33,3	-26,9	
	03:35:42	012:36:02	0,210	10,0	-87,2	0,137	-29,2	-30,0	-30,4	-23,4	-30,9	3,8	-33,3	-26,9	
	03:36:42	012:37:02	0,203	10,0	-92,3	0,131	-29,2	-30,0	-30,4	-23,3	-30,9	3,8	-33,3	-26,9	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 42 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	03:37:40	012:38:00	0,206	10,0	-89,2	0,134	-29,2	-30,0	-30,4	-23,4	-30,9	3,8	-33,3	-26,9
	03:38:43	012:39:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	03:38:43	012:39:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	03:38:43	012:39:03	0,208	10,0	-90,8	0,135	-29,3	-30,0	-30,4	-23,4	-30,9	3,8	-33,3	-27,0
	03:39:43	012:40:03	0,204	10,0	-89,0	0,132	-29,3	-30,0	-30,4	-23,3	-30,9	3,8	-33,3	-27,0
	03:40:44	012:41:04	0,208	10,0	-90,6	0,136	-29,3	-30,0	-30,4	-23,3	-30,9	3,8	-33,3	-27,0
	03:41:40	012:42:00	0,201	10,0	-88,7	0,130	-29,3	-30,0	-30,4	-23,3	-30,9	3,8	-33,3	-27,0
	03:42:43	012:43:03	0,205	10,0	-90,7	0,133	-29,3	-30,0	-30,4	-23,3	-30,9	3,8	-33,4	-27,0
	03:43:43	012:44:03	0,208	10,0	-89,1	0,136	-29,3	-30,1	-30,4	-23,3	-30,9	3,9	-33,4	-27,1
	03:44:44	012:45:04	0,202	10,0	-89,8	0,130	-29,4	-30,1	-30,4	-23,3	-30,9	3,9	-33,4	-27,1
	03:45:40	012:46:00	0,204	10,0	-90,2	0,132	-29,4	-30,1	-30,3	-23,3	-30,9	3,9	-33,4	-27,1
	03:46:40	012:47:00	0,209	10,0	-88,6	0,136	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,1
	03:47:42	012:48:02	0,201	10,0	-92,4	0,130	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,1
	03:48:42	012:49:02	0,205	10,0	-88,3	0,133	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,1
	03:49:43	012:50:03	0,210	10,0	-88,3	0,137	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,2
	03:50:44	012:51:04	0,202	10,0	-87,8	0,131	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,2
	03:51:40	012:52:00	0,206	10,0	-91,2	0,134	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,2
	03:52:41	012:53:01	0,210	10,0	-87,0	0,137	-29,4	-30,1	-30,3	-23,3	-31,0	3,9	-33,5	-27,2
	03:53:42	012:54:02	0,202	10,0	-85,1	0,131	---	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,2
	03:54:42	012:55:02	0,206	10,0	-86,4	0,134	-29,5	-30,1	-30,3	-23,3	-31,0	3,9	-33,5	-27,3
	03:55:44	012:56:04	0,210	10,0	-93,2	0,137	-29,5	-30,1	-30,3	-23,3	-31,0	3,9	-33,4	-27,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 43 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:56:43	012:57:03	0,203	10,0	-88,6	0,131	-29,5	-30,1	-30,3	-23,3	-31,0	3,9	-33,5	-27,3	
	03:57:40	012:58:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:57:40	012:58:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:57:40	012:58:00	0,206	10,0	-91,2	0,134	-29,5	-30,2	-30,3	-23,3	-31,0	4,0	-33,5	-27,3	
	03:58:40	012:59:00	0,210	10,0	-89,8	0,137	-29,5	-30,2	-30,3	-23,3	-31,0	4,0	-33,5	-27,3	
	03:59:40	013:00:00	0,203	10,0	-90,7	0,132	-29,5	-30,2	-30,3	-23,2	-31,0	4,0	-33,5	-27,4	
	04:00:42	013:01:02	0,207	10,0	-89,0	0,134	---	-30,2	-30,3	-23,2	-31,0	4,0	-33,5	-27,4	
	04:01:43	013:02:03	0,200	10,0	-90,8	0,130	-29,6	-30,2	-30,3	-23,2	-31,0	4,0	-33,5	-27,4	
	04:02:43	013:03:03	0,203	10,0	-88,8	0,132	-29,6	-30,2	-30,3	-23,2	-31,0	4,0	-33,5	-27,4	
	04:03:43	013:04:03	0,208	10,0	-90,4	0,135	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,5	-27,4	
	04:04:44	013:05:04	0,205	10,0	-88,9	0,133	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,5	-27,5	
	04:05:40	013:06:00	0,205	10,0	-91,2	0,133	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,5	-27,5	
	04:06:41	013:07:01	0,207	10,0	-89,5	0,135	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,5	-27,5	
	04:07:41	013:08:01	0,200	10,0	-88,5	0,130	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,6	-27,5	
	04:08:42	013:09:02	0,204	10,0	-89,5	0,132	-29,6	-30,2	-30,3	-23,2	-31,1	4,0	-33,6	-27,5	
	04:09:43	013:10:03	0,209	10,0	-88,3	0,136	-29,6	-30,2	-30,3	-23,2	-31,1	4,1	-33,6	-27,5	
	04:10:41	013:11:01	0,201	10,0	-88,9	0,130	-29,7	-30,3	-30,3	-23,2	-31,1	4,1	-33,6	-27,5	
	04:11:42	013:12:02	0,205	10,0	-88,3	0,133	-29,7	-30,3	-30,2	-23,2	-31,1	4,1	-33,6	-27,5	
	04:12:43	013:13:03	0,210	10,0	-92,5	0,137	-29,7	-30,3	-30,2	-23,2	-31,1	4,1	-33,6	-27,6	
	04:13:44	013:14:04	0,203	10,0	-87,7	0,131	-29,7	-30,3	-30,2	-23,2	-31,1	4,1	-33,6	-27,6	
	04:14:44	013:15:04	0,207	10,0	-91,7	0,134	-29,7	-30,3	-30,3	-23,2	-31,1	4,1	-33,6	-27,6	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 44 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	04:15:41	013:16:01	0,210	10,0	-87,5	0,137	-29,7	-30,3	-30,3	-23,1	-31,1	4,1	-33,6	-27,6	
	04:16:43	013:17:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:16:43	013:17:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:16:43	013:17:03	0,202	10,0	-92,8	0,131	-29,7	-30,3	-30,2	-23,1	-31,2	4,1	-33,6	-27,6	
	04:17:40	013:18:00	0,206	10,0	-88,5	0,133	-29,8	-30,3	-30,2	-23,1	-31,2	4,1	-33,6	-27,6	
	04:18:43	013:19:03	0,210	10,0	-90,5	0,137	-29,8	-30,3	-30,2	-23,1	-31,2	4,1	-33,6	-27,7	
	04:19:40	013:20:00	0,201	10,0	-89,4	0,130	-29,7	-30,3	-30,2	-23,1	-31,2	4,1	-33,6	-27,7	
	04:20:40	013:21:00	0,206	10,0	-90,9	0,134	-29,7	-30,3	-30,2	-23,1	-31,2	4,1	-33,6	-27,7	
	04:21:43	013:22:03	0,210	10,0	-89,6	0,137	-29,7	-30,3	-30,2	-23,1	-31,2	4,2	-33,6	-27,7	
	04:22:43	013:23:03	0,204	10,0	-88,3	0,131	-29,8	-30,3	-30,2	-23,1	-31,2	4,2	-33,7	-27,7	
	04:23:44	013:24:04	0,206	10,0	-90,7	0,134	---	-30,3	-30,2	-23,1	-31,2	4,2	-33,6	-27,7	
	04:24:44	013:25:04	0,210	10,0	-86,8	0,137	---	-30,3	-30,2	-23,1	-31,2	4,2	-33,7	-27,7	
	04:25:40	013:26:00	0,202	10,0	-88,7	0,131	-29,7	-30,4	-30,2	-23,1	-31,2	4,2	-33,7	-27,7	
	04:26:40	013:27:00	0,207	10,0	-87,9	0,135	-29,7	-30,4	-30,2	-23,1	-31,2	4,2	-33,7	-27,8	
	04:27:42	013:28:02	0,210	10,0	-90,0	0,138	-29,8	-30,4	-30,2	-23,1	-31,2	4,2	-33,7	-27,8	
	04:28:43	013:29:03	0,203	10,0	-87,4	0,132	-29,8	-30,4	-30,2	-23,1	-31,2	4,2	-33,7	-27,8	
	04:29:40	013:30:00	0,206	10,0	-88,5	0,134	-29,9	-30,4	-30,2	-23,1	-31,3	4,2	-33,7	-27,8	
	04:30:40	013:31:00	0,211	10,0	-88,1	0,138	-29,9	-30,4	-30,2	-23,1	-31,3	4,2	-33,7	-27,8	
	04:31:41	013:32:01	0,203	10,0	-89,9	0,131	-29,9	-30,4	-30,2	-23,1	-31,3	4,2	-33,7	-27,8	
	04:32:42	013:33:02	0,208	10,0	-87,7	0,135	-29,9	-30,4	-30,2	-23,1	-31,3	4,3	-33,7	-27,8	
	04:33:44	013:34:04	0,211	10,0	-92,9	0,138	-29,9	-30,4	-30,2	-23,1	-31,3	4,2	-33,7	-27,9	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 45 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	04:34:44	013:35:04	0,203	10,0	-88,4	0,131	-29,9	-30,4	-30,2	-23,1	-31,3	4,3	-33,7	-27,9	
	04:35:44	013:36:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:35:44	013:36:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:35:44	013:36:04	0,207	10,0	-88,0	0,135	-29,9	-30,4	-30,2	-23,1	-31,3	4,3	-33,7	-27,9	
	04:36:41	013:37:01	0,211	10,0	-88,2	0,138	-29,9	-30,5	-30,2	-23,1	-31,3	4,3	-33,8	-27,9	
	04:37:42	013:38:02	0,203	10,0	-91,4	0,132	-30,0	-30,5	-30,1	-23,1	-31,3	4,3	-33,8	-27,9	
	04:38:44	013:39:04	0,207	10,0	-89,6	0,134	-29,9	-30,5	-30,1	-23,1	-31,3	4,3	-33,8	-27,9	
	04:39:44	013:40:04	0,207	10,0	-90,5	0,133	-29,9	-30,5	-30,1	-23,1	-31,3	4,3	-33,8	-27,9	
	04:40:45	013:41:05	0,204	10,0	-89,0	0,132	-29,9	-30,5	-30,1	-23,0	-31,3	4,3	-33,8	-28,0	
	04:41:41	013:42:01	0,208	10,0	-89,8	0,135	-30,0	-30,5	-30,1	-23,0	-31,3	4,3	-33,8	-28,0	
	04:42:40	013:43:00	0,210	10,0	-89,9	0,137	-29,9	-30,5	-30,1	-23,0	-31,4	4,3	-33,8	-28,0	
	04:43:41	013:44:01	0,203	10,0	-88,8	0,132	-29,9	-30,5	-30,0	-23,1	-31,4	4,3	-33,8	-28,0	
	04:44:42	013:45:02	0,207	10,0	-88,6	0,134	-30,0	-30,5	-30,0	-23,1	-31,4	4,4	-33,8	-28,0	
	04:45:41	013:46:01	0,207	10,0	-89,5	0,134	-30,0	-30,5	-30,0	-23,1	-31,4	4,4	-33,8	-28,0	
	04:46:44	013:47:04	0,203	10,0	-88,4	0,131	-30,0	-30,5	-30,0	-23,0	-31,4	4,4	-33,8	-28,0	
	04:47:45	013:48:05	0,208	10,0	-88,1	0,135	-30,0	-30,5	-29,9	-23,0	-31,4	4,4	-33,8	-28,0	
	04:48:40	013:49:00	0,210	10,0	-90,9	0,137	-30,0	-30,6	-29,9	-23,0	-31,4	4,4	-33,9	-28,0	
	04:49:40	013:50:00	0,202	10,0	-88,0	0,132	-30,0	-30,5	-29,9	-23,0	-31,4	4,4	-33,8	-28,0	
	04:50:44	013:51:04	0,207	10,0	-86,3	0,134	-30,1	-30,6	-29,9	-23,0	-31,4	4,4	-33,8	-28,1	
	04:51:41	013:52:01	0,210	10,0	-88,3	0,136	-30,1	-30,6	-29,8	-23,0	-31,4	4,4	-33,8	-28,1	
	04:52:41	013:53:01	0,202	10,0	-88,1	0,131	-30,0	-30,6	-29,8	-23,0	-31,4	4,4	-33,8	-28,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus			
Document issued by:		Signature:			Page: 46 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	04:53:42	013:54:02	0,207	10,0	-88,2	0,134	-30,0	-30,6	-29,8	-23,0	-31,4	4,4	-33,9	-28,1	
	04:54:44	013:55:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:54:44	013:55:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:54:44	013:55:04	0,210	10,0	-92,9	0,137	-30,0	-30,6	-29,7	-23,0	-31,5	4,4	-33,9	-28,1	
	04:55:40	013:56:00	0,203	10,0	-87,4	0,131	-30,0	-30,6	-29,7	-23,0	-31,4	4,4	-33,9	-28,1	
	04:56:42	013:57:02	0,206	10,0	-85,4	0,134	-30,1	-30,6	-29,6	-23,0	-31,4	4,4	-33,9	-28,1	
	04:57:42	013:58:02	0,210	10,0	-87,4	0,137	-30,1	-30,6	-29,6	-23,0	-31,4	4,4	-33,9	-28,1	
	04:58:40	013:59:00	0,203	10,0	-91,6	0,131	-30,1	-30,6	-29,5	-23,0	-31,4	4,5	-33,9	-28,1	
	04:59:41	014:00:01	0,206	10,0	-90,0	0,133	-30,1	-30,6	-29,4	-23,0	-31,4	4,5	-33,9	-28,2	
	05:00:41	014:01:01	0,210	10,0	-89,8	0,137	-30,1	-30,6	-29,4	-23,0	-31,4	4,5	-33,9	-28,2	
	05:01:42	014:02:02	0,202	10,0	-88,8	0,131	-30,1	-30,6	-29,4	-23,0	-31,5	4,5	-33,9	-28,2	
	05:02:42	014:03:02	0,207	10,0	-87,7	0,134	-30,1	-30,6	-29,3	-23,0	-31,5	4,5	-33,9	-28,2	
	05:03:43	014:04:03	0,210	10,0	-91,7	0,137	-30,1	-30,7	-29,3	-23,0	-31,5	4,5	-33,9	-28,2	
	05:04:40	014:05:00	0,201	10,0	-88,1	0,130	-30,2	-30,7	-29,3	-23,0	-31,5	4,5	-33,9	-28,2	
	05:05:42	014:06:02	0,205	10,0	-87,1	0,133	-30,2	-30,7	-29,3	-23,0	-31,5	4,5	-33,9	-28,2	
	05:06:43	014:07:03	0,209	10,0	-86,6	0,137	-30,2	-30,7	-29,2	-23,0	-31,5	4,5	-33,9	-28,2	
	05:07:43	014:08:03	0,202	10,0	-93,7	0,131	-30,2	-30,7	-29,2	-23,0	-31,5	4,5	-33,9	-28,2	
	05:08:40	014:09:00	0,205	10,0	-87,0	0,133	-30,2	-30,7	-29,1	-23,0	-31,5	4,6	-33,9	-28,2	
	05:09:40	014:10:00	0,209	10,0	-89,8	0,136	-30,2	-30,7	-29,1	-22,9	-31,5	4,5	-33,9	-28,3	
	05:10:40	014:11:00	0,201	10,0	-86,7	0,130	-30,2	-30,7	-29,1	-22,9	-31,5	4,6	-33,9	-28,3	
	05:11:41	014:12:01	0,205	10,0	-87,6	0,133	-30,2	-30,7	-29,1	-22,9	-31,5	4,5	-33,9	-28,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 47 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



05:12:40	014:13:00	0,208	10,0	-87,2	0,136	-30,2	-30,7	-29,1	-22,9	-31,6	4,6	-34,0	-28,3
05:13:42	014:14:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
05:13:42	014:14:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
05:13:42	014:14:02	0,201	10,0	-87,3	0,130	-30,2	-30,7	-29,0	-22,9	-31,5	4,6	-34,0	-28,3
05:14:43	014:15:03	0,205	10,0	-87,3	0,133	-30,2	-30,7	-29,0	-22,9	-31,5	4,6	-34,0	-28,3
05:15:41	014:16:01	0,207	10,0	-92,7	0,135	-30,2	-30,7	-29,0	-22,9	-31,5	4,6	-34,0	-28,3
05:16:44	014:17:04	0,200	10,0	-89,1	0,129	-30,3	-30,7	-29,0	-22,9	-31,5	4,6	-34,0	-28,3
05:17:41	014:18:01	0,204	10,0	-85,9	0,132	-30,2	-30,8	-28,9	-22,9	-31,5	4,6	-34,0	-28,3
05:18:43	014:19:03	0,207	10,0	-88,6	0,135	-30,2	-30,8	-28,9	-22,9	-31,5	4,6	-34,0	-28,3
05:19:44	014:20:04	0,200	10,0	-88,6	0,129	-30,3	-30,8	-28,9	-22,9	-31,6	4,6	-34,0	-28,3
05:20:40	014:21:00	0,203	10,0	-88,1	0,132	-30,3	-30,8	-28,9	-22,9	-31,5	4,6	-34,0	-28,4
05:21:43	014:22:03	0,207	10,0	-90,7	0,135	-30,3	-30,8	-28,9	-22,9	-31,5	4,6	-34,0	-28,4
05:22:43	014:23:03	0,211	10,0	-87,6	0,138	-30,3	-30,8	-28,8	-22,9	-31,6	4,6	-34,0	-28,4
05:23:41	014:24:01	0,203	10,0	-93,3	0,132	-30,3	-30,8	-28,8	-22,9	-31,6	4,7	-34,0	-28,4
05:24:40	014:25:00	0,207	10,0	-88,1	0,134	-30,3	-30,8	-28,8	-22,9	-31,6	4,7	-34,0	-28,4
05:25:41	014:26:01	0,210	10,0	-89,2	0,137	-30,3	-30,8	-28,8	-22,9	-31,6	4,7	-34,0	-28,4
05:26:42	014:27:02	0,202	10,0	-88,2	0,131	-30,3	-30,8	-28,7	-22,9	-31,6	4,7	-34,0	-28,4
05:27:40	014:28:00	0,206	10,0	-91,6	0,134	-30,4	-30,8	-28,7	-22,9	-31,6	4,7	-34,0	-28,4
05:28:41	014:29:01	0,208	10,0	-89,5	0,136	-30,4	-30,9	-28,7	-22,9	-31,6	4,7	-34,0	-28,5
05:29:43	014:30:03	0,202	10,0	-91,0	0,131	-30,3	-30,9	-28,7	-22,9	-31,6	4,7	-34,0	-28,4
05:30:43	014:31:03	0,205	10,0	-89,4	0,133	-30,4	-30,9	-28,6	-22,9	-31,6	4,7	-34,0	-28,5

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 48 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:31:43	014:32:03	0,210	10,0	-89,8	0,137	-30,4	-30,9	-28,6	-22,9	-31,6	4,7	-34,0	-28,5
	05:32:40	014:33:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	05:32:40	014:33:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	05:32:40	014:33:00	0,201	10,0	-89,7	0,130	-30,4	-30,9	-28,6	-22,9	-31,7	4,7	-34,0	-28,4
	05:33:43	014:34:03	0,205	10,0	-88,0	0,133	-30,4	-30,9	-28,6	-22,9	-31,7	4,7	-34,0	-28,4
	05:34:44	014:35:04	0,208	10,0	-87,1	0,136	-30,4	-30,9	-28,5	-22,8	-31,8	4,7	-34,0	-28,5
	05:35:40	014:36:00	0,202	10,0	-88,9	0,130	-30,4	-30,9	-28,5	-22,8	-31,8	4,7	-34,0	-28,5
	05:36:42	014:37:02	0,204	10,0	-88,4	0,132	-30,4	-30,9	-28,5	-22,8	-31,8	4,7	-34,0	-28,5
	05:37:42	014:38:02	0,208	10,0	-87,6	0,135	-30,4	-30,9	-28,5	-22,8	-31,9	4,7	-34,1	-28,5
	05:38:43	014:39:03	0,203	10,0	-88,1	0,131	-30,4	-30,9	-28,4	-22,8	-31,9	4,8	-34,1	-28,5
	05:39:43	014:40:03	0,205	10,0	-87,6	0,132	-30,4	-30,9	-28,4	-22,8	-31,9	4,8	-34,0	-28,5
	05:40:44	014:41:04	0,207	10,0	-86,4	0,135	-30,4	-30,9	-28,4	-22,8	-31,9	4,8	-34,0	-28,5
	05:41:40	014:42:00	0,204	10,0	-88,3	0,130	-30,4	-30,9	-28,3	-22,8	-31,9	4,8	-34,0	-28,5
	05:42:41	014:43:01	0,203	10,0	-91,0	0,131	-30,4	-30,9	-28,3	-22,8	-31,9	4,8	-34,0	-28,5
	05:43:43	014:44:03	0,206	10,0	-87,7	0,134	-30,4	-30,9	-28,2	-22,8	-31,9	4,8	-34,0	-28,5
	05:44:44	014:45:04	0,210	10,0	-88,4	0,137	-30,4	-30,9	-28,2	-22,8	-31,9	4,8	-34,1	-28,5
	05:45:41	014:46:01	0,201	10,0	-88,9	0,130	-30,4	-31,0	-28,2	-22,8	-31,9	4,8	-34,0	-28,5
	05:46:44	014:47:04	0,205	10,0	-89,9	0,133	-30,4	-31,0	-28,1	-22,8	-31,9	4,8	-34,0	-28,6
	05:47:41	014:48:01	0,210	10,0	-89,1	0,136	-30,4	-31,0	-28,1	-22,8	-31,9	4,8	-34,1	-28,6
	05:48:40	014:49:00	0,201	10,0	-85,9	0,130	-30,4	-31,0	-28,0	-22,8	-31,9	4,8	-34,0	-28,6
	05:49:42	014:50:02	0,204	10,0	-89,1	0,132	-30,4	-30,9	-27,9	-22,8	-31,9	4,8	-34,0	-28,5

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 49 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:50:44	014:51:04	0,208	10,0	-90,6	0,135	-30,5	-31,0	-27,8	-22,8	-31,9	4,8	-34,1	-28,6	
	05:51:41	014:52:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:51:41	014:52:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:51:41	014:52:01	0,211	10,0	-89,3	0,138	-30,5	-31,0	-27,8	-22,8	-31,9	4,8	-34,1	-28,6	
	05:52:42	014:53:02	0,204	10,0	-85,7	0,132	-30,5	-31,0	-27,7	-22,8	-31,9	4,9	-34,0	-28,6	
	05:53:42	014:54:02	0,208	10,0	-88,6	0,134	-30,5	-31,0	-27,6	-22,8	-31,9	4,9	-34,0	-28,6	
	05:54:44	014:55:04	0,210	10,0	-92,2	0,138	-30,5	-31,0	-27,5	-22,8	-31,9	4,9	-34,1	-28,6	
	05:55:40	014:56:00	0,202	10,0	-89,1	0,131	-30,5	-31,0	-27,4	-22,8	-31,9	4,9	-34,0	-28,6	
	05:56:43	014:57:03	0,205	10,0	-88,5	0,133	-30,5	-31,0	-27,4	-22,8	-31,9	4,9	-34,1	-28,6	
	05:57:44	014:58:04	0,209	10,0	-92,6	0,136	-30,5	-31,0	-27,3	-22,8	-31,9	4,9	-34,1	-28,6	
	05:58:40	014:59:00	0,200	10,0	-88,8	0,129	-30,5	-31,0	-27,2	-22,7	-31,9	4,9	-34,1	-28,6	
	05:59:40	015:00:00	0,203	10,0	-85,5	0,132	-30,5	-31,0	-27,1	-22,7	-31,9	4,9	-34,1	-28,6	
	06:00:41	015:01:01	0,208	10,0	-87,3	0,135	-30,5	-31,0	-27,0	-22,7	-31,9	4,9	-34,1	-28,6	
	06:01:41	015:02:01	0,208	10,0	-86,5	0,134	-30,5	-31,0	-27,0	-22,7	-31,9	4,9	-34,1	-28,6	
	06:02:42	015:03:02	0,203	10,0	-87,1	0,132	-30,5	-31,0	-26,9	-22,7	-31,9	4,9	-34,0	-28,6	
	06:03:42	015:04:02	0,206	10,0	-91,9	0,134	-30,5	-31,0	-26,8	-22,7	-31,9	4,9	-34,1	-28,7	
	06:04:42	015:05:02	0,209	10,0	-89,5	0,136	-30,5	-31,1	-26,7	-22,7	-31,9	5,0	-34,1	-28,7	
	06:05:40	015:06:00	0,201	10,0	-83,8	0,130	-30,5	-31,1	-26,6	-22,7	-31,9	5,0	-34,0	-28,7	
	06:06:44	015:07:04	0,205	10,0	-89,3	0,133	-30,5	-31,1	-26,5	-22,7	-31,9	5,0	-34,0	-28,7	
	06:07:44	015:08:04	0,210	10,0	-89,6	0,136	-30,5	-31,1	-26,4	-22,7	-31,9	5,0	-34,0	-28,7	
	06:08:42	015:09:02	0,204	10,0	-90,4	0,132	-30,5	-31,1	-26,3	-22,7	-31,9	5,0	-34,0	-28,7	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 50 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:09:44	015:10:04	0,203	10,0	-89,0	0,132	-30,5	-31,1	-26,1	-22,7	-31,9	5,0	-34,0	-28,7	
	06:10:41	015:11:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:10:41	015:11:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:10:41	015:11:01	0,206	10,0	-89,0	0,134	-30,5	-31,1	-26,0	-22,7	-31,9	5,0	-34,0	-28,7	
	06:11:41	015:12:01	0,210	10,0	-88,0	0,137	-30,6	-31,1	-25,9	-22,7	-31,9	5,0	-34,0	-28,7	
	06:12:42	015:13:02	0,202	10,0	-86,9	0,131	-30,6	-31,1	-25,8	-22,7	-31,9	5,0	-34,0	-28,7	
	06:13:40	015:14:00	0,206	10,0	-88,0	0,134	-30,6	-31,1	-25,6	-22,7	-31,9	5,0	-34,0	-28,7	
	06:14:44	015:15:04	0,208	10,0	-88,5	0,136	-30,6	-31,1	-25,5	-22,7	-31,9	5,0	-34,0	-28,7	
	06:15:40	015:16:00	0,201	10,0	-87,7	0,130	-30,6	-31,1	-25,4	-22,7	-31,9	5,0	-34,0	-28,7	
	06:16:42	015:17:02	0,204	10,0	-88,7	0,132	-30,6	-31,1	-25,2	-22,7	-31,9	5,0	-34,0	-28,7	
	06:17:45	015:18:05	0,208	10,0	-87,6	0,135	-30,6	-31,1	-25,1	-22,7	-31,9	5,0	-34,1	-28,8	
	06:18:43	015:19:03	0,211	10,0	-88,1	0,138	-30,6	-31,1	-25,0	-22,7	-32,0	5,1	-34,1	-28,8	
	06:19:40	015:20:00	0,202	10,0	-87,8	0,131	-30,6	-31,1	-24,8	-22,7	-31,9	5,1	-34,0	-28,7	
	06:20:42	015:21:02	0,205	10,0	-87,1	0,133	-30,6	-31,1	-24,7	-22,7	-31,9	5,1	-34,0	-28,7	
	06:21:43	015:22:03	0,210	10,0	-87,0	0,137	-30,6	-31,2	-24,6	-22,7	-32,0	5,1	-34,1	-28,8	
	06:22:44	015:23:04	0,201	10,0	-86,2	0,130	-30,6	-31,2	-24,5	-22,7	-32,0	5,1	-34,0	-28,8	
	06:23:40	015:24:00	0,205	10,0	-88,1	0,133	-30,6	-31,2	-24,3	-22,7	-32,0	5,1	-34,0	-28,8	
	06:24:40	015:25:00	0,208	10,0	-86,8	0,135	-30,6	-31,2	-24,2	-22,7	-32,0	5,1	-34,0	-28,8	
	06:25:43	015:26:03	0,200	10,0	-88,9	0,130	-30,6	-31,2	-24,0	-22,7	-32,0	5,1	-34,0	-28,8	
	06:26:40	015:27:00	0,203	10,0	-88,2	0,131	-30,6	-31,2	-23,9	-22,6	-32,0	5,1	-34,0	-28,8	
	06:27:40	015:28:00	0,207	10,0	-89,5	0,134	-30,6	-31,2	-23,8	-22,6	-32,0	5,1	-34,0	-28,8	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 51 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



06:28:41	015:29:01	0,209	10,0	-89,8	0,136	-30,6	-31,2	-23,6	-22,7	-32,0	5,1	-34,1	-28,8
06:29:42	015:30:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
06:29:42	015:30:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
06:29:42	015:30:02	0,201	10,0	-88,5	0,130	-30,6	-31,2	-23,4	-22,7	-32,0	5,1	-34,0	-28,8
06:30:44	015:31:04	0,205	10,0	-88,5	0,132	-30,6	-31,2	-23,2	-22,6	-32,0	5,1	-34,0	-28,8
06:31:40	015:32:00	0,208	10,0	-88,9	0,135	-30,6	-31,2	-23,1	-22,6	-32,0	5,1	-34,0	-28,8
06:32:42	015:33:02	0,211	10,0	-87,2	0,138	-30,6	-31,2	-22,9	-22,6	-32,0	5,2	-34,0	-28,8
06:33:43	015:34:03	0,203	10,0	-89,0	0,132	-30,6	-31,2	-22,7	-22,6	-32,0	5,2	-34,0	-28,8
06:34:40	015:35:00	0,207	10,0	-88,1	0,133	-30,6	-31,2	-22,4	-22,6	-32,0	5,2	-34,0	-28,8
06:35:41	015:36:01	0,210	10,0	-89,4	0,137	-30,6	-31,2	-22,3	-22,6	-32,0	5,2	-34,0	-28,8
06:36:40	015:37:00	0,201	10,0	-89,0	0,130	-30,7	-31,2	-22,1	-22,6	-32,0	5,2	-34,0	-28,8
06:37:42	015:38:02	0,205	10,0	-88,5	0,132	-30,6	-31,3	-21,9	-22,6	-32,0	5,2	-34,0	-28,8
06:38:43	015:39:03	0,207	10,0	-88,1	0,135	-30,6	-31,3	-21,7	-22,6	-32,0	5,2	-34,0	-28,8
06:39:43	015:40:03	0,208	10,0	-88,1	0,135	-30,6	-31,3	-21,5	-22,6	-32,0	5,2	-34,0	-28,8
06:40:44	015:41:04	0,202	10,0	-87,8	0,131	-30,6	-31,3	-21,3	-22,6	-32,0	5,2	-34,0	-28,8
06:41:43	015:42:03	0,206	10,0	-88,4	0,134	-30,6	-31,3	-21,0	-22,6	-32,1	5,2	-34,0	-28,8
06:42:40	015:43:00	0,208	10,0	-88,4	0,136	-30,6	-31,3	-20,8	-22,6	-32,1	5,2	-34,0	-28,9
06:43:43	015:44:03	0,200	10,0	-87,7	0,130	-30,6	-31,3	-20,6	-22,6	-32,1	5,2	-34,0	-28,8
06:44:41	015:45:01	0,203	10,0	-87,3	0,132	-30,6	-31,3	-20,4	-22,6	-32,1	5,2	-34,0	-28,8
06:45:43	015:46:03	0,207	10,0	-87,6	0,135	-30,6	-31,3	-20,2	-22,6	-32,1	5,2	-34,0	-28,9
06:46:44	015:47:04	0,210	10,0	-86,3	0,137	-30,6	-31,3	-19,9	-22,6	-32,2	5,2	-34,0	-28,9

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 52 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	06:47:40	015:48:00	0,202	10,0	-88,5	0,131	-30,6	-31,3	-19,7	-22,6	-32,2	5,3	-33,9	-28,8	
	06:48:42	015:49:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:48:42	015:49:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:48:42	015:49:02	0,204	10,0	-87,1	0,133	-30,6	-31,3	-19,4	-22,6	-32,2	5,3	-33,9	-28,8	
	06:49:42	015:50:02	0,208	10,0	-89,1	0,136	-30,6	-31,4	-19,2	-22,6	-32,2	5,3	-34,0	-28,9	
	06:50:43	015:51:03	0,212	10,0	-88,6	0,137	-30,6	-31,4	-19,0	-22,6	-32,3	5,3	-33,9	-28,8	
	06:51:40	015:52:00	0,202	10,0	-88,0	0,131	-30,6	-31,4	-18,7	-22,5	-32,3	5,3	-33,9	-28,8	
	06:52:44	015:53:04	0,206	10,0	-93,1	0,134	-30,6	-31,4	-18,4	-22,5	-32,3	5,3	-33,9	-28,8	
	06:53:40	015:54:00	0,209	10,0	-89,2	0,136	-30,6	-31,4	-18,1	-22,5	-32,4	5,3	-33,9	-28,9	
	06:54:42	015:55:02	0,200	10,0	-84,0	0,130	-30,6	-31,4	-17,8	-22,5	-32,4	5,3	-33,9	-28,8	
	06:55:43	015:56:03	0,204	10,0	-86,7	0,132	-30,6	-31,4	-17,5	-22,5	-32,4	5,3	-33,9	-28,8	
	06:56:41	015:57:01	0,207	10,0	-85,1	0,135	-30,6	-31,4	-17,2	-22,5	-32,4	5,3	-33,9	-28,8	
	06:57:41	015:58:01	0,211	10,0	-88,2	0,138	-30,6	-31,4	-16,9	-22,5	-32,5	5,3	-33,8	-28,8	
	06:58:45	015:59:05	0,202	10,0	-86,8	0,131	-30,6	-31,4	-16,6	-22,5	-32,5	5,3	-33,8	-28,8	
	06:59:43	016:00:03	0,206	10,0	-89,0	0,134	-30,6	-31,4	-16,3	-22,5	-32,6	5,3	-33,8	-28,8	
	07:00:41	016:01:01	0,208	10,0	-87,6	0,135	-30,6	-31,5	-16,0	-22,5	-32,6	5,3	-33,8	-28,8	
	07:01:42	016:02:02	0,200	10,0	-89,1	0,129	-30,6	-31,4	-15,7	-22,5	-32,6	5,3	-33,8	-28,8	
	07:02:43	016:03:03	0,202	10,0	-87,9	0,131	-30,6	-31,4	-15,4	-22,5	-32,7	5,3	-33,7	-28,7	
	07:03:43	016:04:03	0,206	10,0	-89,4	0,134	-30,7	-31,5	-15,1	-22,5	-32,8	5,4	-33,7	-28,8	
	07:04:44	016:05:04	0,209	10,0	-88,2	0,136	-30,7	-31,5	-14,8	-22,5	-32,8	5,4	-33,7	-28,8	
	07:05:44	016:06:04	0,201	10,0	-90,0	0,130	-30,7	-31,5	-14,4	-22,5	-32,9	5,4	-33,6	-28,7	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 53 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:06:44	016:07:04	0,203	10,0	-88,9	0,132	-30,7	-31,5	-14,0	-22,5	-32,9	5,4	-33,6	-28,7	
	07:07:40	016:08:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:07:40	016:08:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:07:40	016:08:00	0,345	10,0	-90,8	0,226	-30,7	-31,1	-13,4	-22,2	-32,8	5,4	-33,5	-28,1	
	07:08:40	016:09:00	0,198	10,0	-88,7	0,128	-30,0	-30,2	-12,7	-21,8	-32,0	5,4	-33,6	-27,4	
	07:09:45	016:10:05	0,198	10,0	-88,7	0,128	-29,6	-29,8	-12,0	-21,5	-31,6	5,4	---	-27,2	
	07:10:42	016:11:02	0,199	10,0	-88,5	0,128	-29,8	-30,5	-11,5	-21,8	-32,0	5,4	---	-27,9	
	07:11:42	016:12:02	0,201	10,0	-92,9	0,130	-30,2	-31,0	-11,2	-22,0	-32,5	5,4	---	-28,4	
	07:12:42	016:13:02	0,203	10,0	-88,9	0,131	-30,4	-31,3	-10,9	-22,2	-32,9	5,3	-33,4	-28,7	
	07:13:40	016:14:00	0,206	10,0	-92,2	0,134	-30,6	-31,5	-10,5	-22,3	-33,0	5,3	-33,4	-28,7	
	07:14:42	016:15:02	0,208	10,0	-89,2	0,136	-30,6	-31,6	-10,1	-22,4	-33,1	5,3	-33,4	-28,8	
	07:15:41	016:16:01	0,199	10,0	-88,2	0,129	-30,7	-31,6	-9,7	-22,4	-33,2	5,3	-33,3	-28,7	
	07:16:42	016:17:02	0,202	10,0	-88,9	0,131	-30,7	-31,6	-9,3	-22,4	-33,2	5,3	-33,3	-28,7	
	07:17:44	016:18:04	0,207	10,0	-90,2	0,134	-30,7	-31,6	-8,9	-22,5	-33,2	5,2	-33,2	-28,7	
	07:18:40	016:19:00	0,208	10,0	-89,2	0,136	-30,7	-31,6	-8,5	-22,5	-33,2	5,2	-33,2	-28,7	
	07:19:40	016:20:00	0,200	10,0	-91,0	0,129	-30,7	-31,6	-8,1	-22,5	-33,2	5,2	-33,2	-28,7	
	07:20:41	016:21:01	0,202	10,0	-90,2	0,131	-30,7	-31,6	-7,8	-22,5	-33,2	5,2	-33,1	-28,7	
	07:21:42	016:22:02	0,206	10,0	-88,9	0,134	-30,7	-31,6	-7,4	-22,5	-33,2	5,2	-33,1	-28,7	
	07:22:44	016:23:04	0,208	10,0	-88,5	0,136	-30,7	-31,6	-7,1	-22,5	-33,2	5,2	-33,1	-28,7	
	07:23:40	016:24:00	0,201	10,0	-89,1	0,130	-30,7	-31,6	-6,7	-22,5	-33,2	5,2	-33,0	-28,6	
	07:24:40	016:25:00	0,203	10,0	-91,7	0,132	-30,6	-31,6	-6,4	-22,5	-33,2	5,2	-32,9	-28,6	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	54 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:25:41	016:26:01	0,207	10,0	-87,7	0,134	-30,6	-31,6	-6,1	-22,5	-33,2	5,2	-32,9	-28,6	
	07:26:44	016:27:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:26:44	016:27:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:26:44	016:27:04	0,210	10,0	-87,0	0,137	-30,7	-31,7	-5,8	-22,5	-33,2	5,2	-32,8	-28,6	
	07:27:40	016:28:00	0,201	10,0	-87,6	0,130	-30,6	-31,6	-5,5	-22,5	-33,2	5,2	-32,7	-28,6	
	07:28:41	016:29:01	0,204	10,0	-87,5	0,133	-30,6	-31,6	-5,2	-22,5	-33,2	5,2	-32,6	-28,5	
	07:29:41	016:30:01	0,208	10,0	-87,1	0,136	-30,6	-31,7	-4,9	-22,5	-33,2	5,2	-32,5	-28,5	
	07:30:44	016:31:04	0,204	10,0	-85,1	0,131	-30,6	-31,7	-4,6	-22,5	-33,2	5,3	-32,5	-28,5	
	07:31:44	016:32:04	0,203	10,0	-88,0	0,132	-30,6	-31,6	-4,3	-22,5	-33,2	5,3	-32,4	-28,5	
	07:32:44	016:33:04	0,206	10,0	-87,1	0,134	-30,5	-31,6	-4,0	-22,4	-33,2	5,3	-32,3	-28,4	
	07:33:41	016:34:01	0,210	10,0	-89,9	0,137	-30,5	-31,7	-3,8	-22,4	-33,2	5,3	-32,2	-28,4	
	07:34:46	016:35:06	0,201	10,0	-89,5	0,130	-30,5	-31,7	-3,5	-22,4	-33,2	5,3	-32,1	-28,4	
	07:35:41	016:36:01	0,204	10,0	-89,6	0,132	-30,5	-31,7	-3,3	-22,4	-33,2	5,3	-32,0	-28,3	
	07:36:42	016:37:02	0,207	10,0	-88,1	0,134	-30,5	-31,7	-3,1	-22,4	-33,2	5,3	-32,0	-28,3	
	07:37:40	016:38:00	0,211	10,0	-90,6	0,137	-30,5	-31,7	-2,8	-22,4	-33,2	5,3	-31,9	-28,4	
	07:38:41	016:39:01	0,201	10,0	-88,4	0,131	-30,5	-31,7	-2,5	-22,4	-33,2	5,3	-31,7	-28,3	
	07:39:43	016:40:03	0,205	10,0	-91,7	0,133	-30,5	-31,7	-2,2	-22,4	-33,2	5,3	-31,6	-28,3	
	07:40:41	016:41:01	0,208	10,0	-89,2	0,135	-30,5	-31,7	-2,0	-22,4	-33,3	5,3	-31,6	-28,3	
	07:41:43	016:42:03	0,200	10,0	-90,3	0,129	-30,5	-31,7	-1,7	-22,4	-33,2	5,3	-31,4	-28,3	
	07:42:43	016:43:03	0,202	10,0	-91,2	0,131	-30,4	-31,7	-1,5	-22,4	-33,2	5,3	-31,2	-28,2	
	07:43:44	016:44:04	0,206	10,0	-88,1	0,134	-30,4	-31,7	-1,3	-22,4	-33,2	5,3	-31,1	-28,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 55 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:44:44	016:45:04	0,209	10,0	-87,4	0,136	-30,3	-31,8	-1,1	-22,4	-33,3	5,3	-30,9	-28,2	
	07:45:41	016:46:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:45:41	016:46:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:45:41	016:46:01	0,201	10,0	-88,0	0,130	-30,3	-31,7	-0,9	-22,4	-33,2	5,3	-30,7	-28,1	
	07:46:43	016:47:03	0,204	10,0	-86,3	0,132	-30,3	-31,7	-0,8	-22,4	-33,2	5,3	-30,5	-28,1	
	07:47:43	016:48:03	0,207	10,0	-88,4	0,135	-30,2	-31,7	-0,6	-22,4	-33,2	5,3	-30,2	-28,1	
	07:48:40	016:49:00	0,210	10,0	-87,7	0,137	-30,2	-31,7	-0,5	-22,4	-33,2	5,4	-29,9	-28,1	
	07:49:41	016:50:01	0,202	10,0	-89,1	0,131	-30,1	-31,7	-0,3	-22,4	-33,2	5,4	-29,4	-28,0	
	07:50:43	016:51:03	0,205	10,0	-87,4	0,132	-30,0	-31,7	-0,2	-22,4	-33,2	5,4	-28,5	-28,0	
	07:51:43	016:52:03	0,208	10,0	-88,4	0,135	-30,0	-31,8	-0,1	-22,4	-33,2	5,4	-27,2	-28,0	
	07:52:40	016:53:00	0,212	10,0	-87,3	0,137	-29,9	-31,8	0,1	-22,4	-33,2	5,4	-26,0	-27,9	
	07:53:41	016:54:01	0,203	10,0	-89,6	0,131	-29,8	-31,7	0,2	-22,4	-33,2	5,4	-25,1	-27,9	
	07:54:43	016:55:03	0,205	10,0	-88,7	0,133	-29,7	-31,7	0,3	-22,4	-33,2	5,4	-24,6	-27,9	
	07:55:44	016:56:04	0,209	10,0	-89,3	0,136	-29,6	-31,8	0,4	-22,4	-33,2	5,4	-24,3	-27,9	
	07:56:42	016:57:02	0,207	10,0	-88,6	0,133	-29,4	-31,7	0,5	-22,4	-33,2	5,4	-24,0	-27,8	
	07:57:41	016:58:01	0,203	10,0	-90,7	0,132	-29,2	-31,7	0,6	-22,3	-33,2	5,4	-23,8	-27,8	
	07:58:41	016:59:01	0,206	10,0	-88,9	0,134	-28,9	-31,7	0,7	-22,3	-33,2	5,4	-23,6	-27,7	
	07:59:41	017:00:01	0,211	10,0	-88,1	0,136	-28,6	-31,8	0,8	-22,3	-33,2	5,4	-23,5	-27,7	
	08:00:42	017:01:02	0,201	10,0	-88,7	0,130	-28,2	-31,8	0,9	-22,3	-33,2	5,4	-23,4	-27,6	
	08:01:43	017:02:03	0,204	10,0	-87,9	0,132	-27,6	-31,8	1,1	-22,3	-33,2	5,4	-23,2	-27,5	
	08:02:43	017:03:03	0,207	10,0	-88,3	0,135	-26,9	-31,8	1,2	-22,3	-33,2	5,4	-23,1	-27,5	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 56 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



08:03:43	017:04:03	0,210	10,0	-92,8	0,137	-26,3	-31,8	1,3	-22,3	-33,2	5,4	-23,0	-27,4
08:04:43	017:05:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:04:43	017:05:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:04:43	017:05:03	0,201	10,0	-89,3	0,130	-25,9	-31,8	1,4	-22,3	-33,2	5,4	-22,9	-27,3
08:05:44	017:06:04	0,204	10,0	-88,1	0,132	-25,7	-31,8	1,4	-22,3	-33,2	5,4	-22,8	-27,2
08:06:40	017:07:00	0,208	10,0	-87,9	0,135	-25,6	-31,8	1,5	-22,3	-33,2	5,4	-22,7	-27,1
08:07:41	017:08:01	0,211	10,0	-85,5	0,137	-25,5	-31,8	1,6	-22,3	-33,2	5,4	-22,6	-27,0
08:08:41	017:09:01	0,203	10,0	-86,9	0,131	-25,3	-31,8	1,6	-22,3	-33,2	5,5	-22,5	-26,7
08:09:43	017:10:03	0,205	10,0	-84,8	0,133	-25,2	-31,8	1,7	-22,3	-33,2	5,5	-22,4	-26,5
08:10:41	017:11:01	0,208	10,0	-88,3	0,135	-25,1	-31,8	1,7	-22,3	-33,2	5,5	-22,3	-26,4
08:11:44	017:12:04	0,210	10,0	-87,0	0,138	-25,0	-31,8	1,8	-22,3	-33,2	5,5	-22,2	-26,3
08:12:41	017:13:01	0,202	10,0	-89,7	0,131	-24,8	-31,8	1,8	-22,3	-33,2	5,5	-22,1	-26,1
08:13:42	017:14:02	0,205	10,0	-88,3	0,133	-24,7	-31,8	1,9	-22,3	-33,2	5,5	-22,0	-26,0
08:14:42	017:15:02	0,208	10,0	-89,6	0,136	-24,6	-31,8	1,9	-22,3	-33,2	5,5	-21,9	-25,9
08:15:44	017:16:04	0,211	10,0	-88,4	0,138	-24,5	-31,8	2,0	-22,3	-33,2	5,5	-21,8	-25,8
08:16:44	017:17:04	0,203	10,0	-88,3	0,131	-24,4	-31,8	2,0	-22,3	-33,1	5,5	-21,7	-25,7
08:17:40	017:18:00	0,205	10,0	-88,2	0,133	-24,3	-31,8	2,0	-22,3	-33,1	5,5	-21,6	-25,6
08:18:42	017:19:02	0,209	10,0	-88,0	0,136	-24,2	-31,9	2,1	-22,3	-33,1	5,5	-21,6	-25,4
08:19:44	017:20:04	0,200	10,0	-86,0	0,129	-24,2	-31,8	2,1	-22,3	-33,1	5,5	-21,5	-25,2
08:20:40	017:21:00	0,203	10,0	-89,6	0,132	-24,1	-31,8	2,1	-22,3	-33,1	5,5	-21,4	-25,0
08:21:42	017:22:02	0,206	10,0	-87,3	0,134	-24,0	-31,8	2,2	-22,3	-33,1	5,5	-21,4	-24,8

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature: Page: 57 / 209
Document issued by:

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



08:22:42	017:23:02	0,210	10,0	-90,0	0,137	-23,9	-31,9	2,2	-22,3	-33,1	5,5	-21,3	-24,5
08:23:43	017:24:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:23:43	017:24:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:23:43	017:24:03	0,201	10,0	-88,3	0,130	-23,8	-31,8	2,2	-22,3	-33,1	5,5	-21,2	-24,3
08:24:40	017:25:00	0,204	10,0	-89,8	0,132	-23,7	-31,9	2,3	-22,3	-33,1	5,5	-21,1	-24,1
08:25:41	017:26:01	0,207	10,0	-89,5	0,134	-23,6	-31,9	2,3	-22,3	-33,1	5,5	-21,0	-23,9
08:26:41	017:27:01	0,210	10,0	-89,7	0,137	-23,5	-31,9	2,3	-22,3	-33,1	5,5	-21,0	-23,8
08:27:42	017:28:02	0,201	10,0	-88,6	0,130	-23,4	-31,9	2,4	-22,3	-33,1	5,5	-20,9	-23,7
08:28:45	017:29:05	0,204	10,0	-88,9	0,133	-23,3	-31,9	2,4	-22,3	-33,1	5,5	-20,8	-23,5
08:29:42	017:30:02	0,207	10,0	-88,7	0,134	-23,3	-31,9	2,4	-22,3	-33,1	5,5	-20,7	-23,4
08:30:43	017:31:03	0,211	10,0	-88,5	0,137	-23,2	-31,9	2,4	-22,3	-33,1	5,5	-20,7	-23,3
08:31:40	017:32:00	0,202	10,0	-87,7	0,130	-23,1	-31,9	2,5	-22,3	-33,0	5,6	-20,6	-23,2
08:32:40	017:33:00	0,204	10,0	-87,7	0,133	-23,0	-31,9	2,5	-22,3	-33,0	5,6	-20,5	-23,1
08:33:40	017:34:00	0,207	10,0	-86,6	0,135	-23,0	-31,9	2,5	-22,3	-33,0	5,6	-20,5	-23,0
08:34:44	017:35:04	0,211	10,0	-88,2	0,138	-22,9	-31,9	2,5	-22,3	-33,0	5,6	-20,4	-22,9
08:35:44	017:36:04	0,203	10,0	-87,6	0,131	-22,8	-31,9	2,6	-22,3	-33,0	5,6	-20,3	-22,8
08:36:44	017:37:04	0,205	10,0	-89,1	0,133	-22,7	-31,9	2,6	-22,3	-33,0	5,6	-20,2	-22,7
08:37:40	017:38:00	0,208	10,0	-89,2	0,135	-22,6	-31,9	2,6	-22,3	-33,0	5,6	-20,2	-22,6
08:38:41	017:39:01	0,209	10,0	-88,8	0,136	-22,5	-31,9	2,6	-22,3	-33,0	5,6	-20,1	-22,5
08:39:42	017:40:02	0,202	10,0	-87,6	0,131	-22,5	-31,9	2,6	-22,3	-32,9	5,6	-20,0	-22,4
08:40:42	017:41:02	0,206	10,0	-88,8	0,133	-22,4	-31,9	2,7	-22,3	-32,9	5,6	-20,0	-22,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 58 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



08:41:41	017:42:01	0,208	10,0	-88,1	0,135	-22,4	-31,9	2,7	-22,3	-32,9	5,6	-19,9	-22,2
08:42:42	017:43:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:42:42	017:43:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:42:42	017:43:02	0,200	10,0	-89,6	0,129	-22,3	-31,9	2,7	-22,3	-32,9	5,6	-19,8	-22,1
08:43:42	017:44:02	0,202	10,0	-88,1	0,131	-22,2	-31,9	2,7	-22,3	-32,9	5,6	-19,8	-22,0
08:44:42	017:45:02	0,206	10,0	-90,0	0,134	-22,1	-31,9	2,7	-22,3	-32,9	5,6	-19,7	-21,9
08:45:44	017:46:04	0,209	10,0	-89,1	0,136	-22,1	-32,0	2,8	-22,3	-32,9	5,6	-19,7	-21,8
08:46:40	017:47:00	0,200	10,0	-91,0	0,130	-22,0	-31,9	2,8	-22,3	-32,9	5,6	-19,6	-21,8
08:47:41	017:48:01	0,202	10,0	-89,1	0,131	-21,9	-31,9	2,8	-22,2	-32,8	5,6	-19,5	-21,6
08:48:42	017:49:02	0,206	10,0	-89,6	0,134	-21,8	-31,9	2,8	-22,2	-32,8	5,6	-19,5	-21,5
08:49:43	017:50:03	0,209	10,0	-89,3	0,136	-21,7	-31,9	2,8	-22,2	-32,8	5,6	-19,4	-21,4
08:50:44	017:51:04	0,200	10,0	-89,7	0,129	-21,7	-31,9	2,8	-22,2	-32,8	5,7	-19,3	-21,3
08:51:40	017:52:00	0,202	10,0	-89,0	0,131	-21,6	-31,9	2,9	-22,2	-32,7	5,7	-19,3	-21,2
08:52:40	017:53:00	0,206	10,0	-82,8	0,134	-21,5	-31,9	2,9	-22,2	-32,6	5,7	-19,2	-21,1
08:53:40	017:54:00	0,209	10,0	-87,7	0,136	-21,5	-32,0	2,9	-22,2	-32,6	5,7	-19,1	-21,0
08:54:41	017:55:01	0,200	10,0	-84,4	0,130	-21,4	-31,9	2,9	-22,2	-32,5	5,7	-19,1	-20,9
08:55:41	017:56:01	0,203	10,0	-87,0	0,132	-21,3	-31,9	2,9	-22,2	-32,4	5,7	-19,0	-20,8
08:56:42	017:57:02	0,206	10,0	-85,8	0,134	-21,2	-31,9	2,9	-22,2	-32,3	5,7	-19,0	-20,7
08:57:44	017:58:04	0,210	10,0	-87,3	0,137	-21,2	-32,0	3,0	-22,2	-32,3	5,7	-18,9	-20,6
08:58:41	017:59:01	0,200	10,0	-87,0	0,130	-21,1	-31,9	3,0	-22,2	-32,1	5,7	-18,8	-20,5
08:59:41	018:00:01	0,203	10,0	-87,4	0,132	-21,0	-31,9	3,0	-22,2	-31,9	5,7	-18,8	-20,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 59 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



09:00:43	018:01:03	0,206	10,0	-86,8	0,134	-20,9	-32,0	3,0	-22,2	-31,7	5,7	-18,7	-20,2
09:01:45	018:02:05	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
09:01:45	018:02:05	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
09:01:45	018:02:05	0,210	10,0	-88,5	0,137	-20,9	-32,0	3,0	-22,2	-31,4	5,7	-18,7	-20,1
09:02:44	018:03:04	0,200	10,0	-85,0	0,130	-20,8	-32,0	3,0	-22,2	-31,2	5,7	-18,6	-20,0
09:03:41	018:04:01	0,203	10,0	-87,7	0,132	-20,7	-32,0	3,1	-22,2	-31,0	5,7	-18,6	-20,0
09:04:42	018:05:02	0,206	10,0	-85,4	0,134	-20,6	-32,0	3,1	-22,2	-30,9	5,7	-18,5	-19,8
09:05:45	018:06:05	0,210	10,0	-87,2	0,137	-20,6	-32,0	3,1	-22,2	-30,8	5,8	-18,4	-19,7
09:06:41	018:07:01	0,200	10,0	-86,3	0,130	-20,5	-32,0	3,1	-22,2	-30,8	5,8	-18,4	-19,6
09:07:44	018:08:04	0,204	10,0	-88,6	0,132	-20,4	-32,0	3,1	-22,2	-30,7	5,7	-18,3	-19,5
09:08:40	018:09:00	0,206	10,0	-87,9	0,134	-20,4	-32,0	3,1	-22,2	-30,6	5,8	-18,3	-19,4
09:09:41	018:10:01	0,210	10,0	-89,3	0,137	-20,3	-32,0	3,1	-22,2	-30,6	5,8	-18,2	-19,3
09:10:41	018:11:01	0,200	10,0	-88,5	0,130	-20,2	-32,0	3,2	-22,2	-30,6	5,8	-18,2	-19,2
09:11:43	018:12:03	0,204	10,0	-88,6	0,132	-20,2	-32,0	3,2	-22,2	-30,5	5,8	-18,1	-19,0
09:12:44	018:13:04	0,207	10,0	-86,5	0,135	-20,1	-32,0	3,2	-22,2	-30,5	5,8	-18,1	-18,9
09:13:40	018:14:00	0,210	10,0	-89,6	0,137	-20,0	-32,0	3,2	-22,2	-30,5	5,8	-18,0	-18,7
09:14:40	018:15:00	0,200	10,0	-88,3	0,130	-20,0	-32,0	3,2	-22,2	-30,5	5,8	-18,0	-18,5
09:15:40	018:16:00	0,204	10,0	-92,0	0,133	-19,9	-32,0	3,2	-22,2	-30,4	5,8	-17,9	-18,4
09:16:44	018:17:04	0,207	10,0	-88,7	0,135	-19,9	-32,0	3,2	-22,2	-30,4	5,8	-17,9	-18,2
09:17:40	018:18:00	0,210	10,0	-91,8	0,137	-19,8	-32,0	3,2	-22,2	-30,3	5,8	-17,8	-18,1
09:18:42	018:19:02	0,201	10,0	-88,4	0,131	-19,7	-32,0	3,3	-22,2	-30,3	5,8	-17,8	-17,9

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 60 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	09:19:44	018:20:04	0,204	10,0	-90,0	0,133	-19,7	-32,0	3,3	-22,2	-30,3	5,8	-17,7	-17,8
	09:20:44	018:21:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	09:20:44	018:21:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	09:20:44	018:21:04	0,209	10,0	-88,4	0,135	-19,6	-32,1	3,3	-22,2	-30,3	5,8	-17,7	-17,7
	09:21:44	018:22:04	0,206	10,0	-90,3	0,134	-19,6	-32,1	3,3	-22,2	-30,2	5,8	-17,6	-17,6
	09:22:44	018:23:04	0,201	10,0	-88,1	0,131	-19,5	-32,0	3,3	-22,2	-30,2	5,8	-17,6	-17,5
	09:23:44	018:24:04	0,204	10,0	-84,1	0,133	-19,4	-32,0	3,3	-22,2	-30,2	5,8	-17,5	-17,4
	09:24:41	018:25:01	0,208	10,0	-86,7	0,136	-19,4	-32,1	3,3	-22,2	-30,2	5,8	-17,5	-17,3
	09:25:41	018:26:01	0,201	10,0	-85,1	0,130	-19,3	-32,1	3,3	-22,2	-30,1	5,8	-17,4	-17,2
	09:26:42	018:27:02	0,202	10,0	-87,8	0,132	-19,2	-32,0	3,3	-22,2	-30,1	5,9	-17,4	-17,1
	09:27:42	018:28:02	0,205	10,0	-88,0	0,133	-19,2	-32,0	3,4	-22,2	-30,1	5,9	-17,3	-17,0
	09:28:44	018:29:04	0,208	10,0	-87,5	0,136	-19,1	-32,1	3,4	-22,2	-30,1	5,9	-17,3	-16,9
	09:29:40	018:30:00	0,206	10,0	-87,7	0,133	-19,1	-32,1	3,4	-22,2	-30,0	5,9	-17,2	-16,8
	09:30:40	018:31:00	0,202	10,0	-87,7	0,131	-19,0	-32,0	3,4	-22,2	-30,0	5,9	-17,2	-16,7
	09:31:41	018:32:01	0,205	10,0	-86,7	0,133	-18,9	-32,1	3,4	-22,2	-29,9	5,9	-17,1	-16,6
	09:32:42	018:33:02	0,209	10,0	-87,3	0,136	-18,9	-32,1	3,4	-22,2	-29,9	5,9	-17,1	-16,5
	09:33:43	018:34:03	0,199	10,0	-85,5	0,129	-18,8	-32,1	3,4	-22,2	-29,9	5,9	-17,0	-16,4
	09:34:43	018:35:03	0,202	10,0	-87,1	0,132	-18,7	-32,1	3,4	-22,2	-29,8	5,9	-17,0	-16,3
	09:35:40	018:36:00	0,205	10,0	-84,9	0,133	-18,7	-32,1	3,4	-22,1	-29,8	5,9	-16,9	-16,3
	09:36:41	018:37:01	0,208	10,0	-87,7	0,136	-18,6	-32,1	3,4	-22,2	-29,8	5,9	-16,9	-16,2
	09:37:44	018:38:04	0,199	10,0	-90,8	0,129	-18,6	-32,1	3,5	-22,1	-29,8	5,9	-16,9	-16,1

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 61 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	09:38:44	018:39:04	0,202	10,0	-88,8	0,131	-18,5	-32,1	3,5	-22,1	-29,7	5,9	-16,8	-16,0	
	09:39:41	018:40:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	09:39:41	018:40:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	09:39:41	018:40:01	0,206	10,0	-85,3	0,133	-18,5	-32,1	3,5	-22,1	-29,7	5,9	-16,8	-16,0	
	09:40:42	018:41:02	0,208	10,0	-87,3	0,136	-18,4	-32,1	3,5	-22,2	-29,7	5,9	-16,7	-15,9	
	09:41:44	018:42:04	0,200	10,0	-84,2	0,129	-18,3	-32,1	3,5	-22,1	-29,6	5,9	-16,7	-15,8	
	09:42:42	018:43:02	0,202	10,0	-88,2	0,131	-18,3	-32,1	3,5	-22,1	-29,6	5,9	-16,6	-15,7	
	09:43:42	018:44:02	0,205	10,0	-86,2	0,133	-18,2	-32,1	3,5	-22,1	-29,6	5,9	-16,6	-15,7	
	09:44:43	018:45:03	0,208	10,0	-88,5	0,136	-18,1	-32,1	3,5	-22,1	-29,5	5,9	-16,5	-15,6	
	09:45:43	018:46:03	0,206	10,0	-88,6	0,133	-18,1	-32,1	3,5	-22,1	-29,5	5,9	-16,5	-15,5	
	09:46:44	018:47:04	0,203	10,0	-88,2	0,131	-18,0	-32,1	3,5	-22,1	-29,5	6,0	-16,4	-15,4	
	09:47:41	018:48:01	0,204	10,0	-86,9	0,133	-18,0	-32,1	3,6	-22,1	-29,4	6,0	-16,4	-15,3	
	09:48:42	018:49:02	0,208	10,0	-89,3	0,135	-17,9	-32,2	3,6	-22,1	-29,4	6,0	-16,4	-15,3	
	09:49:43	018:50:03	0,211	10,0	-87,8	0,138	-17,9	-32,1	3,6	-22,1	-29,4	6,0	-16,3	-15,2	
	09:50:44	018:51:04	0,203	10,0	-89,0	0,131	-17,8	-32,1	3,6	-22,1	-29,3	6,0	-16,3	-15,1	
	09:51:45	018:52:05	0,204	10,0	-88,6	0,133	-17,8	-32,1	3,6	-22,1	-29,3	6,0	-16,2	-15,0	
	09:52:41	018:53:01	0,208	10,0	-89,5	0,135	-17,7	-32,2	3,6	-22,1	-29,3	6,0	-16,2	-15,0	
	09:53:41	018:54:01	0,210	10,0	-88,2	0,137	-17,7	-32,2	3,6	-22,1	-29,2	6,0	-16,1	-14,9	
	09:54:43	018:55:03	0,203	10,0	-90,2	0,131	-17,6	-32,2	3,6	-22,1	-29,2	6,0	-16,1	-14,8	
	09:55:41	018:56:01	0,204	10,0	-89,5	0,132	-17,6	-32,2	3,6	-22,1	-29,1	6,0	-16,0	-14,7	
	09:56:41	018:57:01	0,207	10,0	-90,9	0,135	-17,5	-32,2	3,6	-22,1	-29,1	6,0	-16,0	-14,7	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature: Page: 62 / 209
Document issued by:

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	09:57:40	018:58:00	0,210	10,0	-89,4	0,137	-17,5	-32,3	3,7	-22,1	-29,1	6,0	-15,9	-14,6	
	09:58:42	018:59:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	09:58:42	018:59:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	09:58:42	018:59:02	0,201	10,0	-85,0	0,130	-17,5	-32,2	3,7	-22,1	-29,1	6,0	-15,9	-14,5	
	09:59:43	019:00:03	0,205	10,0	-88,9	0,132	-17,4	-32,3	3,7	-22,1	-29,0	6,0	-15,8	-14,4	
	10:00:43	019:01:03	0,207	10,0	-89,2	0,134	-17,3	-32,3	3,7	-22,1	-29,0	6,1	-15,8	-14,4	
	10:01:40	019:02:00	0,210	10,0	-88,0	0,137	-17,3	-32,3	3,7	-22,1	-29,0	6,1	-15,8	-14,3	
	10:02:42	019:03:02	0,201	10,0	-84,3	0,130	-17,2	-32,3	3,7	-22,1	-28,9	6,1	-15,7	-14,2	
	10:03:40	019:04:00	0,203	10,0	-88,2	0,132	-17,2	-32,4	3,7	-22,1	-28,9	6,1	-15,7	-14,1	
	10:04:40	019:05:00	0,206	10,0	-84,6	0,134	-17,1	-32,4	3,7	-22,1	-28,8	6,1	-15,6	-14,0	
	10:05:44	019:06:04	0,209	10,0	-88,6	0,136	-17,1	-32,5	3,8	-22,1	-28,8	6,1	-15,6	-14,0	
	10:06:40	019:07:00	0,200	10,0	-85,2	0,129	-17,0	-32,5	3,8	-22,1	-28,8	6,1	-15,6	-13,9	
	10:07:41	019:08:01	0,202	10,0	-89,0	0,131	-17,0	-32,5	3,8	-22,1	-28,7	6,1	-15,5	-13,8	
	10:08:43	019:09:03	0,204	10,0	-83,3	0,133	-16,9	-32,6	3,8	-22,1	-28,7	6,1	-15,5	-13,7	
	10:09:44	019:10:04	0,209	10,0	-88,4	0,135	-16,9	-32,6	3,8	-22,1	-28,7	6,1	-15,4	-13,7	
	10:10:44	019:11:04	0,212	10,0	-88,7	0,138	-16,8	-32,7	3,8	-22,1	-28,7	6,1	-15,4	-13,6	
	10:11:40	019:12:00	0,202	10,0	-88,2	0,130	-16,7	-32,7	3,8	-22,1	-28,6	6,1	-15,4	-13,5	
	10:12:41	019:13:01	0,204	10,0	-84,3	0,132	-16,7	-32,7	3,8	-22,1	-28,6	6,1	-15,3	-13,4	
	10:13:42	019:14:02	0,208	10,0	-88,3	0,135	-16,6	-32,9	3,8	-22,1	-28,6	6,1	-15,3	-13,4	
	10:14:43	019:15:03	0,210	10,0	-86,9	0,137	-16,6	-32,9	3,9	-22,1	-28,5	6,1	-15,2	-13,3	
	10:15:43	019:16:03	0,202	10,0	-89,7	0,130	-16,5	-32,9	3,9	-22,1	-28,5	6,1	-15,2	-13,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 63 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



10:16:41	019:17:01	0,203	10,0	-88,9	0,131	-16,5	-32,9	3,9	-22,0	-28,4	6,1	-15,2	-13,2
10:17:42	019:18:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
10:17:42	019:18:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
10:17:42	019:18:02	0,207	10,0	-90,1	0,134	-16,4	-33,0	3,9	-22,1	-28,4	6,2	-15,1	-13,1
10:18:42	019:19:02	0,208	10,0	-87,5	0,136	-16,4	-33,0	3,9	-22,1	-28,4	6,2	-15,1	-13,1
10:19:44	019:20:04	0,200	10,0	-90,4	0,130	-16,3	-33,0	3,9	-22,1	-28,3	6,2	-15,0	-13,0
10:20:42	019:21:02	0,202	10,0	-89,9	0,131	-16,3	-33,0	3,9	-22,0	-28,3	6,2	-15,0	-12,9
10:21:44	019:22:04	0,207	10,0	-89,8	0,133	-16,2	-33,1	3,9	-22,0	-28,2	6,2	-14,9	-12,8
10:22:44	019:23:04	0,207	10,0	-91,9	0,135	-16,2	-33,1	3,9	-22,1	-28,2	6,2	-14,9	-12,7
10:23:44	019:24:04	0,210	10,0	-88,1	0,137	-16,1	-33,1	3,9	-22,0	-28,1	6,2	-14,9	-12,6
10:24:40	019:25:00	0,202	10,0	-86,8	0,130	-16,1	-33,1	3,9	-22,0	-28,1	6,2	-14,8	-12,5
10:25:44	019:26:04	0,205	10,0	-89,4	0,132	-16,0	-33,1	4,0	-22,0	-28,1	6,2	-14,8	-12,5
10:26:40	019:27:00	0,207	10,0	-87,8	0,134	-16,0	-33,2	4,0	-22,1	-28,0	6,2	-14,8	-12,4
10:27:42	019:28:02	0,210	10,0	-92,7	0,137	-15,9	-33,2	4,0	-22,1	-28,0	6,2	-14,7	-12,3
10:28:43	019:29:03	0,201	10,0	-88,9	0,130	-15,9	-33,2	4,0	-22,0	-27,9	6,2	-14,7	-12,2
10:29:44	019:30:04	0,203	10,0	-84,5	0,132	-15,8	-33,2	4,0	-22,0	-27,9	6,2	-14,7	-12,1
10:30:41	019:31:01	0,206	10,0	-87,2	0,134	-15,8	-33,2	4,0	-22,1	-27,8	6,2	-14,6	-12,0
10:31:44	019:32:04	0,208	10,0	-84,0	0,136	-15,7	-33,3	4,0	-22,1	-27,8	6,2	-14,6	-12,0
10:32:42	019:33:02	0,200	10,0	-88,2	0,129	-15,7	-33,2	4,0	-22,0	-27,7	6,3	-14,6	-11,9
10:33:40	019:34:00	0,201	10,0	-87,0	0,130	-15,6	-33,2	4,0	-22,0	-27,7	6,2	-14,5	-11,8
10:34:42	019:35:02	0,205	10,0	-89,6	0,133	-15,6	-33,3	4,0	-22,0	-27,6	6,2	-14,5	-11,8

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 64 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



10:35:43	019:36:03	0,208	10,0	-88,1	0,135	-15,5	-33,3	4,0	-22,1	-27,6	6,3	-14,4	-11,7
10:36:43	019:37:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
10:36:43	019:37:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
10:36:43	019:37:03	0,210	10,0	-91,1	0,137	-15,5	-33,3	4,1	-22,0	-27,6	6,3	-14,4	-11,6
10:37:40	019:38:00	0,201	10,0	-89,6	0,130	-15,4	-33,3	4,1	-22,0	-27,5	6,3	-14,4	-11,5
10:38:42	019:39:02	0,204	10,0	-87,8	0,132	-15,4	-33,3	4,1	-22,0	-27,5	6,3	-14,3	-11,5
10:39:42	019:40:02	0,206	10,0	-89,5	0,134	-15,3	-33,4	4,1	-22,0	-27,4	6,3	-14,3	-11,4
10:40:44	019:41:04	0,210	10,0	-89,7	0,136	-15,3	-33,4	4,1	-22,0	-27,4	6,3	-14,3	-11,3
10:41:44	019:42:04	0,200	10,0	-89,4	0,129	-15,2	-33,3	4,1	-22,0	-27,4	6,3	-14,2	-11,3
10:42:40	019:43:00	0,202	10,0	-84,6	0,131	-15,2	-33,3	4,1	-22,0	-27,3	6,3	-14,2	-11,2
10:43:40	019:44:00	0,204	10,0	-88,3	0,133	-15,1	-33,4	4,1	-22,0	-27,2	6,3	-14,1	-11,1
10:44:43	019:45:03	0,207	10,0	-84,8	0,135	-15,1	-33,4	4,1	-22,0	-27,2	6,3	-14,1	-11,1
10:45:40	019:46:00	0,211	10,0	-87,4	0,138	-15,0	-33,4	4,2	-22,0	-27,2	6,3	-14,1	-11,0
10:46:40	019:47:00	0,201	10,0	-85,3	0,130	-15,0	-33,4	4,2	-22,0	-27,1	6,3	-14,0	-10,9
10:47:42	019:48:02	0,204	10,0	-88,6	0,133	-14,9	-33,4	4,2	-22,0	-27,1	6,3	-14,0	-10,9
10:48:43	019:49:03	0,207	10,0	-86,9	0,135	-14,9	-33,4	4,2	-22,0	-27,0	6,3	-14,0	-10,8
10:49:44	019:50:04	0,210	10,0	-89,7	0,137	-14,8	-33,4	4,2	-22,0	-27,0	6,3	-13,9	-10,8
10:50:40	019:51:00	0,201	10,0	-88,9	0,130	-14,8	-33,4	4,2	-22,0	-26,9	6,3	-13,9	-10,7
10:51:43	019:52:03	0,204	10,0	-90,5	0,132	-14,7	-33,4	4,2	-22,0	-26,9	6,3	-13,8	-10,6
10:52:44	019:53:04	0,205	10,0	-87,8	0,134	-14,7	-33,5	4,2	-22,0	-26,8	6,3	-13,8	-10,6
10:53:44	019:54:04	0,209	10,0	-90,6	0,136	-14,6	-33,5	4,2	-22,0	-26,8	6,3	-13,8	-10,5

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 65 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	10:54:40	019:55:00	0,203	10,0	-89,1	0,130	-14,6	-33,4	4,2	-22,0	-26,8	6,3	-13,8	-10,5
	10:55:41	019:56:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
	10:55:41	019:56:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
	10:55:41	019:56:01	0,203	10,0	-82,5	0,131	-14,5	-33,4	4,2	-22,0	-26,7	6,3	-13,7	-10,4
	10:56:42	019:57:02	0,205	10,0	-88,9	0,133	-14,5	-33,5	4,3	-22,0	-26,7	6,3	-13,7	-10,3
	10:57:43	019:58:03	0,207	10,0	-85,8	0,135	-14,4	-33,5	4,3	-22,0	-26,6	6,4	-13,6	-10,3
	10:58:45	019:59:05	0,210	10,0	-89,2	0,137	-14,4	-33,5	4,3	-22,0	-26,6	6,4	-13,6	-10,2
	10:59:42	020:00:02	0,201	10,0	-83,2	0,130	-14,3	-33,4	4,3	-22,0	-26,5	6,4	-13,6	-10,2
	11:00:40	020:01:00	0,204	10,0	-87,3	0,132	-14,3	-33,5	4,3	-22,0	-26,5	6,4	-13,5	-10,1
	11:01:40	020:02:00	0,207	10,0	-84,4	0,134	-14,2	-33,5	4,3	-22,0	-26,4	6,4	-13,5	-10,1
	11:02:43	020:03:03	0,210	10,0	-87,9	0,137	-14,2	-33,5	4,3	-22,0	-26,4	6,4	-13,4	-10,0
	11:03:41	020:04:01	0,201	10,0	-85,9	0,130	-14,1	-33,5	4,3	-21,9	-26,4	6,4	-13,4	-9,9
	11:04:44	020:05:04	0,203	10,0	-88,1	0,132	-14,1	-33,5	4,3	-21,9	-26,3	6,4	-13,4	-9,9
	11:05:45	020:06:05	0,205	10,0	-86,2	0,134	-14,0	-33,5	4,3	-21,9	-26,3	6,4	-13,3	-9,8
	11:06:41	020:07:01	0,210	10,0	-89,5	0,136	-14,0	-33,5	4,4	-22,0	-26,2	6,4	-13,3	-9,8
	11:07:42	020:08:02	0,200	10,0	-87,1	0,129	-13,9	-33,5	4,4	-21,9	-26,2	6,4	-13,3	-9,7
	11:08:43	020:09:03	0,203	10,0	-89,9	0,131	-13,9	-33,5	4,4	-21,9	-26,2	6,4	-13,2	-9,6
	11:09:44	020:10:04	0,204	10,0	-89,1	0,133	-13,9	-33,5	4,4	-21,9	-26,1	6,4	-13,2	-9,6
	11:10:43	020:11:03	0,208	10,0	-90,1	0,135	-13,8	-33,6	4,4	-21,9	-26,1	6,4	-13,2	-9,6
	11:11:44	020:12:04	0,210	10,0	-88,0	0,137	-13,8	-33,5	4,4	-21,9	-26,0	6,4	-13,1	-9,5
	11:12:40	020:13:00	0,202	10,0	-89,4	0,130	-13,7	-33,5	4,4	-21,9	-26,0	6,4	-13,1	-9,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 66 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



11:13:42	020:14:02	0,204	10,0	-88,6	0,132	-13,7	-33,5	4,4	-21,9	-25,9	6,5	-13,0	-9,4
11:14:44	020:15:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
11:14:44	020:15:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
11:14:44	020:15:04	0,207	10,0	-88,6	0,134	-13,6	-33,6	4,4	-21,9	-25,9	6,5	-13,0	-9,3
11:15:40	020:16:00	0,208	10,0	-87,4	0,136	-13,6	-33,6	4,4	-21,9	-25,8	6,5	-13,0	-9,3
11:16:42	020:17:02	0,200	10,0	-80,5	0,130	-13,5	-33,5	4,4	-21,9	-25,8	6,5	-12,9	-9,3
11:17:43	020:18:03	0,203	10,0	-88,7	0,131	-13,5	-33,5	4,5	-21,9	-25,7	6,5	-12,9	-9,2
11:18:40	020:19:00	0,207	10,0	-82,4	0,133	-13,4	-33,6	4,5	-21,9	-25,7	6,5	-12,9	-9,2
11:19:41	020:20:01	0,208	10,0	-87,2	0,136	-13,4	-33,6	4,5	-21,9	-25,7	6,5	-12,8	-9,1
11:20:41	020:21:01	0,210	10,0	-85,9	0,138	-13,4	-33,6	4,5	-21,9	-25,6	6,5	-12,8	-9,1
11:21:44	020:22:04	0,202	10,0	-87,7	0,131	-13,3	-33,5	4,5	-21,9	-25,6	6,5	-12,7	-9,1
11:22:40	020:23:00	0,204	10,0	-90,4	0,132	-13,3	-33,6	4,5	-21,9	-25,5	6,5	-12,7	-9,0
11:23:41	020:24:01	0,207	10,0	-88,5	0,134	-13,2	-33,6	4,5	-21,9	-25,5	6,5	-12,7	-8,9
11:24:44	020:25:04	0,209	10,0	-84,4	0,136	-13,2	-33,6	4,5	-21,9	-25,4	6,5	-12,7	-8,9
11:25:41	020:26:01	0,200	10,0	-88,2	0,129	-13,1	-33,6	4,5	-21,9	-25,4	6,5	-12,6	-8,8
11:26:42	020:27:02	0,202	10,0	-86,5	0,131	-13,1	-33,6	4,5	-21,9	-25,3	6,5	-12,6	-8,8
11:27:43	020:28:03	0,205	10,0	-89,8	0,133	-13,1	-33,6	4,5	-21,9	-25,3	6,5	-12,6	-8,8
11:28:44	020:29:04	0,207	10,0	-88,1	0,135	-13,0	-33,6	4,6	-21,9	-25,3	6,5	-12,6	-8,7
11:29:40	020:30:00	0,210	10,0	-90,1	0,137	-13,0	-33,6	4,6	-21,9	-25,2	6,5	-12,5	-8,7
11:30:41	020:31:01	0,201	10,0	-88,7	0,130	-12,9	-33,5	4,6	-21,9	-25,2	6,5	-12,5	-8,6
11:31:43	020:32:03	0,203	10,0	-81,4	0,132	-12,9	-33,5	4,6	-21,9	-25,1	6,5	-12,5	-8,6

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 67 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



11:32:40	020:33:00	0,205	10,0	-89,4	0,134	-12,8	-33,6	4,6	-21,9	-25,1	6,5	-12,4	-8,5
11:33:40	020:34:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
11:33:40	020:34:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
11:33:40	020:34:00	0,208	10,0	-85,2	0,136	-12,8	-33,5	4,6	-21,9	-25,1	6,6	-12,4	-8,5
11:34:40	020:35:00	0,200	10,0	-88,5	0,129	-12,7	-33,4	4,6	-21,9	-25,0	6,6	-12,4	-8,4
11:35:42	020:36:02	0,202	10,0	-84,2	0,131	-12,7	-33,4	4,6	-21,9	-25,0	6,6	-12,4	-8,4
11:36:44	020:37:04	0,206	10,0	-88,0	0,133	-12,6	-33,4	4,6	-21,9	-24,9	6,6	-12,3	-8,3
11:37:42	020:38:02	0,207	10,0	-86,7	0,135	-12,6	-33,4	4,6	-21,9	-24,9	6,6	-12,3	-8,3
11:38:42	020:39:02	0,212	10,0	-90,0	0,137	-12,5	-33,4	4,7	-21,9	-24,8	6,6	-12,2	-8,2
11:39:45	020:40:05	0,200	10,0	-89,5	0,129	-12,5	-33,3	4,7	-21,8	-24,8	6,6	-12,2	-8,1
11:40:40	020:41:00	0,203	10,0	-90,1	0,131	-12,4	-33,3	4,7	-21,8	-24,7	6,6	-12,2	-8,1
11:41:40	020:42:00	0,204	10,0	-87,8	0,133	-12,4	-33,3	4,7	-21,9	-24,7	6,6	-12,1	-8,0
11:42:42	020:43:02	0,208	10,0	-79,7	0,136	-12,4	-33,3	4,7	-21,9	-24,7	6,6	-12,1	-8,0
11:43:44	020:44:04	0,211	10,0	-88,9	0,138	-12,3	-33,2	4,7	-21,9	-24,6	6,6	-12,1	-7,9
11:44:41	020:45:01	0,201	10,0	-92,0	0,130	-12,3	-33,2	4,7	-21,8	-24,5	6,6	-12,0	-7,9
11:45:41	020:46:01	0,203	10,0	-89,3	0,132	-12,2	-33,2	4,7	-21,8	-24,5	6,6	-12,0	-7,8
11:46:42	020:47:02	0,206	10,0	-84,7	0,134	-12,2	-33,2	4,7	-21,9	-24,5	6,6	-12,0	-7,8
11:47:43	020:48:03	0,208	10,0	-87,8	0,136	-12,1	-33,2	4,7	-21,9	-24,4	6,6	-12,0	-7,7
11:48:40	020:49:00	0,199	10,0	-84,6	0,129	-12,1	-33,1	4,7	-21,8	-24,4	6,6	-11,9	-7,7
11:49:40	020:50:00	0,202	10,0	-89,1	0,131	-12,0	-33,0	4,8	-21,8	-24,3	6,6	-11,9	-7,7
11:50:40	020:51:00	0,204	10,0	-87,9	0,132	-12,0	-33,0	4,8	-21,8	-24,2	6,6	-11,8	-7,6

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 68 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



11:51:41	020:52:01	0,207	10,0	-89,0	0,135	-12,0	-33,0	4,8	-21,8	-24,2	6,6	-11,8	-7,5
11:52:42	020:53:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
11:52:42	020:53:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
11:52:42	020:53:02	0,210	10,0	-87,0	0,137	-11,9	-33,0	4,8	-21,8	-24,2	6,7	-11,8	-7,5
11:53:40	020:54:00	0,200	10,0	-80,1	0,130	-11,9	-32,9	4,8	-21,8	-24,1	6,7	-11,7	-7,4
11:54:40	020:55:00	0,202	10,0	-88,6	0,131	-11,8	-32,8	4,8	-21,8	-24,1	6,7	-11,7	-7,4
11:55:42	020:56:02	0,205	10,0	-83,1	0,133	-11,8	-32,8	4,8	-21,8	-24,0	6,7	-11,6	-7,3
11:56:42	020:57:02	0,207	10,0	-88,5	0,135	-11,8	-32,8	4,8	-21,8	-24,0	6,7	-11,6	-7,3
11:57:42	020:58:02	0,210	10,0	-84,7	0,137	-11,7	-32,7	4,8	-21,8	-24,0	6,7	-11,6	-7,3
11:58:40	020:59:00	0,201	10,0	-89,0	0,130	-11,7	-32,6	4,8	-21,8	-23,9	6,7	-11,5	-7,2
11:59:42	021:00:02	0,203	10,0	-85,9	0,132	-11,6	-32,5	4,8	-21,8	-23,8	6,7	-11,5	-7,2
12:00:42	021:01:02	0,206	10,0	-89,3	0,134	-11,6	-32,5	4,9	-21,8	-23,8	6,7	-11,5	-7,1
12:01:43	021:02:03	0,208	10,0	-82,8	0,136	-11,6	-32,4	4,9	-21,8	-23,8	6,7	-11,4	-7,1
12:02:44	021:03:04	0,199	10,0	-87,9	0,129	-11,5	-32,3	4,9	-21,8	-23,7	6,7	-11,4	-7,0
12:03:44	021:04:04	0,201	10,0	-85,5	0,130	-11,5	-32,2	4,9	-21,8	-23,7	6,7	-11,4	-7,0
12:04:40	021:05:00	0,204	10,0	-91,1	0,132	-11,4	-32,0	4,9	-21,8	-23,6	6,7	-11,3	-6,9
12:05:40	021:06:00	0,205	10,0	-88,9	0,134	-11,4	-31,9	4,9	-21,8	-23,6	6,7	-11,3	-6,9
12:06:40	021:07:00	0,208	10,0	-83,7	0,136	-11,4	-31,8	4,9	-21,8	-23,6	6,7	-11,3	-6,9
12:07:42	021:08:02	0,200	10,0	-89,2	0,129	-11,3	-31,6	4,9	-21,8	-23,5	6,7	-11,2	-6,8
12:08:43	021:09:03	0,202	10,0	-88,1	0,131	-11,3	-31,5	4,9	-21,8	-23,5	6,7	-11,2	-6,8
12:09:43	021:10:03	0,204	10,0	-88,7	0,132	-11,2	-31,4	4,9	-21,8	-23,4	6,7	-11,2	-6,7

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 69 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



12:10:44	021:11:04	0,207	10,0	-84,4	0,134	-11,2	-31,4	4,9	-21,8	-23,4	6,7	-11,1	-6,7
12:11:40	021:12:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
12:11:40	021:12:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
12:11:40	021:12:00	0,209	10,0	-88,8	0,137	-11,2	-31,3	4,9	-21,8	-23,3	6,7	-11,1	-6,6
12:12:40	021:13:00	0,200	10,0	-86,8	0,129	-11,1	-31,2	5,0	-21,8	-23,3	6,8	-11,0	-6,6
12:13:41	021:14:01	0,204	10,0	-90,6	0,131	-11,0	-31,2	5,0	-21,8	-23,2	6,8	-11,0	-6,5
12:14:42	021:15:02	0,206	10,0	-88,3	0,133	-11,0	-31,1	5,0	-21,8	-23,2	6,8	-11,0	-6,5
12:15:43	021:16:03	0,207	10,0	-86,6	0,135	-11,0	-31,1	5,0	-21,8	-23,1	6,8	-10,9	-6,5
12:16:43	021:17:03	0,210	10,0	-88,9	0,137	-10,9	-31,1	5,0	-21,8	-23,1	6,8	-10,9	-6,4
12:17:43	021:18:03	0,200	10,0	-83,1	0,130	-10,9	-31,0	5,0	-21,8	-23,0	6,8	-10,8	-6,4
12:18:44	021:19:04	0,203	10,0	-87,4	0,132	-10,8	-30,9	5,0	-21,7	-23,0	6,8	-10,8	-6,3
12:19:40	021:20:00	0,205	10,0	-85,8	0,133	-10,8	-30,9	5,0	-21,8	-23,0	6,8	-10,8	-6,3
12:20:43	021:21:03	0,208	10,0	-88,7	0,135	-10,8	-30,9	5,0	-21,8	-22,9	6,8	-10,7	-6,3
12:21:42	021:22:02	0,210	10,0	-86,9	0,137	-10,7	-30,8	5,0	-21,8	-22,9	6,8	-10,7	-6,2
12:22:42	021:23:02	0,201	10,0	-87,1	0,130	-10,7	-30,7	5,1	-21,7	-22,8	6,8	-10,7	-6,2
12:23:43	021:24:03	0,203	10,0	-88,5	0,132	-10,6	-30,6	5,1	-21,8	-22,8	6,8	-10,6	-6,1
12:24:45	021:25:05	0,206	10,0	-86,8	0,134	-10,6	-30,6	5,1	-21,8	-22,7	6,8	-10,6	-6,1
12:25:40	021:26:00	0,209	10,0	-87,1	0,136	-10,6	-30,5	5,1	-21,8	-22,7	6,8	-10,6	-6,0
12:26:42	021:27:02	0,211	10,0	-83,8	0,136	-10,5	-30,4	5,1	-21,8	-22,7	6,8	-10,5	-6,0
12:27:42	021:28:02	0,202	10,0	-89,8	0,131	-10,5	-30,3	5,1	-21,7	-22,6	6,8	-10,5	-5,9
12:28:42	021:29:02	0,203	10,0	-87,2	0,132	-10,4	-30,3	5,1	-21,7	-22,6	6,9	-10,5	-5,9

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 70 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	12:29:42	021:30:02	0,207	10,0	-92,3	0,134	-10,4	-30,3	5,1	-21,8	-22,5	6,9	-10,5	-5,8	
	12:30:43	021:31:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	12:30:43	021:31:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	12:30:43	021:31:03	0,208	10,0	-88,9	0,136	-10,4	-30,2	5,1	-21,8	-22,5	6,9	-10,4	-5,8	
	12:31:44	021:32:04	0,198	10,0	-91,1	0,129	-10,3	-30,2	5,1	-21,7	-22,5	6,9	-10,4	-5,7	
	12:32:44	021:33:04	0,201	10,0	-89,6	0,130	-10,3	-30,1	5,2	-21,7	-22,4	6,9	-10,3	-5,6	
	12:33:40	021:34:00	0,203	10,0	-83,8	0,132	-10,3	-30,0	5,2	-21,7	-22,4	6,9	-10,3	-5,6	
	12:34:40	021:35:00	0,206	10,0	-89,3	0,134	-10,2	-30,0	5,2	-21,8	-22,3	6,9	-10,3	-5,6	
	12:35:40	021:36:00	0,209	10,0	-87,0	0,136	-10,2	-30,0	5,2	-21,8	-22,3	6,9	-10,3	-5,5	
	12:36:41	021:37:01	0,201	10,0	-88,4	0,129	-10,1	-29,9	5,2	-21,7	-22,3	6,9	-10,2	-5,5	
	12:37:42	021:38:02	0,201	10,0	-83,6	0,130	-10,1	-29,8	5,2	-21,7	-22,2	6,9	-10,2	-5,4	
	12:38:42	021:39:02	0,204	10,0	-89,0	0,132	-10,1	-29,8	5,2	-21,7	-22,1	6,9	-10,2	-5,4	
	12:39:42	021:40:02	0,207	10,0	-87,5	0,134	-10,0	-29,8	5,2	-21,8	-22,1	6,9	-10,1	-5,3	
	12:40:42	021:41:02	0,210	10,0	-90,9	0,137	-10,0	-29,7	5,2	-21,8	-22,1	6,9	-10,1	-5,3	
	12:41:44	021:42:04	0,199	10,0	-88,8	0,129	-9,9	-29,6	5,2	-21,7	-22,0	6,9	-10,1	-5,3	
	12:42:42	021:43:02	0,202	10,0	-83,2	0,131	-9,9	-29,6	5,2	-21,7	-21,9	6,9	-10,1	-5,2	
	12:43:42	021:44:02	0,204	10,0	-89,2	0,132	-9,9	-29,6	5,2	-21,7	-21,9	6,9	-10,0	-5,2	
	12:44:44	021:45:04	0,207	10,0	-82,1	0,134	-9,8	-29,5	5,2	-21,7	-21,9	7,0	-10,0	-5,1	
	12:45:44	021:46:04	0,209	10,0	-87,5	0,137	-9,8	-29,5	5,3	-21,7	-21,8	7,0	-10,0	-5,1	
	12:46:40	021:47:00	0,199	10,0	-85,4	0,129	-9,8	-29,4	5,3	-21,7	-21,8	7,0	-9,9	-5,0	
	12:47:41	021:48:01	0,202	10,0	-88,2	0,131	-9,7	-29,3	5,3	-21,7	-21,7	7,0	-9,9	-5,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 71 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	12:48:42	021:49:02	0,204	10,0	-86,4	0,132	-9,7	-29,3	5,3	-21,7	-21,6	7,0	-9,9	-5,0	
	12:49:43	021:50:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	12:49:43	021:50:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	12:49:43	021:50:03	0,207	10,0	-89,1	0,135	-9,7	-29,3	5,3	-21,7	-21,6	7,0	-9,8	-4,9	
	12:50:43	021:51:03	0,209	10,0	-87,7	0,136	-9,6	-29,2	5,3	-21,7	-21,6	7,0	-9,8	-4,9	
	12:51:43	021:52:03	0,200	10,0	-90,2	0,129	-9,6	-29,2	5,3	-21,7	-21,5	7,0	-9,8	-4,8	
	12:52:44	021:53:04	0,201	10,0	-88,0	0,130	-9,5	-29,1	5,3	-21,7	-21,4	7,0	-9,7	-4,8	
	12:53:44	021:54:04	0,204	10,0	-80,9	0,132	-9,5	-29,1	5,3	-21,7	-21,4	7,0	-9,7	-4,7	
	12:54:44	021:55:04	0,207	10,0	-88,8	0,135	-9,5	-29,1	5,3	-21,7	-21,3	7,0	-9,6	-4,7	
	12:55:40	021:56:00	0,209	10,0	-82,2	0,136	-9,4	-29,0	5,3	-21,7	-21,3	7,0	-9,6	-4,7	
	12:56:40	021:57:00	0,200	10,0	-89,0	0,130	-9,4	-29,0	5,3	-21,7	-21,2	7,0	-9,6	-4,6	
	12:57:42	021:58:02	0,202	10,0	-86,6	0,131	-9,4	-28,9	5,4	-21,7	-21,1	7,0	-9,6	-4,6	
	12:58:44	021:59:04	0,205	10,0	-91,1	0,133	-9,3	-28,9	5,4	-21,7	-21,1	7,0	-9,5	-4,5	
	12:59:40	022:00:00	0,208	10,0	-88,8	0,135	-9,3	-28,9	5,4	-21,7	-21,0	7,1	-9,5	-4,5	
	13:00:43	022:01:03	0,210	10,0	-81,6	0,137	-9,3	-28,8	5,4	-21,7	-20,9	7,1	-9,5	-4,4	
	13:01:45	022:02:05	0,200	10,0	-88,0	0,129	-9,2	-28,7	5,4	-21,7	-20,9	7,1	-9,5	-4,4	
	13:02:44	022:03:04	0,201	10,0	-84,4	0,130	-9,2	-28,7	5,4	-21,6	-20,8	7,1	-9,4	-4,3	
	13:03:45	022:04:05	0,204	10,0	-88,1	0,132	-9,1	-28,7	5,4	-21,7	-20,7	7,1	-9,4	-4,3	
	13:04:42	022:05:02	0,205	10,0	-86,1	0,134	-9,1	-28,6	5,4	-21,7	-20,6	7,1	-9,4	-4,2	
	13:05:41	022:06:01	0,209	10,0	-89,5	0,136	-9,1	-28,6	5,4	-21,7	-20,5	7,1	-9,4	-4,2	
	13:06:43	022:07:03	0,200	10,0	-88,8	0,129	-9,1	-28,5	5,4	-21,6	-20,5	7,1	-9,3	-4,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 72 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



13:07:44	022:08:04	0,202	10,0	-81,3	0,131	-9,0	-28,4	5,5	-21,6	-20,4	7,1	-9,3	-4,1
13:08:40	022:09:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
13:08:40	022:09:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
13:08:40	022:09:00	0,204	10,0	-88,8	0,132	-9,0	-28,4	5,5	-21,6	-20,3	7,1	-9,3	-4,1
13:09:41	022:10:01	0,206	10,0	-83,0	0,134	-8,9	-28,4	5,5	-21,7	-20,3	7,1	-9,2	-4,0
13:10:45	022:11:05	0,210	10,0	-88,2	0,137	-8,9	-28,3	5,5	-21,7	-20,2	7,1	-9,2	-4,0
13:11:40	022:12:00	0,199	10,0	-85,2	0,129	-8,9	-28,3	5,5	-21,6	-20,1	7,1	-9,2	-4,0
13:12:43	022:13:03	0,202	10,0	-88,7	0,131	-8,9	-28,2	5,5	-21,6	-20,1	7,1	-9,1	-3,9
13:13:45	022:14:05	0,204	10,0	-86,6	0,132	-8,8	-28,2	5,5	-21,6	-20,0	7,1	-9,1	-3,9
13:14:42	022:15:02	0,207	10,0	-89,6	0,134	-8,8	-28,2	5,5	-21,6	-20,0	7,1	-9,1	-3,8
13:15:42	022:16:02	0,209	10,0	-88,4	0,136	-8,8	-28,1	5,5	-21,6	-19,9	7,1	-9,1	-3,8
13:16:43	022:17:03	0,199	10,0	-89,0	0,129	-8,7	-28,1	5,5	-21,6	-19,9	7,2	-9,1	-3,7
13:17:40	022:18:00	0,201	10,0	-87,1	0,130	-8,7	-28,0	5,5	-21,6	-19,8	7,2	-9,0	-3,7
13:18:43	022:19:03	0,204	10,0	-89,2	0,132	-8,6	-28,0	5,5	-21,6	-19,7	7,2	-9,0	-3,7
13:19:43	022:20:03	0,206	10,0	-87,5	0,134	-8,6	-28,0	5,6	-21,6	-19,7	7,2	-9,0	-3,6
13:20:44	022:21:04	0,208	10,0	-90,5	0,136	-8,6	-27,9	5,6	-21,6	-19,6	7,2	-9,0	-3,6
13:21:40	022:22:00	0,210	10,0	-88,5	0,138	-8,5	-27,9	5,6	-21,6	-19,6	7,2	-8,9	-3,5
13:22:41	022:23:01	0,202	10,0	-81,3	0,130	-8,5	-27,8	5,6	-21,6	-19,5	7,2	-8,9	-3,5
13:23:42	022:24:02	0,204	10,0	-88,6	0,132	-8,4	-27,8	5,6	-21,6	-19,5	7,2	-8,9	-3,4
13:24:44	022:25:04	0,206	10,0	-82,1	0,134	-8,4	-27,8	5,6	-21,6	-19,4	7,2	-8,8	-3,4
13:25:40	022:26:00	0,208	10,0	-89,3	0,136	-8,4	-27,7	5,6	-21,6	-19,4	7,2	-8,8	-3,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 73 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	13:26:41	022:27:01	0,210	10,0	-86,8	0,138	-8,3	-27,6	5,6	-21,5	-19,3	7,2	-8,8	-3,3	
	13:27:44	022:28:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	13:27:44	022:28:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	13:27:44	022:28:04	0,201	10,0	-89,7	0,130	-8,3	-27,6	5,6	-21,5	-19,2	7,2	-8,7	-3,3	
	13:28:40	022:29:00	0,203	10,0	-87,5	0,132	-8,3	-27,5	5,6	-21,5	-19,2	7,2	-8,7	-3,2	
	13:29:41	022:30:01	0,206	10,0	-90,9	0,134	-8,2	-27,5	5,6	-21,6	-19,1	7,2	-8,7	-3,2	
	13:30:42	022:31:02	0,207	10,0	-90,0	0,135	-8,2	-27,5	5,6	-21,6	-19,1	7,2	-8,7	-3,2	
	13:31:44	022:32:04	0,211	10,0	-90,5	0,138	-8,1	-27,4	5,7	-21,6	-19,0	7,2	-8,6	-3,1	
	13:32:42	022:33:02	0,201	10,0	-89,0	0,130	-8,1	-27,3	5,7	-21,5	-19,0	7,3	-8,6	-3,1	
	13:33:41	022:34:01	0,203	10,0	-85,0	0,131	-8,1	-27,3	5,7	-21,5	-18,9	7,3	-8,6	-3,0	
	13:34:42	022:35:02	0,205	10,0	-89,2	0,133	-8,1	-27,3	5,7	-21,5	-18,9	7,3	-8,6	-3,0	
	13:35:42	022:36:02	0,208	10,0	-87,5	0,135	-8,0	-27,3	5,7	-21,5	-18,9	7,3	-8,5	-3,0	
	13:36:44	022:37:04	0,211	10,0	-89,3	0,137	-8,0	-27,2	5,7	-21,5	-18,8	7,3	-8,5	-2,9	
	13:37:44	022:38:04	0,200	10,0	-84,0	0,130	-7,9	-27,1	5,7	-21,5	-18,7	7,3	-8,5	-2,9	
	13:38:41	022:39:01	0,202	10,0	-88,7	0,131	-7,9	-27,1	5,7	-21,5	-18,7	7,3	-8,4	-2,8	
	13:39:43	022:40:03	0,204	10,0	-84,2	0,133	-7,9	-27,1	5,7	-21,5	-18,6	7,3	-8,4	-2,8	
	13:40:43	022:41:03	0,208	10,0	-88,1	0,135	-7,8	-27,1	5,7	-21,5	-18,6	7,3	-8,4	-2,8	
	13:41:40	022:42:00	0,210	10,0	-86,0	0,137	-7,8	-27,0	5,7	-21,5	-18,5	7,3	-8,3	-2,7	
	13:42:43	022:43:03	0,200	10,0	-89,5	0,129	-7,8	-26,9	5,8	-21,5	-18,5	7,3	-8,3	-2,7	
	13:43:43	022:44:03	0,202	10,0	-88,7	0,131	-7,7	-26,9	5,8	-21,5	-18,4	7,3	-8,3	-2,7	
	13:44:40	022:45:00	0,204	10,0	-81,4	0,132	-7,7	-26,9	5,8	-21,5	-18,4	7,3	-8,3	-2,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 74 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	13:45:43	022:46:03	0,207	10,0	-88,4	0,134	-7,7	-26,9	5,8	-21,5	-18,3	7,3	-8,2	-2,6	
	13:46:44	022:47:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	13:46:44	022:47:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	13:46:44	022:47:04	0,208	10,0	-83,8	0,136	-7,6	-26,8	5,8	-21,5	-18,3	7,4	-8,2	-2,6	
	13:47:40	022:48:00	0,199	10,0	-89,7	0,129	-7,6	-26,7	5,8	-21,4	-18,2	7,4	-8,2	-2,5	
	13:48:41	022:49:01	0,201	10,0	-87,1	0,130	-7,6	-26,7	5,8	-21,4	-18,2	7,4	-8,1	-2,5	
	13:49:43	022:50:03	0,203	10,0	-92,9	0,132	-7,5	-26,6	5,8	-21,5	-18,1	7,4	-8,1	-2,4	
	13:50:41	022:51:01	0,206	10,0	-89,2	0,134	-7,5	-26,6	5,8	-21,5	-18,1	7,4	-8,1	-2,4	
	13:51:40	022:52:00	0,208	10,0	-81,6	0,136	-7,5	-26,6	5,8	-21,5	-18,0	7,4	-8,1	-2,4	
	13:52:41	022:53:01	0,210	10,0	-87,3	0,135	-7,4	-26,5	5,8	-21,5	-18,0	7,4	-8,0	-2,3	
	13:53:42	022:54:02	0,201	10,0	-84,4	0,130	-7,4	-26,5	5,8	-21,4	-17,9	7,4	-8,0	-2,3	
	13:54:41	022:55:01	0,203	10,0	-88,8	0,131	-7,4	-26,4	5,9	-21,4	-17,8	7,4	-8,0	-2,2	
	13:55:43	022:56:03	0,205	10,0	-86,8	0,133	-7,3	-26,4	5,9	-21,4	-17,8	7,4	-8,0	-2,2	
	13:56:40	022:57:00	0,207	10,0	-90,9	0,135	-7,3	-26,4	5,9	-21,4	-17,7	7,4	-7,9	-2,1	
	13:57:41	022:58:01	0,210	10,0	-88,2	0,137	-7,3	-26,3	5,9	-21,4	-17,7	7,4	-7,9	-2,1	
	13:58:43	022:59:03	0,201	10,0	-88,8	0,130	-7,2	-26,3	5,9	-21,4	-17,6	7,4	-7,9	-2,1	
	13:59:44	023:00:04	0,201	10,0	-88,8	0,131	-7,2	-26,2	5,9	-21,4	-17,6	7,4	-7,8	-2,0	
	14:00:40	023:01:00	0,204	10,0	-83,3	0,133	-7,2	-26,2	5,9	-21,4	-17,5	7,4	-7,8	-2,0	
	14:01:42	023:02:02	0,207	10,0	-88,6	0,135	-7,1	-26,2	5,9	-21,4	-17,5	7,4	-7,8	-1,9	
	14:02:40	023:03:00	0,208	10,0	-86,9	0,137	-7,1	-26,1	5,9	-21,4	-17,4	7,4	-7,8	-1,9	
	14:03:42	023:04:02	0,200	10,0	-82,8	0,130	-7,1	-26,1	5,9	-21,4	-17,4	7,5	-7,7	-1,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 75 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	14:04:42	023:05:02	0,201	10,0	-88,9	0,131	-7,0	-26,0	5,9	-21,4	-17,3	7,5	-7,7	-1,9	
	14:05:42	023:06:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	14:05:42	023:06:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	14:05:42	023:06:02	0,204	10,0	-83,4	0,132	-7,0	-26,0	6,0	-21,4	-17,3	7,5	-7,7	-1,8	
	14:06:44	023:07:04	0,206	10,0	-89,0	0,134	-7,0	-26,0	6,0	-21,4	-17,2	7,5	-7,6	-1,8	
	14:07:40	023:08:00	0,208	10,0	-83,9	0,136	-7,0	-25,9	6,0	-21,4	-17,2	7,5	-7,6	-1,8	
	14:08:41	023:09:01	0,199	10,0	-89,4	0,129	-6,9	-25,9	6,0	-21,4	-17,1	7,5	-7,6	-1,7	
	14:09:43	023:10:03	0,200	10,0	-86,8	0,130	-6,9	-25,8	6,0	-21,3	-17,1	7,5	-7,6	-1,7	
	14:10:40	023:11:00	0,203	10,0	-91,5	0,132	-6,9	-25,8	6,0	-21,3	-17,0	7,5	-7,5	-1,7	
	14:11:41	023:12:01	0,204	10,0	-88,9	0,133	-6,8	-25,8	6,0	-21,4	-17,0	7,5	-7,5	-1,6	
	14:12:41	023:13:01	0,208	10,0	-81,8	0,135	-6,8	-25,7	6,0	-21,4	-16,9	7,5	-7,5	-1,6	
	14:13:42	023:14:02	0,211	10,0	-88,5	0,138	-6,8	-25,7	6,0	-21,3	-16,9	7,5	-7,5	-1,6	
	14:14:44	023:15:04	0,200	10,0	-85,3	0,130	-6,7	-25,6	6,0	-21,3	-16,8	7,5	-7,4	-1,5	
	14:15:40	023:16:00	0,203	10,0	-91,7	0,131	-6,7	-25,6	6,0	-21,3	-16,8	7,5	-7,4	-1,5	
	14:16:42	023:17:02	0,205	10,0	-88,8	0,133	-6,7	-25,6	6,0	-21,3	-16,7	7,5	-7,4	-1,4	
	14:17:43	023:18:03	0,207	10,0	-79,6	0,135	-6,7	-25,5	6,0	-21,3	-16,7	7,5	-7,4	-1,4	
	14:18:40	023:19:00	0,210	10,0	-88,8	0,137	-6,6	-25,5	6,0	-21,3	-16,6	7,5	-7,3	-1,4	
	14:19:42	023:20:02	0,199	10,0	-82,2	0,129	-6,6	-25,4	6,0	-21,3	-16,6	7,5	-7,3	-1,3	
	14:20:43	023:21:03	0,201	10,0	-88,7	0,131	-6,6	-25,4	6,1	-21,3	-16,5	7,5	-7,3	-1,3	
	14:21:41	023:22:01	0,203	10,0	-86,5	0,132	-6,5	-25,4	6,1	-21,3	-16,5	7,5	-7,2	-1,3	
	14:22:42	023:23:02	0,206	10,0	-87,6	0,134	-6,5	-25,4	6,1	-21,3	-16,4	7,6	-7,2	-1,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 76 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



14:23:43	023:24:03	0,208	10,0	-89,0	0,136	-6,5	-25,3	6,1	-21,3	-16,4	7,6	-7,2	-1,2
14:24:44	023:25:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
14:24:44	023:25:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
14:24:44	023:25:04	0,210	10,0	-89,1	0,138	-6,4	-25,3	6,1	-21,3	-16,3	7,6	-7,2	-1,2
14:25:40	023:26:00	0,200	10,0	-89,1	0,130	-6,4	-25,2	6,1	-21,2	-16,2	7,6	-7,1	-1,1
14:26:40	023:27:00	0,202	10,0	-82,1	0,131	-6,3	-25,1	6,1	-21,2	-16,2	7,6	-7,1	-1,1
14:27:42	023:28:02	0,204	10,0	-88,5	0,133	-6,3	-25,1	6,1	-21,2	-16,1	7,6	-7,1	-1,0
14:28:43	023:29:03	0,207	10,0	-83,6	0,135	-6,3	-25,1	6,1	-21,3	-16,1	7,6	-7,1	-1,0
14:29:45	023:30:05	0,210	10,0	-88,8	0,137	-6,3	-25,1	6,2	-21,2	-16,1	7,6	-7,0	-1,0
14:30:40	023:31:00	0,199	10,0	-87,2	0,129	-6,2	-25,0	6,2	-21,2	-16,0	7,6	-7,0	-0,9
14:31:43	023:32:03	0,201	10,0	-90,3	0,130	-6,2	-25,0	6,2	-21,2	-15,9	7,6	-7,0	-0,9
14:32:44	023:33:04	0,204	10,0	-88,8	0,132	-6,1	-24,9	6,2	-21,2	-15,9	7,6	-7,0	-0,8
14:33:44	023:34:04	0,207	10,0	-81,5	0,134	-6,1	-24,9	6,2	-21,2	-15,9	7,6	-6,9	-0,8
14:34:43	023:35:03	0,208	10,0	-87,4	0,136	-6,1	-24,9	6,2	-21,2	-15,8	7,6	-6,9	-0,7
14:35:43	023:36:03	0,210	10,0	-84,8	0,137	-6,0	-24,8	6,2	-21,2	-15,7	7,6	-6,9	-0,7
14:36:44	023:37:04	0,200	10,0	-90,7	0,130	-6,0	-24,8	6,2	-21,1	-15,7	7,6	-6,9	-0,7
14:37:40	023:38:00	0,202	10,0	-88,7	0,131	-6,0	-24,7	6,2	-21,1	-15,6	7,7	-6,8	-0,6
14:38:42	023:39:02	0,204	10,0	-87,6	0,133	-5,9	-24,7	6,2	-21,1	-15,5	7,7	-6,8	-0,5
14:39:43	023:40:03	0,207	10,0	-89,3	0,134	-5,9	-24,7	6,2	-21,1	-15,5	7,7	-6,8	-0,5
14:40:40	023:41:00	0,208	10,0	-83,0	0,136	-5,9	-24,6	6,2	-21,1	-15,4	7,7	-6,8	-0,5
14:41:42	023:42:02	0,199	10,0	-87,6	0,129	-5,9	-24,6	6,3	-21,1	-15,4	7,7	-6,7	-0,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 77 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	14:42:44	023:43:04	0,200	10,0	-85,2	0,130	-5,8	-24,5	6,3	-21,1	-15,3	7,7	-6,7	-0,4	
	14:43:44	023:44:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	14:43:44	023:44:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	14:43:44	023:44:04	0,203	10,0	-89,9	0,132	-5,8	-24,5	6,3	-21,1	-15,2	7,7	-6,7	-0,3	
	14:44:41	023:45:01	0,204	10,0	-87,1	0,133	-5,8	-24,5	6,3	-21,1	-15,2	7,7	-6,7	-0,3	
	14:45:41	023:46:01	0,208	10,0	-83,3	0,135	-5,7	-24,4	6,3	-21,1	-15,1	7,7	-6,6	-0,3	
	14:46:44	023:47:04	0,211	10,0	-88,8	0,137	-5,7	-24,4	6,3	-21,1	-15,1	7,7	-6,6	-0,2	
	14:47:44	023:48:04	0,199	10,0	-82,6	0,129	-5,7	-24,3	6,3	-21,0	-15,0	7,7	-6,6	-0,2	
	14:48:40	023:49:00	0,201	10,0	-89,4	0,131	-5,6	-24,3	6,3	-21,0	-15,0	7,7	-6,5	-0,1	
	14:49:42	023:50:02	0,203	10,0	-86,8	0,132	-5,6	-24,2	6,3	-21,1	-14,9	7,7	-6,5	-0,1	
	14:50:40	023:51:00	0,205	10,0	-92,9	0,134	-5,6	-24,2	6,3	-21,1	-14,9	7,7	-6,5	-0,3	
	14:51:43	023:52:03	0,207	10,0	-89,1	0,136	-5,6	-24,2	6,3	-21,1	-14,8	7,7	-6,5	-0,5	
	14:52:44	023:53:04	0,210	10,0	-81,4	0,138	-5,5	-24,1	6,3	-21,0	-14,8	7,7	-6,4	-0,6	
	14:53:40	023:54:00	0,200	10,0	-87,7	0,130	-5,5	-24,1	6,4	-21,0	-14,7	7,7	-6,4	-0,5	
	14:54:41	023:55:01	0,201	10,0	-83,8	0,131	-5,5	-24,0	6,4	-21,0	-14,7	7,7	-6,4	-0,5	
	14:55:41	023:56:01	0,204	10,0	-89,6	0,133	-5,4	-24,0	6,4	-21,0	-14,6	7,7	-6,3	-0,4	
	14:56:43	023:57:03	0,206	10,0	-86,6	0,134	-5,4	-24,0	6,4	-21,0	-14,6	7,7	-6,3	-0,3	
	14:57:44	023:58:04	0,208	10,0	-92,6	0,136	-5,4	-23,9	6,4	-21,0	-14,5	7,8	-6,3	-0,3	
	14:58:41	023:59:01	0,200	10,0	-89,1	0,129	-5,3	-23,9	6,4	-21,0	-14,5	7,8	-6,3	-0,2	
	14:59:42	024:00:02	0,200	10,0	-85,4	0,130	-5,3	-23,8	6,4	-21,0	-14,4	7,8	-6,2	-0,2	
	15:00:42	024:01:02	0,202	10,0	-88,6	0,131	-5,3	-23,8	6,4	-21,0	-14,4	7,8	-6,2	-0,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus Signature: Page: 78 / 209
 Document issued by:

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



15:01:40	024:02:00	0,204	10,0	-83,9	0,133	-5,3	-23,8	6,4	-21,0	-14,3	7,8	-6,2	0,0	
15:02:41	024:03:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
15:02:41	024:03:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
15:02:41	024:03:01	0,208	10,0	-88,9	0,135	-5,2	-23,7	6,4	-21,0	-14,3	7,8	-6,2	0,1	
15:03:41	024:04:01	0,208	10,0	-86,8	0,136	-5,2	-23,7	6,4	-20,9	-14,2	7,8	-6,2	0,1	
15:04:42	024:05:02	0,199	10,0	-90,9	0,129	-5,2	-23,6	6,4	-20,9	-14,2	7,8	-6,1	0,2	
15:05:43	024:06:03	0,200	10,0	-87,9	0,130	-5,1	-23,6	6,5	-20,9	-14,1	7,8	-6,1	0,3	
15:06:44	024:07:04	0,203	10,0	-84,0	0,132	-5,1	-23,5	6,5	-20,9	-14,0	7,8	-6,1	0,5	
15:07:44	024:08:04	0,205	10,0	-89,2	0,133	-5,1	-23,5	6,5	-20,9	-14,0	7,8	-6,1	0,5	
15:08:44	024:09:04	0,207	10,0	-82,9	0,135	-5,0	-23,5	6,5	-20,9	-14,0	7,8	-6,0	0,6	
15:09:44	024:10:04	0,210	10,0	-88,2	0,137	-5,0	-23,4	6,5	-20,9	-13,9	7,8	-6,0	0,7	
15:10:41	024:11:01	0,199	10,0	-85,9	0,129	-5,0	-23,3	6,5	-20,8	-13,9	7,8	-6,0	0,8	
15:11:41	024:12:01	0,201	10,0	-91,9	0,131	-4,9	-23,3	6,5	-20,8	-13,8	7,8	-6,0	0,9	
15:12:44	024:13:04	0,203	10,0	-89,3	0,132	-4,9	-23,2	6,5	-20,8	-13,7	7,8	-5,9	0,9	
15:13:42	024:14:02	0,205	10,0	-86,3	0,134	-4,9	-23,2	6,5	-20,9	-13,7	7,8	-5,9	1,0	
15:14:43	024:15:03	0,208	10,0	-89,0	0,136	-4,9	-23,2	6,5	-20,8	-13,7	7,8	-5,9	1,2	
15:15:44	024:16:04	0,210	10,0	-82,0	0,137	-4,8	-23,1	6,5	-20,8	-13,6	7,8	-5,9	1,3	
15:16:41	024:17:01	0,199	10,0	-90,0	0,129	-4,8	-23,0	6,6	-20,8	-13,5	7,8	-5,8	1,5	
15:17:42	024:18:02	0,201	10,0	-86,1	0,130	-4,7	-23,0	6,6	-20,8	-13,5	7,9	-5,8	1,7	
15:18:42	024:19:02	0,203	10,0	-92,7	0,132	-4,7	-23,0	6,6	-20,8	-13,4	7,9	-5,8	1,8	
15:19:40	024:20:00	0,204	10,0	-88,8	0,134	-4,7	-23,0	6,6	-20,8	-13,4	7,9	-5,8	1,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 79 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



15:20:41	024:21:01	0,208	10,0	-85,8	0,136	-4,7	-22,9	6,6	-20,7	-13,3	7,9	-5,7	2,1
15:21:44	024:22:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:21:44	024:22:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:21:44	024:22:04	0,210	10,0	-89,6	0,138	-4,6	-22,9	6,6	-20,7	-13,3	7,9	-5,7	2,2
15:22:44	024:23:04	0,199	10,0	-87,6	0,130	-4,6	-22,8	6,6	-20,7	-13,2	7,9	-5,7	2,3
15:23:41	024:24:01	0,201	10,0	-89,6	0,131	-4,5	-22,8	6,6	-20,7	-13,2	7,9	-5,6	2,4
15:24:43	024:25:03	0,203	10,0	-82,5	0,132	-4,5	-22,7	6,6	-20,7	-13,1	7,9	-5,6	2,5
15:25:43	024:26:03	0,206	10,0	-87,8	0,134	-4,5	-22,7	6,6	-20,7	-13,1	7,9	-5,6	2,6
15:26:43	024:27:03	0,208	10,0	-85,1	0,136	-4,5	-22,7	6,6	-20,7	-13,0	7,9	-5,6	2,8
15:27:41	024:28:01	0,201	10,0	-89,5	0,130	-4,4	-22,6	6,6	-20,7	-13,0	7,9	-5,5	2,9
15:28:41	024:29:01	0,199	10,0	-87,2	0,129	-4,4	-22,5	6,6	-20,6	-12,9	7,9	-5,5	3,0
15:29:43	024:30:03	0,202	10,0	-91,2	0,131	-4,4	-22,5	6,6	-20,6	-12,9	7,9	-5,5	3,2
15:30:44	024:31:04	0,205	10,0	-89,0	0,133	-4,4	-22,5	6,7	-20,6	-12,8	7,9	-5,5	3,3
15:31:40	024:32:00	0,206	10,0	-81,4	0,134	-4,3	-22,5	6,7	-20,6	-12,8	7,9	-5,5	3,4
15:32:43	024:33:03	0,208	10,0	-87,6	0,137	-4,3	-22,4	6,7	-20,6	-12,8	7,9	-5,4	3,6
15:33:40	024:34:00	0,200	10,0	-85,0	0,129	-4,3	-22,4	6,7	-20,6	-12,7	7,9	-5,4	3,7
15:34:43	024:35:03	0,200	10,0	-90,5	0,130	-4,2	-22,3	6,7	-20,5	-12,6	7,9	-5,3	3,8
15:35:43	024:36:03	0,202	10,0	-87,4	0,131	-4,2	-22,2	6,7	-20,5	-12,6	7,9	-5,3	3,9
15:36:42	024:37:02	0,204	10,0	-84,8	0,133	-4,2	-22,2	6,7	-20,6	-12,6	7,9	-5,3	4,0
15:37:42	024:38:02	0,207	10,0	-89,0	0,135	-4,2	-22,2	6,7	-20,5	-12,5	7,9	-5,3	4,1
15:38:43	024:39:03	0,209	10,0	-82,7	0,136	-4,1	-22,1	6,7	-20,5	-12,5	7,9	-5,2	4,3

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 80 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



15:39:44	024:40:04	0,198	10,0	-87,9	0,129	-4,1	-22,1	6,7	-20,5	-12,4	7,9	-5,2	4,4	
15:40:40	024:41:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
15:40:40	024:41:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
15:40:40	024:41:00	0,200	10,0	-85,5	0,130	-4,1	-22,0	6,7	-20,5	-12,4	7,9	-5,2	4,5	
15:41:41	024:42:01	0,202	10,0	-90,0	0,131	-4,0	-22,0	6,8	-20,5	-12,3	8,0	-5,2	4,7	
15:42:44	024:43:04	0,204	10,0	-87,5	0,133	-4,0	-21,9	6,8	-20,5	-12,3	8,0	-5,1	4,8	
15:43:41	024:44:01	0,207	10,0	-82,6	0,135	-4,0	-21,9	6,8	-20,5	-12,3	8,0	-5,1	5,0	
15:44:42	024:45:02	0,208	10,0	-88,8	0,137	-3,9	-21,8	6,8	-20,5	-12,2	8,0	-5,1	5,1	
15:45:43	024:46:03	0,199	10,0	-80,5	0,129	-3,9	-21,8	6,8	-20,4	-12,2	8,0	-5,1	5,2	
15:46:44	024:47:04	0,200	10,0	-87,4	0,130	-3,9	-21,7	6,8	-20,4	-12,1	8,0	-5,0	5,3	
15:47:44	024:48:04	0,202	10,0	-83,8	0,131	-3,8	-21,7	6,8	-20,4	-12,1	8,0	-5,0	5,4	
15:48:44	024:49:04	0,204	10,0	-89,6	0,133	-3,8	-21,6	6,8	-20,4	-12,0	8,0	-5,0	5,5	
15:49:44	024:50:04	0,207	10,0	-87,9	0,135	-3,8	-21,6	6,8	-20,4	-12,0	8,0	-5,0	5,6	
15:50:40	024:51:00	0,209	10,0	-80,6	0,137	-3,8	-21,5	6,8	-20,4	-11,9	8,0	-4,9	5,7	
15:51:40	024:52:00	0,199	10,0	-88,3	0,129	-3,7	-21,4	6,8	-20,4	-11,9	8,0	-4,9	5,8	
15:52:42	024:53:02	0,200	10,0	-82,9	0,130	-3,7	-21,4	6,8	-20,4	-11,8	8,0	-4,9	5,9	
15:53:42	024:54:02	0,202	10,0	-90,1	0,132	-3,7	-21,3	6,8	-20,4	-11,8	8,0	-4,9	6,0	
15:54:44	024:55:04	0,204	10,0	-87,2	0,133	-3,7	-21,3	6,8	-20,4	-11,8	8,0	-4,8	6,2	
15:55:40	024:56:00	0,207	10,0	-89,0	0,135	-3,7	-21,2	6,9	-20,4	-11,7	8,0	-4,8	6,3	
15:56:43	024:57:03	0,209	10,0	-88,9	0,137	-3,6	-21,2	6,9	-20,4	-11,7	8,0	-4,8	6,5	
15:57:45	024:58:05	0,199	10,0	-80,9	0,129	-3,6	-21,1	6,9	-20,4	-11,6	8,0	-4,7	6,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 81 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	15:58:41	024:59:01	0,200	10,0	-87,8	0,130	-3,5	-21,1	6,9	-20,3	-11,6	8,1	-4,7	6,8	
	15:59:41	025:00:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	15:59:41	025:00:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	15:59:41	025:00:01	0,202	10,0	-84,8	0,131	-3,5	-21,0	6,9	-20,3	-11,5	8,1	-4,7	6,9	
	16:00:42	025:01:02	0,204	10,0	-90,3	0,133	-3,5	-21,0	6,9	-20,4	-11,5	8,1	-4,7	7,1	
	16:01:42	025:02:02	0,207	10,0	-86,9	0,135	-3,5	-21,0	6,9	-20,4	-11,4	8,1	-4,6	7,2	
	16:02:43	025:03:03	0,208	10,0	-91,0	0,137	-3,4	-20,9	6,9	-20,3	-11,4	8,1	-4,6	7,4	
	16:03:43	025:04:03	0,198	10,0	-89,3	0,129	-3,4	-20,8	6,9	-20,3	-11,4	8,1	-4,6	7,5	
	16:04:43	025:05:03	0,199	10,0	-82,1	0,130	-3,3	-20,8	6,9	-20,3	-11,3	8,1	-4,6	7,6	
	16:05:41	025:06:01	0,201	10,0	-88,3	0,131	-3,3	-20,7	6,9	-20,3	-11,3	8,1	-4,5	7,7	
	16:06:42	025:07:02	0,203	10,0	-83,3	0,132	-3,3	-20,7	7,0	-20,3	-11,2	8,1	-4,5	7,8	
	16:07:43	025:08:03	0,205	10,0	-89,6	0,134	-3,3	-20,7	7,0	-20,3	-11,2	8,1	-4,5	7,9	
	16:08:44	025:09:04	0,208	10,0	-87,0	0,136	-3,3	-20,6	7,0	-20,3	-11,2	8,1	-4,5	7,9	
	16:09:40	025:10:00	0,209	10,0	-83,0	0,138	-3,3	-20,6	7,0	-20,3	-11,2	8,1	-4,4	8,0	
	16:10:42	025:11:02	0,199	10,0	-89,0	0,129	-3,3	-20,5	7,0	-20,2	-11,1	8,1	-4,4	8,1	
	16:11:44	025:12:04	0,201	10,0	-80,6	0,131	-3,2	-20,4	7,0	-20,2	-11,0	8,1	-4,4	8,2	
	16:12:42	025:13:02	0,203	10,0	-88,1	0,132	-3,2	-20,4	7,0	-20,2	-11,0	8,1	-4,4	8,3	
	16:13:42	025:14:02	0,204	10,0	-85,0	0,134	-3,1	-20,4	7,0	-20,2	-11,0	8,1	-4,3	8,3	
	16:14:40	025:15:00	0,208	10,0	-85,3	0,136	-3,1	-20,3	7,0	-20,2	-10,9	8,1	-4,3	8,3	
	16:15:42	025:16:02	0,210	10,0	-88,3	0,138	-3,1	-20,3	7,0	-20,2	-10,9	8,1	-4,3	8,4	
	16:16:44	025:17:04	0,199	10,0	-82,2	0,129	-3,1	-20,2	7,0	-20,2	-10,8	8,1	-4,3	8,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 82 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



16:17:40	025:18:00	0,200	10,0	-88,4	0,131	-3,0	-20,1	7,0	-20,2	-10,8	8,1	-4,2	8,5
16:18:42	025:19:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
16:18:42	025:19:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
16:18:42	025:19:02	0,203	10,0	-82,6	0,132	-3,0	-20,1	7,0	-20,2	-10,8	8,1	-4,2	8,5
16:19:43	025:20:03	0,205	10,0	-89,7	0,134	-3,0	-20,1	7,1	-20,2	-10,7	8,1	-4,2	8,6
16:20:43	025:21:03	0,206	10,0	-87,6	0,135	-3,0	-20,0	7,1	-20,1	-10,7	8,1	-4,2	8,6
16:21:44	025:22:04	0,209	10,0	-91,8	0,137	-2,9	-20,0	7,1	-20,1	-10,6	8,1	-4,1	8,6
16:22:42	025:23:02	0,198	10,0	-89,1	0,129	-2,9	-19,9	7,1	-20,1	-10,6	8,2	-4,1	8,7
16:23:42	025:24:02	0,199	10,0	-80,9	0,130	-2,9	-19,9	7,1	-20,1	-10,5	8,2	-4,1	8,7
16:24:43	025:25:03	0,201	10,0	-87,4	0,132	-2,9	-19,8	7,1	-20,1	-10,5	8,2	-4,1	8,7
16:25:43	025:26:03	0,203	10,0	-84,0	0,133	-2,8	-19,8	7,1	-20,1	-10,5	8,2	-4,1	8,7
16:26:43	025:27:03	0,205	10,0	-90,6	0,135	-2,8	-19,8	7,1	-20,1	-10,5	8,2	-4,0	8,8
16:27:40	025:28:00	0,207	10,0	-89,0	0,136	-2,8	-19,7	7,1	-20,1	-10,4	8,2	-4,0	8,8
16:28:40	025:29:00	0,198	10,0	-90,2	0,129	-2,8	-19,7	7,1	-20,0	-10,4	8,2	-4,0	8,8
16:29:42	025:30:02	0,198	10,0	-89,7	0,129	-2,7	-19,6	7,1	-20,0	-10,3	8,2	-4,0	8,9
16:30:42	025:31:02	0,200	10,0	-92,5	0,131	-2,7	-19,5	7,1	-20,0	-10,3	8,2	-3,9	8,9
16:31:44	025:32:04	0,202	10,0	-89,3	0,132	-2,7	-19,5	7,1	-20,0	-10,3	8,2	-3,9	8,9
16:32:43	025:33:03	0,203	10,0	-81,6	0,134	-2,7	-19,5	7,2	-20,0	-10,2	8,2	-3,9	8,9
16:33:43	025:34:03	0,207	10,0	-88,3	0,135	-2,7	-19,4	7,2	-20,0	-10,2	8,2	-3,9	8,9
16:34:40	025:35:00	0,210	10,0	-83,6	0,137	-2,6	-19,4	7,2	-20,0	-10,2	8,2	-3,9	9,0
16:35:40	025:36:00	0,200	10,0	-90,8	0,129	-2,6	-19,3	7,2	-19,9	-10,1	8,2	-3,8	9,0

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 83 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



16:36:41	025:37:01	0,199	10,0	-87,1	0,130	-2,6	-19,2	7,2	-19,9	-10,0	8,2	-3,8	9,0
16:37:41	025:38:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
16:37:41	025:38:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
16:37:41	025:38:01	0,201	10,0	-79,2	0,132	-2,5	-19,2	7,2	-19,9	-10,0	8,2	-3,8	9,0
16:38:41	025:39:01	0,204	10,0	-87,5	0,133	-2,5	-19,1	7,2	-19,9	-10,0	8,2	-3,8	9,1
16:39:41	025:40:01	0,206	10,0	-83,5	0,135	-2,5	-19,1	7,2	-19,9	-9,9	8,2	-3,8	9,1
16:40:41	025:41:01	0,207	10,0	-90,6	0,137	-2,5	-19,0	7,2	-19,9	-9,9	8,2	-3,7	9,1
16:41:41	025:42:01	0,198	10,0	-87,8	0,129	-2,4	-18,9	7,2	-19,9	-9,8	8,2	-3,7	9,1
16:42:40	025:43:00	0,199	10,0	-80,7	0,130	-2,4	-18,9	7,2	-19,8	-9,8	8,2	-3,7	9,1
16:43:41	025:44:01	0,202	10,0	-88,9	0,131	-2,4	-18,8	7,2	-19,8	-9,8	8,2	-3,6	9,2
16:44:43	025:45:03	0,202	10,0	-81,9	0,132	-2,4	-18,8	7,3	-19,8	-9,7	8,2	-3,6	9,2
16:45:44	025:46:04	0,205	10,0	-89,1	0,134	-2,4	-18,8	7,3	-19,8	-9,7	8,2	-3,6	9,2
16:46:44	025:47:04	0,207	10,0	-86,1	0,136	-2,3	-18,7	7,3	-19,8	-9,7	8,2	-3,6	9,2
16:47:42	025:48:02	0,209	10,0	-82,7	0,138	-2,3	-18,6	7,3	-19,8	-9,6	8,2	-3,6	9,2
16:48:42	025:49:02	0,198	10,0	-89,0	0,129	-2,3	-18,5	7,3	-19,8	-9,6	8,3	-3,5	9,2
16:49:42	025:50:02	0,198	10,0	-81,7	0,130	-2,2	-18,5	7,3	-19,7	-9,5	8,3	-3,5	9,3
16:50:44	025:51:04	0,201	10,0	-88,7	0,132	-2,2	-18,5	7,3	-19,8	-9,5	8,3	-3,5	9,3
16:51:40	025:52:00	0,203	10,0	-86,4	0,133	-2,2	-18,4	7,3	-19,8	-9,5	8,3	-3,5	9,3
16:52:42	025:53:02	0,206	10,0	-91,2	0,135	-2,2	-18,4	7,3	-19,7	-9,5	8,3	-3,5	9,3
16:53:42	025:54:02	0,208	10,0	-88,2	0,136	-2,2	-18,3	7,3	-19,7	-9,4	8,3	-3,4	9,3
16:54:41	025:55:01	0,198	10,0	-80,9	0,129	-2,1	-18,3	7,3	-19,7	-9,4	8,3	-3,4	9,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 84 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	16:55:42	025:56:02	0,199	10,0	-89,0	0,130	-2,1	-18,2	7,3	-19,6	-9,3	8,3	-3,4	9,3	
	16:56:42	025:57:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	16:56:42	025:57:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	16:56:42	025:57:02	0,200	10,0	-82,4	0,131	-2,1	-18,2	7,3	-19,6	-9,3	8,3	-3,4	9,4	
	16:57:43	025:58:03	0,202	10,0	-89,3	0,132	-2,1	-18,1	7,3	-19,7	-9,3	8,3	-3,4	9,4	
	16:58:45	025:59:05	0,204	10,0	-86,6	0,134	-2,1	-18,1	7,3	-19,6	-9,2	8,3	-3,3	9,4	
	16:59:42	026:00:02	0,207	10,0	-92,0	0,136	-2,0	-18,0	7,4	-19,6	-9,2	8,3	-3,3	9,4	
	17:00:42	026:01:02	0,208	10,0	-89,3	0,137	-2,0	-18,0	7,4	-19,6	-9,1	8,3	-3,3	9,4	
	17:01:42	026:02:02	0,198	10,0	-81,0	0,129	-2,0	-17,9	7,4	-19,6	-9,1	8,3	-3,3	9,4	
	17:02:43	026:03:03	0,199	10,0	-88,0	0,130	-2,0	-17,8	7,4	-19,5	-9,1	8,3	-3,2	9,4	
	17:03:44	026:04:04	0,202	10,0	-85,2	0,131	-1,9	-17,8	7,4	-19,5	-9,0	8,3	-3,2	9,5	
	17:04:41	026:05:01	0,203	10,0	-91,9	0,133	-1,9	-17,8	7,4	-19,6	-9,0	8,3	-3,2	9,5	
	17:05:41	026:06:01	0,204	10,0	-88,7	0,134	-1,9	-17,7	7,4	-19,5	-9,0	8,3	-3,2	9,5	
	17:06:42	026:07:02	0,208	10,0	-79,7	0,136	-1,9	-17,7	7,4	-19,5	-8,9	8,3	-3,1	9,5	
	17:07:44	026:08:04	0,198	10,0	-88,0	0,129	-1,9	-17,6	7,4	-19,5	-8,9	8,3	-3,1	9,5	
	17:08:41	026:09:01	0,198	10,0	-83,8	0,129	-1,8	-17,5	7,4	-19,5	-8,8	8,3	-3,1	9,5	
	17:09:42	026:10:02	0,199	10,0	-89,1	0,131	-1,8	-17,5	7,5	-19,5	-8,8	8,3	-3,1	9,5	
	17:10:40	026:11:00	0,201	10,0	-86,9	0,132	-1,8	-17,4	7,5	-19,5	-8,8	8,3	-3,1	9,6	
	17:11:40	026:12:00	0,203	10,0	-81,3	0,133	-1,8	-17,4	7,5	-19,5	-8,8	8,3	-3,0	9,6	
	17:12:41	026:13:01	0,206	10,0	-88,8	0,135	-1,8	-17,3	7,5	-19,4	-8,7	8,3	-3,0	9,6	
	17:13:43	026:14:03	0,208	10,0	-80,7	0,137	-1,8	-17,2	7,5	-19,4	-8,7	8,3	-3,0	9,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 85 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



17:14:44	026:15:04	0,197	10,0	-88,7	0,129	-1,7	-17,1	7,5	-19,4	-8,6	8,3	-3,0	9,6	
17:15:41	026:16:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
17:15:41	026:16:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
17:15:41	026:16:01	0,198	10,0	-84,5	0,129	-1,7	-17,0	7,5	-19,3	-8,6	8,4	-2,9	9,6	
17:16:41	026:17:01	0,200	10,0	-90,0	0,131	-1,7	-16,9	7,5	-19,3	-8,6	8,4	-2,9	9,6	
17:17:42	026:18:02	0,201	10,0	-87,1	0,132	-1,7	-16,8	7,5	-19,3	-8,5	8,4	-2,9	9,6	
17:18:43	026:19:03	0,204	10,0	-78,3	0,134	-1,7	-16,8	7,5	-19,3	-8,5	8,4	-2,9	9,6	
17:19:44	026:20:04	0,206	10,0	-88,8	0,136	-1,6	-16,7	7,5	-19,3	-8,5	8,4	-2,9	9,7	
17:20:44	026:21:04	0,208	10,0	-81,2	0,137	-1,6	-16,6	7,5	-19,3	-8,4	8,4	-2,9	9,7	
17:21:42	026:22:02	0,198	10,0	-89,2	0,129	-1,6	-16,5	7,5	-19,2	-8,4	8,4	-2,8	9,7	
17:22:43	026:23:03	0,198	10,0	-85,5	0,130	-1,6	-16,5	7,5	-19,2	-8,4	8,4	-2,8	9,7	
17:23:44	026:24:04	0,201	10,0	-85,1	0,131	-1,5	-16,4	7,5	-19,2	-8,3	8,4	-2,8	9,7	
17:24:45	026:25:05	0,202	10,0	-88,6	0,132	-1,5	-16,4	7,5	-19,2	-8,3	8,4	-2,8	9,7	
17:25:41	026:26:01	0,204	10,0	-80,2	0,134	-1,5	-16,3	7,6	-19,1	-8,3	8,4	-2,8	9,7	
17:26:41	026:27:01	0,205	10,0	-87,8	0,136	-1,5	-16,2	7,6	-19,1	-8,3	8,4	-2,7	9,7	
17:27:41	026:28:01	0,208	10,0	-85,0	0,137	-1,5	-16,1	7,6	-19,1	-8,2	8,4	-2,7	9,7	
17:28:41	026:29:01	0,198	10,0	-90,6	0,129	-1,5	-16,0	7,6	-19,0	-8,1	8,4	-2,7	9,7	
17:29:42	026:30:02	0,198	10,0	-86,4	0,130	-1,5	-16,0	7,6	-19,0	-8,1	8,4	-2,6	9,8	
17:30:43	026:31:03	0,201	10,0	-78,9	0,131	-1,4	-15,9	7,6	-19,0	-8,1	8,4	-2,6	9,8	
17:31:44	026:32:04	0,202	10,0	-88,1	0,133	-1,5	-15,9	7,6	-19,0	-8,1	8,4	-2,6	9,8	
17:32:40	026:33:00	0,204	10,0	-83,8	0,134	-1,5	-15,8	7,6	-19,0	-8,0	8,4	-2,6	9,8	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 86 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



17:33:41	026:34:01	0,206	10,0	-89,8	0,136	-1,4	-15,7	7,6	-18,9	-8,0	8,4	-2,6	9,8	
17:34:42	026:35:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
17:34:42	026:35:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
17:34:42	026:35:02	0,207	10,0	-87,4	0,137	-1,4	-15,7	7,6	-18,9	-7,9	8,4	-2,6	9,8	
17:35:44	026:36:04	0,197	10,0	-82,7	0,129	-1,4	-15,6	7,6	-18,8	-7,9	8,4	-2,5	9,8	
17:36:44	026:37:04	0,198	10,0	-88,9	0,130	-1,4	-15,5	7,6	-18,8	-7,9	8,4	-2,5	9,8	
17:37:42	026:38:02	0,199	10,0	-80,8	0,131	-1,4	-15,4	7,6	-18,8	-7,8	8,4	-2,5	9,8	
17:38:43	026:39:03	0,202	10,0	-87,7	0,133	-1,3	-15,4	7,6	-18,8	-7,8	8,4	-2,5	9,8	
17:39:40	026:40:00	0,204	10,0	-84,6	0,134	-1,3	-15,4	7,7	-18,8	-7,8	8,4	-2,5	9,9	
17:40:43	026:41:03	0,205	10,0	-89,3	0,136	-1,3	-15,3	7,7	-18,7	-7,7	8,4	-2,4	9,9	
17:41:40	026:42:00	0,208	10,0	-87,1	0,137	-1,3	-15,2	7,7	-18,7	-7,7	8,4	-2,4	9,9	
17:42:40	026:43:00	0,197	10,0	-79,5	0,129	-1,3	-15,1	7,7	-18,6	-7,6	8,4	-2,4	9,9	
17:43:42	026:44:02	0,197	10,0	-89,0	0,130	-1,2	-15,0	7,7	-18,6	-7,6	8,4	-2,3	9,9	
17:44:43	026:45:03	0,199	10,0	-80,5	0,131	-1,2	-15,0	7,7	-18,6	-7,5	8,4	-2,3	9,9	
17:45:44	026:46:04	0,202	10,0	-87,8	0,132	-1,2	-14,9	7,7	-18,6	-7,5	8,4	-2,3	9,9	
17:46:41	026:47:01	0,203	10,0	-85,5	0,133	-1,2	-14,8	7,7	-18,6	-7,5	8,5	-2,3	9,9	
17:47:43	026:48:03	0,205	10,0	-81,9	0,135	-1,2	-14,7	7,7	-18,5	-7,5	8,5	-2,3	9,9	
17:48:44	026:49:04	0,207	10,0	-88,8	0,137	-1,1	-14,6	7,7	-18,5	-7,4	8,5	-2,3	9,9	
17:49:40	026:50:00	0,197	10,0	-80,6	0,129	-1,1	-14,5	7,7	-18,5	-7,4	8,5	-2,2	9,9	
17:50:41	026:51:01	0,198	10,0	-88,1	0,130	-1,1	-14,4	7,7	-18,4	-7,3	8,5	-2,2	9,9	
17:51:43	026:52:03	0,199	10,0	-85,1	0,130	-1,0	-14,4	7,7	-18,4	-7,3	8,5	-2,2	10,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 87 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	17:52:40	026:53:00	0,200	10,0	-90,9	0,132	-1,0	-14,3	7,8	-18,4	-7,3	8,5	-2,2	9,9	
	17:53:41	026:54:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	17:53:41	026:54:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	17:53:41	026:54:01	0,202	10,0	-87,6	0,133	-1,0	-14,2	7,8	-18,4	-7,2	8,5	-2,1	10,0	
	17:54:41	026:55:01	0,205	10,0	-79,3	0,135	-1,0	-14,2	7,8	-18,4	-7,2	8,5	-2,1	10,0	
	17:55:44	026:56:04	0,207	10,0	-88,7	0,137	-1,0	-14,1	7,8	-18,3	-7,2	8,5	-2,1	10,0	
	17:56:44	026:57:04	0,197	10,0	-82,9	0,129	-1,0	-14,0	7,8	-18,3	-7,1	8,5	-2,1	10,0	
	17:57:40	026:58:00	0,197	10,0	-90,4	0,130	-0,9	-13,9	7,8	-18,2	-7,1	8,5	-2,1	10,0	
	17:58:44	026:59:04	0,198	10,0	-87,0	0,130	-0,9	-13,8	7,8	-18,2	-7,0	8,5	-2,0	10,0	
	17:59:42	027:00:02	0,201	10,0	-81,6	0,132	-0,9	-13,7	7,8	-18,2	-7,0	8,5	-2,0	10,0	
	18:00:41	027:01:01	0,202	10,0	-89,1	0,133	-0,9	-13,7	7,8	-18,2	-7,0	8,5	-2,0	10,0	
	18:01:44	027:02:04	0,204	10,0	-81,0	0,135	-0,9	-13,6	7,8	-18,1	-6,9	8,5	-2,0	10,0	
	18:02:41	027:03:01	0,206	10,0	-89,1	0,136	-0,9	-13,4	7,8	-18,1	-6,9	8,5	-2,0	10,0	
	18:03:42	027:04:02	0,208	10,0	-86,3	0,138	-0,8	-13,3	7,8	-18,0	-6,8	8,5	-2,0	10,0	
	18:04:43	027:05:03	0,197	10,0	-90,0	0,129	-0,8	-13,2	7,8	-18,0	-6,8	8,5	-1,9	10,1	
	18:05:40	027:06:00	0,197	10,0	-87,4	0,130	-0,8	-13,1	7,9	-18,0	-6,8	8,5	-1,9	10,1	
	18:06:41	027:07:01	0,199	10,0	-79,6	0,131	-0,8	-13,0	7,9	-17,9	-6,7	8,5	-1,9	10,1	
	18:07:44	027:08:04	0,201	10,0	-89,7	0,133	-0,7	-12,9	7,9	-17,9	-6,7	8,5	-1,9	10,1	
	18:08:40	027:09:00	0,203	10,0	-83,2	0,134	-0,7	-12,7	7,9	-17,9	-6,7	8,5	-1,9	10,1	
	18:09:42	027:10:02	0,204	10,0	-87,9	0,136	-0,7	-12,6	7,9	-17,9	-6,6	8,5	-1,8	10,1	
	18:10:40	027:11:00	0,207	10,0	-85,0	0,137	-0,7	-12,5	7,9	-17,8	-6,6	8,5	-1,8	10,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 88 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



18:11:40	027:12:00	0,197	10,0	-89,2	0,129	-0,7	-12,3	7,9	-17,7	-6,6	8,5	-1,8	10,1	
18:12:42	027:13:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
18:12:42	027:13:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
18:12:42	027:13:02	0,197	10,0	-88,7	0,130	-0,6	-12,1	7,9	-17,7	-6,5	8,5	-1,8	10,1	
18:13:41	027:14:01	0,198	10,0	-80,1	0,131	-0,6	-12,0	7,9	-17,7	-6,5	8,5	-1,8	10,1	
18:14:42	027:15:02	0,200	10,0	-88,8	0,132	-0,6	-11,9	7,9	-17,7	-6,5	8,5	-1,8	10,1	
18:15:43	027:16:03	0,202	10,0	-84,5	0,134	-0,6	-11,8	7,9	-17,6	-6,4	8,6	-1,7	10,1	
18:16:41	027:17:01	0,204	10,0	-88,6	0,135	-0,6	-11,6	7,9	-17,6	-6,4	8,6	-1,7	10,1	
18:17:42	027:18:02	0,205	10,0	-85,8	0,137	-0,6	-11,5	7,9	-17,5	-6,4	8,6	-1,7	10,1	
18:18:40	027:19:00	0,196	10,0	-89,2	0,129	-0,5	-11,4	7,9	-17,5	-6,3	8,6	-1,7	10,2	
18:19:42	027:20:02	0,196	10,0	-88,2	0,129	-0,5	-11,2	7,9	-17,4	-6,3	8,6	-1,6	10,2	
18:20:42	027:21:02	0,197	10,0	-79,1	0,130	-0,5	-11,1	8,0	-17,4	-6,2	8,6	-1,6	10,2	
18:21:42	027:22:02	0,199	10,0	-87,8	0,132	-0,4	-11,0	8,0	-17,3	-6,2	8,6	-1,6	10,2	
18:22:42	027:23:02	0,201	10,0	-84,7	0,133	-0,4	-10,9	8,0	-17,3	-6,2	8,6	-1,6	10,2	
18:23:40	027:24:00	0,203	10,0	-83,1	0,135	-0,4	-10,8	8,0	-17,2	-6,2	8,6	-1,6	10,2	
18:24:41	027:25:01	0,204	10,0	-88,7	0,136	-0,4	-10,7	8,0	-17,2	-6,1	8,6	-1,6	10,2	
18:25:44	027:26:04	0,207	10,0	-79,6	0,138	-0,4	-10,6	8,0	-17,1	-6,1	8,6	-1,5	10,2	
18:26:44	027:27:04	0,197	10,0	-87,8	0,129	-0,3	-10,4	8,0	-17,1	-6,0	8,6	-1,5	10,2	
18:27:43	027:28:03	0,196	10,0	-85,1	0,130	-0,3	-10,3	8,0	-17,0	-6,0	8,6	-1,5	10,2	
18:28:44	027:29:04	0,198	10,0	-88,3	0,131	-0,3	-10,3	8,0	-17,0	-6,0	8,6	-1,5	10,2	
18:29:40	027:30:00	0,199	10,0	-88,7	0,132	-0,3	-10,2	8,0	-17,0	-5,9	8,6	-1,5	10,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 89 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



18:30:42	027:31:02	0,201	10,0	-84,8	0,134	-0,3	-10,1	8,0	-16,9	-5,9	8,6	-1,5	10,3
18:31:40	027:32:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:31:40	027:32:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:31:40	027:32:00	0,203	10,0	-89,1	0,135	-0,3	-10,0	8,0	-16,8	-5,9	8,6	-1,4	10,3
18:32:43	027:33:03	0,205	10,0	-83,6	0,136	-0,2	-9,9	8,0	-16,8	-5,8	8,6	-1,4	10,3
18:33:41	027:34:01	0,196	10,0	-88,3	0,129	-0,2	-9,7	8,0	-16,7	-5,8	8,6	-1,4	10,3
18:34:41	027:35:01	0,196	10,0	-85,3	0,129	-0,2	-9,6	8,0	-16,6	-5,7	8,6	-1,4	10,3
18:35:41	027:36:01	0,197	10,0	-78,7	0,130	-0,1	-9,6	8,1	-16,6	-5,7	8,6	-1,4	10,3
18:36:41	027:37:01	0,198	10,0	-88,4	0,131	-0,1	-9,5	8,1	-16,5	-5,7	8,6	-1,3	10,3
18:37:41	027:38:01	0,201	10,0	-79,3	0,133	-0,1	-9,4	8,1	-16,5	-5,7	8,6	-1,3	10,3
18:38:43	027:39:03	0,202	10,0	-88,0	0,134	-0,1	-9,3	8,1	-16,4	-5,6	8,6	-1,3	10,3
18:39:40	027:40:00	0,204	10,0	-85,8	0,135	-0,1	-9,2	8,1	-16,4	-5,6	8,6	-1,3	10,3
18:40:42	027:41:02	0,206	10,0	-89,4	0,137	-0,1	-9,1	8,1	-16,3	-5,6	8,6	-1,3	10,3
18:41:40	027:42:00	0,195	10,0	-88,4	0,129	0,0	-9,0	8,1	-16,2	-5,5	8,6	-1,3	10,3
18:42:41	027:43:01	0,196	10,0	-79,7	0,129	0,0	-8,9	8,1	-16,1	-5,5	8,6	-1,2	10,3
18:43:45	027:44:05	0,197	10,0	-88,6	0,131	0,0	-8,8	8,1	-16,0	-5,4	8,6	-1,2	10,3
18:44:40	027:45:00	0,198	10,0	-84,0	0,132	0,0	-8,7	8,1	-16,0	-5,4	8,6	-1,2	10,3
18:45:42	027:46:02	0,202	10,0	-90,3	0,133	0,0	-8,6	8,1	-15,9	-5,4	8,7	-1,2	10,3
18:46:42	027:47:02	0,202	10,0	-86,6	0,135	0,0	-8,5	8,1	-15,9	-5,4	8,7	-1,2	10,3
18:47:40	027:48:00	0,204	10,0	-78,7	0,136	0,0	-8,4	8,1	-15,8	-5,3	8,6	-1,2	10,4
18:48:42	027:49:02	0,206	10,0	-88,0	0,138	0,1	-8,3	8,1	-15,7	-5,3	8,7	-1,1	10,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 90 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



18:49:43	027:50:03	0,195	10,0	-83,8	0,129	0,1	-8,2	8,1	-15,6	-5,2	8,7	-1,1	10,4
18:50:44	027:51:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:50:44	027:51:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:50:44	027:51:04	0,196	10,0	-90,2	0,130	0,2	-8,1	8,1	-15,6	-5,2	8,7	-1,1	10,4
18:51:44	027:52:04	0,197	10,0	-87,4	0,131	0,2	-8,0	8,1	-15,5	-5,2	8,7	-1,1	10,4
18:52:40	027:53:00	0,199	10,0	-85,5	0,132	0,2	-8,0	8,1	-15,5	-5,1	8,7	-1,1	10,4
18:53:42	027:54:02	0,201	10,0	-89,5	0,134	0,2	-7,8	8,2	-15,4	-5,1	8,7	-1,1	10,4
18:54:42	027:55:02	0,203	10,0	-83,7	0,135	0,2	-7,7	8,2	-15,3	-5,1	8,7	-1,0	10,4
18:55:43	027:56:03	0,204	10,0	-89,4	0,137	0,2	-7,7	8,2	-15,2	-5,1	8,7	-1,0	10,4
18:56:41	027:57:01	0,206	10,0	-83,1	0,138	0,2	-7,5	8,2	-15,2	-5,0	8,7	-1,0	10,4
18:57:41	027:58:01	0,195	10,0	-89,7	0,129	0,3	-7,4	8,2	-15,1	-5,0	8,7	-1,0	10,4
18:58:44	027:59:04	0,196	10,0	-86,3	0,130	0,3	-7,3	8,2	-15,0	-5,0	8,7	-1,0	10,4
18:59:40	028:00:00	0,197	10,0	-83,1	0,131	0,3	-7,2	8,2	-14,9	-4,9	8,7	-0,9	10,4
19:00:40	028:01:00	0,199	10,0	-89,0	0,132	0,3	-7,1	8,2	-14,8	-4,9	8,7	-0,9	10,4
19:01:41	028:02:01	0,200	10,0	-82,6	0,134	0,3	-7,0	8,2	-14,7	-4,9	8,7	-0,9	10,4
19:02:41	028:03:01	0,202	10,0	-90,5	0,135	0,4	-6,9	8,2	-14,6	-4,8	8,7	-0,9	10,4
19:03:44	028:04:04	0,204	10,0	-86,5	0,137	0,4	-6,8	8,2	-14,5	-4,8	8,7	-0,9	10,4
19:04:40	028:05:00	0,195	10,0	-79,1	0,129	0,4	-6,7	8,2	-14,5	-4,8	8,7	-0,9	10,5
19:05:41	028:06:01	0,194	10,0	-88,1	0,129	0,4	-6,6	8,2	-14,3	-4,7	8,7	-0,8	10,5
19:06:44	028:07:04	0,196	10,0	-82,1	0,130	0,4	-6,5	8,2	-14,2	-4,7	8,7	-0,8	10,5
19:07:40	028:08:00	0,197	10,0	-90,6	0,131	0,5	-6,4	8,2	-14,2	-4,7	8,7	-0,8	10,5

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 91 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



19:08:41	028:09:01	0,198	10,0	-86,4	0,132	0,5	-6,3	8,2	-14,0	-4,6	8,7	-0,8	10,5
19:09:44	028:10:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
19:09:44	028:10:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
19:09:44	028:10:04	0,200	10,0	-79,9	0,134	0,5	-6,2	8,3	-13,9	-4,6	8,7	-0,8	10,5
19:10:41	028:11:01	0,202	10,0	-88,6	0,135	0,5	-6,0	8,3	-13,8	-4,6	8,7	-0,8	10,5
19:11:42	028:12:02	0,204	10,0	-82,4	0,136	0,5	-5,9	8,3	-13,7	-4,5	8,7	-0,7	10,5
19:12:44	028:13:04	0,195	10,0	-90,0	0,129	0,5	-5,8	8,3	-13,5	-4,5	8,7	-0,7	10,5
19:13:41	028:14:01	0,193	10,0	-86,4	0,129	0,6	-5,7	8,3	-13,4	-4,5	8,7	-0,7	10,5
19:14:42	028:15:02	0,194	10,0	-78,1	0,129	0,6	-5,6	8,3	-13,3	-4,4	8,7	-0,7	10,5
19:15:42	028:16:02	0,197	10,0	-87,7	0,131	0,6	-5,5	8,3	-13,1	-4,4	8,7	-0,7	10,5
19:16:42	028:17:02	0,197	10,0	-83,0	0,132	0,6	-5,4	8,3	-13,0	-4,4	8,7	-0,6	10,5
19:17:43	028:18:03	0,200	10,0	-90,9	0,133	0,6	-5,2	8,3	-12,9	-4,3	8,7	-0,6	10,5
19:18:43	028:19:03	0,201	10,0	-87,8	0,135	0,6	-5,1	8,3	-12,8	-4,3	8,7	-0,6	10,5
19:19:44	028:20:04	0,203	10,0	-78,9	0,136	0,7	-4,9	8,3	-12,6	-4,3	8,7	-0,6	10,5
19:20:40	028:21:00	0,205	10,0	-87,8	0,138	0,7	-4,8	8,3	-12,5	-4,2	8,7	-0,6	10,5
19:21:43	028:22:03	0,194	10,0	-84,1	0,129	0,7	-4,6	8,3	-12,3	-4,2	8,7	-0,6	10,5
19:22:45	028:23:05	0,194	10,0	-87,8	0,130	0,7	-4,4	8,3	-12,2	-4,2	8,7	-0,5	10,6
19:23:42	028:24:02	0,195	10,0	-88,2	0,130	0,8	-4,3	8,3	-12,1	-4,1	8,7	-0,5	10,6
19:24:43	028:25:03	0,198	10,0	-78,8	0,132	0,8	-4,1	8,4	-11,9	-4,1	8,8	-0,5	10,6
19:25:42	028:26:02	0,200	10,0	-88,3	0,133	0,8	-4,0	8,4	-11,8	-4,1	8,8	-0,5	10,6
19:26:43	028:27:03	0,200	10,0	-85,2	0,134	0,8	-3,8	8,4	-11,7	-4,0	8,8	-0,5	10,6

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 92 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



19:27:41	028:28:01	0,202	10,0	-92,6	0,136	0,8	-3,6	8,4	-11,5	-4,0	8,8	-0,5	10,6
19:28:43	028:29:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
19:28:43	028:29:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
19:28:43	028:29:03	0,204	10,0	-88,7	0,137	0,8	-3,4	8,4	-11,4	-4,0	8,8	-0,4	10,6
19:29:44	028:30:04	0,194	10,0	-80,2	0,129	0,9	-3,2	8,4	-11,3	-3,9	8,8	-0,4	10,6
19:30:40	028:31:00	0,194	10,0	-88,3	0,129	0,9	-3,1	8,4	-11,2	-3,9	8,8	-0,4	10,6
19:31:40	028:32:00	0,195	10,0	-83,7	0,130	0,9	-2,9	8,4	-11,1	-3,8	8,8	-0,4	10,6
19:32:40	028:33:00	0,196	10,0	-90,2	0,131	0,9	-2,6	8,4	-11,0	-3,8	8,8	-0,4	10,6
19:33:42	028:34:02	0,197	10,0	-85,4	0,132	0,9	-2,3	8,4	-10,9	-3,8	8,8	-0,4	10,6
19:34:42	028:35:02	0,199	10,0	-78,4	0,134	0,9	-2,0	8,4	-10,8	-3,8	8,8	-0,4	10,6
19:35:40	028:36:00	0,201	10,0	-88,8	0,135	0,9	-1,7	8,4	-10,7	-3,7	8,8	-0,3	10,6
19:36:41	028:37:01	0,203	10,0	-82,1	0,137	1,0	-1,4	8,4	-10,5	-3,7	8,8	-0,3	10,6
19:37:43	028:38:03	0,194	10,0	-89,9	0,129	1,0	-1,1	8,4	-10,4	-3,6	8,8	-0,3	10,6
19:38:43	028:39:03	0,194	10,0	-85,8	0,129	1,0	-0,9	8,4	-10,3	-3,6	8,8	-0,3	10,7
19:39:43	028:40:03	0,194	10,0	-81,0	0,130	1,0	-0,6	8,4	-10,2	-3,6	8,8	-0,3	10,7
19:40:43	028:41:03	0,195	10,0	-89,3	0,131	1,0	-0,4	8,4	-10,1	-3,5	8,8	-0,3	10,7
19:41:41	028:42:01	0,196	10,0	-81,8	0,132	1,0	-0,3	8,4	-10,0	-3,5	8,8	-0,3	10,7
19:42:42	028:43:02	0,198	10,0	-88,0	0,133	1,1	-0,1	8,4	-9,8	-3,5	8,8	-0,2	10,7
19:43:43	028:44:03	0,200	10,0	-84,5	0,135	1,1	0,0	8,4	-9,7	-3,5	8,8	-0,2	10,7
19:44:43	028:45:03	0,202	10,0	-90,3	0,136	1,1	0,2	8,4	-9,6	-3,4	8,8	-0,2	10,7
19:45:44	028:46:04	0,203	10,0	-86,3	0,137	1,1	0,3	8,5	-9,5	-3,4	8,8	-0,2	10,7

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 93 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	19:46:41	028:47:01	0,193	10,0	-78,9	0,129	1,1	0,4	8,5	-9,4	-3,4	8,8	-0,2	10,7	
	19:47:41	028:48:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	19:47:41	028:48:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	19:47:41	028:48:01	0,193	10,0	-87,9	0,129	1,1	0,6	8,5	-9,3	-3,3	8,8	-0,1	10,7	
	19:48:42	028:49:02	0,194	10,0	-83,3	0,130	1,2	0,7	8,5	-9,2	-3,3	8,8	-0,1	10,7	
	19:49:44	028:50:04	0,195	10,0	-90,4	0,131	1,2	0,8	8,5	-9,1	-3,3	8,8	-0,1	10,7	
	19:50:40	028:51:00	0,198	10,0	-87,9	0,132	1,2	0,9	8,5	-9,0	-3,2	8,8	-0,1	10,7	
	19:51:41	028:52:01	0,199	10,0	-79,1	0,134	1,2	1,1	8,5	-8,9	-3,2	8,8	-0,1	10,7	
	19:52:42	028:53:02	0,200	10,0	-87,7	0,135	1,2	1,3	8,5	-8,7	-3,2	8,8	-0,1	10,7	
	19:53:42	028:54:02	0,202	10,0	-84,0	0,136	1,2	1,4	8,5	-8,6	-3,1	8,8	-0,1	10,7	
	19:54:44	028:55:04	0,193	10,0	-90,6	0,129	1,3	1,6	8,5	-8,5	-3,1	8,8	-0,1	10,7	
	19:55:41	028:56:01	0,192	10,0	-87,6	0,129	1,3	1,8	8,5	-8,4	-3,1	8,8	-0,1	10,7	
	19:56:42	028:57:02	0,192	10,0	-79,7	0,129	1,3	2,0	8,5	-8,2	-3,0	8,8	-0,1	10,7	
	19:57:43	028:58:03	0,193	10,0	-89,2	0,130	1,3	2,1	8,5	-8,1	-3,0	8,8	-0,1	10,7	
	19:58:40	028:59:00	0,195	10,0	-82,7	0,131	1,3	2,2	8,5	-8,1	-3,0	8,8	-0,1	10,7	
	19:59:41	029:00:01	0,197	10,0	-90,3	0,133	1,3	2,4	8,5	-7,9	-2,9	8,9	-0,1	10,7	
	20:00:42	029:01:02	0,198	10,0	-86,4	0,134	1,3	2,5	8,5	-7,8	-2,9	8,9	-0,1	10,7	
	20:01:44	029:02:04	0,200	10,0	-79,3	0,135	1,3	2,6	8,5	-7,7	-2,9	8,9	-0,1	10,8	
	20:02:45	029:03:05	0,202	10,0	-89,2	0,137	1,4	2,7	8,5	-7,6	-2,9	8,8	-0,1	10,8	
	20:03:41	029:04:01	0,197	10,0	-82,2	0,132	1,4	2,8	8,6	-7,5	-2,8	8,9	0,0	10,8	
	20:04:42	029:05:02	0,192	10,0	-89,5	0,129	1,4	2,9	8,6	-7,4	-2,8	8,9	0,0	10,8	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	94 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:05:44	029:06:04	0,192	10,0	-86,2	0,130	1,5	3,1	8,6	-7,3	-2,8	8,9	0,0	10,8	
	20:06:40	029:07:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:06:40	029:07:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:06:40	029:07:00	0,193	10,0	-80,7	0,131	1,5	3,1	8,6	-7,2	-2,7	8,9	0,0	10,8	
	20:07:42	029:08:02	0,195	10,0	-89,0	0,132	1,5	3,2	8,6	-7,1	-2,7	8,9	0,0	10,8	
	20:08:42	029:09:02	0,196	10,0	-81,1	0,133	1,5	3,3	8,6	-7,0	-2,7	8,8	0,0	10,8	
	20:09:42	029:10:02	0,198	10,0	-89,4	0,135	1,5	3,4	8,6	-6,9	-2,6	8,9	0,1	10,8	
	20:10:43	029:11:03	0,200	10,0	-85,8	0,136	1,5	3,5	8,6	-6,9	-2,6	8,9	0,1	10,8	
	20:11:44	029:12:04	0,202	10,0	-90,9	0,137	1,5	3,6	8,6	-6,8	-2,6	8,9	0,1	10,8	
	20:12:41	029:13:01	0,192	10,0	-89,2	0,129	1,6	3,6	8,6	-6,7	-2,5	8,9	0,1	10,8	
	20:13:42	029:14:02	0,191	10,0	-81,4	0,129	1,6	3,7	8,6	-6,6	-2,5	8,9	0,1	10,8	
	20:14:43	029:15:03	0,192	10,0	-88,7	0,130	1,6	3,8	8,6	-6,5	-2,5	8,9	0,1	10,8	
	20:15:40	029:16:00	0,193	10,0	-86,2	0,131	1,6	3,8	8,6	-6,5	-2,5	8,9	0,1	10,8	
	20:16:42	029:17:02	0,195	10,0	-78,3	0,132	1,7	3,9	8,6	-6,4	-2,4	8,9	0,1	10,8	
	20:17:43	029:18:03	0,196	10,0	-88,6	0,133	1,7	4,0	8,6	-6,3	-2,4	8,9	0,2	10,8	
	20:18:43	029:19:03	0,198	10,0	-79,8	0,134	1,7	4,0	8,6	-6,2	-2,4	8,9	0,2	10,8	
	20:19:44	029:20:04	0,199	10,0	-87,9	0,136	1,7	4,1	8,6	-6,1	-2,3	8,9	0,2	10,8	
	20:20:40	029:21:00	0,202	10,0	-85,4	0,137	1,7	4,2	8,6	-6,1	-2,3	8,9	0,2	10,8	
	20:21:42	029:22:02	0,192	10,0	-90,9	0,129	1,7	4,3	8,7	-5,9	-2,3	8,9	0,2	10,8	
	20:22:43	029:23:03	0,190	10,0	-87,4	0,129	1,8	4,4	8,7	-5,9	-2,2	8,9	0,2	10,9	
	20:23:44	029:24:04	0,191	10,0	-80,2	0,130	1,8	4,5	8,7	-5,8	-2,2	8,9	0,2	10,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 95 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:24:40	029:25:00	0,192	10,0	-89,5	0,131	1,8	4,5	8,7	-5,7	-2,2	8,9	0,2	10,9	
	20:25:41	029:26:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:25:41	029:26:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:25:41	029:26:01	0,193	10,0	-82,9	0,132	1,8	4,6	8,7	-5,6	-2,2	8,9	0,2	10,9	
	20:26:41	029:27:01	0,195	10,0	-88,7	0,133	1,8	4,7	8,7	-5,6	-2,1	8,9	0,2	10,9	
	20:27:41	029:28:01	0,197	10,0	-85,1	0,134	1,8	4,8	8,7	-5,5	-2,1	8,9	0,2	10,9	
	20:28:42	029:29:02	0,200	10,0	-76,6	0,135	1,8	4,8	8,7	-5,4	-2,1	8,9	0,2	10,9	
	20:29:42	029:30:02	0,200	10,0	-88,8	0,137	1,9	4,9	8,7	-5,3	-2,0	8,9	0,2	10,9	
	20:30:43	029:31:03	0,192	10,0	-80,8	0,130	1,9	5,0	8,7	-5,2	-2,0	8,9	0,3	10,9	
	20:31:43	029:32:03	0,190	10,0	-89,3	0,129	1,9	5,0	8,7	-5,1	-2,0	8,9	0,3	10,9	
	20:32:44	029:33:04	0,190	10,0	-86,1	0,130	1,9	5,1	8,7	-5,0	-1,9	8,9	0,4	10,9	
	20:33:40	029:34:00	0,192	10,0	-79,7	0,131	2,0	5,1	8,7	-5,0	-1,9	8,9	0,4	10,9	
	20:34:40	029:35:00	0,193	10,0	-88,9	0,132	2,0	5,2	8,7	-4,9	-1,9	8,9	0,4	10,9	
	20:35:42	029:36:02	0,194	10,0	-79,8	0,133	2,0	5,2	8,7	-4,8	-1,9	8,9	0,5	10,9	
	20:36:40	029:37:00	0,196	10,0	-88,6	0,134	2,0	5,2	8,7	-4,8	-1,9	8,9	0,5	10,9	
	20:37:42	029:38:02	0,198	10,0	-84,9	0,135	2,0	5,3	8,7	-4,7	-1,8	8,9	0,5	10,9	
	20:38:42	029:39:02	0,199	10,0	-90,5	0,137	2,0	5,3	8,8	-4,6	-1,8	8,9	0,5	10,9	
	20:39:43	029:40:03	0,201	10,0	-86,5	0,138	2,0	5,3	8,8	-4,5	-1,8	8,9	0,6	10,9	
	20:40:44	029:41:04	0,190	10,0	-79,0	0,129	2,1	5,4	8,8	-4,4	-1,7	8,9	0,6	10,9	
	20:41:44	029:42:04	0,190	10,0	-89,2	0,130	2,1	5,4	8,8	-4,4	-1,7	8,9	0,7	10,9	
	20:42:40	029:43:00	0,191	10,0	-82,4	0,130	2,1	5,5	8,8	-4,3	-1,7	8,9	0,7	10,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 96 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:43:43	029:44:03	0,192	10,0	-90,1	0,132	2,1	5,5	8,8	-4,2	-1,6	8,9	0,8	10,9	
	20:44:40	029:45:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:44:40	029:45:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:44:40	029:45:00	0,194	10,0	-88,0	0,133	2,1	5,5	8,8	-4,2	-1,6	8,9	0,8	10,9	
	20:45:40	029:46:00	0,196	10,0	-78,6	0,134	2,1	5,5	8,8	-4,1	-1,6	8,9	0,8	10,9	
	20:46:40	029:47:00	0,197	10,0	-88,2	0,135	2,2	5,6	8,8	-4,0	-1,6	8,9	0,8	10,9	
	20:47:42	029:48:02	0,198	10,0	-84,4	0,136	2,2	5,6	8,8	-4,0	-1,6	9,0	0,9	11,0	
	20:48:44	029:49:04	0,200	10,0	-80,9	0,138	2,2	5,6	8,8	-3,9	-1,5	9,0	0,9	11,0	
	20:49:43	029:50:03	0,189	10,0	-89,2	0,129	2,2	5,6	8,8	-3,8	-1,5	9,0	0,9	11,0	
	20:50:40	029:51:00	0,189	10,0	-81,8	0,129	2,3	5,7	8,8	-3,7	-1,5	9,0	0,9	11,0	
	20:51:42	029:52:02	0,190	10,0	-87,9	0,130	2,3	5,7	8,8	-3,6	-1,4	9,0	1,0	11,0	
	20:52:44	029:53:04	0,191	10,0	-83,8	0,131	2,3	5,7	8,8	-3,6	-1,4	9,0	1,0	11,0	
	20:53:41	029:54:01	0,192	10,0	-82,3	0,132	2,3	5,7	8,8	-3,5	-1,4	9,0	1,0	11,0	
	20:54:41	029:55:01	0,194	10,0	-88,9	0,133	2,3	5,7	8,8	-3,4	-1,4	9,0	1,0	11,0	
	20:55:42	029:56:02	0,195	10,0	-79,9	0,135	2,3	5,8	8,8	-3,4	-1,3	9,0	1,0	11,0	
	20:56:40	029:57:00	0,197	10,0	-88,8	0,136	2,3	5,8	8,8	-3,3	-1,3	9,0	1,0	11,0	
	20:57:40	029:58:00	0,199	10,0	-85,5	0,137	2,4	5,8	8,8	-3,2	-1,3	9,0	1,1	11,0	
	20:58:41	029:59:01	0,190	10,0	-89,2	0,130	2,4	5,8	8,8	-3,1	-1,2	9,0	1,1	11,0	
	20:59:44	030:00:04	0,188	10,0	-89,0	0,129	2,4	5,9	8,8	-3,0	-1,2	9,0	1,1	11,0	
	21:00:44	030:01:04	0,189	10,0	-80,4	0,130	2,4	5,9	8,8	-3,0	-1,2	9,0	1,1	11,0	
	21:01:40	030:02:00	0,191	10,0	-89,5	0,131	2,4	5,9	8,8	-2,9	-1,1	9,0	1,1	11,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 97 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:02:40	030:03:00	0,192	10,0	-87,0	0,131	2,4	5,9	8,9	-2,9	-1,1	9,0	1,1	11,1	
	21:03:43	030:04:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:03:43	030:04:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:03:43	030:04:03	0,193	10,0	-90,3	0,133	2,5	5,9	8,9	-2,8	-1,1	9,0	1,1	11,0	
	21:04:43	030:05:03	0,194	10,0	-86,1	0,134	2,5	5,9	8,9	-2,7	-1,1	9,0	1,2	11,1	
	21:05:41	030:06:01	0,195	10,0	-78,2	0,135	2,5	6,0	8,9	-2,6	-1,1	9,0	1,2	11,1	
	21:06:42	030:07:02	0,197	10,0	-87,8	0,136	2,5	6,0	8,9	-2,6	-1,0	9,0	1,2	11,1	
	21:07:43	030:08:03	0,198	10,0	-83,3	0,137	2,5	6,0	8,9	-2,5	-1,0	9,0	1,2	11,1	
	21:08:43	030:09:03	0,188	10,0	-87,4	0,129	2,6	6,0	8,9	-2,4	-1,0	9,0	1,3	11,1	
	21:09:44	030:10:04	0,187	10,0	-88,1	0,129	2,6	6,0	8,9	-2,3	-1,0	9,0	1,3	11,1	
	21:10:44	030:11:04	0,189	10,0	-80,7	0,130	2,6	6,1	8,9	-2,3	-0,9	9,0	1,3	11,1	
	21:11:41	030:12:01	0,189	10,0	-89,2	0,131	2,6	6,1	8,9	-2,2	-0,9	9,0	1,3	11,1	
	21:12:44	030:13:04	0,190	10,0	-81,7	0,132	2,6	6,1	8,9	-2,1	-0,9	9,0	1,3	11,1	
	21:13:44	030:14:04	0,192	10,0	-88,8	0,133	2,6	6,1	8,9	-2,1	-0,8	9,0	1,3	11,1	
	21:14:41	030:15:01	0,193	10,0	-86,4	0,134	2,6	6,1	8,9	-2,0	-0,8	9,0	1,3	11,1	
	21:15:41	030:16:01	0,195	10,0	-90,8	0,135	2,6	6,1	8,9	-1,9	-0,8	9,0	1,4	11,1	
	21:16:42	030:17:02	0,196	10,0	-88,8	0,136	2,7	6,2	8,9	-1,8	-0,8	9,0	1,4	11,1	
	21:17:43	030:18:03	0,198	10,0	-82,7	0,137	2,7	6,2	8,9	-1,8	-0,7	9,0	1,4	11,1	
	21:18:43	030:19:03	0,187	10,0	-89,3	0,129	2,7	6,2	8,9	-1,7	-0,7	9,0	1,4	11,1	
	21:19:43	030:20:03	0,187	10,0	-81,4	0,129	2,7	6,2	9,0	-1,6	-0,7	9,1	1,4	11,1	
	21:20:44	030:21:04	0,188	10,0	-89,8	0,130	2,8	6,2	9,0	-1,5	-0,6	9,0	1,5	11,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 98 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:21:44	030:22:04	0,188	10,0	-85,8	0,131	2,8	6,2	9,0	-1,5	-0,6	9,0	1,5	11,1	
	21:22:44	030:23:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:22:44	030:23:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:22:44	030:23:04	0,190	10,0	-78,6	0,132	2,8	6,2	9,0	-1,4	-0,6	9,0	1,5	11,1	
	21:23:40	030:24:00	0,191	10,0	-88,5	0,133	2,8	6,3	9,0	-1,3	-0,6	9,0	1,5	11,1	
	21:24:43	030:25:03	0,193	10,0	-82,5	0,134	2,8	6,3	9,0	-1,2	-0,6	9,1	1,5	11,1	
	21:25:44	030:26:04	0,194	10,0	-90,1	0,135	2,8	6,3	9,0	-1,2	-0,6	9,0	1,5	11,1	
	21:26:42	030:27:02	0,196	10,0	-86,2	0,136	2,8	6,3	9,0	-1,1	-0,5	9,1	1,5	11,2	
	21:27:44	030:28:04	0,197	10,0	-78,5	0,138	2,9	6,3	9,0	-1,0	-0,5	9,1	1,6	11,1	
	21:28:40	030:29:00	0,186	10,0	-88,3	0,129	2,9	6,3	9,0	-0,9	-0,5	9,1	1,6	11,1	
	21:29:43	030:30:03	0,186	10,0	-82,4	0,129	2,9	6,4	9,0	-0,8	-0,5	9,1	1,6	11,1	
	21:30:41	030:31:01	0,187	10,0	-90,4	0,130	2,9	6,4	9,0	-0,8	-0,5	9,1	1,6	11,1	
	21:31:41	030:32:01	0,187	10,0	-86,7	0,131	2,9	6,4	9,0	-0,7	-0,4	9,1	1,6	11,2	
	21:32:44	030:33:04	0,190	10,0	-78,9	0,132	2,9	6,4	9,0	-0,6	-0,4	9,1	1,6	11,2	
	21:33:45	030:34:05	0,190	10,0	-89,1	0,133	2,9	6,4	9,0	-0,6	-0,4	9,1	1,6	11,2	
	21:34:41	030:35:01	0,191	10,0	-82,4	0,134	2,9	6,4	9,0	-0,5	-0,3	9,1	1,6	11,2	
	21:35:43	030:36:03	0,193	10,0	-90,4	0,135	2,9	6,4	9,0	-0,4	-0,3	9,1	1,6	11,2	
	21:36:45	030:37:05	0,194	10,0	-86,0	0,136	3,0	6,5	9,0	-0,4	-0,3	9,1	1,7	11,2	
	21:37:41	030:38:01	0,196	10,0	-78,7	0,137	3,0	6,5	9,0	-0,3	-0,3	9,1	1,7	11,2	
	21:38:42	030:39:02	0,186	10,0	-88,0	0,129	3,0	6,5	9,0	-0,2	-0,2	9,1	1,7	11,2	
	21:39:43	030:40:03	0,185	10,0	-83,0	0,129	3,0	6,5	9,0	-0,1	-0,2	9,1	1,7	11,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 99 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	21:40:41	030:41:01	0,185	10,0	-91,1	0,130	3,0	6,5	9,0	0,0	-0,2	9,1	1,7	11,2	
	21:41:43	030:42:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:41:43	030:42:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:41:43	030:42:03	0,187	10,0	-86,4	0,130	3,0	6,5	9,0	0,1	-0,2	9,1	1,7	11,2	
	21:42:44	030:43:04	0,187	10,0	-78,5	0,131	3,1	6,5	9,0	0,1	-0,2	9,1	1,7	11,2	
	21:43:41	030:44:01	0,189	10,0	-87,5	0,132	3,1	6,5	9,1	0,2	-0,2	9,1	1,7	11,2	
	21:44:41	030:45:01	0,190	10,0	-83,8	0,133	3,1	6,6	9,1	0,3	-0,2	9,1	1,8	11,2	
	21:45:44	030:46:04	0,191	10,0	-83,8	0,135	3,1	6,6	9,1	0,4	-0,2	9,1	1,8	11,2	
	21:46:41	030:47:01	0,193	10,0	-88,8	0,136	3,1	6,6	9,1	0,5	-0,1	9,1	1,8	11,2	
	21:47:42	030:48:02	0,195	10,0	-79,4	0,137	3,1	6,6	9,1	0,6	-0,1	9,1	1,8	11,2	
	21:48:43	030:49:03	0,192	10,0	-88,6	0,134	3,1	6,6	9,1	0,7	-0,1	9,1	1,8	11,2	
	21:49:41	030:50:01	0,184	10,0	-83,9	0,129	3,2	6,6	9,1	0,8	-0,1	9,1	1,9	11,2	
	21:50:44	030:51:04	0,184	10,0	-90,6	0,130	3,2	6,6	9,1	0,9	-0,1	9,1	1,9	11,2	
	21:51:44	030:52:04	0,185	10,0	-86,7	0,130	3,2	6,6	9,1	1,0	-0,1	9,1	1,9	11,2	
	21:52:41	030:53:01	0,186	10,0	-79,3	0,131	3,2	6,7	9,1	1,0	-0,1	9,1	1,9	11,2	
	21:53:41	030:54:01	0,187	10,0	-89,3	0,132	3,2	6,7	9,1	1,2	0,0	9,1	1,9	11,2	
	21:54:42	030:55:02	0,188	10,0	-82,1	0,133	3,2	6,7	9,1	1,3	0,0	9,1	1,9	11,2	
	21:55:43	030:56:03	0,191	10,0	-88,8	0,134	3,3	6,7	9,1	1,4	0,0	9,1	1,9	11,2	
	21:56:43	030:57:03	0,191	10,0	-86,1	0,135	3,3	6,7	9,1	1,5	-0,1	9,1	1,9	11,2	
	21:57:43	030:58:03	0,193	10,0	-89,4	0,136	3,3	6,7	9,1	1,6	0,0	9,1	2,0	11,3	
	21:58:43	030:59:03	0,194	10,0	-88,3	0,137	3,3	6,7	9,1	1,7	0,0	9,1	2,0	11,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 100 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:59:45	031:00:05	0,184	10,0	-79,2	0,130	3,3	6,8	9,1	1,8	0,1	9,1	2,0	11,3	
	22:00:40	031:01:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:00:40	031:01:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:00:40	031:01:00	0,184	10,0	-88,0	0,129	3,3	6,8	9,1	1,8	0,1	9,1	2,0	11,3	
	22:01:40	031:02:00	0,183	10,0	-84,4	0,130	3,4	6,8	9,1	1,9	0,0	9,1	2,0	11,3	
	22:02:41	031:03:01	0,184	10,0	-76,8	0,131	3,4	6,8	9,1	2,0	0,0	9,1	2,1	11,3	
	22:03:42	031:04:02	0,185	10,0	-88,8	0,131	3,4	6,8	9,1	2,1	0,1	9,1	2,1	11,3	
	22:04:43	031:05:03	0,186	10,0	-79,9	0,132	3,4	6,8	9,1	2,2	0,1	9,1	2,1	11,3	
	22:05:41	031:06:01	0,188	10,0	-89,0	0,134	3,4	6,8	9,1	2,3	0,2	9,1	2,1	11,3	
	22:06:42	031:07:02	0,189	10,0	-85,4	0,134	3,4	6,8	9,2	2,4	0,4	9,1	2,1	11,3	
	22:07:40	031:08:00	0,190	10,0	-88,6	0,135	3,4	6,9	9,2	2,5	0,5	9,1	2,1	11,3	
	22:08:44	031:09:04	0,191	10,0	-89,2	0,137	3,5	6,9	9,2	2,6	0,5	9,1	2,1	11,3	
	22:09:40	031:10:00	0,193	10,0	-81,4	0,137	3,5	6,9	9,2	2,7	0,5	9,1	2,1	11,3	
	22:10:41	031:11:01	0,183	10,0	-88,9	0,129	3,5	6,9	9,2	2,8	0,5	9,1	2,1	11,3	
	22:11:44	031:12:04	0,182	10,0	-84,9	0,129	3,5	6,9	9,2	3,0	0,6	9,1	2,2	11,3	
	22:12:41	031:13:01	0,183	10,0	-83,1	0,130	3,5	6,9	9,2	3,1	0,8	9,1	2,2	11,3	
	22:13:45	031:14:05	0,183	10,0	-88,8	0,131	3,5	6,9	9,2	3,2	1,1	9,1	2,2	11,3	
	22:14:41	031:15:01	0,184	10,0	-82,1	0,131	3,6	6,9	9,2	3,3	1,2	9,1	2,2	11,3	
	22:15:43	031:16:03	0,185	10,0	-90,9	0,133	3,6	7,0	9,2	3,4	1,3	9,1	2,2	11,3	
	22:16:41	031:17:01	0,186	10,0	-87,7	0,133	3,6	7,0	9,2	3,5	1,4	9,1	2,2	11,3	
	22:17:44	031:18:04	0,188	10,0	-78,6	0,134	3,6	7,0	9,2	3,7	1,5	9,1	2,3	11,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 101 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	22:18:40	031:19:00	0,189	10,0	-87,8	0,135	3,6	7,0	9,2	3,8	1,6	9,1	2,3	11,4	
	22:19:40	031:20:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:19:40	031:20:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:19:40	031:20:00	0,190	10,0	-84,0	0,136	3,6	7,0	9,2	3,9	1,7	9,1	2,3	11,4	
	22:20:44	031:21:04	0,192	10,0	-89,0	0,137	3,7	7,0	9,2	4,1	1,8	9,1	2,3	11,4	
	22:21:44	031:22:04	0,181	10,0	-87,8	0,129	3,7	7,0	9,2	4,3	1,8	9,1	2,3	11,4	
	22:22:43	031:23:03	0,180	10,0	-78,5	0,129	3,7	7,0	9,2	4,5	1,8	9,1	2,4	11,4	
	22:23:45	031:24:05	0,181	10,0	-88,0	0,130	3,7	7,1	9,2	4,8	1,9	9,1	2,4	11,4	
	22:24:40	031:25:00	0,181	10,0	-85,8	0,130	3,7	7,1	9,2	5,1	1,9	9,1	2,4	11,4	
	22:25:40	031:26:00	0,182	10,0	-88,2	0,131	3,7	7,1	9,3	5,4	1,9	9,1	2,4	11,4	
	22:26:40	031:27:00	0,183	10,0	-88,8	0,132	3,7	7,1	9,3	5,6	2,0	9,1	2,4	11,4	
	22:27:41	031:28:01	0,184	10,0	-79,8	0,133	3,8	7,1	9,3	6,0	2,0	9,1	2,4	11,4	
	22:28:42	031:29:02	0,186	10,0	-89,8	0,134	3,8	7,1	9,3	6,2	2,1	9,1	2,4	11,4	
	22:29:40	031:30:00	0,187	10,0	-85,8	0,135	3,8	7,1	9,3	6,4	2,1	9,1	2,4	11,4	
	22:30:42	031:31:02	0,188	10,0	-77,8	0,136	3,8	7,1	9,3	6,6	2,1	9,1	2,5	11,4	
	22:31:42	031:32:02	0,189	10,0	-88,8	0,137	3,8	7,1	9,3	6,7	2,1	9,1	2,5	11,4	
	22:32:40	031:33:00	0,190	10,0	-83,3	0,138	3,8	7,1	9,3	6,9	2,1	9,1	2,5	11,4	
	22:33:41	031:34:01	0,181	10,0	-90,4	0,130	3,8	7,1	9,3	7,0	2,2	9,1	2,5	11,4	
	22:34:43	031:35:03	0,179	10,0	-87,5	0,130	3,8	7,2	9,3	7,2	2,2	9,1	2,5	11,4	
	22:35:43	031:36:03	0,180	10,0	-81,2	0,130	3,9	7,2	9,3	7,3	2,2	9,1	2,6	11,4	
	22:36:45	031:37:05	0,181	10,0	-89,3	0,131	3,9	7,2	9,3	7,4	2,2	9,1	2,6	11,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	102 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



22:37:42	031:38:02	0,182	10,0	-83,2	0,132	3,9	7,2	9,3	7,5	2,3	9,1	2,6	11,4
22:38:43	031:39:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
22:38:43	031:39:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
22:38:43	031:39:03	0,183	10,0	-90,6	0,133	3,9	7,2	9,3	7,6	2,3	9,1	2,6	11,4
22:39:43	031:40:03	0,184	10,0	-87,3	0,134	3,9	7,2	9,3	7,7	2,3	9,1	2,6	11,4
22:40:41	031:41:01	0,185	10,0	-78,7	0,135	3,9	7,2	9,3	7,7	2,3	9,1	2,6	11,4
22:41:41	031:42:01	0,187	10,0	-88,1	0,136	3,9	7,2	9,3	7,8	2,4	9,2	2,7	11,4
22:42:42	031:43:02	0,188	10,0	-84,3	0,137	3,9	7,2	9,3	7,9	2,4	9,2	2,7	11,4
22:43:40	031:44:00	0,189	10,0	-81,0	0,138	3,9	7,2	9,3	8,0	2,4	9,2	2,7	11,4
22:44:43	031:45:03	0,179	10,0	-88,7	0,129	4,0	7,3	9,3	8,2	2,4	9,2	2,7	11,5
22:45:44	031:46:04	0,178	10,0	-81,0	0,129	4,0	7,3	9,3	8,3	2,5	9,2	2,7	11,5
22:46:41	031:47:01	0,179	10,0	-89,6	0,130	4,0	7,3	9,3	8,4	2,5	9,2	2,8	11,5
22:47:42	031:48:02	0,180	10,0	-87,5	0,131	4,0	7,3	9,3	8,5	2,5	9,2	2,8	11,5
22:48:43	031:49:03	0,180	10,0	-77,9	0,132	4,0	7,3	9,3	8,6	2,5	9,2	2,8	11,5
22:49:44	031:50:04	0,182	10,0	-88,0	0,133	4,0	7,3	9,4	8,7	2,6	9,2	2,8	11,5
22:50:40	031:51:00	0,183	10,0	-85,1	0,133	4,0	7,3	9,4	8,7	2,6	9,2	2,8	11,5
22:51:40	031:52:00	0,185	10,0	-88,1	0,135	4,1	7,3	9,4	8,8	2,6	9,2	2,8	11,5
22:52:41	031:53:01	0,186	10,0	-89,1	0,135	4,1	7,3	9,4	8,9	2,6	9,2	2,8	11,5
22:53:41	031:54:01	0,186	10,0	-80,7	0,136	4,1	7,3	9,4	8,9	2,6	9,2	2,9	11,5
22:54:44	031:55:04	0,188	10,0	-88,1	0,137	4,1	7,4	9,4	9,0	2,7	9,2	2,9	11,5
22:55:44	031:56:04	0,178	10,0	-84,3	0,129	4,1	7,4	9,4	9,0	2,7	9,2	2,9	11,5

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	103 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	22:56:41	031:57:01	0,177	10,0	-90,5	0,129	4,2	7,4	9,4	9,1	2,7	9,2	2,9	11,5	
	22:57:45	031:58:05	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:57:45	031:58:05	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:57:45	031:58:05	0,177	10,0	-85,6	0,130	4,2	7,4	9,4	9,1	2,8	9,2	2,9	11,5	
	22:58:41	031:59:01	0,178	10,0	-77,6	0,131	4,2	7,4	9,4	9,1	2,8	9,2	3,0	11,5	
	22:59:43	032:00:03	0,179	10,0	-88,0	0,131	4,2	7,4	9,4	9,1	2,8	9,2	3,0	11,5	
	23:00:44	032:01:04	0,180	10,0	-83,4	0,132	4,2	7,4	9,4	9,1	2,8	9,2	3,0	11,5	
	23:01:42	032:02:02	0,181	10,0	-88,3	0,133	4,2	7,4	9,4	9,2	2,8	9,2	3,0	11,5	
	23:02:41	032:03:01	0,182	10,0	-88,3	0,134	4,2	7,4	9,4	9,2	2,9	9,2	3,0	11,5	
	23:03:41	032:04:01	0,183	10,0	-78,8	0,135	4,2	7,4	9,4	9,2	2,9	9,2	3,0	11,5	
	23:04:41	032:05:01	0,184	10,0	-89,0	0,136	4,3	7,4	9,4	9,2	2,9	9,2	3,0	11,5	
	23:05:41	032:06:01	0,185	10,0	-83,0	0,137	4,3	7,5	9,4	9,3	2,9	9,2	3,1	11,5	
	23:06:43	032:07:03	0,187	10,0	-90,2	0,138	4,3	7,5	9,4	9,3	3,0	9,2	3,1	11,5	
	23:07:43	032:08:03	0,176	10,0	-86,1	0,129	4,3	7,5	9,4	9,3	3,0	9,2	3,1	11,5	
	23:08:41	032:09:01	0,176	10,0	-78,3	0,130	4,3	7,5	9,4	9,3	3,0	9,2	3,1	11,5	
	23:09:42	032:10:02	0,176	10,0	-88,1	0,130	4,3	7,5	9,4	9,3	3,0	9,2	3,2	11,5	
	23:10:43	032:11:03	0,177	10,0	-83,1	0,131	4,4	7,5	9,4	9,3	3,1	9,2	3,2	11,6	
	23:11:41	032:12:01	0,178	10,0	-79,0	0,132	4,4	7,5	9,4	9,3	3,1	9,2	3,2	11,6	
	23:12:43	032:13:03	0,179	10,0	-89,0	0,133	4,4	7,5	9,4	9,4	3,1	9,2	3,2	11,6	
	23:13:44	032:14:04	0,181	10,0	-79,6	0,134	4,4	7,5	9,5	9,4	3,1	9,2	3,2	11,6	
	23:14:41	032:15:01	0,182	10,0	-88,9	0,135	4,4	7,5	9,5	9,4	3,1	9,2	3,2	11,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 104 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:15:42	032:16:02	0,184	10,0	-85,3	0,135	4,4	7,6	9,5	9,4	3,1	9,2	3,2	11,6	
	23:16:43	032:17:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:16:43	032:17:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:16:43	032:17:03	0,184	10,0	-90,3	0,137	4,4	7,6	9,5	9,4	3,1	9,2	3,3	11,6	
	23:17:40	032:18:00	0,185	10,0	-89,2	0,137	4,5	7,6	9,5	9,4	3,2	9,2	3,3	11,6	
	23:18:40	032:19:00	0,176	10,0	-79,7	0,130	4,5	7,6	9,5	9,5	3,2	9,2	3,3	11,6	
	23:19:41	032:20:01	0,175	10,0	-88,8	0,130	4,5	7,6	9,5	9,5	3,2	9,2	3,3	11,6	
	23:20:42	032:21:02	0,176	10,0	-83,2	0,130	4,5	7,6	9,5	9,5	3,3	9,2	3,3	11,6	
	23:21:44	032:22:04	0,176	10,0	-88,7	0,131	4,5	7,6	9,5	9,5	3,3	9,2	3,3	11,6	
	23:22:40	032:23:00	0,177	10,0	-88,6	0,132	4,5	7,6	9,5	9,5	3,3	9,2	3,4	11,6	
	23:23:41	032:24:01	0,178	10,0	-84,4	0,132	4,5	7,6	9,5	9,5	3,3	9,2	3,4	11,6	
	23:24:41	032:25:01	0,179	10,0	-89,6	0,133	4,5	7,6	9,5	9,5	3,4	9,2	3,4	11,6	
	23:25:42	032:26:02	0,181	10,0	-82,2	0,134	4,6	7,6	9,5	9,6	3,4	9,2	3,4	11,6	
	23:26:40	032:27:00	0,181	10,0	-90,6	0,135	4,6	7,6	9,5	9,6	3,4	9,2	3,4	11,6	
	23:27:40	032:28:00	0,182	10,0	-86,2	0,136	4,6	7,7	9,5	9,6	3,5	9,2	3,4	11,6	
	23:28:42	032:29:02	0,183	10,0	-77,1	0,137	4,6	7,7	9,5	9,6	3,5	9,2	3,5	11,6	
	23:29:43	032:30:03	0,175	10,0	-88,2	0,130	4,6	7,7	9,5	9,6	3,5	9,2	3,5	11,6	
	23:30:40	032:31:00	0,174	10,0	-84,7	0,129	4,6	7,7	9,5	9,6	3,5	9,2	3,5	11,6	
	23:31:42	032:32:02	0,174	10,0	-88,4	0,130	4,7	7,7	9,5	9,6	3,5	9,2	3,5	11,6	
	23:32:43	032:33:03	0,175	10,0	-88,8	0,131	4,7	7,7	9,5	9,6	3,6	9,2	3,5	11,6	
	23:33:44	032:34:04	0,176	10,0	-79,3	0,131	4,7	7,7	9,5	9,7	3,6	9,2	3,6	11,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 105 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:34:40	032:35:00	0,176	10,0	-89,0	0,132	4,7	7,7	9,5	9,7	3,6	9,2	3,6	11,6	
	23:35:42	032:36:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:35:42	032:36:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:35:42	032:36:02	0,177	10,0	-85,1	0,133	4,7	7,7	9,5	9,7	3,6	9,2	3,6	11,6	
	23:36:44	032:37:04	0,179	10,0	-76,5	0,134	4,7	7,7	9,5	9,7	3,6	9,2	3,6	11,6	
	23:37:45	032:38:05	0,180	10,0	-88,8	0,135	4,7	7,7	9,5	9,7	3,7	9,2	3,6	11,6	
	23:38:43	032:39:03	0,181	10,0	-80,7	0,136	4,8	7,8	9,6	9,7	3,7	9,2	3,6	11,6	
	23:39:41	032:40:01	0,183	10,0	-89,5	0,137	4,8	7,8	9,6	9,7	3,7	9,2	3,7	11,6	
	23:40:44	032:41:04	0,183	10,0	-85,4	0,138	4,8	7,8	9,6	9,8	3,7	9,2	3,7	11,6	
	23:41:40	032:42:00	0,173	10,0	-78,7	0,130	4,8	7,8	9,6	9,8	3,7	9,2	3,7	11,6	
	23:42:40	032:43:00	0,172	10,0	-89,1	0,130	4,8	7,8	9,6	9,8	3,8	9,2	3,7	11,7	
	23:43:43	032:44:03	0,173	10,0	-80,5	0,130	4,8	7,8	9,6	9,8	3,8	9,2	3,7	11,6	
	23:44:44	032:45:04	0,174	10,0	-89,7	0,131	4,8	7,8	9,6	9,8	3,8	9,2	3,7	11,7	
	23:45:41	032:46:01	0,175	10,0	-85,7	0,132	4,8	7,8	9,6	9,8	3,8	9,2	3,7	11,7	
	23:46:43	032:47:03	0,176	10,0	-78,1	0,133	4,9	7,8	9,6	9,8	3,8	9,2	3,8	11,7	
	23:47:40	032:48:00	0,177	10,0	-88,2	0,134	4,9	7,8	9,6	9,8	3,8	9,2	3,8	11,7	
	23:48:41	032:49:01	0,178	10,0	-83,6	0,135	4,9	7,8	9,6	9,8	3,8	9,2	3,8	11,7	
	23:49:43	032:50:03	0,179	10,0	-88,4	0,136	4,9	7,8	9,6	9,8	3,9	9,2	3,8	11,7	
	23:50:43	032:51:03	0,180	10,0	-87,7	0,136	4,9	7,9	9,6	9,8	3,9	9,2	3,8	11,7	
	23:51:44	032:52:04	0,181	10,0	-78,5	0,137	4,9	7,9	9,6	9,8	3,9	9,2	3,8	11,7	
	23:52:40	032:53:00	0,173	10,0	-89,5	0,130	4,9	7,9	9,6	9,8	3,9	9,2	3,9	11,7	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 106 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	23:53:40	032:54:00	0,171	10,0	-85,6	0,129	4,9	7,9	9,6	9,8	4,0	9,2	3,9	11,7	
	23:54:41	032:55:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:54:41	032:55:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:54:41	032:55:01	0,171	10,0	-89,7	0,130	4,9	7,9	9,6	9,9	4,0	9,2	3,9	11,7	
	23:55:42	032:56:02	0,172	10,0	-89,2	0,131	5,0	7,9	9,6	9,9	4,0	9,2	3,9	11,7	
	23:56:44	032:57:04	0,173	10,0	-80,9	0,131	5,0	7,9	9,6	9,9	4,1	9,2	3,9	11,7	
	23:57:41	032:58:01	0,174	10,0	-90,1	0,132	5,0	7,9	9,6	9,9	4,1	9,2	3,9	11,7	
	23:58:41	032:59:01	0,175	10,0	-86,8	0,133	5,0	7,9	9,6	9,9	4,1	9,2	3,9	11,7	
	23:59:42	033:00:02	0,176	10,0	-80,6	0,134	5,0	7,9	9,6	9,9	4,1	9,2	4,0	11,7	
2021-11-25	00:00:43	033:01:03	0,177	10,0	-89,1	0,135	5,0	7,9	9,6	9,9	4,2	9,2	4,0	11,7	
	00:01:43	033:02:03	0,178	10,0	-80,1	0,136	5,1	7,9	9,6	9,9	4,2	9,2	4,0	11,7	
	00:02:44	033:03:04	0,180	10,0	-89,8	0,137	5,1	8,0	9,6	9,9	4,2	9,2	4,0	11,7	
	00:03:40	033:04:00	0,180	10,0	-85,9	0,138	5,1	8,0	9,6	9,9	4,2	9,2	4,0	11,7	
	00:04:40	033:05:00	0,171	10,0	-77,2	0,130	5,1	8,0	9,7	9,9	4,3	9,2	4,1	11,7	
	00:05:41	033:06:01	0,170	10,0	-87,8	0,130	5,1	8,0	9,7	10,0	4,3	9,3	4,1	11,7	
	00:06:41	033:07:01	0,171	10,0	-83,8	0,130	5,1	8,0	9,7	10,0	4,3	9,2	4,1	11,7	
	00:07:42	033:08:02	0,173	10,0	-91,4	0,131	5,1	8,0	9,6	10,0	4,3	9,3	4,1	11,7	
	00:08:40	033:09:00	0,173	10,0	-89,2	0,132	5,2	8,0	9,7	10,0	4,3	9,3	4,1	11,7	
	00:09:40	033:10:00	0,174	10,0	-80,7	0,133	5,2	8,0	9,7	10,0	4,3	9,3	4,1	11,7	
	00:10:40	033:11:00	0,175	10,0	-89,3	0,134	5,2	8,0	9,7	10,0	4,4	9,2	4,1	11,7	
	00:11:44	033:12:04	0,175	10,0	-85,4	0,135	5,2	8,0	9,7	10,0	4,4	9,2	4,2	11,7	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	107 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	00:12:40	033:13:00	0,177	10,0	-90,5	0,136	5,2	8,0	9,7	10,0	4,4	9,2	4,2	11,7	
	00:13:43	033:14:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:13:43	033:14:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:13:43	033:14:03	0,179	10,0	-89,4	0,136	5,2	8,0	9,7	10,0	4,4	9,3	4,2	11,7	
	00:14:40	033:15:00	0,180	10,0	-81,3	0,137	5,2	8,0	9,7	10,0	4,5	9,3	4,2	11,7	
	00:15:40	033:16:00	0,170	10,0	-88,3	0,130	5,3	8,0	9,7	10,0	4,5	9,3	4,2	11,7	
	00:16:41	033:17:01	0,169	10,0	-84,3	0,129	5,3	8,0	9,7	10,0	4,5	9,3	4,2	11,7	
	00:17:43	033:18:03	0,169	10,0	-79,3	0,130	5,3	8,0	9,7	10,0	4,5	9,2	4,3	11,7	
	00:18:44	033:19:04	0,169	10,0	-88,5	0,131	5,3	8,1	9,7	10,0	4,5	9,2	4,3	11,7	
	00:19:40	033:20:00	0,170	10,0	-79,8	0,131	5,3	8,1	9,7	10,0	4,6	9,2	4,3	11,7	
	00:20:40	033:21:00	0,171	10,0	-89,8	0,132	5,3	8,1	9,7	10,0	4,6	9,3	4,3	11,8	
	00:21:40	033:22:00	0,173	10,0	-84,9	0,133	5,3	8,1	9,7	10,1	4,6	9,3	4,3	11,8	
	00:22:42	033:23:02	0,174	10,0	-76,2	0,134	5,3	8,1	9,7	10,1	4,6	9,3	4,3	11,8	
	00:23:44	033:24:04	0,174	10,0	-88,2	0,135	5,4	8,1	9,7	10,1	4,7	9,3	4,3	11,8	
	00:24:40	033:25:00	0,175	10,0	-84,0	0,136	5,4	8,1	9,7	10,1	4,7	9,3	4,3	11,8	
	00:25:42	033:26:02	0,176	10,0	-91,4	0,137	5,4	8,1	9,7	10,1	4,7	9,3	4,3	11,8	
	00:26:42	033:27:02	0,178	10,0	-87,9	0,137	5,4	8,1	9,7	10,1	4,7	9,3	4,4	11,8	
	00:27:43	033:28:03	0,169	10,0	-78,4	0,130	5,4	8,1	9,7	10,1	4,7	9,3	4,4	11,8	
	00:28:40	033:29:00	0,168	10,0	-89,2	0,130	5,4	8,1	9,7	10,1	4,7	9,3	4,4	11,8	
	00:29:42	033:30:02	0,169	10,0	-85,4	0,130	5,5	8,1	9,7	10,1	4,8	9,3	4,4	11,8	
	00:30:44	033:31:04	0,168	10,0	-87,4	0,131	5,5	8,2	9,7	10,1	4,8	9,3	4,4	11,8	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	108 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	00:31:40	033:32:00	0,169	10,0	-89,1	0,132	5,5	8,1	9,7	10,1	4,8	9,3	4,5	11,8	
	00:32:41	033:33:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:32:41	033:33:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:32:41	033:33:01	0,170	10,0	-81,7	0,132	5,5	8,2	9,7	10,1	4,8	9,3	4,5	11,8	
	00:33:42	033:34:02	0,172	10,0	-90,6	0,133	5,5	8,2	9,7	10,1	4,9	9,3	4,5	11,8	
	00:34:42	033:35:02	0,172	10,0	-86,5	0,134	5,5	8,2	9,8	10,1	4,9	9,3	4,5	11,8	
	00:35:44	033:36:04	0,173	10,0	-77,7	0,135	5,5	8,2	9,7	10,1	4,9	9,3	4,5	11,8	
	00:36:40	033:37:00	0,174	10,0	-87,9	0,136	5,6	8,2	9,7	10,1	4,9	9,3	4,5	11,8	
	00:37:43	033:38:03	0,175	10,0	-83,6	0,136	5,6	8,2	9,7	10,1	4,9	9,3	4,5	11,8	
	00:38:41	033:39:01	0,177	10,0	-79,3	0,137	5,6	8,2	9,7	10,1	5,0	9,3	4,6	11,8	
	00:39:41	033:40:01	0,166	10,0	-88,8	0,130	5,6	8,2	9,7	10,2	5,0	9,3	4,6	11,8	
	00:40:41	033:41:01	0,166	10,0	-79,4	0,130	5,6	8,2	9,8	10,2	5,0	9,3	4,6	11,8	
	00:41:41	033:42:01	0,166	10,0	-89,3	0,130	5,6	8,2	9,8	10,1	5,0	9,3	4,6	11,8	
	00:42:42	033:43:02	0,167	10,0	-85,2	0,131	5,7	8,2	9,8	10,1	5,1	9,3	4,6	11,8	
	00:43:43	033:44:03	0,169	10,0	-81,9	0,132	5,7	8,3	9,8	10,2	5,1	9,3	4,6	11,8	
	00:44:43	033:45:03	0,169	10,0	-89,2	0,133	5,7	8,3	9,8	10,2	5,1	9,3	4,7	11,8	
	00:45:45	033:46:05	0,170	10,0	-81,3	0,133	5,7	8,2	9,8	10,2	5,1	9,3	4,7	11,8	
	00:46:40	033:47:00	0,171	10,0	-90,5	0,134	5,7	8,2	9,8	10,2	5,1	9,3	4,7	11,8	
	00:47:40	033:48:00	0,172	10,0	-86,6	0,135	5,7	8,3	9,8	10,2	5,1	9,3	4,7	11,8	
	00:48:42	033:49:02	0,173	10,0	-79,1	0,136	5,7	8,3	9,8	10,2	5,1	9,3	4,7	11,8	
	00:49:43	033:50:03	0,174	10,0	-89,3	0,137	5,7	8,3	9,8	10,2	5,2	9,3	4,7	11,8	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 109 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	00:50:44	033:51:04	0,175	10,0	-81,2	0,137	5,8	8,3	9,8	10,2	5,2	9,3	4,7	11,8	
	00:51:41	033:52:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:51:41	033:52:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:51:41	033:52:01	0,166	10,0	-90,3	0,130	5,8	8,3	9,8	10,2	5,2	9,3	4,7	11,8	
	00:52:42	033:53:02	0,165	10,0	-85,8	0,130	5,8	8,3	9,8	10,2	5,3	9,3	4,8	11,8	
	00:53:42	033:54:02	0,166	10,0	-76,7	0,130	5,8	8,3	9,8	10,2	5,3	9,3	4,8	11,8	
	00:54:43	033:55:03	0,166	10,0	-88,0	0,131	5,8	8,3	9,8	10,2	5,3	9,3	4,8	11,8	
	00:55:44	033:56:04	0,167	10,0	-84,2	0,132	5,8	8,3	9,8	10,3	5,3	9,3	4,8	11,8	
	00:56:41	033:57:01	0,168	10,0	-92,2	0,133	5,8	8,3	9,8	10,3	5,3	9,3	4,8	11,8	
	00:57:43	033:58:03	0,169	10,0	-88,4	0,133	5,8	8,3	9,8	10,3	5,3	9,3	4,8	11,8	
	00:58:44	033:59:04	0,170	10,0	-78,8	0,134	5,8	8,3	9,8	10,3	5,4	9,3	4,9	11,8	
	00:59:40	034:00:00	0,171	10,0	-89,7	0,135	5,9	8,3	9,8	10,3	5,4	9,3	4,9	11,8	
	01:00:40	034:01:00	0,172	10,0	-85,5	0,136	5,9	8,3	9,8	10,3	5,4	9,3	4,9	11,9	
	01:01:41	034:02:01	0,173	10,0	-76,8	0,137	5,9	8,3	9,8	10,3	5,4	9,3	4,9	11,8	
	01:02:41	034:03:01	0,174	10,0	-88,6	0,138	5,9	8,3	9,8	10,3	5,5	9,3	4,9	11,8	
	01:03:43	034:04:03	0,164	10,0	-82,4	0,130	5,9	8,3	9,8	10,3	5,5	9,3	4,9	11,9	
	01:04:44	034:05:04	0,164	10,0	-90,2	0,130	5,9	8,4	9,8	10,3	5,5	9,3	5,0	11,9	
	01:05:44	034:06:04	0,164	10,0	-88,0	0,131	5,9	8,4	9,8	10,3	5,5	9,3	5,0	11,9	
	01:06:44	034:07:04	0,165	10,0	-78,6	0,131	6,0	8,4	9,8	10,3	5,5	9,3	5,0	11,8	
	01:07:41	034:08:01	0,166	10,0	-89,6	0,132	6,0	8,4	9,8	10,3	5,5	9,3	5,0	11,9	
	01:08:42	034:09:02	0,167	10,0	-84,8	0,133	6,0	8,4	9,9	10,3	5,5	9,3	5,0	11,9	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	110 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:09:43	034:10:03	0,168	10,0	-75,7	0,134	6,0	8,4	9,9	10,3	5,6	9,3	5,0	11,9	
	01:10:40	034:11:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:10:40	034:11:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:10:40	034:11:00	0,169	10,0	-87,7	0,135	6,0	8,4	9,8	10,3	5,6	9,3	5,0	11,9	
	01:11:40	034:12:00	0,170	10,0	-83,8	0,136	6,0	8,4	9,9	10,3	5,6	9,3	5,0	11,9	
	01:12:41	034:13:01	0,171	10,0	-76,8	0,137	6,0	8,4	9,9	10,3	5,7	9,3	5,1	11,9	
	01:13:41	034:14:01	0,172	10,0	-88,9	0,137	6,0	8,4	9,9	10,3	5,7	9,3	5,1	11,9	
	01:14:44	034:15:04	0,163	10,0	-80,5	0,130	6,1	8,4	9,9	10,4	5,7	9,3	5,1	11,9	
	01:15:40	034:16:00	0,163	10,0	-90,1	0,130	6,1	8,4	9,9	10,4	5,8	9,3	5,1	11,9	
	01:16:41	034:17:01	0,164	10,0	-85,5	0,130	6,1	8,4	9,9	10,4	5,8	9,3	5,1	11,9	
	01:17:42	034:18:02	0,164	10,0	-77,5	0,131	6,1	8,4	9,8	10,4	5,8	9,3	5,1	11,9	
	01:18:42	034:19:02	0,165	10,0	-88,8	0,132	6,1	8,4	9,8	10,4	5,8	9,3	5,2	11,9	
	01:19:44	034:20:04	0,167	10,0	-82,6	0,133	6,1	8,4	9,9	10,4	5,9	9,3	5,2	11,9	
	01:20:40	034:21:00	0,167	10,0	-76,9	0,134	6,2	8,4	9,9	10,4	5,9	9,3	5,2	11,9	
	01:21:42	034:22:02	0,168	10,0	-88,4	0,134	6,2	8,4	9,9	10,4	5,9	9,3	5,2	11,9	
	01:22:43	034:23:03	0,169	10,0	-79,2	0,135	6,2	8,4	9,9	10,4	6,0	9,3	5,2	11,9	
	01:23:43	034:24:03	0,171	10,0	-89,4	0,136	6,2	8,4	9,9	10,4	6,0	9,3	5,2	11,9	
	01:24:45	034:25:05	0,171	10,0	-85,5	0,137	6,2	8,4	9,9	10,4	6,0	9,3	5,2	11,9	
	01:25:41	034:26:01	0,172	10,0	-82,4	0,138	6,2	8,5	9,9	10,4	6,1	9,3	5,3	11,9	
	01:26:42	034:27:02	0,163	10,0	-89,2	0,130	6,2	8,5	9,9	10,4	6,1	9,3	5,3	11,9	
	01:27:43	034:28:03	0,161	10,0	-81,2	0,130	6,3	8,5	9,9	10,4	6,1	9,3	5,3	11,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature:
Document issued by: Page: 111 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:28:44	034:29:04	0,162	10,0	-89,0	0,131	6,3	8,5	9,9	10,4	6,2	9,3	5,2	11,9	
	01:29:44	034:30:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:29:44	034:30:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:29:44	034:30:04	0,164	10,0	-85,1	0,131	6,3	8,5	9,9	10,4	6,2	9,3	5,2	11,9	
	01:30:41	034:31:01	0,164	10,0	-76,5	0,132	6,3	8,5	9,9	10,4	6,2	9,3	5,2	11,9	
	01:31:42	034:32:02	0,165	10,0	-87,7	0,133	6,3	8,5	9,9	10,4	6,2	9,3	5,3	11,9	
	01:32:42	034:33:02	0,166	10,0	-83,2	0,134	6,3	8,5	9,9	10,4	6,3	9,3	5,3	11,9	
	01:33:43	034:34:03	0,167	10,0	-88,8	0,135	6,3	8,5	9,9	10,4	6,3	9,3	5,3	11,9	
	01:34:44	034:35:04	0,168	10,0	-88,7	0,135	6,3	8,5	9,9	10,4	6,4	9,3	5,3	11,9	
	01:35:40	034:36:00	0,169	10,0	-79,7	0,136	6,4	8,5	9,9	10,4	6,4	9,3	5,3	11,9	
	01:36:40	034:37:00	0,170	10,0	-88,7	0,137	6,4	8,5	9,9	10,4	6,4	9,3	5,4	11,9	
	01:37:41	034:38:01	0,165	10,0	-85,0	0,131	6,4	8,5	9,9	10,4	6,4	9,3	5,4	11,9	
	01:38:41	034:39:01	0,161	10,0	-78,2	0,130	6,4	8,5	9,9	10,4	6,5	9,3	5,4	11,9	
	01:39:40	034:40:00	0,161	10,0	-89,2	0,130	6,4	8,5	9,9	10,4	6,5	9,3	5,4	11,9	
	01:40:40	034:41:00	0,162	10,0	-81,4	0,131	6,4	8,5	9,9	10,4	6,6	9,3	5,4	11,9	
	01:41:40	034:42:00	0,163	10,0	-90,3	0,132	6,5	8,5	9,9	10,4	6,6	9,3	5,4	11,9	
	01:42:41	034:43:01	0,163	10,0	-85,9	0,132	6,5	8,5	9,9	10,4	6,6	9,3	5,5	11,9	
	01:43:42	034:44:02	0,165	10,0	-76,7	0,133	6,5	8,5	9,9	10,4	6,7	9,3	5,5	11,9	
	01:44:42	034:45:02	0,166	10,0	-87,7	0,134	6,5	8,5	9,9	10,4	6,7	9,3	5,5	11,9	
	01:45:43	034:46:03	0,166	10,0	-83,7	0,135	6,5	8,5	9,9	10,4	6,7	9,3	5,5	11,9	
	01:46:40	034:47:00	0,169	10,0	-87,4	0,136	6,5	8,5	9,9	10,4	6,7	9,3	5,5	11,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 112 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:47:42	034:48:02	0,168	10,0	-88,6	0,137	6,5	8,5	9,9	10,5	6,8	9,3	5,5	11,9	
	01:48:44	034:49:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:48:44	034:49:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:48:44	034:49:04	0,170	10,0	-79,0	0,138	6,6	8,6	9,9	10,5	6,8	9,3	5,6	12,0	
	01:49:44	034:50:04	0,161	10,0	-88,8	0,130	6,6	8,6	9,9	10,5	6,8	9,3	5,6	11,9	
	01:50:44	034:51:04	0,161	10,0	-84,7	0,130	6,6	8,6	9,9	10,5	6,9	9,3	5,6	11,9	
	01:51:42	034:52:02	0,160	10,0	-79,2	0,131	6,6	8,6	9,9	10,5	6,9	9,3	5,6	11,9	
	01:52:44	034:53:04	0,161	10,0	-89,2	0,132	6,6	8,6	9,9	10,5	6,9	9,3	5,6	11,9	
	01:53:41	034:54:01	0,162	10,0	-81,9	0,132	6,6	8,6	9,9	10,5	7,0	9,3	5,7	11,9	
	01:54:44	034:55:04	0,163	10,0	-90,1	0,133	6,7	8,6	9,9	10,5	7,0	9,3	5,7	11,9	
	01:55:40	034:56:00	0,164	10,0	-87,3	0,134	6,7	8,6	9,9	10,5	7,0	9,3	5,7	11,9	
	01:56:42	034:57:02	0,165	10,0	-77,8	0,135	6,7	8,6	9,9	10,5	7,1	9,3	5,7	11,9	
	01:57:42	034:58:02	0,166	10,0	-88,7	0,136	6,7	8,6	9,9	10,5	7,1	9,3	5,7	11,9	
	01:58:42	034:59:02	0,167	10,0	-85,2	0,136	6,7	8,6	9,9	10,5	7,1	9,3	5,7	11,9	
	01:59:44	035:00:04	0,168	10,0	-92,3	0,137	6,7	8,6	9,9	10,5	7,2	9,3	5,8	12,0	
	02:00:44	035:01:04	0,160	10,0	-88,9	0,130	6,7	8,6	9,9	10,5	7,2	9,3	5,8	11,9	
	02:01:41	035:02:01	0,160	10,0	-80,1	0,130	6,8	8,6	9,9	10,5	7,2	9,3	5,8	11,9	
	02:02:42	035:03:02	0,160	10,0	-89,9	0,131	6,8	8,6	9,9	10,5	7,2	9,3	5,8	11,9	
	02:03:43	035:04:03	0,161	10,0	-86,1	0,131	6,8	8,6	9,9	10,5	7,3	9,3	5,8	11,9	
	02:04:40	035:05:00	0,161	10,0	-77,8	0,132	6,8	8,6	9,9	10,5	7,3	9,3	5,8	12,0	
	02:05:42	035:06:02	0,162	10,0	-88,2	0,133	6,8	8,6	9,9	10,5	7,3	9,3	5,8	11,9	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 113 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:06:44	035:07:04	0,163	10,0	-83,2	0,133	6,8	8,6	9,9	10,5	7,3	9,3	5,9	11,9	
	02:07:44	035:08:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:07:44	035:08:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:07:44	035:08:04	0,164	10,0	-90,8	0,134	6,9	8,6	9,9	10,5	7,4	9,3	5,9	11,9	
	02:08:41	035:09:01	0,166	10,0	-87,6	0,135	6,9	8,6	9,9	10,5	7,4	9,3	5,9	11,9	
	02:09:40	035:10:00	0,167	10,0	-78,6	0,136	6,9	8,6	9,9	10,5	7,4	9,3	5,9	11,9	
	02:10:44	035:11:04	0,167	10,0	-88,6	0,137	6,9	8,6	9,9	10,5	7,4	9,3	5,9	11,9	
	02:11:41	035:12:01	0,168	10,0	-85,9	0,138	6,9	8,6	9,9	10,5	7,5	9,3	6,0	11,9	
	02:12:44	035:13:04	0,159	10,0	-91,8	0,130	7,0	8,6	9,9	10,5	7,5	9,3	6,0	12,0	
	02:13:45	035:14:05	0,159	10,0	-88,0	0,130	7,0	8,6	9,9	10,5	7,5	9,4	6,0	12,0	
	02:14:42	035:15:02	0,159	10,0	-79,4	0,131	7,0	8,7	9,9	10,5	7,6	9,3	6,0	12,0	
	02:15:44	035:16:04	0,159	10,0	-89,0	0,132	7,0	8,7	9,9	10,5	7,6	9,3	6,0	12,0	
	02:16:44	035:17:04	0,160	10,0	-85,6	0,132	7,0	8,7	9,9	10,5	7,6	9,4	6,1	11,9	
	02:17:40	035:18:00	0,161	10,0	-92,0	0,133	7,0	8,7	9,9	10,5	7,7	9,4	6,1	11,9	
	02:18:41	035:19:01	0,162	10,0	-89,5	0,134	7,0	8,7	10,0	10,5	7,7	9,3	6,1	11,9	
	02:19:42	035:20:02	0,163	10,0	-81,1	0,135	7,1	8,7	10,0	10,5	7,8	9,3	6,1	12,0	
	02:20:43	035:21:03	0,164	10,0	-89,5	0,136	7,1	8,7	10,0	10,5	7,8	9,3	6,1	12,0	
	02:21:44	035:22:04	0,165	10,0	-85,6	0,136	7,1	8,7	10,0	10,5	7,9	9,4	6,1	12,0	
	02:22:46	035:23:06	0,166	10,0	-75,2	0,137	7,1	8,7	9,9	10,5	7,9	9,3	6,2	11,9	
	02:23:43	035:24:03	0,158	10,0	-88,9	0,130	7,2	8,7	9,9	10,5	8,0	9,3	6,2	12,0	
	02:24:43	035:25:03	0,158	10,0	-81,8	0,130	7,2	8,7	9,9	10,5	8,0	9,3	6,2	12,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	114 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:25:41	035:26:01	0,158	10,0	-91,2	0,130	7,2	8,7	9,9	10,5	8,0	9,4	6,2	12,0	
	02:26:41	035:27:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:26:41	035:27:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:26:41	035:27:01	0,159	10,0	-87,0	0,131	7,2	8,7	10,0	10,5	8,0	9,4	6,2	12,0	
	02:27:41	035:28:01	0,159	10,0	-78,0	0,132	7,2	8,7	10,0	10,5	8,0	9,3	6,3	12,0	
	02:28:44	035:29:04	0,160	10,0	-87,8	0,133	7,3	8,7	10,0	10,5	8,0	9,3	6,3	12,0	
	02:29:41	035:30:01	0,161	10,0	-83,7	0,133	7,3	8,7	10,0	10,5	8,1	9,3	6,3	12,0	
	02:30:43	035:31:03	0,162	10,0	-86,7	0,134	7,3	8,7	10,0	10,5	8,1	9,4	6,3	12,0	
	02:31:43	035:32:03	0,164	10,0	-89,1	0,135	7,3	8,7	10,0	10,5	8,1	9,4	6,4	12,0	
	02:32:42	035:33:02	0,164	10,0	-81,1	0,136	7,4	8,7	10,0	10,5	8,1	9,4	6,4	12,0	
	02:33:42	035:34:02	0,165	10,0	-90,4	0,137	7,4	8,7	10,0	10,6	8,1	9,3	6,4	12,0	
	02:34:44	035:35:04	0,166	10,0	-85,9	0,138	7,4	8,7	10,0	10,5	8,1	9,4	6,4	12,0	
	02:35:44	035:36:04	0,157	10,0	-77,4	0,130	7,4	8,7	10,0	10,5	8,1	9,4	6,4	12,0	
	02:36:43	035:37:03	0,157	10,0	-88,0	0,130	7,4	8,7	10,0	10,5	8,1	9,4	6,5	12,0	
	02:37:44	035:38:04	0,158	10,0	-83,3	0,131	7,5	8,7	10,0	10,5	8,2	9,4	6,5	12,0	
	02:38:43	035:39:03	0,158	10,0	-88,2	0,132	7,5	8,7	10,0	10,5	8,2	9,4	6,5	12,0	
	02:39:40	035:40:00	0,159	10,0	-89,0	0,132	7,5	8,7	10,0	10,5	8,2	9,3	6,5	12,0	
	02:40:42	035:41:02	0,160	10,0	-79,3	0,133	7,5	8,7	10,0	10,5	8,1	9,3	6,6	12,0	
	02:41:44	035:42:04	0,161	10,0	-88,4	0,134	7,5	8,7	10,0	10,5	8,2	9,3	6,6	12,0	
	02:42:44	035:43:04	0,162	10,0	-84,8	0,135	7,6	8,7	10,0	10,5	8,2	9,3	6,6	12,0	
	02:43:40	035:44:00	0,162	10,0	-84,3	0,135	7,6	8,7	10,0	10,5	8,2	9,4	6,6	12,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 115 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:44:42	035:45:02	0,164	10,0	-89,4	0,136	7,6	8,7	10,0	10,5	8,2	9,4	6,6	12,0	
	02:45:40	035:46:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:45:40	035:46:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:45:40	035:46:00	0,164	10,0	-81,5	0,137	7,6	8,7	10,0	10,5	8,2	9,4	6,7	12,0	
	02:46:42	035:47:02	0,157	10,0	-89,7	0,130	7,6	8,7	10,0	10,5	8,2	9,4	6,7	12,0	
	02:47:43	035:48:03	0,156	10,0	-84,4	0,130	7,7	8,7	10,0	10,5	8,2	9,4	6,7	12,0	
	02:48:45	035:49:05	0,156	10,0	-76,3	0,130	7,7	8,7	10,0	10,5	8,2	9,4	6,7	12,0	
	02:49:43	035:50:03	0,157	10,0	-88,5	0,131	7,7	8,7	10,0	10,5	8,2	9,4	6,7	12,0	
	02:50:40	035:51:00	0,158	10,0	-83,3	0,132	7,7	8,7	10,0	10,5	8,2	9,4	6,8	12,0	
	02:51:41	035:52:01	0,159	10,0	-83,0	0,132	7,8	8,7	10,0	10,5	8,2	9,4	6,8	12,0	
	02:52:44	035:53:04	0,159	10,0	-88,7	0,133	7,8	8,7	10,0	10,5	8,2	9,4	6,8	12,0	
	02:53:44	035:54:04	0,160	10,0	-79,5	0,134	7,8	8,7	10,0	10,5	8,3	9,4	6,8	12,0	
	02:54:40	035:55:00	0,161	10,0	-89,7	0,135	7,9	8,7	10,0	10,5	8,3	9,4	6,8	12,0	
	02:55:43	035:56:03	0,161	10,0	-84,6	0,136	7,9	8,7	10,0	10,5	8,3	9,4	6,9	12,0	
	02:56:40	035:57:00	0,163	10,0	-77,0	0,136	7,9	8,8	10,0	10,5	8,3	9,4	6,9	12,0	
	02:57:42	035:58:02	0,163	10,0	-88,8	0,137	7,9	8,8	10,0	10,5	8,3	9,4	6,9	12,0	
	02:58:42	035:59:02	0,157	10,0	-82,4	0,130	8,0	8,8	10,0	10,5	8,3	9,4	7,0	12,0	
	02:59:43	036:00:03	0,155	10,0	-91,7	0,130	8,0	8,8	10,0	10,5	8,3	9,4	7,0	12,0	
	03:00:44	036:01:04	0,155	10,0	-88,5	0,130	8,0	8,8	10,0	10,5	8,3	9,4	7,0	12,0	
	03:01:41	036:02:01	0,156	10,0	-79,3	0,131	8,0	8,8	10,0	10,5	8,3	9,4	7,0	12,0	
	03:02:43	036:03:03	0,157	10,0	-89,3	0,132	8,1	8,8	10,0	10,5	8,3	9,4	7,0	12,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 116 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:03:44	036:04:04	0,158	10,0	-84,6	0,132	8,1	8,8	10,0	10,5	8,3	9,4	7,1	12,0	
	03:04:44	036:05:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:04:44	036:05:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:04:44	036:05:04	0,159	10,0	-75,6	0,133	8,2	8,8	10,0	10,5	8,3	9,4	7,1	12,0	
	03:05:44	036:06:04	0,160	10,0	-87,8	0,134	8,2	8,8	10,0	10,5	8,3	9,4	7,1	12,0	
	03:06:41	036:07:01	0,160	10,0	-83,6	0,135	8,2	8,8	10,0	10,5	8,3	9,4	7,1	12,0	
	03:07:42	036:08:02	0,161	10,0	-89,0	0,136	8,3	8,8	10,0	10,5	8,4	9,4	7,1	12,0	
	03:08:42	036:09:02	0,161	10,0	-89,0	0,137	8,3	8,8	10,0	10,5	8,4	9,4	7,1	12,0	
	03:09:40	036:10:00	0,162	10,0	-80,3	0,137	8,4	8,8	10,0	10,5	8,4	9,4	7,2	12,0	
	03:10:40	036:11:00	0,154	10,0	-89,3	0,130	8,4	8,8	10,0	10,5	8,4	9,4	7,2	12,0	
	03:11:43	036:12:03	0,154	10,0	-84,7	0,130	8,5	8,8	10,0	10,5	8,4	9,4	7,2	12,0	
	03:12:43	036:13:03	0,154	10,0	-80,4	0,131	8,5	8,8	10,0	10,5	8,4	9,4	7,2	12,0	
	03:13:40	036:14:00	0,155	10,0	-89,2	0,131	8,5	8,8	10,0	10,5	8,4	9,4	7,2	12,0	
	03:14:41	036:15:01	0,156	10,0	-82,3	0,132	8,6	8,8	10,0	10,6	8,4	9,4	7,3	12,0	
	03:15:41	036:16:01	0,157	10,0	-90,6	0,133	8,6	8,8	10,0	10,5	8,4	9,4	7,3	12,0	
	03:16:45	036:17:05	0,158	10,0	-86,2	0,133	8,6	8,8	10,0	10,5	8,4	9,4	7,3	12,0	
	03:17:42	036:18:02	0,159	10,0	-78,2	0,134	8,6	8,8	10,0	10,5	8,4	9,4	7,4	12,0	
	03:18:44	036:19:04	0,160	10,0	-88,4	0,135	8,7	8,8	10,0	10,5	8,4	9,4	7,4	12,1	
	03:19:43	036:20:03	0,160	10,0	-84,4	0,136	8,7	8,8	10,0	10,6	8,5	9,4	7,4	12,0	
	03:20:44	036:21:04	0,161	10,0	-88,7	0,137	8,7	8,8	10,0	10,6	8,5	9,4	7,4	12,0	
	03:21:40	036:22:00	0,162	10,0	-89,4	0,137	8,7	8,8	10,0	10,6	8,5	9,4	7,4	12,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 117 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	03:22:40	036:23:00	0,153	10,0	-81,9	0,130	8,7	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:23:40	036:24:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:23:40	036:24:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:23:40	036:24:00	0,153	10,0	-88,9	0,130	8,7	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:24:42	036:25:02	0,153	10,0	-85,1	0,131	8,7	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:25:40	036:26:00	0,154	10,0	-76,8	0,131	8,7	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:26:42	036:27:02	0,155	10,0	-88,8	0,132	8,8	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:27:43	036:28:03	0,156	10,0	-80,7	0,133	8,8	8,8	10,0	10,5	8,5	9,4	7,5	12,0	
	03:28:43	036:29:03	0,157	10,0	-90,5	0,134	8,8	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:29:44	036:30:04	0,158	10,0	-86,7	0,134	8,8	8,8	10,0	10,6	8,5	9,4	7,5	12,0	
	03:30:40	036:31:00	0,159	10,0	-78,7	0,135	8,8	8,8	10,0	10,6	8,5	9,4	7,6	12,0	
	03:31:40	036:32:00	0,160	10,0	-88,5	0,136	8,8	8,8	10,0	10,6	8,5	9,4	7,6	12,0	
	03:32:44	036:33:04	0,160	10,0	-83,8	0,137	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:33:40	036:34:00	0,161	10,0	-76,9	0,137	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:34:42	036:35:02	0,152	10,0	-88,8	0,130	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:35:42	036:36:02	0,152	10,0	-81,5	0,130	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:36:44	036:37:04	0,153	10,0	-79,2	0,131	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:37:40	036:38:00	0,153	10,0	-88,3	0,131	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:38:43	036:39:03	0,154	10,0	-78,0	0,132	8,8	8,8	10,0	10,6	8,6	9,4	7,6	12,0	
	03:39:43	036:40:03	0,155	10,0	-88,8	0,133	8,8	8,9	10,0	10,6	8,6	9,4	7,6	12,0	
	03:40:43	036:41:03	0,156	10,0	-84,4	0,134	8,8	8,9	10,0	10,6	8,6	9,4	7,7	12,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 118 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:41:44	036:42:04	0,157	10,0	-84,3	0,135	8,9	8,9	10,1	10,6	8,6	9,4	7,7	12,0	
	03:42:42	036:43:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:42:42	036:43:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:42:42	036:43:02	0,158	10,0	-89,0	0,135	8,9	8,9	10,0	10,6	8,6	9,4	7,7	12,0	
	03:43:42	036:44:02	0,159	10,0	-81,7	0,136	8,9	8,9	10,0	10,6	8,6	9,4	7,7	12,1	
	03:44:44	036:45:04	0,159	10,0	-90,1	0,137	8,9	8,9	10,0	10,6	8,6	9,4	7,7	12,1	
	03:45:41	036:46:01	0,160	10,0	-88,1	0,137	8,9	8,9	10,0	10,6	8,7	9,4	7,7	12,1	
	03:46:42	036:47:02	0,151	10,0	-77,9	0,130	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,1	
	03:47:44	036:48:04	0,151	10,0	-88,0	0,130	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,1	
	03:48:41	036:49:01	0,151	10,0	-84,7	0,131	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,0	
	03:49:44	036:50:04	0,152	10,0	-90,5	0,131	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,0	
	03:50:40	036:51:00	0,153	10,0	-87,4	0,132	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,0	
	03:51:40	036:52:00	0,154	10,0	-78,4	0,133	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,1	
	03:52:43	036:53:03	0,154	10,0	-88,4	0,134	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,1	
	03:53:44	036:54:04	0,156	10,0	-82,5	0,134	8,9	8,9	10,0	10,6	8,7	9,4	7,8	12,0	
	03:54:40	036:55:00	0,157	10,0	-75,7	0,135	8,9	8,9	10,0	10,6	8,7	9,4	7,9	12,0	
	03:55:41	036:56:01	0,158	10,0	-88,6	0,136	8,9	8,9	10,0	10,6	8,7	9,4	7,9	12,0	
	03:56:42	036:57:02	0,158	10,0	-79,6	0,137	8,9	8,9	10,1	10,6	8,7	9,4	7,9	12,0	
	03:57:42	036:58:02	0,160	10,0	-89,5	0,138	8,9	8,9	10,1	10,6	8,8	9,4	7,9	12,0	
	03:58:44	036:59:04	0,150	10,0	-84,8	0,130	8,9	8,9	10,0	10,6	8,7	9,4	7,9	12,0	
	03:59:40	037:00:00	0,150	10,0	-77,1	0,130	8,9	8,9	10,0	10,6	8,8	9,4	7,9	12,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 119 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	04:00:41	037:01:01	0,151	10,0	-88,0	0,131	8,9	8,9	10,1	10,6	8,8	9,4	8,0	12,1	
	04:01:41	037:02:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:01:41	037:02:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:01:41	037:02:01	0,151	10,0	-82,8	0,131	8,9	8,9	10,1	10,6	8,8	9,4	8,0	12,0	
	04:02:43	037:03:03	0,152	10,0	-91,2	0,132	8,9	8,9	10,0	10,6	8,8	9,4	8,0	12,0	
	04:03:43	037:04:03	0,153	10,0	-88,6	0,133	8,9	8,9	10,0	10,6	8,8	9,4	8,0	12,0	
	04:04:40	037:05:00	0,154	10,0	-79,9	0,134	8,9	8,9	10,1	10,6	8,8	9,4	8,0	12,0	
	04:05:40	037:06:00	0,155	10,0	-90,2	0,135	8,9	8,9	10,1	10,6	8,8	9,4	8,0	12,0	
	04:06:42	037:07:02	0,156	10,0	-86,4	0,135	8,9	8,9	10,1	10,6	8,8	9,4	8,1	12,1	
	04:07:44	037:08:04	0,157	10,0	-77,7	0,136	8,9	8,9	10,1	10,6	8,8	9,4	8,1	12,1	
	04:08:41	037:09:01	0,158	10,0	-88,3	0,137	8,9	8,9	10,1	10,6	8,8	9,4	8,1	12,1	
	04:09:42	037:10:02	0,157	10,0	-83,2	0,137	8,9	8,9	10,1	10,6	8,8	9,4	8,1	12,0	
	04:10:44	037:11:04	0,149	10,0	-90,7	0,130	8,9	8,9	10,1	10,6	8,8	9,4	8,1	12,1	
	04:11:41	037:12:01	0,149	10,0	-88,0	0,130	8,9	8,9	10,0	10,6	8,8	9,4	8,1	12,1	
	04:12:42	037:13:02	0,149	10,0	-78,9	0,131	8,9	8,9	10,1	10,6	8,8	9,4	8,2	12,0	
	04:13:43	037:14:03	0,150	10,0	-88,6	0,132	9,0	8,9	10,1	10,6	8,8	9,4	8,2	12,0	
	04:14:43	037:15:03	0,151	10,0	-83,1	0,132	8,9	8,9	10,1	10,6	8,8	9,4	8,2	12,1	
	04:15:43	037:16:03	0,152	10,0	-77,8	0,133	8,9	8,9	10,1	10,6	8,9	9,4	8,2	12,0	
	04:16:44	037:17:04	0,153	10,0	-88,6	0,134	8,9	8,9	10,1	10,6	8,9	9,4	8,2	12,1	
	04:17:41	037:18:01	0,154	10,0	-79,3	0,134	9,0	8,9	10,1	10,6	8,9	9,4	8,2	12,1	
	04:18:42	037:19:02	0,155	10,0	-89,2	0,135	9,0	8,9	10,1	10,6	8,9	9,4	8,2	12,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 120 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	04:19:43	037:20:03	0,155	10,0	-85,2	0,136	9,0	8,9	10,1	10,6	8,9	9,4	8,3	12,1	
	04:20:43	037:21:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:20:43	037:21:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:20:43	037:21:03	0,156	10,0	-77,4	0,137	8,9	8,9	10,1	10,6	8,9	9,4	8,3	12,1	
	04:21:44	037:22:04	0,157	10,0	-89,0	0,137	9,0	9,0	10,1	10,6	8,9	9,4	8,3	12,1	
	04:22:44	037:23:04	0,148	10,0	-80,9	0,130	9,0	9,0	10,1	10,6	8,9	9,4	8,3	12,1	
	04:23:44	037:24:04	0,148	10,0	-90,4	0,130	9,0	9,0	10,1	10,6	8,9	9,4	8,3	12,1	
	04:24:40	037:25:00	0,148	10,0	-88,6	0,131	9,0	9,0	10,1	10,6	8,9	9,4	8,3	12,1	
	04:25:41	037:26:01	0,149	10,0	-78,8	0,131	9,0	9,0	10,1	10,6	8,9	9,4	8,4	12,1	
	04:26:44	037:27:04	0,150	10,0	-89,1	0,132	9,0	9,0	10,1	10,6	9,0	9,4	8,4	12,1	
	04:27:40	037:28:00	0,150	10,0	-85,1	0,133	9,0	9,0	10,1	10,6	8,9	9,4	8,4	12,1	
	04:28:41	037:29:01	0,151	10,0	-78,1	0,134	9,0	9,0	10,1	10,6	9,0	9,4	8,4	12,1	
	04:29:42	037:30:02	0,152	10,0	-89,2	0,134	9,0	9,0	10,1	10,6	9,0	9,4	8,4	12,1	
	04:30:41	037:31:01	0,153	10,0	-82,1	0,135	9,0	9,0	10,1	10,6	9,0	9,4	8,4	12,1	
	04:31:41	037:32:01	0,154	10,0	-87,1	0,136	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:32:41	037:33:01	0,155	10,0	-88,2	0,136	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:33:44	037:34:04	0,156	10,0	-77,9	0,137	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:34:41	037:35:01	0,147	10,0	-88,9	0,130	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:35:42	037:36:02	0,147	10,0	-85,1	0,130	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:36:42	037:37:02	0,147	10,0	-88,0	0,131	9,0	9,0	10,1	10,6	9,0	9,4	8,5	12,1	
	04:37:42	037:38:02	0,148	10,0	-89,2	0,131	9,0	9,0	10,1	10,6	9,1	9,4	8,5	12,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 121 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	04:38:42	037:39:02	0,148	10,0	-81,5	0,132	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:39:42	037:40:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:39:42	037:40:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:39:42	037:40:02	0,149	10,0	-90,1	0,133	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:40:42	037:41:02	0,150	10,0	-86,7	0,133	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:41:41	037:42:01	0,151	10,0	-77,9	0,134	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:42:41	037:43:01	0,152	10,0	-87,9	0,135	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:43:43	037:44:03	0,152	10,0	-83,6	0,136	9,0	9,0	10,1	10,6	9,1	9,4	8,6	12,1	
	04:44:44	037:45:04	0,153	10,0	-89,4	0,136	9,0	9,0	10,1	10,6	9,1	9,4	8,7	12,1	
	04:45:40	037:46:00	0,154	10,0	-89,4	0,137	9,0	9,0	10,1	10,6	9,1	9,4	8,7	12,1	
	04:46:43	037:47:03	0,155	10,0	-83,4	0,138	9,0	9,0	10,1	10,6	9,1	9,4	8,7	12,1	
	04:47:40	037:48:00	0,145	10,0	-89,8	0,130	9,0	9,0	10,1	10,6	9,1	9,4	8,7	12,1	
	04:48:41	037:49:01	0,145	10,0	-83,0	0,130	9,0	9,0	10,1	10,6	9,1	9,4	8,7	12,1	
	04:49:41	037:50:01	0,146	10,0	-89,9	0,131	9,0	9,0	10,1	10,7	9,1	9,4	8,7	12,1	
	04:50:42	037:51:02	0,147	10,0	-85,1	0,132	9,0	9,0	10,1	10,7	9,1	9,4	8,7	12,1	
	04:51:42	037:52:02	0,148	10,0	-77,2	0,133	9,0	9,0	10,1	10,7	9,1	9,4	8,7	12,1	
	04:52:44	037:53:04	0,149	10,0	-89,2	0,133	9,0	9,0	10,1	10,6	9,2	9,4	8,8	12,1	
	04:53:40	037:54:00	0,149	10,0	-82,6	0,134	9,0	9,0	10,1	10,6	9,2	9,4	8,8	12,1	
	04:54:41	037:55:01	0,150	10,0	-89,6	0,135	9,0	9,0	10,1	10,6	9,1	9,4	8,8	12,1	
	04:55:42	037:56:02	0,151	10,0	-88,2	0,135	9,0	9,0	10,1	10,6	9,1	9,4	8,8	12,1	
	04:56:42	037:57:02	0,152	10,0	-78,8	0,136	9,0	9,0	10,1	10,7	9,2	9,4	8,8	12,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	122 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	04:57:43	037:58:03	0,153	10,0	-88,8	0,137	9,0	9,0	10,1	10,7	9,2	9,4	8,8	12,1	
	04:58:40	037:59:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:58:40	037:59:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:58:40	037:59:00	0,154	10,0	-84,5	0,138	9,0	9,0	10,1	10,7	9,2	9,4	8,8	12,1	
	04:59:40	038:00:00	0,145	10,0	-77,0	0,130	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:00:40	038:01:00	0,145	10,0	-89,1	0,130	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:01:42	038:02:02	0,145	10,0	-82,5	0,131	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:02:44	038:03:04	0,146	10,0	-90,5	0,132	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:03:44	038:04:04	0,147	10,0	-87,4	0,133	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:04:41	038:05:01	0,148	10,0	-78,7	0,133	9,0	9,0	10,1	10,7	9,2	9,4	8,9	12,1	
	05:05:45	038:06:05	0,149	10,0	-88,3	0,134	9,0	9,0	10,1	10,7	9,3	9,4	8,9	12,1	
	05:06:40	038:07:00	0,149	10,0	-86,1	0,135	9,0	9,0	10,1	10,7	9,3	9,4	8,9	12,1	
	05:07:40	038:08:00	0,150	10,0	-91,1	0,136	9,0	9,0	10,1	10,7	9,3	9,4	9,0	12,1	
	05:08:43	038:09:03	0,151	10,0	-87,3	0,136	9,0	9,0	10,1	10,7	9,3	9,4	9,0	12,1	
	05:09:41	038:10:01	0,152	10,0	-78,6	0,137	9,0	9,0	10,1	10,7	9,3	9,4	9,0	12,1	
	05:10:42	038:11:02	0,153	10,0	-88,1	0,138	9,1	9,1	10,1	10,7	9,3	9,4	9,0	12,1	
	05:11:44	038:12:04	0,144	10,0	-83,8	0,130	9,1	9,1	10,1	10,7	9,3	9,4	9,0	12,1	
	05:12:40	038:13:00	0,144	10,0	-87,9	0,131	9,0	9,1	10,1	10,7	9,3	9,4	9,0	12,1	
	05:13:42	038:14:02	0,145	10,0	-89,1	0,131	9,1	9,0	10,1	10,7	9,3	9,4	9,0	12,1	
	05:14:43	038:15:03	0,145	10,0	-80,4	0,132	9,1	9,1	10,1	10,7	9,3	9,4	9,0	12,1	
	05:15:40	038:16:00	0,146	10,0	-90,6	0,133	9,1	9,1	10,1	10,7	9,3	9,4	9,1	12,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	123 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:16:41	038:17:01	0,146	10,0	-86,4	0,133	9,1	9,1	10,1	10,7	9,3	9,4	9,1	12,1	
	05:17:41	038:18:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:17:41	038:18:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:17:41	038:18:01	0,147	10,0	-77,3	0,134	9,1	9,1	10,1	10,7	9,3	9,4	9,1	12,1	
	05:18:45	038:19:05	0,149	10,0	-88,9	0,135	9,1	9,1	10,1	10,7	9,3	9,4	9,1	12,1	
	05:19:43	038:20:03	0,149	10,0	-82,9	0,135	9,1	9,1	10,1	10,7	9,4	9,4	9,1	12,1	
	05:20:43	038:21:03	0,150	10,0	-83,0	0,136	9,1	9,1	10,1	10,7	9,4	9,4	9,1	12,1	
	05:21:40	038:22:00	0,150	10,0	-88,9	0,137	9,1	9,1	10,1	10,7	9,4	9,4	9,1	12,1	
	05:22:43	038:23:03	0,151	10,0	-79,7	0,138	9,1	9,1	10,1	10,7	9,4	9,4	9,2	12,1	
	05:23:43	038:24:03	0,143	10,0	-89,6	0,130	9,1	9,1	10,1	10,7	9,4	9,4	9,2	12,1	
	05:24:44	038:25:04	0,143	10,0	-85,1	0,130	9,1	9,1	10,1	10,7	9,4	9,4	9,2	12,1	
	05:25:41	038:26:01	0,143	10,0	-80,2	0,131	9,0	9,1	10,1	10,7	9,4	9,4	9,2	12,1	
	05:26:42	038:27:02	0,144	10,0	-89,2	0,132	9,1	9,1	10,1	10,7	9,4	9,4	9,2	12,1	
	05:27:43	038:28:03	0,145	10,0	-81,6	0,132	9,1	9,1	10,1	10,7	9,4	9,4	9,3	12,1	
	05:28:44	038:29:04	0,145	10,0	-89,8	0,133	9,1	9,1	10,1	10,7	9,4	9,4	9,3	12,1	
	05:29:45	038:30:05	0,146	10,0	-86,9	0,134	9,0	9,1	10,1	10,7	9,4	9,4	9,3	12,1	
	05:30:43	038:31:03	0,147	10,0	-78,3	0,135	9,0	9,1	10,1	10,7	9,4	9,4	9,3	12,1	
	05:31:43	038:32:03	0,148	10,0	-88,5	0,135	9,0	9,1	10,1	10,7	9,4	9,4	9,3	12,1	
	05:32:43	038:33:03	0,148	10,0	-84,9	0,136	9,1	9,1	10,1	10,7	9,4	9,4	9,4	12,1	
	05:33:43	038:34:03	0,149	10,0	-92,4	0,137	9,1	9,1	10,1	10,7	9,5	9,4	9,4	12,1	
	05:34:44	038:35:04	0,150	10,0	-89,2	0,137	9,0	9,1	10,1	10,7	9,5	9,4	9,4	12,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	124 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:35:44	038:36:04	0,142	10,0	-80,7	0,130	9,1	9,1	10,1	10,7	9,5	9,4	9,4	12,1	
	05:36:40	038:37:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:36:40	038:37:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:36:40	038:37:00	0,142	10,0	-90,1	0,130	9,0	9,1	10,1	10,7	9,5	9,4	9,4	12,1	
	05:37:43	038:38:03	0,142	10,0	-85,1	0,131	9,1	9,1	10,1	10,7	9,5	9,4	9,4	12,1	
	05:38:40	038:39:00	0,143	10,0	-77,8	0,132	9,1	9,1	10,1	10,7	9,5	9,4	9,5	12,1	
	05:39:41	038:40:01	0,144	10,0	-87,9	0,132	9,1	9,1	10,1	10,7	9,5	9,4	9,5	12,1	
	05:40:40	038:41:00	0,144	10,0	-83,9	0,133	9,1	9,1	10,1	10,7	9,5	9,4	9,5	12,1	
	05:41:41	038:42:01	0,145	10,0	-90,0	0,134	9,1	9,1	10,1	10,7	9,5	9,4	9,5	12,1	
	05:42:42	038:43:02	0,146	10,0	-89,1	0,134	9,1	9,1	10,1	10,7	9,5	9,4	9,5	12,1	
	05:43:42	038:44:02	0,146	10,0	-80,8	0,135	9,1	9,1	10,1	10,7	9,5	9,4	9,6	12,1	
	05:44:44	038:45:04	0,147	10,0	-89,6	0,136	9,1	9,1	10,1	10,7	9,5	9,4	9,6	12,1	
	05:45:42	038:46:02	0,147	10,0	-85,2	0,136	9,1	9,1	10,1	10,7	9,5	9,4	9,6	12,1	
	05:46:42	038:47:02	0,148	10,0	-77,3	0,137	9,1	9,1	10,1	10,7	9,5	9,4	9,6	12,1	
	05:47:45	038:48:05	0,141	10,0	-88,8	0,130	9,1	9,1	10,1	10,7	9,6	9,4	9,6	12,1	
	05:48:40	038:49:00	0,141	10,0	-83,0	0,130	9,1	9,1	10,1	10,7	9,6	9,4	9,7	12,1	
	05:49:42	038:50:02	0,141	10,0	-90,9	0,131	9,1	9,1	10,1	10,7	9,6	9,4	9,7	12,1	
	05:50:43	038:51:03	0,142	10,0	-86,5	0,131	9,1	9,1	10,1	10,7	9,6	9,4	9,7	12,1	
	05:51:43	038:52:03	0,142	10,0	-78,2	0,132	9,1	9,1	10,1	10,7	9,6	9,5	9,7	12,1	
	05:52:40	038:53:00	0,143	10,0	-89,1	0,133	9,1	9,1	10,1	10,7	9,6	9,5	9,8	12,1	
	05:53:41	038:54:01	0,143	10,0	-82,0	0,133	9,1	9,1	10,2	10,7	9,6	9,4	9,8	12,1	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	125 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	05:54:44	038:55:04	0,144	10,0	-90,7	0,134	9,1	9,1	10,1	10,7	9,6	9,4	9,8	12,1	
	05:55:44	038:56:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:55:44	038:56:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:55:44	038:56:04	0,145	10,0	-86,5	0,135	9,1	9,1	10,1	10,7	9,6	9,4	9,8	12,1	
	05:56:44	038:57:04	0,145	10,0	-77,3	0,136	9,1	9,1	10,1	10,7	9,6	9,4	9,9	12,1	
	05:57:42	038:58:02	0,146	10,0	-88,4	0,136	9,1	9,1	10,1	10,7	9,6	9,4	9,9	12,1	
	05:58:44	038:59:04	0,146	10,0	-83,9	0,137	9,1	9,1	10,1	10,7	9,7	9,4	9,9	12,1	
	05:59:40	039:00:00	0,147	10,0	-78,9	0,138	9,1	9,1	10,1	10,7	9,7	9,4	9,9	12,1	
	06:00:41	039:01:01	0,140	10,0	-89,1	0,130	9,1	9,1	10,2	10,7	9,7	9,4	10,0	12,1	
	06:01:43	039:02:03	0,140	10,0	-81,5	0,130	9,1	9,1	10,2	10,7	9,7	9,4	10,0	12,1	
	06:02:40	039:03:00	0,141	10,0	-90,4	0,131	9,1	9,1	10,1	10,7	9,7	9,4	10,0	12,1	
	06:03:42	039:04:02	0,141	10,0	-86,9	0,132	9,1	9,1	10,1	10,7	9,7	9,4	10,0	12,1	
	06:04:42	039:05:02	0,142	10,0	-78,2	0,132	9,1	9,1	10,2	10,7	9,7	9,4	10,1	12,1	
	06:05:44	039:06:04	0,142	10,0	-89,4	0,133	9,1	9,1	10,2	10,7	9,7	9,4	10,1	12,1	
	06:06:40	039:07:00	0,143	10,0	-83,8	0,134	9,1	9,1	10,2	10,7	9,7	9,4	10,1	12,1	
	06:07:42	039:08:02	0,144	10,0	-90,5	0,135	9,1	9,1	10,1	10,7	9,7	9,4	10,1	12,1	
	06:08:40	039:09:00	0,144	10,0	-86,7	0,135	9,1	9,1	10,1	10,7	9,8	9,4	10,2	12,2	
	06:09:44	039:10:04	0,145	10,0	-77,4	0,136	9,1	9,1	10,2	10,7	9,8	9,4	10,2	12,2	
	06:10:40	039:11:00	0,145	10,0	-88,4	0,137	9,1	9,1	10,1	10,7	9,8	9,4	10,2	12,1	
	06:11:40	039:12:00	0,146	10,0	-82,6	0,137	9,1	9,1	10,2	10,7	9,8	9,4	10,3	12,2	
	06:12:41	039:13:01	0,139	10,0	-82,7	0,130	9,1	9,1	10,2	10,7	9,8	9,4	10,3	12,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 126 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	06:13:45	039:14:05	0,139	10,0	-88,3	0,130	9,1	9,1	10,2	10,7	9,8	9,4	10,3	12,2	
	06:14:42	039:15:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:14:42	039:15:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:14:42	039:15:02	0,139	10,0	-78,7	0,131	9,1	9,1	10,2	10,7	9,9	9,4	10,3	12,2	
	06:15:44	039:16:04	0,140	10,0	-89,5	0,132	9,1	9,1	10,2	10,7	9,9	9,4	10,4	12,2	
	06:16:44	039:17:04	0,141	10,0	-85,0	0,132	9,1	9,1	10,2	10,7	9,9	9,4	10,4	12,2	
	06:17:41	039:18:01	0,141	10,0	-77,4	0,133	9,1	9,1	10,2	10,7	9,9	9,4	10,4	12,2	
	06:18:41	039:19:01	0,142	10,0	-88,4	0,134	9,1	9,1	10,2	10,7	9,9	9,4	10,4	12,2	
	06:19:41	039:20:01	0,143	10,0	-82,5	0,134	9,1	9,1	10,2	10,7	9,9	9,4	10,4	12,2	
	06:20:42	039:21:02	0,143	10,0	-90,6	0,135	9,1	9,1	10,2	10,7	9,9	9,4	10,5	12,2	
	06:21:43	039:22:03	0,143	10,0	-86,5	0,135	9,1	9,1	10,2	10,7	9,9	9,4	10,5	12,2	
	06:22:40	039:23:00	0,144	10,0	-78,7	0,136	9,1	9,1	10,2	10,7	9,9	9,4	10,5	12,2	
	06:23:41	039:24:01	0,145	10,0	-88,1	0,137	9,1	9,1	10,2	10,7	9,9	9,4	10,5	12,2	
	06:24:43	039:25:03	0,145	10,0	-83,0	0,138	9,1	9,1	10,2	10,7	9,9	9,4	10,5	12,2	
	06:25:43	039:26:03	0,138	10,0	-77,0	0,130	9,1	9,1	10,2	10,7	9,9	9,4	10,6	12,2	
	06:26:43	039:27:03	0,138	10,0	-89,1	0,131	9,1	9,1	10,2	10,7	10,0	9,4	10,6	12,2	
	06:27:42	039:28:02	0,139	10,0	-80,2	0,131	9,1	9,1	10,2	10,7	10,0	9,4	10,6	12,2	
	06:28:43	039:29:03	0,139	10,0	-89,1	0,132	9,1	9,1	10,2	10,7	10,0	9,4	10,6	12,2	
	06:29:41	039:30:01	0,140	10,0	-84,9	0,132	9,1	9,1	10,2	10,7	10,0	9,4	10,7	12,2	
	06:30:43	039:31:03	0,141	10,0	-77,6	0,133	9,1	9,1	10,2	10,7	10,0	9,4	10,7	12,2	
	06:31:43	039:32:03	0,141	10,0	-89,1	0,134	9,1	9,1	10,2	10,7	10,0	9,4	10,7	12,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 127 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:32:40	039:33:00	0,142	10,0	-82,2	0,135	9,1	9,1	10,2	10,7	10,0	9,4	10,7	12,2	
	06:33:42	039:34:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:33:42	039:34:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:33:42	039:34:02	0,142	10,0	-81,3	0,135	9,1	9,1	10,2	10,7	10,0	9,4	10,8	12,2	
	06:34:44	039:35:04	0,143	10,0	-88,7	0,136	9,1	9,2	10,2	10,8	10,1	9,4	10,8	12,2	
	06:35:40	039:36:00	0,143	10,0	-79,6	0,137	9,1	9,2	10,2	10,8	10,1	9,4	10,9	12,2	
	06:36:41	039:37:01	0,144	10,0	-88,9	0,137	9,1	9,2	10,2	10,8	10,1	9,4	10,9	12,2	
	06:37:43	039:38:03	0,137	10,0	-84,9	0,130	9,1	9,2	10,2	10,8	10,1	9,4	11,0	12,2	
	06:38:45	039:39:05	0,137	10,0	-76,9	0,131	9,1	9,2	10,2	10,8	10,1	9,4	11,0	12,2	
	06:39:40	039:40:00	0,138	10,0	-89,1	0,131	9,1	9,2	10,2	10,8	10,1	9,4	11,0	12,2	
	06:40:41	039:41:01	0,138	10,0	-80,8	0,132	9,1	9,2	10,2	10,8	10,1	9,4	11,1	12,2	
	06:41:42	039:42:02	0,139	10,0	-89,7	0,132	9,1	9,2	10,2	10,8	10,1	9,4	11,1	12,2	
	06:42:45	039:43:05	0,140	10,0	-86,1	0,133	9,1	9,1	10,2	10,8	10,1	9,4	11,1	12,2	
	06:43:42	039:44:02	0,140	10,0	-80,4	0,134	9,1	9,1	10,2	10,8	10,1	9,4	11,1	12,2	
	06:44:40	039:45:00	0,141	10,0	-89,5	0,134	9,1	9,2	10,2	10,8	10,2	9,4	11,2	12,2	
	06:45:40	039:46:00	0,141	10,0	-80,2	0,135	9,1	9,2	10,2	10,8	10,2	9,4	11,2	12,2	
	06:46:40	039:47:00	0,142	10,0	-89,0	0,136	9,1	9,2	10,2	10,8	10,2	9,4	11,2	12,2	
	06:47:42	039:48:02	0,143	10,0	-85,3	0,137	9,1	9,1	10,2	10,8	10,2	9,4	11,2	12,2	
	06:48:42	039:49:02	0,143	10,0	-78,0	0,137	9,1	9,1	10,2	10,8	10,2	9,4	11,2	12,2	
	06:49:40	039:50:00	0,136	10,0	-89,3	0,130	9,1	9,1	10,2	10,8	10,2	9,4	11,3	12,2	
	06:50:40	039:51:00	0,136	10,0	-81,1	0,130	9,1	9,2	10,2	10,8	10,2	9,4	11,3	12,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 128 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:51:40	039:52:00	0,137	10,0	-89,9	0,131	9,1	9,2	10,2	10,8	10,2	9,4	11,3	12,2	
	06:52:41	039:53:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:52:41	039:53:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:52:41	039:53:01	0,137	10,0	-85,3	0,132	9,1	9,2	10,2	10,8	10,2	9,4	11,3	12,2	
	06:53:42	039:54:02	0,138	10,0	-81,1	0,132	9,1	9,2	10,2	10,8	10,2	9,4	11,3	12,2	
	06:54:40	039:55:00	0,139	10,0	-89,3	0,133	9,1	9,2	10,2	10,8	10,2	9,4	11,3	12,2	
	06:55:42	039:56:02	0,139	10,0	-81,8	0,134	9,1	9,1	10,2	10,8	10,3	9,4	11,3	12,2	
	06:56:45	039:57:05	0,140	10,0	-82,6	0,135	9,1	9,2	10,2	10,8	10,3	9,4	11,3	12,2	
	06:57:41	039:58:01	0,141	10,0	-88,7	0,135	9,1	9,2	10,2	10,8	10,3	9,4	11,3	12,2	
	06:58:44	039:59:04	0,141	10,0	-78,7	0,136	9,1	9,2	10,2	10,8	10,3	9,4	11,3	12,2	
	06:59:44	040:00:04	0,142	10,0	-89,5	0,137	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:00:44	040:01:04	0,143	10,0	-85,6	0,137	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:01:42	040:02:02	0,136	10,0	-79,0	0,130	9,1	9,2	10,2	10,8	10,3	9,5	11,4	12,2	
	07:02:43	040:03:03	0,136	10,0	-89,3	0,131	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:03:40	040:04:00	0,136	10,0	-82,3	0,131	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:04:43	040:05:03	0,137	10,0	-75,8	0,132	9,1	9,1	10,2	10,8	10,3	9,5	11,4	12,2	
	07:05:42	040:06:02	0,138	10,0	-88,6	0,133	9,1	9,1	10,2	10,8	10,3	9,5	11,4	12,2	
	07:06:44	040:07:04	0,138	10,0	-78,8	0,133	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:07:40	040:08:00	0,139	10,0	-89,6	0,134	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:08:42	040:09:02	0,140	10,0	-84,5	0,135	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:09:41	040:10:01	0,140	10,0	-77,4	0,135	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	129 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:10:43	040:11:03	0,141	10,0	-89,2	0,136	9,1	9,2	10,2	10,8	10,3	9,4	11,4	12,2	
	07:11:44	040:12:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:11:44	040:12:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:11:44	040:12:04	0,142	10,0	-81,8	0,137	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:12:40	040:13:00	0,142	10,0	-79,3	0,138	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:13:42	040:14:02	0,135	10,0	-88,2	0,130	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:14:44	040:15:04	0,136	10,0	-77,9	0,131	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:15:40	040:16:00	0,136	10,0	-89,1	0,132	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:16:41	040:17:01	0,137	10,0	-84,0	0,132	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:17:42	040:18:02	0,137	10,0	-92,4	0,133	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:18:43	040:19:03	0,138	10,0	-89,6	0,133	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:19:43	040:20:03	0,139	10,0	-80,4	0,134	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:20:41	040:21:01	0,139	10,0	-88,6	0,135	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:21:43	040:22:03	0,140	10,0	-83,6	0,135	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:22:45	040:23:05	0,141	10,0	-78,3	0,136	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:23:42	040:24:02	0,141	10,0	-89,1	0,137	9,1	9,2	10,2	10,8	10,4	9,5	11,4	12,2	
	07:24:42	040:25:02	0,142	10,0	-81,9	0,138	9,1	9,2	10,2	10,8	10,4	9,4	11,4	12,2	
	07:25:40	040:26:00	0,135	10,0	-81,7	0,130	9,1	9,2	10,2	10,8	10,4	9,4	11,5	12,2	
	07:26:42	040:27:02	0,135	10,0	-88,8	0,131	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:27:44	040:28:04	0,136	10,0	-79,1	0,131	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:28:40	040:29:00	0,136	10,0	-89,7	0,132	9,1	9,2	10,2	10,8	10,4	9,4	11,5	12,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 130 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:29:42	040:30:02	0,137	10,0	-85,2	0,133	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:30:43	040:31:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:30:43	040:31:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:30:43	040:31:03	0,138	10,0	-86,0	0,134	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:31:40	040:32:00	0,138	10,0	-89,6	0,134	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:32:40	040:33:00	0,139	10,0	-81,9	0,135	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:33:41	040:34:01	0,140	10,0	-89,6	0,136	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:34:42	040:35:02	0,140	10,0	-84,5	0,136	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:35:40	040:36:00	0,141	10,0	-77,3	0,137	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:36:41	040:37:01	0,141	10,0	-89,3	0,138	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:37:43	040:38:03	0,134	10,0	-81,7	0,130	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:38:40	040:39:00	0,135	10,0	-90,3	0,131	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:39:41	040:40:01	0,135	10,0	-85,5	0,131	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,3	
	07:40:44	040:41:04	0,136	10,0	-77,5	0,132	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:41:40	040:42:00	0,136	10,0	-87,8	0,133	9,1	9,2	10,2	10,8	10,5	9,4	11,5	12,2	
	07:42:41	040:43:01	0,137	10,0	-83,2	0,133	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:43:41	040:44:01	0,138	10,0	-82,0	0,134	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:44:43	040:45:03	0,138	10,0	-88,8	0,135	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:45:41	040:46:01	0,139	10,0	-80,2	0,136	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:46:42	040:47:02	0,140	10,0	-90,3	0,136	9,1	9,2	10,2	10,8	10,5	9,5	11,5	12,2	
	07:47:44	040:48:04	0,140	10,0	-85,5	0,137	9,1	9,2	10,2	10,8	10,6	9,4	11,5	12,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 131 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:48:43	040:49:03	0,141	10,0	-76,7	0,138	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,3	
	07:49:44	040:50:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:49:44	040:50:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:49:44	040:50:04	0,134	10,0	-88,0	0,130	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:50:44	040:51:04	0,134	10,0	-83,5	0,131	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:51:40	040:52:00	0,135	10,0	-75,0	0,132	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:52:42	040:53:02	0,135	10,0	-88,6	0,132	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,2	
	07:53:43	040:54:03	0,136	10,0	-81,7	0,133	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	
	07:54:44	040:55:04	0,137	10,0	-90,3	0,134	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:55:44	040:56:04	0,137	10,0	-88,8	0,134	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:56:44	040:57:04	0,138	10,0	-79,3	0,135	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:57:41	040:58:01	0,138	10,0	-89,2	0,136	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,3	
	07:58:41	040:59:01	0,139	10,0	-83,9	0,136	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	07:59:41	041:00:01	0,140	10,0	-75,8	0,137	9,1	9,2	10,2	10,8	10,6	9,5	11,5	12,2	
	08:00:42	041:01:02	0,137	10,0	-88,1	0,132	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	
	08:01:42	041:02:02	0,133	10,0	-82,0	0,130	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	
	08:02:44	041:03:04	0,134	10,0	-91,8	0,131	9,1	9,2	10,2	10,9	10,6	9,4	11,5	12,3	
	08:03:44	041:04:04	0,134	10,0	-88,7	0,132	9,1	9,2	10,2	10,9	10,6	9,4	11,5	12,3	
	08:04:45	041:05:05	0,135	10,0	-78,6	0,132	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	
	08:05:41	041:06:01	0,136	10,0	-89,4	0,133	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	
	08:06:41	041:07:01	0,136	10,0	-84,3	0,134	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus			
Document issued by:			Signature:	Page: 132 / 209	

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



08:07:44	041:08:04	0,137	10,0	-78,1	0,134	9,1	9,2	10,2	10,9	10,6	9,5	11,6	12,3
08:08:41	041:09:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:08:41	041:09:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:08:41	041:09:01	0,138	10,0	-88,9	0,135	9,1	9,2	10,2	10,9	10,6	9,5	11,6	12,3
08:09:42	041:10:02	0,138	10,0	-82,3	0,136	9,1	9,2	10,2	10,9	10,6	9,5	11,5	12,3
08:10:44	041:11:04	0,140	10,0	-90,5	0,137	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:11:40	041:12:00	0,139	10,0	-87,8	0,137	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:12:41	041:13:01	0,140	10,0	-79,6	0,138	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:13:43	041:14:03	0,133	10,0	-89,3	0,130	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:14:40	041:15:00	0,133	10,0	-82,6	0,131	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:15:42	041:16:02	0,134	10,0	-76,2	0,132	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:16:44	041:17:04	0,134	10,0	-88,7	0,132	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:17:40	041:18:00	0,136	10,0	-79,7	0,133	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:18:40	041:19:00	0,135	10,0	-89,8	0,134	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:19:42	041:20:02	0,136	10,0	-84,4	0,134	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:20:43	041:21:03	0,137	10,0	-76,5	0,135	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:21:43	041:22:03	0,137	10,0	-88,9	0,136	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:22:43	041:23:03	0,138	10,0	-82,0	0,136	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:23:40	041:24:00	0,139	10,0	-76,7	0,137	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:24:40	041:25:00	0,139	10,0	-88,8	0,138	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:25:41	041:26:01	0,132	10,0	-80,8	0,130	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 133 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



08:26:42	041:27:02	0,132	10,0	-90,7	0,131	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:27:43	041:28:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:27:43	041:28:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:27:43	041:28:03	0,133	10,0	-86,7	0,131	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:28:43	041:29:03	0,134	10,0	-77,6	0,132	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:29:40	041:30:00	0,134	10,0	-88,7	0,133	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:30:41	041:31:01	0,135	10,0	-84,6	0,134	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:31:41	041:32:01	0,136	10,0	-90,7	0,134	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:32:43	041:33:03	0,136	10,0	-89,1	0,135	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:33:44	041:34:04	0,137	10,0	-80,3	0,136	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:34:43	041:35:03	0,137	10,0	-90,1	0,136	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:35:44	041:36:04	0,138	10,0	-85,5	0,137	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:36:40	041:37:00	0,139	10,0	-77,1	0,138	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:37:41	041:38:01	0,131	10,0	-87,9	0,130	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:38:41	041:39:01	0,132	10,0	-83,4	0,131	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:39:42	041:40:02	0,133	10,0	-76,6	0,132	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:40:45	041:41:05	0,133	10,0	-88,9	0,132	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:41:40	041:42:00	0,134	10,0	-81,3	0,133	9,1	9,2	10,2	10,9	10,7	9,5	11,6	12,3
08:42:41	041:43:01	0,135	10,0	-91,1	0,134	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:43:41	041:44:01	0,136	10,0	-87,3	0,134	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:44:43	041:45:03	0,136	10,0	-78,6	0,135	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 134 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



08:45:44	041:46:04	0,136	10,0	-89,2	0,136	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:46:44	041:47:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
08:46:44	041:47:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
08:46:44	041:47:04	0,137	10,0	-82,8	0,136	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:47:42	041:48:02	0,138	10,0	-90,0	0,137	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:48:43	041:49:03	0,138	10,0	-84,9	0,138	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:49:43	041:50:03	0,131	10,0	-77,2	0,130	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:50:44	041:51:04	0,132	10,0	-89,2	0,131	9,1	9,2	10,3	10,9	10,7	9,5	11,6	12,3
08:51:44	041:52:04	0,133	10,0	-82,4	0,132	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:52:44	041:53:04	0,133	10,0	-78,2	0,133	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:53:40	041:54:00	0,134	10,0	-88,9	0,133	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:54:42	041:55:02	0,134	10,0	-79,1	0,134	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:55:44	041:56:04	0,135	10,0	-89,2	0,135	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:56:43	041:57:03	0,135	10,0	-84,1	0,135	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:57:43	041:58:03	0,136	10,0	-75,8	0,136	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:58:43	041:59:03	0,137	10,0	-87,8	0,137	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
08:59:44	042:00:04	0,137	10,0	-83,1	0,137	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
09:00:44	042:01:04	0,130	10,0	-90,1	0,130	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
09:01:40	042:02:00	0,131	10,0	-89,1	0,131	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
09:02:41	042:03:01	0,132	10,0	-80,6	0,131	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3
09:03:42	042:04:02	0,132	10,0	-90,1	0,132	9,2	9,2	10,3	10,9	10,8	9,5	11,6	12,3

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 135 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	09:04:42	042:05:02	0,133	10,0	-85,7	0,133	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:05:43	042:06:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	09:05:43	042:06:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	09:05:43	042:06:03	0,133	10,0	-77,3	0,133	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:06:44	042:07:04	0,134	10,0	-88,7	0,134	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:07:40	042:08:00	0,134	10,0	-83,0	0,135	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:08:42	042:09:02	0,135	10,0	-89,8	0,136	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:09:44	042:10:04	0,136	10,0	-88,7	0,136	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3	
	09:10:41	042:11:01	0,136	10,0	-79,5	0,137	9,1	9,2	10,3	10,9	10,8	9,5	11,6	12,3	
	09:11:42	042:12:02	0,137	10,0	-89,5	0,138	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:12:43	042:13:03	0,130	10,0	-84,4	0,130	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:13:40	042:14:00	0,131	10,0	-76,7	0,131	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:14:41	042:15:01	0,131	10,0	-89,1	0,132	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:15:41	042:16:01	0,132	10,0	-82,3	0,132	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:16:41	042:17:01	0,133	10,0	-89,5	0,133	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:17:41	042:18:01	0,133	10,0	-88,2	0,133	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:18:43	042:19:03	0,134	10,0	-78,2	0,134	9,1	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:19:40	042:20:00	0,135	10,0	-88,3	0,135	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:20:41	042:21:01	0,135	10,0	-84,4	0,136	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:21:41	042:22:01	0,136	10,0	-86,3	0,136	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:22:42	042:23:02	0,136	10,0	-89,3	0,137	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 136 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	09:23:43	042:24:03	0,137	10,0	-79,8	0,138	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,3	
	09:24:40	042:25:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	09:24:40	042:25:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	09:24:40	042:25:00	0,130	10,0	-88,8	0,130	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:25:42	042:26:02	0,130	10,0	-84,0	0,131	9,2	9,2	10,3	10,9	10,8	9,5	11,7	12,4	
	09:26:43	042:27:03	0,131	10,0	-74,3	0,132	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,4	
	09:27:43	042:28:03	0,132	10,0	-89,0	0,132	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,3	
	09:28:43	042:29:03	0,132	10,0	-80,5	0,133	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,4	
	09:29:40	042:30:00	0,133	10,0	-90,2	0,134	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,4	
	09:30:41	042:31:01	0,134	10,0	-85,0	0,134	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,4	
	09:31:42	042:32:02	0,134	10,0	-77,3	0,135	9,2	9,2	10,3	11,0	10,8	9,5	11,7	12,4	
	09:32:44	042:33:04	0,135	10,0	-89,4	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:33:44	042:34:04	0,135	10,0	-81,8	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:34:41	042:35:01	0,136	10,0	-90,5	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:35:40	042:36:00	0,137	10,0	-88,7	0,138	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:36:44	042:37:04	0,130	10,0	-79,1	0,130	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:37:40	042:38:00	0,130	10,0	-89,9	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:38:40	042:39:00	0,131	10,0	-85,4	0,132	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:39:42	042:40:02	0,132	10,0	-76,5	0,133	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:40:42	042:41:02	0,132	10,0	-88,7	0,133	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	09:41:45	042:42:05	0,133	10,0	-82,6	0,134	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 137 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



09:42:41	042:43:01	0,134	10,0	-86,4	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:43:43	042:44:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
09:43:43	042:44:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
09:43:43	042:44:03	0,134	10,0	-88,3	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:44:40	042:45:00	0,135	10,0	-79,4	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:45:40	042:46:00	0,136	10,0	-88,8	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:46:43	042:47:03	0,136	10,0	-84,0	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:47:40	042:48:00	0,129	10,0	-76,3	0,130	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:48:41	042:49:01	0,130	10,0	-89,1	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:49:44	042:50:04	0,130	10,0	-81,8	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:50:40	042:51:00	0,132	10,0	-87,5	0,132	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:51:42	042:52:02	0,132	10,0	-88,2	0,133	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:52:42	042:53:02	0,132	10,0	-78,2	0,133	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:53:40	042:54:00	0,133	10,0	-89,9	0,134	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:54:43	042:55:03	0,134	10,0	-84,9	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:55:44	042:56:04	0,134	10,0	-75,8	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:56:40	042:57:00	0,135	10,0	-87,9	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:57:43	042:58:03	0,136	10,0	-82,8	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:58:44	042:59:04	0,136	10,0	-78,0	0,138	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
09:59:41	043:00:01	0,129	10,0	-88,8	0,130	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:00:44	043:01:04	0,130	10,0	-78,6	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 138 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



10:01:42	043:02:02	0,131	10,0	-89,2	0,132	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:02:43	043:03:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
10:02:43	043:03:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
10:02:43	043:03:03	0,132	10,0	-84,2	0,132	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:03:41	043:04:01	0,132	10,0	-74,9	0,133	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:04:43	043:05:03	0,133	10,0	-89,1	0,134	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:05:43	043:06:03	0,133	10,0	-81,6	0,134	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:06:40	043:07:00	0,134	10,0	-90,7	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:07:44	043:08:04	0,134	10,0	-85,3	0,136	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:08:40	043:09:00	0,135	10,0	-78,7	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:09:40	043:10:00	0,136	10,0	-87,7	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:10:43	043:11:03	0,131	10,0	-82,5	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:11:43	043:12:03	0,129	10,0	-82,2	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:12:43	043:13:03	0,130	10,0	-88,9	0,131	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:13:44	043:14:04	0,131	10,0	-80,5	0,132	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:14:40	043:15:00	0,131	10,0	-90,5	0,133	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4
10:15:43	043:16:03	0,132	10,0	-86,3	0,133	9,2	9,3	10,3	11,0	10,9	9,5	11,8	12,4
10:16:40	043:17:00	0,133	10,0	-77,4	0,134	9,2	9,3	10,3	11,0	10,9	9,5	11,8	12,4
10:17:42	043:18:02	0,133	10,0	-87,8	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,8	12,4
10:18:43	043:19:03	0,134	10,0	-83,7	0,135	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4
10:19:44	043:20:04	0,135	10,0	-80,1	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 139 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	10:20:41	043:21:01	0,135	10,0	-89,3	0,137	9,2	9,2	10,3	11,0	10,9	9,5	11,7	12,4	
	10:21:43	043:22:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	10:21:43	043:22:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	10:21:43	043:22:03	0,136	10,0	-79,9	0,138	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	10:22:45	043:23:05	0,129	10,0	-88,9	0,130	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:23:42	043:24:02	0,129	10,0	-84,2	0,131	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	10:24:43	043:25:03	0,130	10,0	-85,9	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:25:40	043:26:00	0,131	10,0	-89,2	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:26:43	043:27:03	0,131	10,0	-83,0	0,133	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	10:27:42	043:28:02	0,132	10,0	-90,3	0,134	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:28:42	043:29:02	0,133	10,0	-86,7	0,134	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4	
	10:29:43	043:30:03	0,133	10,0	-77,9	0,135	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4	
	10:30:43	043:31:03	0,134	10,0	-88,7	0,136	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4	
	10:31:40	043:32:00	0,136	10,0	-83,5	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:32:40	043:33:00	0,135	10,0	-77,8	0,137	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:33:42	043:34:02	0,128	10,0	-89,2	0,130	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:34:43	043:35:03	0,129	10,0	-81,2	0,130	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:35:45	043:36:05	0,130	10,0	-86,7	0,131	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	10:36:42	043:37:02	0,130	10,0	-88,5	0,132	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4	
	10:37:43	043:38:03	0,132	10,0	-78,4	0,132	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4	
	10:38:43	043:39:03	0,132	10,0	-88,5	0,133	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 140 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



10:39:40	043:40:00	0,132	10,0	-83,7	0,134	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:40:40	043:41:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
10:40:40	043:41:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
10:40:40	043:41:00	0,133	10,0	-74,9	0,135	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
10:41:41	043:42:01	0,135	10,0	-88,5	0,135	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
10:42:43	043:43:03	0,134	10,0	-82,1	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
10:43:40	043:44:00	0,135	10,0	-90,2	0,137	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
10:44:42	043:45:02	0,135	10,0	-88,5	0,137	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
10:45:42	043:46:02	0,128	10,0	-78,4	0,130	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:46:44	043:47:04	0,129	10,0	-88,7	0,131	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:47:40	043:48:00	0,130	10,0	-84,4	0,131	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:48:42	043:49:02	0,130	10,0	-79,3	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
10:49:44	043:50:04	0,132	10,0	-89,2	0,133	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
10:50:41	043:51:01	0,131	10,0	-81,4	0,133	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4
10:51:42	043:52:02	0,132	10,0	-90,4	0,134	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4
10:52:45	043:53:05	0,133	10,0	-86,2	0,135	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
10:53:42	043:54:02	0,134	10,0	-77,5	0,135	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
10:54:40	043:55:00	0,134	10,0	-89,1	0,136	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
10:55:40	043:56:00	0,135	10,0	-83,7	0,137	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:56:40	043:57:00	0,135	10,0	-77,7	0,138	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4
10:57:41	043:58:01	0,128	10,0	-89,1	0,130	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 141 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



10:58:42	043:59:02	0,129	10,0	-82,0	0,131	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
10:59:41	044:00:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
10:59:41	044:00:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
10:59:41	044:00:01	0,129	10,0	-81,7	0,131	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:00:44	044:01:04	0,130	10,0	-88,8	0,132	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:01:41	044:02:01	0,131	10,0	-79,4	0,133	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:02:43	044:03:03	0,131	10,0	-89,1	0,134	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:03:40	044:04:00	0,132	10,0	-85,0	0,134	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:04:41	044:05:01	0,133	10,0	-78,2	0,135	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:05:43	044:06:03	0,133	10,0	-89,3	0,135	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4
11:06:40	044:07:00	0,134	10,0	-81,8	0,136	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:07:40	044:08:00	0,135	10,0	-90,5	0,137	9,2	9,3	10,3	11,0	11,0	9,5	11,8	12,4
11:08:43	044:09:03	0,135	10,0	-86,8	0,137	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4
11:09:40	044:10:00	0,128	10,0	-80,4	0,130	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:10:40	044:11:00	0,129	10,0	-88,9	0,131	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:11:41	044:12:01	0,129	10,0	-83,7	0,131	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:12:41	044:13:01	0,130	10,0	-76,7	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4
11:13:42	044:14:02	0,131	10,0	-89,1	0,133	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:14:43	044:15:03	0,131	10,0	-81,4	0,134	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:15:44	044:16:04	0,132	10,0	-91,6	0,134	9,1	9,2	10,3	11,0	11,0	9,5	11,8	12,4
11:16:42	044:17:02	0,133	10,0	-88,8	0,135	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 142 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	11:17:43	044:18:03	0,133	10,0	-78,6	0,136	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,5	
	11:18:43	044:19:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	11:18:43	044:19:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	11:18:43	044:19:03	0,134	10,0	-88,7	0,137	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,5	
	11:19:40	044:20:00	0,135	10,0	-84,0	0,137	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:20:42	044:21:02	0,135	10,0	-75,9	0,137	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:21:43	044:22:03	0,128	10,0	-88,9	0,131	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,5	
	11:22:44	044:23:04	0,129	10,0	-81,4	0,131	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,5	
	11:23:41	044:24:01	0,131	10,0	-76,4	0,132	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,5	
	11:24:41	044:25:01	0,130	10,0	-88,4	0,132	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	11:25:44	044:26:04	0,131	10,0	-78,2	0,133	9,2	9,2	10,4	11,0	11,0	9,5	11,8	12,4	
	11:26:40	044:27:00	0,132	10,0	-89,8	0,134	9,2	9,2	10,4	11,0	11,0	9,5	11,8	12,4	
	11:27:43	044:28:03	0,132	10,0	-84,6	0,135	9,2	9,2	10,4	11,0	11,0	9,6	11,8	12,4	
	11:28:43	044:29:03	0,133	10,0	-76,8	0,135	9,2	9,2	10,4	11,0	11,0	9,5	11,8	12,4	
	11:29:41	044:30:01	0,134	10,0	-89,1	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	11:30:41	044:31:01	0,134	10,0	-82,0	0,137	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	11:31:41	044:32:01	0,135	10,0	-89,9	0,137	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:32:41	044:33:01	0,128	10,0	-87,8	0,130	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:33:42	044:34:02	0,129	10,0	-77,6	0,131	9,2	9,2	10,3	11,1	11,0	9,5	11,8	12,4	
	11:34:44	044:35:04	0,129	10,0	-87,9	0,131	9,1	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:35:45	044:36:05	0,130	10,0	-83,2	0,132	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 143 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	11:36:41	044:37:01	0,130	10,0	-84,0	0,133	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	11:37:42	044:38:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	11:37:42	044:38:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	11:37:42	044:38:02	0,132	10,0	-89,0	0,133	9,2	9,2	10,3	11,0	11,0	9,6	11,7	12,4	
	11:38:43	044:39:03	0,132	10,0	-79,1	0,134	9,1	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	11:39:44	044:40:04	0,133	10,0	-89,1	0,135	9,1	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	11:40:44	044:41:04	0,133	10,0	-84,1	0,135	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	11:41:41	044:42:01	0,134	10,0	-76,3	0,136	9,1	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	11:42:42	044:43:02	0,134	10,0	-88,5	0,137	9,1	9,2	10,3	11,0	11,1	9,5	11,7	12,4	
	11:43:42	044:44:02	0,135	10,0	-81,7	0,138	9,1	9,3	10,3	11,0	11,0	9,5	11,7	12,4	
	11:44:44	044:45:04	0,128	10,0	-92,0	0,130	9,2	9,2	10,4	11,0	11,0	9,5	11,7	12,4	
	11:45:40	044:46:00	0,130	10,0	-89,0	0,131	9,1	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	11:46:41	044:47:01	0,130	10,0	-79,7	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,8	12,4	
	11:47:44	044:48:04	0,130	10,0	-88,5	0,132	9,2	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	11:48:44	044:49:04	0,131	10,0	-83,4	0,133	9,2	9,3	10,4	11,1	11,0	9,6	11,7	12,4	
	11:49:43	044:50:03	0,131	10,0	-76,5	0,134	9,2	9,2	10,4	11,0	11,0	9,6	11,7	12,4	
	11:50:43	044:51:03	0,132	10,0	-88,9	0,134	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	11:51:43	044:52:03	0,133	10,0	-81,4	0,135	9,2	9,2	10,3	11,0	11,1	9,6	11,8	12,4	
	11:52:40	044:53:00	0,133	10,0	-84,2	0,136	9,2	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	11:53:40	044:54:00	0,134	10,0	-88,9	0,136	9,2	9,2	10,4	11,0	11,0	9,5	11,7	12,4	
	11:54:43	044:55:03	0,136	10,0	-79,3	0,137	9,2	9,2	10,4	11,1	11,0	9,5	11,8	12,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 144 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	11:55:41	044:56:01	0,128	10,0	-89,8	0,130	9,2	9,3	10,4	11,0	11,0	9,6	11,7	12,4	
	11:56:42	044:57:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	11:56:42	044:57:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	11:56:42	044:57:02	0,129	10,0	-84,8	0,130	9,2	9,3	10,4	11,1	11,0	9,5	11,7	12,4	
	11:57:44	044:58:04	0,129	10,0	-79,8	0,131	9,1	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	11:58:41	044:59:01	0,130	10,0	-89,3	0,132	9,2	9,2	10,3	11,0	11,1	9,5	11,7	12,4	
	11:59:42	045:00:02	0,130	10,0	-81,9	0,132	9,2	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	12:00:44	045:01:04	0,131	10,0	-90,0	0,133	9,2	9,2	10,3	11,1	11,1	9,5	11,7	12,4	
	12:01:41	045:02:01	0,131	10,0	-86,2	0,134	9,1	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	12:02:41	045:03:01	0,133	10,0	-78,0	0,135	9,2	9,2	10,3	11,0	11,0	9,6	11,7	12,4	
	12:03:42	045:04:02	0,133	10,0	-89,0	0,135	9,2	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	12:04:43	045:05:03	0,133	10,0	-82,2	0,136	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	12:05:40	045:06:00	0,134	10,0	-75,9	0,137	9,2	9,2	10,3	11,1	11,0	9,6	11,7	12,4	
	12:06:41	045:07:01	0,135	10,0	-88,7	0,137	9,2	9,2	10,3	11,1	11,0	9,6	11,7	12,4	
	12:07:42	045:08:02	0,128	10,0	-79,3	0,130	9,2	9,2	10,3	11,1	11,0	9,6	11,7	12,4	
	12:08:42	045:09:02	0,128	10,0	-89,9	0,131	9,2	9,2	10,3	11,1	11,1	9,6	11,7	12,4	
	12:09:43	045:10:03	0,129	10,0	-85,2	0,131	9,2	9,3	10,3	11,0	11,1	9,6	11,7	12,4	
	12:10:44	045:11:04	0,130	10,0	-76,3	0,132	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	12:11:44	045:12:04	0,130	10,0	-88,2	0,133	9,2	9,3	10,3	11,0	11,0	9,5	11,7	12,4	
	12:12:44	045:13:04	0,132	10,0	-82,8	0,133	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	12:13:46	045:14:06	0,132	10,0	-84,8	0,134	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 145 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	12:14:42	045:15:02	0,132	10,0	-88,9	0,135	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,5	
	12:15:45	045:16:05	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	12:15:45	045:16:05	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	12:15:45	045:16:05	0,133	10,0	-78,7	0,136	9,2	9,3	10,3	11,0	11,0	9,5	11,7	12,5	
	12:16:42	045:17:02	0,134	10,0	-89,0	0,136	9,2	9,3	10,3	11,0	11,0	9,5	11,7	12,4	
	12:17:43	045:18:03	0,134	10,0	-84,4	0,137	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	12:18:43	045:19:03	0,135	10,0	-76,0	0,138	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	12:19:44	045:20:04	0,128	10,0	-89,0	0,130	9,2	9,3	10,3	11,1	11,0	9,5	11,7	12,5	
	12:20:43	045:21:03	0,129	10,0	-80,5	0,131	9,2	9,3	10,3	11,1	11,1	9,5	11,7	12,5	
	12:21:43	045:22:03	0,129	10,0	-90,1	0,132	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,5	
	12:22:40	045:23:00	0,130	10,0	-86,4	0,132	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	12:23:41	045:24:01	0,131	10,0	-77,0	0,133	9,2	9,2	10,3	11,0	11,1	9,5	11,8	12,4	
	12:24:44	045:25:04	0,131	10,0	-88,0	0,134	9,2	9,2	10,3	11,0	11,0	9,6	11,8	12,4	
	12:25:40	045:26:00	0,132	10,0	-83,3	0,134	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	12:26:42	045:27:02	0,133	10,0	-80,0	0,135	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,5	
	12:27:43	045:28:03	0,133	10,0	-88,9	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,5	
	12:28:40	045:29:00	0,134	10,0	-79,8	0,136	9,2	9,2	10,3	11,0	11,0	9,5	11,7	12,4	
	12:29:40	045:30:00	0,135	10,0	-89,4	0,137	9,2	9,2	10,3	11,0	11,0	9,6	11,7	12,4	
	12:30:41	045:31:01	0,134	10,0	-84,0	0,136	9,1	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	12:31:41	045:32:01	0,128	10,0	-76,1	0,131	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	12:32:44	045:33:04	0,129	10,0	-88,7	0,131	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 146 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	12:33:44	045:34:04	0,130	10,0	-81,4	0,132	9,1	9,3	10,3	11,1	11,0	9,5	11,8	12,4	
	12:34:40	045:35:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	12:34:40	045:35:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	12:34:40	045:35:00	0,130	10,0	-89,0	0,133	9,1	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	12:35:42	045:36:02	0,131	10,0	-88,7	0,133	9,1	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	12:36:43	045:37:03	0,131	10,0	-79,7	0,134	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	12:37:43	045:38:03	0,132	10,0	-89,7	0,135	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	12:38:40	045:39:00	0,133	10,0	-85,8	0,135	9,2	9,2	10,3	11,1	11,0	9,6	11,7	12,4	
	12:39:41	045:40:01	0,134	10,0	-76,5	0,136	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4	
	12:40:43	045:41:03	0,134	10,0	-87,6	0,137	9,2	9,2	10,3	11,1	11,0	9,5	11,7	12,4	
	12:41:43	045:42:03	0,135	10,0	-83,1	0,137	9,2	9,3	10,3	11,1	11,1	9,5	11,7	12,4	
	12:42:43	045:43:03	0,128	10,0	-90,0	0,130	9,2	9,3	10,3	11,1	11,0	9,5	11,7	12,4	
	12:43:41	045:44:01	0,128	10,0	-88,9	0,131	9,2	9,3	10,3	11,1	11,0	9,5	11,7	12,4	
	12:44:41	045:45:01	0,129	10,0	-79,8	0,132	9,2	9,2	10,3	11,1	11,0	9,5	11,7	12,4	
	12:45:42	045:46:02	0,130	10,0	-89,4	0,132	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4	
	12:46:42	045:47:02	0,130	10,0	-84,3	0,133	9,2	9,3	10,3	11,1	11,1	9,5	11,8	12,4	
	12:47:41	045:48:01	0,131	10,0	-75,6	0,134	9,2	9,3	10,3	11,1	11,1	9,5	11,7	12,4	
	12:48:42	045:49:02	0,132	10,0	-87,7	0,134	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	12:49:43	045:50:03	0,132	10,0	-82,4	0,135	9,2	9,2	10,3	11,0	11,0	9,6	11,7	12,4	
	12:50:44	045:51:04	0,133	10,0	-90,4	0,136	9,2	9,2	10,3	11,0	11,1	9,6	11,7	12,4	
	12:51:40	045:52:00	0,134	10,0	-89,4	0,136	9,2	9,3	10,3	11,0	11,1	9,6	11,7	12,4	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 147 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	12:52:40	045:53:00	0,134	10,0	-80,6	0,137	9,2	9,3	10,3	11,0	11,0	9,5	11,7	12,4	
	12:53:41	045:54:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	12:53:41	045:54:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	12:53:41	045:54:01	0,135	10,0	-88,3	0,138	9,1	9,3	10,4	11,0	11,0	9,5	11,7	12,4	
	12:54:41	045:55:01	0,128	10,0	-83,1	0,130	9,2	9,3	10,4	11,1	11,0	9,5	11,7	12,5	
	12:55:42	045:56:02	0,129	10,0	-77,6	0,131	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,5	
	12:56:45	045:57:05	0,129	10,0	-88,7	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	12:57:44	045:58:04	0,130	10,0	-78,8	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,4	
	12:58:45	045:59:05	0,131	10,0	-89,2	0,133	9,2	9,3	10,3	11,0	11,1	9,6	11,8	12,4	
	12:59:43	046:00:03	0,131	10,0	-84,4	0,134	9,2	9,3	10,3	11,0	11,1	9,6	11,7	12,4	
	13:00:43	046:01:03	0,132	10,0	-76,7	0,135	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4	
	13:01:43	046:02:03	0,133	10,0	-88,9	0,135	9,2	9,3	10,3	11,1	11,0	9,5	11,8	12,4	
	13:02:44	046:03:04	0,133	10,0	-81,8	0,136	9,2	9,3	10,3	11,0	11,1	9,5	11,8	12,4	
	13:03:44	046:04:04	0,134	10,0	-89,4	0,137	9,2	9,3	10,3	11,0	11,1	9,5	11,8	12,4	
	13:04:44	046:05:04	0,135	10,0	-87,9	0,137	9,2	9,3	10,4	11,0	11,1	9,5	11,8	12,4	
	13:05:44	046:06:04	0,128	10,0	-78,0	0,130	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	13:06:44	046:07:04	0,128	10,0	-88,2	0,131	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4	
	13:07:40	046:08:00	0,129	10,0	-83,5	0,131	9,2	9,3	10,3	11,0	11,1	9,6	11,8	12,4	
	13:08:44	046:09:04	0,130	10,0	-76,9	0,132	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	
	13:09:40	046:10:00	0,130	10,0	-88,9	0,133	9,2	9,3	10,3	11,0	11,1	9,6	11,7	12,4	
	13:10:40	046:11:00	0,131	10,0	-79,9	0,133	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 148 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	13:11:41	046:12:01	0,132	10,0	-89,7	0,134	9,2	9,3	10,4	11,0	11,0	9,6	11,7	12,4	
	13:12:42	046:13:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	13:12:42	046:13:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	13:12:42	046:13:02	0,132	10,0	-84,3	0,135	9,2	9,3	10,4	11,0	11,0	9,6	11,8	12,4	
	13:13:44	046:14:04	0,133	10,0	-76,4	0,136	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	13:14:40	046:15:00	0,134	10,0	-88,3	0,136	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,5	
	13:15:41	046:16:01	0,134	10,0	-82,3	0,137	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	13:16:42	046:17:02	0,135	10,0	-80,9	0,138	9,1	9,3	10,3	11,0	11,1	9,6	11,8	12,4	
	13:17:43	046:18:03	0,128	10,0	-88,5	0,130	9,2	9,3	10,3	11,0	11,1	9,5	11,8	12,4	
	13:18:43	046:19:03	0,128	10,0	-78,2	0,131	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	
	13:19:43	046:20:03	0,129	10,0	-88,5	0,132	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	
	13:20:41	046:21:01	0,131	10,0	-83,1	0,132	9,2	9,3	10,3	11,0	11,0	9,5	11,7	12,4	
	13:21:41	046:22:01	0,130	10,0	-76,6	0,133	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	13:22:40	046:23:00	0,131	10,0	-88,4	0,134	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4	
	13:23:43	046:24:03	0,133	10,0	-81,8	0,134	9,2	9,3	10,3	11,1	11,1	9,5	11,7	12,4	
	13:24:40	046:25:00	0,133	10,0	-91,1	0,135	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	
	13:25:40	046:26:00	0,134	10,0	-88,5	0,136	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	13:26:40	046:27:00	0,134	10,0	-78,9	0,136	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	13:27:40	046:28:00	0,135	10,0	-89,5	0,137	9,2	9,3	10,4	11,1	11,0	9,6	11,7	12,4	
	13:28:43	046:29:03	0,135	10,0	-85,1	0,138	9,2	9,3	10,4	11,1	11,0	9,6	11,7	12,4	
	13:29:43	046:30:03	0,128	10,0	-91,3	0,131	9,2	9,3	10,4	11,0	11,0	9,6	11,7	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 149 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	13:30:40	046:31:00	0,129	10,0	-89,4	0,131	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	13:31:41	046:32:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	13:31:41	046:32:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	13:31:41	046:32:01	0,129	10,0	-80,2	0,132	9,2	9,3	10,3	11,0	11,1	9,5	11,7	12,4	
	13:32:43	046:33:03	0,131	10,0	-89,0	0,133	9,2	9,3	10,3	11,0	11,1	9,6	11,7	12,4	
	13:33:44	046:34:04	0,131	10,0	-83,3	0,133	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	13:34:45	046:35:05	0,131	10,0	-78,0	0,134	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	13:35:42	046:36:02	0,132	10,0	-88,5	0,135	9,2	9,3	10,3	11,0	11,0	9,6	11,8	12,4	
	13:36:44	046:37:04	0,133	10,0	-82,1	0,135	9,2	9,3	10,3	11,0	11,0	9,6	11,7	12,4	
	13:37:40	046:38:00	0,133	10,0	-90,2	0,136	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	13:38:40	046:39:00	0,134	10,0	-87,5	0,137	9,2	9,3	10,3	11,1	11,0	9,5	11,7	12,4	
	13:39:40	046:40:00	0,135	10,0	-77,6	0,137	9,2	9,2	10,4	11,1	11,0	9,5	11,7	12,4	
	13:40:41	046:41:01	0,128	10,0	-88,1	0,130	9,2	9,3	10,4	11,1	11,1	9,6	11,7	12,4	
	13:41:43	046:42:03	0,128	10,0	-83,3	0,131	9,2	9,3	10,4	11,0	11,1	9,6	11,7	12,4	
	13:42:44	046:43:04	0,129	10,0	-81,3	0,132	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	13:43:41	046:44:01	0,130	10,0	-88,9	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,4	
	13:44:41	046:45:01	0,130	10,0	-81,2	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:45:41	046:46:01	0,131	10,0	-88,4	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:46:42	046:47:02	0,132	10,0	-87,3	0,134	9,2	9,3	10,3	11,1	11,0	9,6	11,7	12,4	
	13:47:43	046:48:03	0,132	10,0	-77,0	0,135	9,2	9,3	10,4	11,1	11,1	9,6	11,7	12,4	
	13:48:44	046:49:04	0,133	10,0	-87,7	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,5	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 150 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	13:49:40	046:50:00	0,133	10,0	-83,5	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:50:40	046:51:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	13:50:40	046:51:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	13:50:40	046:51:00	0,134	10,0	-77,3	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,4	
	13:51:41	046:52:01	0,135	10,0	-89,0	0,138	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:52:43	046:53:03	0,128	10,0	-79,9	0,130	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4	
	13:53:41	046:54:01	0,130	10,0	-89,7	0,131	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4	
	13:54:41	046:55:01	0,130	10,0	-85,2	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:55:41	046:56:01	0,130	10,0	-76,6	0,132	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	13:56:41	046:57:01	0,131	10,0	-88,4	0,133	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	13:57:44	046:58:04	0,131	10,0	-81,9	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	13:58:41	046:59:01	0,132	10,0	-88,0	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	13:59:41	047:00:01	0,132	10,0	-88,1	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:00:43	047:01:03	0,133	10,0	-78,0	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:01:40	047:02:00	0,134	10,0	-89,1	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:02:41	047:03:01	0,134	10,0	-84,0	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:03:43	047:04:03	0,127	10,0	-75,5	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,4	
	14:04:43	047:05:03	0,128	10,0	-88,5	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:05:43	047:06:03	0,129	10,0	-81,9	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:06:40	047:07:00	0,129	10,0	-89,3	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:07:41	047:08:01	0,130	10,0	-88,8	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	151 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



14:08:43	047:09:03	0,131	10,0	-78,3	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:09:40	047:10:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
14:09:40	047:10:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
14:09:40	047:10:00	0,131	10,0	-88,7	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:10:45	047:11:05	0,132	10,0	-83,4	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
14:11:41	047:12:01	0,133	10,0	-77,6	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:12:43	047:13:03	0,134	10,0	-89,0	0,136	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:13:40	047:14:00	0,134	10,0	-81,6	0,137	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:14:40	047:15:00	0,135	10,0	-85,6	0,138	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:15:42	047:16:02	0,128	10,0	-88,0	0,130	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:16:41	047:17:01	0,128	10,0	-78,6	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:17:43	047:18:03	0,129	10,0	-88,7	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:18:40	047:19:00	0,130	10,0	-84,0	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:19:42	047:20:02	0,130	10,0	-76,7	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:20:40	047:21:00	0,131	10,0	-89,0	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:21:41	047:22:01	0,132	10,0	-81,6	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:22:42	047:23:02	0,132	10,0	-88,3	0,135	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
14:23:40	047:24:00	0,133	10,0	-88,0	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
14:24:41	047:25:01	0,134	10,0	-77,5	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
14:25:42	047:26:02	0,134	10,0	-87,8	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
14:26:44	047:27:04	0,135	10,0	-82,5	0,138	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 152 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	14:27:40	047:28:00	0,128	10,0	-78,8	0,130	9,2	9,3	10,4	11,0	11,1	9,6	11,8	12,4	
	14:28:41	047:29:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	14:28:41	047:29:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	14:28:41	047:29:01	0,128	10,0	-89,5	0,131	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:29:43	047:30:03	0,129	10,0	-80,7	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:30:44	047:31:04	0,130	10,0	-88,3	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:31:40	047:32:00	0,130	10,0	-82,2	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:32:42	047:33:02	0,131	10,0	-79,3	0,133	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:33:44	047:34:04	0,132	10,0	-89,2	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:34:40	047:35:00	0,132	10,0	-81,9	0,135	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:35:41	047:36:01	0,133	10,0	-89,7	0,136	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:36:43	047:37:03	0,134	10,0	-86,6	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:37:44	047:38:04	0,134	10,0	-76,5	0,137	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:38:41	047:39:01	0,135	10,0	-88,4	0,138	9,2	9,3	10,3	11,0	11,1	9,6	11,8	12,4	
	14:39:42	047:40:02	0,129	10,0	-84,0	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
	14:40:42	047:41:02	0,129	10,0	-90,3	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:41:42	047:42:02	0,129	10,0	-89,1	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:42:42	047:43:02	0,131	10,0	-80,2	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:43:42	047:44:02	0,131	10,0	-89,8	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:44:42	047:45:02	0,131	10,0	-85,3	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:45:43	047:46:03	0,133	10,0	-77,6	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 153 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	14:46:40	047:47:00	0,132	10,0	-88,8	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:47:44	047:48:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	14:47:44	047:48:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	14:47:44	047:48:04	0,133	10,0	-82,8	0,136	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:48:40	047:49:00	0,134	10,0	-77,5	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:49:41	047:50:01	0,134	10,0	-88,8	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:50:43	047:51:03	0,135	10,0	-79,4	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:51:40	047:52:00	0,128	10,0	-89,5	0,131	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:52:42	047:53:02	0,129	10,0	-83,9	0,131	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:53:41	047:54:01	0,129	10,0	-78,2	0,132	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
	14:54:40	047:55:00	0,130	10,0	-89,4	0,132	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:55:41	047:56:01	0,130	10,0	-81,4	0,133	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:56:44	047:57:04	0,131	10,0	-90,1	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
	14:57:40	047:58:00	0,132	10,0	-86,9	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:58:41	047:59:01	0,133	10,0	-77,5	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	14:59:41	048:00:01	0,133	10,0	-88,2	0,136	9,2	9,3	10,3	11,0	11,1	9,6	11,8	12,5	
	15:00:43	048:01:03	0,134	10,0	-83,8	0,137	9,2	9,3	10,3	11,0	11,1	9,6	11,8	12,5	
	15:01:40	048:02:00	0,135	10,0	-77,4	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	15:02:41	048:03:01	0,128	10,0	-88,9	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	15:03:41	048:04:01	0,129	10,0	-83,6	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
	15:04:44	048:05:04	0,129	10,0	-89,3	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	154 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



15:05:44	048:06:04	0,129	10,0	-84,3	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:06:41	048:07:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:06:41	048:07:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:06:41	048:07:01	0,130	10,0	-77,7	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4
15:07:43	048:08:03	0,131	10,0	-89,3	0,133	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5
15:08:44	048:09:04	0,131	10,0	-82,3	0,134	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5
15:09:42	048:10:02	0,132	10,0	-88,9	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:10:40	048:11:00	0,133	10,0	-88,6	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:11:41	048:12:01	0,133	10,0	-79,0	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:12:43	048:13:03	0,135	10,0	-89,1	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:13:41	048:14:01	0,135	10,0	-83,3	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:14:44	048:15:04	0,128	10,0	-90,6	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:15:40	048:16:00	0,129	10,0	-88,6	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:16:43	048:17:03	0,129	10,0	-78,8	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:17:44	048:18:04	0,129	10,0	-89,4	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:18:42	048:19:02	0,130	10,0	-83,7	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:19:42	048:20:02	0,131	10,0	-84,7	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:20:43	048:21:03	0,131	10,0	-89,0	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:21:45	048:22:05	0,132	10,0	-78,7	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:22:41	048:23:01	0,133	10,0	-88,3	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:23:41	048:24:01	0,133	10,0	-83,0	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 155 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



15:24:43	048:25:03	0,134	10,0	-76,1	0,137	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:25:43	048:26:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
15:25:43	048:26:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
15:25:43	048:26:03	0,135	10,0	-89,0	0,138	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:26:40	048:27:00	0,128	10,0	-81,8	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:27:44	048:28:04	0,128	10,0	-89,0	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:28:42	048:29:02	0,129	10,0	-88,5	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:29:43	048:30:03	0,130	10,0	-78,3	0,132	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5
15:30:43	048:31:03	0,130	10,0	-88,9	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:31:43	048:32:03	0,131	10,0	-84,8	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4
15:32:43	048:33:03	0,132	10,0	-92,3	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:33:43	048:34:03	0,132	10,0	-89,4	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:34:40	048:35:00	0,133	10,0	-83,7	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4
15:35:42	048:36:02	0,134	10,0	-89,1	0,136	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4
15:36:43	048:37:03	0,134	10,0	-87,9	0,137	9,2	9,3	10,3	11,1	11,0	9,6	11,8	12,4
15:37:40	048:38:00	0,135	10,0	-79,1	0,138	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:38:42	048:39:02	0,128	10,0	-89,3	0,130	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:39:43	048:40:03	0,128	10,0	-82,9	0,131	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:40:44	048:41:04	0,129	10,0	-77,8	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,7	12,5
15:41:41	048:42:01	0,130	10,0	-88,8	0,132	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5
15:42:44	048:43:04	0,130	10,0	-79,0	0,133	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 156 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



15:43:41	048:44:01	0,131	10,0	-89,0	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,5	
15:44:43	048:45:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
15:44:43	048:45:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
15:44:43	048:45:03	0,132	10,2	-83,2	0,134	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
15:45:42	048:46:02	0,132	10,4	-75,5	0,135	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
15:46:43	048:47:03	0,133	10,6	-88,7	0,136	9,2	9,3	10,3	11,1	11,1	9,6	11,8	12,4	
15:47:44	048:48:04	0,135	10,8	-82,4	0,137	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,4	
15:48:41	048:49:01	0,134	11,0	-88,7	0,137	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
15:49:44	048:50:04	0,128	11,2	-88,8	0,130	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
15:50:41	048:51:01	0,128	11,4	-80,5	0,131	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
15:51:43	048:52:03	0,129	11,6	-89,7	0,131	9,2	9,3	10,4	11,1	11,1	9,6	11,8	12,5	
15:52:41	048:53:01	0,129	11,8	-85,4	0,132	9,3	9,3	10,5	11,1	11,1	9,6	11,8	12,5	
15:53:42	048:54:02	0,130	12,0	-76,4	0,133	9,3	9,4	10,5	11,1	11,2	9,7	11,8	12,5	
15:54:44	048:55:04	0,131	12,2	-87,6	0,134	9,3	9,4	10,6	11,2	11,2	9,7	11,9	12,6	
15:55:40	048:56:00	0,131	12,4	-82,7	0,134	9,3	9,4	10,6	11,2	11,2	9,7	11,9	12,6	
15:56:43	048:57:03	0,132	12,6	-79,9	0,135	9,4	9,5	10,7	11,2	11,3	9,8	11,9	12,7	
15:57:45	048:58:05	0,133	12,8	-88,8	0,136	9,4	9,5	10,7	11,3	11,3	9,8	12,0	12,8	
15:58:42	048:59:02	0,134	13,0	-80,7	0,136	9,5	9,6	10,8	11,3	11,4	9,9	12,0	12,8	
15:59:43	049:00:03	0,134	13,2	-90,2	0,137	9,5	9,6	10,9	11,4	11,4	9,9	12,1	12,9	
16:00:40	049:01:00	0,135	13,4	-86,7	0,138	9,6	9,7	11,0	11,4	11,5	10,0	12,2	12,9	
16:01:43	049:02:03	0,129	13,6	-76,8	0,130	9,6	9,7	11,1	11,5	11,6	10,1	12,2	13,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 157 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



16:02:41	049:03:01	0,129	13,8	-87,9	0,131	9,7	9,8	11,1	11,6	11,6	10,2	12,3	13,1	
16:03:42	049:04:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
16:03:42	049:04:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
16:03:42	049:04:02	0,129	14,0	-82,3	0,132	9,8	9,9	11,2	11,6	11,7	10,3	12,4	13,2	
16:04:44	049:05:04	0,131	14,2	-91,3	0,133	9,9	10,0	11,3	11,7	11,8	10,3	12,4	13,3	
16:05:40	049:06:00	0,132	14,4	-89,2	0,133	10,0	10,0	11,4	11,8	11,9	10,4	12,5	13,3	
16:06:41	049:07:01	0,131	14,6	-81,4	0,134	10,1	10,1	11,5	11,9	12,0	10,5	12,6	13,4	
16:07:42	049:08:02	0,132	14,8	-90,0	0,135	10,1	10,2	11,6	12,0	12,0	10,6	12,7	13,5	
16:08:42	049:09:02	0,133	15,0	-85,9	0,135	10,2	10,3	11,7	12,1	12,1	10,7	12,8	13,6	
16:09:40	049:10:00	0,133	15,2	-77,0	0,136	10,3	10,4	11,8	12,1	12,2	10,8	12,9	13,7	
16:10:44	049:11:04	0,134	15,4	-87,9	0,137	10,4	10,5	12,0	12,3	12,4	10,9	13,0	13,9	
16:11:40	049:12:00	0,136	15,6	-82,3	0,137	10,5	10,6	12,1	12,3	12,4	11,0	13,0	13,9	
16:12:44	049:13:04	0,128	15,8	-79,3	0,130	10,6	10,7	12,2	12,4	12,6	11,1	13,2	14,1	
16:13:44	049:14:04	0,129	16,0	-88,9	0,131	10,7	10,8	12,3	12,5	12,7	11,2	13,3	14,2	
16:14:40	049:15:00	0,129	16,2	-81,6	0,132	10,8	10,9	12,4	12,6	12,8	11,3	13,4	14,3	
16:15:43	049:16:03	0,130	16,4	-90,8	0,132	10,9	11,0	12,5	12,7	12,9	11,4	13,5	14,4	
16:16:45	049:17:05	0,131	16,6	-86,0	0,133	11,0	11,1	12,7	12,9	13,0	11,6	13,6	14,5	
16:17:42	049:18:02	0,131	16,8	-83,6	0,134	11,1	11,2	12,8	13,0	13,1	11,7	13,7	14,6	
16:18:43	049:19:03	0,132	17,0	-88,7	0,135	11,2	11,3	12,9	13,1	13,2	11,8	13,8	14,7	
16:19:43	049:20:03	0,133	17,2	-83,3	0,135	11,3	11,4	13,0	13,2	13,3	11,9	13,9	14,8	
16:20:44	049:21:04	0,134	17,4	-76,3	0,136	11,4	11,5	13,2	13,3	13,4	12,0	14,0	15,0	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 158 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



16:21:44	049:22:04	0,134	17,6	-89,0	0,137	11,5	11,7	13,3	13,4	13,6	12,1	14,2	15,1	
16:22:40	049:23:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
16:22:40	049:23:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
16:22:40	049:23:00	0,135	17,8	-82,8	0,137	11,6	11,8	13,4	13,5	13,7	12,3	14,3	15,2	
16:23:43	049:24:03	0,128	18,0	-84,4	0,130	11,8	11,9	13,6	13,7	13,8	12,4	14,4	15,3	
16:24:43	049:25:03	0,130	18,2	-88,3	0,131	11,9	12,0	13,7	13,8	13,9	12,5	14,5	15,5	
16:25:41	049:26:01	0,129	18,4	-78,5	0,132	12,0	12,2	13,9	13,9	14,1	12,7	14,6	15,6	
16:26:43	049:27:03	0,130	18,6	-88,7	0,133	12,2	12,3	14,0	14,1	14,2	12,8	14,8	15,7	
16:27:44	049:28:04	0,131	18,8	-83,0	0,133	12,3	12,4	14,1	14,2	14,3	12,9	14,9	15,8	
16:28:41	049:29:01	0,132	19,0	-77,1	0,134	12,4	12,5	14,3	14,3	14,5	13,0	15,0	16,0	
16:29:42	049:30:02	0,133	19,2	-88,9	0,135	12,6	12,7	14,4	14,4	14,6	13,2	15,2	16,1	
16:30:44	049:31:04	0,133	19,4	-82,1	0,135	12,7	12,8	14,5	14,6	14,8	13,3	15,3	16,2	
16:31:42	049:32:02	0,134	19,6	-79,4	0,136	12,8	12,9	14,7	14,7	14,9	13,4	15,4	16,4	
16:32:40	049:33:00	0,135	19,8	-89,2	0,137	12,9	13,1	14,8	14,8	15,0	13,6	15,6	16,5	
16:33:40	049:34:00	0,135	20,0	-79,5	0,138	13,1	13,2	15,0	15,0	15,2	13,7	15,7	16,7	
16:34:41	049:35:01	0,128	20,2	-88,4	0,130	13,2	13,3	15,1	15,1	15,3	13,9	15,8	16,8	
16:35:43	049:36:03	0,129	20,4	-83,2	0,131	13,3	13,5	15,3	15,3	15,4	14,0	16,0	16,9	
16:36:43	049:37:03	0,130	20,6	-75,2	0,132	13,5	13,6	15,4	15,4	15,6	14,1	16,1	17,1	
16:37:44	049:38:04	0,130	20,8	-88,7	0,133	13,6	13,7	15,6	15,5	15,7	14,3	16,2	17,2	
16:38:44	049:39:04	0,132	21,0	-79,9	0,133	13,7	13,9	15,7	15,7	15,9	14,4	16,4	17,4	
16:39:41	049:40:01	0,132	21,2	-90,1	0,134	13,9	14,0	15,9	15,8	16,0	14,5	16,5	17,5	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	159 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



16:40:42	049:41:02	0,133	21,4	-85,8	0,135	14,0	14,2	16,0	16,0	16,2	14,7	16,6	17,6	
16:41:42	049:42:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
16:41:42	049:42:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
16:41:42	049:42:02	0,134	21,6	-76,3	0,136	14,1	14,3	16,2	16,1	16,3	14,8	16,8	17,8	
16:42:42	049:43:02	0,134	21,8	-87,6	0,137	14,3	14,4	16,4	16,3	16,5	15,0	16,9	17,9	
16:43:43	049:44:03	0,135	22,0	-82,9	0,137	14,4	14,6	16,5	16,4	16,6	15,1	17,1	18,1	
16:44:44	049:45:04	0,128	22,2	-89,2	0,130	14,6	14,7	16,7	16,6	16,8	15,3	17,2	18,2	
16:45:41	049:46:01	0,129	22,4	-89,3	0,131	14,7	14,9	16,8	16,7	16,9	15,4	17,4	18,4	
16:46:41	049:47:01	0,129	22,6	-80,7	0,131	14,8	15,0	17,0	16,8	17,1	15,6	17,5	18,5	
16:47:41	049:48:01	0,130	22,8	-88,9	0,132	15,0	15,2	17,1	17,0	17,2	15,7	17,6	18,6	
16:48:42	049:49:02	0,131	23,0	-83,8	0,133	15,1	15,3	17,3	17,1	17,4	15,8	17,8	18,8	
16:49:43	049:50:03	0,132	23,2	-79,8	0,134	15,2	15,5	17,5	17,3	17,5	16,0	17,9	18,9	
16:50:43	049:51:03	0,133	23,4	-89,1	0,135	15,3	15,6	17,6	17,5	17,7	16,1	18,1	19,1	
16:51:43	049:52:03	0,133	23,6	-82,2	0,135	15,6	15,8	17,8	17,6	17,8	16,3	18,2	19,2	
16:52:44	049:53:04	0,134	23,8	-78,7	0,136	15,7	15,9	18,0	17,8	18,0	16,4	18,4	19,4	
16:53:40	049:54:00	0,135	24,0	-88,4	0,137	15,9	16,1	18,1	17,9	18,2	16,6	18,6	19,5	
16:54:41	049:55:01	0,136	24,2	-78,7	0,138	16,1	16,3	18,3	18,1	18,3	16,7	18,7	19,7	
16:55:42	049:56:02	0,130	24,4	-88,9	0,130	16,2	16,4	18,5	18,2	18,5	16,9	18,8	19,8	
16:56:43	049:57:03	0,129	24,6	-83,6	0,131	16,3	16,6	18,6	18,4	18,6	17,1	19,0	20,0	
16:57:44	049:58:04	0,130	24,8	-78,2	0,132	16,4	16,7	18,8	18,6	18,8	17,2	19,1	20,1	
16:58:42	049:59:02	0,131	25,0	-89,1	0,133	16,6	16,9	18,9	18,7	18,9	17,4	19,3	20,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature:
Document issued by: Page: 160 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	16:59:42	050:00:02	0,132	25,0	-81,6	0,133	16,7	17,0	19,1	18,9	19,1	17,5	19,5	20,4	
	17:00:40	050:01:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	17:00:40	050:01:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	17:00:40	050:01:00	0,133	25,0	-88,3	0,134	16,9	17,2	19,3	19,1	19,3	17,7	19,6	20,6	
	17:01:41	050:02:01	0,133	25,0	-87,9	0,135	17,1	17,3	19,5	19,2	19,4	17,8	19,8	20,7	
	17:02:42	050:03:02	0,134	25,0	-78,1	0,136	17,2	17,5	19,6	19,4	19,6	18,0	19,9	20,9	
	17:03:42	050:04:02	0,135	25,0	-88,8	0,137	17,3	17,6	19,8	19,5	19,8	18,1	20,1	21,0	
	17:04:42	050:05:02	0,136	25,0	-84,0	0,137	17,5	17,8	19,9	19,7	19,9	18,3	20,2	21,2	
	17:05:42	050:06:02	0,129	25,0	-90,5	0,130	17,6	18,0	20,1	19,8	20,1	18,4	20,3	21,3	
	17:06:44	050:07:04	0,129	25,0	-89,2	0,131	17,8	18,1	20,2	20,0	20,2	18,6	20,5	21,4	
	17:07:44	050:08:04	0,130	25,0	-81,5	0,131	17,9	18,2	20,4	20,1	20,4	18,7	20,6	21,6	
	17:08:40	050:09:00	0,131	25,0	-89,9	0,132	18,1	18,4	20,5	20,3	20,5	18,8	20,8	21,7	
	17:09:41	050:10:01	0,132	25,0	-85,5	0,133	18,2	18,5	20,6	20,4	20,7	18,9	20,9	21,8	
	17:10:41	050:11:01	0,132	25,0	-76,6	0,134	18,3	18,6	20,7	20,6	20,8	19,0	21,0	21,9	
	17:11:43	050:12:03	0,133	25,0	-88,7	0,134	18,4	18,7	20,8	20,7	20,9	19,1	21,1	22,0	
	17:12:44	050:13:04	0,134	25,0	-82,7	0,135	18,5	18,8	20,9	20,8	21,0	19,2	21,3	22,1	
	17:13:44	050:14:04	0,135	25,0	-90,1	0,136	18,6	19,0	21,0	20,9	21,1	19,3	21,4	22,2	
	17:14:44	050:15:04	0,135	25,0	-85,6	0,137	18,7	19,1	21,1	21,0	21,3	19,4	21,5	22,3	
	17:15:41	050:16:01	0,136	25,0	-77,0	0,137	18,8	19,2	21,2	21,1	21,3	19,5	21,5	22,4	
	17:16:41	050:17:01	0,129	25,0	-88,6	0,130	18,9	19,3	21,3	21,3	21,4	19,6	21,6	22,5	
	17:17:44	050:18:04	0,130	25,0	-82,3	0,131	19,0	19,4	21,4	21,4	21,6	19,7	21,7	22,6	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 161 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



17:18:44	050:19:04	0,131	25,0	-90,5	0,132	19,1	19,5	21,5	21,5	21,7	19,8	21,8	22,6
17:19:42	050:20:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
17:19:42	050:20:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
17:19:42	050:20:02	0,131	25,0	-87,7	0,132	19,2	19,6	21,6	21,6	21,8	19,8	21,9	22,7
17:20:42	050:21:02	0,132	25,0	-78,6	0,133	19,3	19,6	21,7	21,6	21,8	19,9	22,0	22,8
17:21:42	050:22:02	0,133	25,0	-88,2	0,134	19,3	19,7	21,7	21,7	21,9	20,0	22,0	22,8
17:22:43	050:23:03	0,133	25,0	-83,8	0,134	19,4	19,8	21,8	21,8	22,0	20,0	22,1	22,9
17:23:45	050:24:05	0,134	25,0	-91,6	0,135	19,5	19,9	21,9	21,9	22,1	20,1	22,2	23,0
17:24:41	050:25:01	0,135	25,0	-88,9	0,136	19,5	20,0	22,0	22,0	22,2	20,2	22,3	23,0
17:25:43	050:26:03	0,136	25,0	-80,7	0,137	19,6	20,0	22,0	22,1	22,2	20,2	22,3	23,1
17:26:45	050:27:05	0,136	25,0	-89,8	0,138	19,7	20,1	22,1	22,1	22,3	20,3	22,4	23,1
17:27:40	050:28:00	0,129	25,0	-87,4	0,130	19,7	20,2	22,2	22,2	22,4	20,3	22,5	23,2
17:28:40	050:29:00	0,131	25,0	-77,0	0,131	19,8	20,2	22,2	22,3	22,4	20,4	22,5	23,2
17:29:44	050:30:04	0,131	25,0	-87,9	0,132	19,8	20,3	22,3	22,3	22,5	20,4	22,6	23,3
17:30:40	050:31:00	0,132	25,0	-84,3	0,132	19,9	20,3	22,3	22,4	22,6	20,5	22,6	23,3
17:31:41	050:32:01	0,132	25,0	-88,4	0,133	19,9	20,4	22,4	22,5	22,6	20,5	22,7	23,4
17:32:43	050:33:03	0,134	25,0	-87,9	0,134	20,0	20,5	22,5	22,5	22,7	20,6	22,7	23,4
17:33:43	050:34:03	0,135	25,0	-78,6	0,135	20,0	20,5	22,5	22,6	22,7	20,6	22,8	23,5
17:34:43	050:35:03	0,135	25,0	-88,4	0,135	20,1	20,6	22,6	22,7	22,8	20,7	22,8	23,5
17:35:44	050:36:04	0,135	25,0	-84,0	0,136	20,1	20,6	22,6	22,7	22,9	20,7	22,9	23,5
17:36:40	050:37:00	0,136	25,0	-90,6	0,137	20,2	20,7	22,7	22,8	22,9	20,7	22,9	23,6

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 162 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	17:37:40	050:38:00	0,136	25,0	-89,4	0,138	20,2	20,7	22,7	22,8	22,9	20,8	22,9	23,6	
	17:38:40	050:39:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	17:38:40	050:39:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	17:38:40	050:39:00	0,129	25,0	-79,8	0,130	20,3	20,8	22,8	22,9	23,0	20,8	23,0	23,6	
	17:39:42	050:40:02	0,130	25,0	-88,0	0,131	20,3	20,8	22,8	22,9	23,1	20,9	23,0	23,7	
	17:40:44	050:41:04	0,131	25,0	-83,2	0,131	20,3	20,9	22,9	23,0	23,1	20,9	23,1	23,7	
	17:41:44	050:42:04	0,132	25,0	-84,6	0,132	20,4	20,9	22,9	23,0	23,1	20,9	23,1	23,7	
	17:42:41	050:43:01	0,132	25,0	-88,7	0,133	20,4	21,0	22,9	23,1	23,2	21,0	23,1	23,8	
	17:43:43	050:44:03	0,133	25,0	-78,7	0,134	20,4	21,0	23,0	23,1	23,2	21,0	23,2	23,8	
	17:44:44	050:45:04	0,134	25,0	-88,9	0,135	20,5	21,0	23,0	23,2	23,3	21,0	23,2	23,8	
	17:45:41	050:46:01	0,134	25,0	-84,2	0,135	20,5	21,1	23,1	23,2	23,3	21,0	23,2	23,9	
	17:46:41	050:47:01	0,135	25,0	-76,4	0,136	20,5	21,1	23,1	23,2	23,3	21,1	23,3	23,9	
	17:47:42	050:48:02	0,136	25,0	-88,6	0,137	20,6	21,2	23,1	23,3	23,4	21,1	23,3	23,9	
	17:48:42	050:49:02	0,136	25,0	-82,0	0,137	20,6	21,2	23,2	23,3	23,4	21,1	23,3	23,9	
	17:49:42	050:50:02	0,129	25,0	-81,1	0,130	20,7	21,2	23,2	23,4	23,5	21,2	23,3	24,0	
	17:50:42	050:51:02	0,130	25,0	-88,7	0,131	20,7	21,3	23,2	23,4	23,5	21,2	23,4	24,0	
	17:51:43	050:52:03	0,132	25,0	-79,1	0,131	20,7	21,3	23,3	23,4	23,5	21,2	23,4	24,0	
	17:52:43	050:53:03	0,132	25,0	-89,2	0,132	20,8	21,3	23,3	23,5	23,5	21,2	23,4	24,0	
	17:53:43	050:54:03	0,132	25,0	-84,3	0,133	20,8	21,4	23,3	23,5	23,6	21,3	23,5	24,1	
	17:54:40	050:55:00	0,133	25,0	-75,9	0,134	20,8	21,4	23,4	23,5	23,6	21,3	23,5	24,1	
	17:55:41	050:56:01	0,134	25,0	-87,6	0,134	20,8	21,4	23,4	23,6	23,6	21,3	23,5	24,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 163 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	17:56:43	050:57:03	0,134	25,0	-82,4	0,135	20,8	21,4	23,4	23,6	23,7	21,3	23,5	24,1	
	17:57:41	050:58:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	17:57:41	050:58:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	17:57:41	050:58:01	0,135	25,0	-86,0	0,136	20,8	21,5	23,4	23,6	23,7	21,3	23,6	24,1	
	17:58:41	050:59:01	0,137	25,0	-88,6	0,137	20,9	21,5	23,5	23,7	23,7	21,4	23,6	24,2	
	17:59:42	051:00:02	0,136	25,0	-79,8	0,137	20,9	21,5	23,5	23,7	23,7	21,4	23,6	24,2	
	18:00:40	051:01:00	0,129	25,0	-89,5	0,130	20,9	21,6	23,5	23,7	23,8	21,4	23,6	24,2	
	18:01:41	051:02:01	0,131	25,0	-84,8	0,131	21,0	21,6	23,5	23,7	23,8	21,4	23,6	24,2	
	18:02:43	051:03:03	0,131	25,0	-81,3	0,131	21,0	21,6	23,6	23,8	23,8	21,5	23,7	24,2	
	18:03:44	051:04:04	0,131	25,0	-89,1	0,132	21,0	21,6	23,6	23,8	23,8	21,5	23,7	24,3	
	18:04:40	051:05:00	0,132	25,0	-81,1	0,133	21,0	21,7	23,6	23,8	23,9	21,5	23,7	24,3	
	18:05:42	051:06:02	0,133	25,0	-89,5	0,134	21,1	21,7	23,6	23,8	23,9	21,5	23,7	24,3	
	18:06:43	051:07:03	0,134	25,0	-84,7	0,134	21,1	21,7	23,7	23,9	23,9	21,5	23,7	24,3	
	18:07:44	051:08:04	0,135	25,0	-82,4	0,135	21,1	21,7	23,7	23,9	23,9	21,5	23,8	24,3	
	18:08:44	051:09:04	0,135	25,0	-89,1	0,136	21,1	21,8	23,7	23,9	23,9	21,6	23,8	24,3	
	18:09:44	051:10:04	0,136	25,0	-81,1	0,137	21,1	21,8	23,7	24,0	24,0	21,6	23,8	24,3	
	18:10:40	051:11:00	0,136	25,0	-89,7	0,137	21,2	21,8	23,7	24,0	24,0	21,6	23,8	24,4	
	18:11:40	051:12:00	0,129	25,0	-85,3	0,130	21,2	21,8	23,8	24,0	24,0	21,6	23,8	24,4	
	18:12:42	051:13:02	0,131	25,0	-82,0	0,131	21,2	21,8	23,8	24,0	24,0	21,6	23,9	24,4	
	18:13:43	051:14:03	0,131	25,0	-89,1	0,131	21,2	21,9	23,8	24,0	24,1	21,6	23,9	24,4	
	18:14:43	051:15:03	0,131	25,0	-79,3	0,132	21,2	21,9	23,8	24,0	24,1	21,6	23,9	24,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	164 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



18:15:40	051:16:00	0,132	25,0	-88,5	0,133	21,2	21,9	23,8	24,1	24,1	21,6	23,9	24,4	
18:16:41	051:17:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
18:16:41	051:17:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
18:16:41	051:17:01	0,133	25,0	-84,3	0,133	21,2	21,9	23,8	24,1	24,1	21,7	23,9	24,4	
18:17:45	051:18:05	0,134	25,0	-82,0	0,134	21,3	21,9	23,9	24,1	24,1	21,7	23,9	24,5	
18:18:41	051:19:01	0,134	25,0	-88,9	0,135	21,3	22,0	23,9	24,1	24,1	21,7	23,9	24,5	
18:19:41	051:20:01	0,135	25,0	-80,1	0,136	21,3	22,0	23,9	24,1	24,1	21,7	23,9	24,5	
18:20:41	051:21:01	0,136	25,0	-89,6	0,137	21,3	22,0	23,9	24,2	24,2	21,7	23,9	24,5	
18:21:42	051:22:02	0,136	25,0	-84,9	0,137	21,3	22,0	23,9	24,2	24,2	21,7	24,0	24,5	
18:22:40	051:23:00	0,130	25,0	-76,6	0,130	21,3	22,0	23,9	24,2	24,2	21,7	24,0	24,5	
18:23:40	051:24:00	0,130	25,0	-87,6	0,131	21,4	22,0	23,9	24,2	24,2	21,8	24,0	24,5	
18:24:42	051:25:02	0,130	25,0	-82,8	0,131	21,4	22,0	24,0	24,2	24,2	21,8	24,0	24,5	
18:25:43	051:26:03	0,131	25,0	-91,1	0,132	21,4	22,1	24,0	24,2	24,2	21,8	24,0	24,5	
18:26:43	051:27:03	0,132	25,0	-88,5	0,133	21,4	22,1	24,0	24,2	24,3	21,8	24,0	24,5	
18:27:43	051:28:03	0,133	25,0	-78,3	0,133	21,4	22,1	24,0	24,3	24,3	21,8	24,0	24,6	
18:28:40	051:29:00	0,133	25,0	-88,6	0,134	21,4	22,1	24,0	24,3	24,3	21,8	24,1	24,6	
18:29:41	051:30:01	0,134	25,0	-84,2	0,135	21,4	22,1	24,0	24,3	24,3	21,8	24,1	24,6	
18:30:41	051:31:01	0,135	25,0	-84,0	0,136	21,4	22,1	24,0	24,3	24,3	21,8	24,1	24,6	
18:31:44	051:32:04	0,135	25,0	-87,8	0,136	21,4	22,1	24,0	24,3	24,3	21,8	24,1	24,6	
18:32:44	051:33:04	0,136	25,0	-78,1	0,137	21,5	22,2	24,1	24,3	24,3	21,9	24,1	24,6	
18:33:44	051:34:04	0,129	25,0	-88,4	0,130	21,5	22,2	24,1	24,3	24,3	21,9	24,1	24,6	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	165 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



18:34:41	051:35:01	0,129	25,0	-84,3	0,130	21,5	22,2	24,1	24,3	24,3	21,9	24,1	24,6
18:35:41	051:36:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8
18:35:41	051:36:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C
18:35:41	051:36:01	0,130	25,0	-77,5	0,131	21,5	22,2	24,1	24,3	24,3	21,9	24,1	24,6
18:36:40	051:37:00	0,131	25,0	-88,9	0,132	21,5	22,2	24,1	24,3	24,4	21,9	24,1	24,6
18:37:41	051:38:01	0,132	25,0	-80,4	0,133	21,5	22,2	24,1	24,4	24,4	21,9	24,1	24,7
18:38:41	051:39:01	0,132	25,0	-88,2	0,133	21,5	22,2	24,1	24,4	24,4	21,9	24,1	24,7
18:39:42	051:40:02	0,134	25,0	-83,8	0,134	21,5	22,2	24,1	24,4	24,4	21,9	24,2	24,7
18:40:40	051:41:00	0,134	25,0	-75,2	0,135	21,6	22,2	24,1	24,4	24,4	21,9	24,2	24,7
18:41:40	051:42:00	0,134	25,0	-88,7	0,135	21,6	22,3	24,1	24,4	24,4	21,9	24,2	24,7
18:42:41	051:43:01	0,135	25,0	-79,2	0,136	21,6	22,3	24,1	24,4	24,4	21,9	24,2	24,7
18:43:42	051:44:02	0,136	25,0	-88,9	0,137	21,6	22,3	24,2	24,4	24,4	21,9	24,2	24,7
18:44:43	051:45:03	0,136	25,0	-83,9	0,138	21,6	22,3	24,2	24,4	24,4	22,0	24,2	24,7
18:45:40	051:46:00	0,129	25,0	-75,7	0,130	21,6	22,3	24,2	24,4	24,4	22,0	24,2	24,7
18:46:40	051:47:00	0,130	25,0	-88,5	0,131	21,6	22,3	24,2	24,4	24,4	22,0	24,2	24,7
18:47:41	051:48:01	0,131	25,0	-81,7	0,132	21,6	22,3	24,2	24,4	24,4	22,0	24,2	24,7
18:48:41	051:49:01	0,131	25,0	-88,7	0,133	21,6	22,3	24,2	24,5	24,5	22,0	24,2	24,7
18:49:44	051:50:04	0,132	25,0	-86,5	0,133	21,6	22,3	24,2	24,5	24,5	22,0	24,2	24,7
18:50:40	051:51:00	0,133	25,0	-80,9	0,134	21,6	22,3	24,2	24,5	24,5	22,0	24,2	24,8
18:51:42	051:52:02	0,133	25,0	-89,0	0,135	21,6	22,3	24,2	24,5	24,5	22,0	24,2	24,8
18:52:44	051:53:04	0,134	25,0	-81,3	0,135	21,6	22,3	24,2	24,5	24,5	22,0	24,3	24,7

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 166 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	18:53:40	051:54:00	0,135	25,0	-85,4	0,136	21,6	22,3	24,2	24,5	24,5	22,0	24,3	24,8	
	18:54:40	051:55:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	18:54:40	051:55:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	18:54:40	051:55:00	0,135	25,0	-87,8	0,137	21,6	22,4	24,2	24,5	24,5	22,0	24,3	24,8	
	18:55:43	051:56:03	0,136	25,0	-77,1	0,138	21,7	22,4	24,2	24,5	24,5	22,0	24,3	24,8	
	18:56:42	051:57:02	0,129	25,0	-87,9	0,130	21,7	22,4	24,3	24,5	24,5	22,0	24,3	24,8	
	18:57:44	051:58:04	0,130	25,0	-83,1	0,131	21,7	22,4	24,3	24,5	24,5	22,0	24,3	24,8	
	18:58:44	051:59:04	0,130	25,0	-84,2	0,132	21,7	22,4	24,3	24,5	24,6	22,0	24,3	24,8	
	18:59:43	052:00:03	0,132	25,0	-88,6	0,132	21,7	22,4	24,3	24,5	24,5	22,0	24,3	24,8	
	19:00:44	052:01:04	0,132	25,0	-80,1	0,133	21,7	22,4	24,3	24,5	24,6	22,1	24,3	24,8	
	19:01:41	052:02:01	0,132	25,0	-90,0	0,134	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:02:43	052:03:03	0,133	25,0	-85,2	0,134	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:03:44	052:04:04	0,135	25,0	-76,9	0,135	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:04:44	052:05:04	0,135	25,0	-88,6	0,136	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:05:41	052:06:01	0,135	25,0	-82,3	0,136	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:06:42	052:07:02	0,136	25,0	-91,3	0,137	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:07:43	052:08:03	0,129	25,0	-88,1	0,130	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:08:44	052:09:04	0,129	25,0	-77,6	0,130	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,8	
	19:09:41	052:10:01	0,130	25,0	-88,4	0,131	21,7	22,4	24,3	24,6	24,6	22,1	24,3	24,9	
	19:10:40	052:11:00	0,131	25,0	-85,3	0,132	21,7	22,5	24,3	24,6	24,6	22,1	24,3	24,9	
	19:11:41	052:12:01	0,131	25,0	-92,2	0,133	21,7	22,5	24,3	24,6	24,6	22,1	24,4	24,9	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	167 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



19:12:41	052:13:01	0,132	25,0	-88,9	0,133	21,7	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:13:42	052:14:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
19:13:42	052:14:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
19:13:42	052:14:02	0,134	25,0	-79,6	0,134	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:14:43	052:15:03	0,133	25,0	-88,8	0,135	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:15:40	052:16:00	0,134	25,0	-84,0	0,135	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:16:42	052:17:02	0,135	25,0	-76,4	0,136	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:17:44	052:18:04	0,135	25,0	-88,8	0,137	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:18:40	052:19:00	0,136	25,0	-81,7	0,138	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:19:41	052:20:01	0,129	25,0	-90,0	0,130	21,8	22,5	24,3	24,6	24,6	22,1	24,4	24,9	
19:20:40	052:21:00	0,131	25,0	-87,2	0,131	21,8	22,5	24,3	24,7	24,6	22,1	24,4	24,9	
19:21:41	052:22:01	0,130	25,0	-77,9	0,132	21,8	22,5	24,4	24,7	24,6	22,1	24,4	24,9	
19:22:43	052:23:03	0,131	25,0	-88,9	0,133	21,8	22,5	24,4	24,7	24,6	22,1	24,4	24,9	
19:23:40	052:24:00	0,132	25,0	-82,9	0,133	21,8	22,5	24,4	24,7	24,6	22,2	24,4	24,9	
19:24:42	052:25:02	0,133	25,0	-79,2	0,134	21,8	22,5	24,4	24,7	24,6	22,2	24,4	24,9	
19:25:40	052:26:00	0,133	25,0	-88,4	0,135	21,8	22,5	24,4	24,7	24,6	22,2	24,4	24,9	
19:26:41	052:27:01	0,135	25,0	-78,5	0,135	21,8	22,5	24,4	24,7	24,7	22,2	24,4	24,9	
19:27:42	052:28:02	0,136	25,0	-88,4	0,136	21,8	22,5	24,4	24,7	24,7	22,2	24,4	24,9	
19:28:43	052:29:03	0,135	25,0	-84,4	0,137	21,8	22,5	24,4	24,7	24,7	22,2	24,4	24,9	
19:29:44	052:30:04	0,136	25,0	-78,5	0,138	21,8	22,6	24,4	24,7	24,7	22,2	24,4	24,9	
19:30:40	052:31:00	0,129	25,0	-88,9	0,130	21,8	22,6	24,4	24,7	24,7	22,2	24,4	24,9	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 168 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



19:31:41	052:32:01	0,129	25,0	-78,9	0,131	21,8	22,6	24,4	24,7	24,7	22,2	24,4	24,9	
19:32:44	052:33:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
19:32:44	052:33:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
19:32:44	052:33:04	0,130	25,0	-87,7	0,132	21,8	22,6	24,4	24,7	24,7	22,2	24,4	24,9	
19:33:40	052:34:00	0,131	25,0	-83,5	0,132	21,9	22,6	24,4	24,7	24,7	22,2	24,4	24,9	
19:34:41	052:35:01	0,132	25,0	-75,9	0,133	21,8	22,6	24,4	24,7	24,7	22,2	24,5	24,9	
19:35:41	052:36:01	0,132	25,0	-88,7	0,134	21,8	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:36:44	052:37:04	0,133	25,0	-78,8	0,134	21,8	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:37:44	052:38:04	0,134	25,0	-88,4	0,135	21,8	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:38:41	052:39:01	0,134	25,0	-83,8	0,136	21,8	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:39:44	052:40:04	0,135	25,0	-81,1	0,137	21,8	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:40:40	052:41:00	0,136	25,0	-88,7	0,137	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:41:40	052:42:00	0,128	25,0	-80,8	0,130	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:42:40	052:43:00	0,129	25,0	-89,9	0,131	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:43:42	052:44:02	0,130	25,0	-85,3	0,131	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:44:43	052:45:03	0,131	25,0	-76,3	0,132	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:45:40	052:46:00	0,131	25,0	-87,3	0,133	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:46:41	052:47:01	0,133	25,0	-83,4	0,133	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:47:44	052:48:04	0,133	25,0	-90,1	0,134	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:48:44	052:49:04	0,133	25,0	-85,8	0,135	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
19:49:46	052:50:06	0,134	25,0	-76,3	0,136	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	169 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	19:50:42	052:51:02	0,135	25,0	-87,7	0,137	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
	19:51:41	052:52:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	19:51:41	052:52:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	19:51:41	052:52:01	0,135	25,0	-83,9	0,137	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
	19:52:44	052:53:04	0,129	25,0	-89,7	0,130	21,9	22,6	24,4	24,7	24,7	22,2	24,5	25,0	
	19:53:40	052:54:00	0,129	25,0	-88,7	0,130	21,9	22,7	24,4	24,7	24,7	22,3	24,5	25,0	
	19:54:40	052:55:00	0,131	25,0	-79,9	0,131	21,9	22,7	24,4	24,7	24,7	22,3	24,5	25,0	
	19:55:40	052:56:00	0,130	25,0	-89,3	0,132	21,9	22,6	24,4	24,7	24,7	22,3	24,5	25,0	
	19:56:40	052:57:00	0,131	25,0	-84,5	0,133	21,9	22,6	24,4	24,7	24,8	22,3	24,5	25,0	
	19:57:41	052:58:01	0,132	25,0	-75,3	0,133	21,9	22,7	24,4	24,7	24,8	22,3	24,5	25,0	
	19:58:43	052:59:03	0,133	25,0	-88,3	0,134	21,9	22,7	24,4	24,8	24,7	22,3	24,5	25,0	
	19:59:44	053:00:04	0,133	25,0	-81,9	0,135	21,9	22,6	24,4	24,8	24,7	22,3	24,5	25,0	
	20:00:41	053:01:01	0,134	25,0	-90,9	0,136	21,9	22,6	24,4	24,8	24,7	22,3	24,5	25,0	
	20:01:43	053:02:03	0,134	25,0	-86,1	0,136	21,9	22,6	24,4	24,8	24,7	22,3	24,5	25,0	
	20:02:44	053:03:04	0,135	25,0	-76,6	0,137	21,9	22,6	24,5	24,8	24,7	22,3	24,5	25,0	
	20:03:40	053:04:00	0,136	25,0	-87,4	0,137	21,9	22,7	24,5	24,8	24,8	22,3	24,5	25,0	
	20:04:41	053:05:01	0,129	25,0	-83,3	0,130	21,9	22,6	24,5	24,8	24,8	22,3	24,5	25,0	
	20:05:44	053:06:04	0,130	25,0	-89,9	0,131	21,9	22,6	24,5	24,8	24,8	22,3	24,6	25,0	
	20:06:41	053:07:01	0,131	25,0	-86,3	0,132	21,9	22,6	24,5	24,8	24,8	22,3	24,5	25,0	
	20:07:43	053:08:03	0,131	25,0	-77,6	0,132	21,9	22,7	24,5	24,8	24,8	22,3	24,5	25,0	
	20:08:43	053:09:03	0,132	25,0	-88,7	0,133	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,0	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	170 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:09:44	053:10:04	0,132	25,0	-81,5	0,134	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,0	
	20:10:41	053:11:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:10:41	053:11:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:10:41	053:11:01	0,133	25,0	-87,1	0,135	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:11:42	053:12:02	0,134	25,0	-88,1	0,135	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,0	
	20:12:41	053:13:01	0,134	25,0	-78,0	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,0	
	20:13:41	053:14:01	0,135	25,0	-88,0	0,137	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,0	
	20:14:41	053:15:01	0,136	25,0	-83,9	0,137	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:15:41	053:16:01	0,129	25,0	-82,3	0,130	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:16:41	053:17:01	0,129	25,0	-88,6	0,131	21,9	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:17:41	053:18:01	0,131	25,0	-80,2	0,131	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,0	
	20:18:43	053:19:03	0,131	25,0	-89,8	0,132	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:19:40	053:20:00	0,132	25,0	-86,5	0,133	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:20:41	053:21:01	0,132	25,0	-78,7	0,133	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:21:43	053:22:03	0,133	25,0	-87,5	0,134	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:22:40	053:23:00	0,133	25,0	-83,4	0,135	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:23:42	053:24:02	0,135	25,0	-77,1	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:24:43	053:25:03	0,135	25,0	-88,6	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:25:43	053:26:03	0,135	25,0	-80,5	0,137	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:26:44	053:27:04	0,128	25,0	-90,0	0,130	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:27:41	053:28:01	0,129	25,0	-86,9	0,130	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 171 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:28:41	053:29:01	0,130	25,0	-77,8	0,131	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:29:42	053:30:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:29:42	053:30:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:29:42	053:30:02	0,131	25,0	-89,0	0,132	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:30:43	053:31:03	0,131	25,0	-82,3	0,133	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:31:44	053:32:04	0,132	25,0	-89,8	0,133	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:32:40	053:33:00	0,133	25,0	-86,5	0,134	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:33:41	053:34:01	0,133	25,0	-76,9	0,135	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:34:42	053:35:02	0,134	25,0	-87,4	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:35:43	053:36:03	0,134	25,0	-83,0	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:36:40	053:37:00	0,136	25,0	-90,7	0,137	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:37:40	053:38:00	0,135	25,0	-87,0	0,138	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:38:40	053:39:00	0,129	25,0	-77,4	0,130	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:39:41	053:40:01	0,129	25,0	-88,2	0,131	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:40:43	053:41:03	0,130	25,0	-82,4	0,132	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:41:43	053:42:03	0,132	25,0	-82,8	0,132	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:42:43	053:43:03	0,131	25,0	-88,3	0,133	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:43:44	053:44:04	0,132	25,0	-77,9	0,134	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:44:44	053:45:04	0,133	25,0	-88,4	0,135	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:45:41	053:46:01	0,134	25,0	-83,9	0,135	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:46:43	053:47:03	0,134	25,0	-88,6	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 172 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	20:47:40	053:48:00	0,135	25,0	-88,9	0,137	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:48:41	053:49:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	20:48:41	053:49:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	20:48:41	053:49:01	0,135	25,0	-81,0	0,137	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:49:43	053:50:03	0,129	25,0	-89,3	0,130	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:50:43	053:51:03	0,130	25,0	-84,4	0,131	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:51:43	053:52:03	0,130	25,0	-77,6	0,131	22,0	22,7	24,5	24,8	24,8	22,4	24,6	25,1	
	20:52:45	053:53:05	0,131	25,0	-89,1	0,132	22,0	22,7	24,5	24,8	24,8	22,4	24,6	25,1	
	20:53:40	053:54:00	0,131	25,0	-82,5	0,133	22,0	22,8	24,5	24,8	24,8	22,4	24,6	25,1	
	20:54:42	053:55:02	0,132	25,0	-88,4	0,134	22,0	22,8	24,5	24,8	24,8	22,4	24,6	25,1	
	20:55:42	053:56:02	0,133	25,0	-87,2	0,134	22,0	22,8	24,5	24,8	24,8	22,3	24,6	25,1	
	20:56:43	053:57:03	0,133	25,0	-76,7	0,135	22,0	22,8	24,5	24,8	24,8	22,3	24,6	25,1	
	20:57:44	053:58:04	0,134	25,0	-87,7	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,6	25,1	
	20:58:43	053:59:03	0,134	25,0	-82,6	0,136	22,0	22,7	24,5	24,8	24,8	22,3	24,7	25,1	
	20:59:44	054:00:04	0,136	25,0	-76,0	0,137	22,0	22,7	24,5	24,8	24,9	22,3	24,6	25,1	
	21:00:44	054:01:04	0,129	25,0	-88,4	0,130	22,0	22,8	24,5	24,8	24,9	22,3	24,6	25,1	
	21:01:44	054:02:04	0,129	25,0	-79,0	0,131	22,0	22,8	24,5	24,8	24,9	22,3	24,6	25,1	
	21:02:41	054:03:01	0,129	25,0	-88,9	0,131	22,0	22,8	24,5	24,8	24,9	22,3	24,6	25,2	
	21:03:43	054:04:03	0,130	25,0	-84,0	0,132	22,0	22,8	24,5	24,8	24,9	22,3	24,7	25,2	
	21:04:44	054:05:04	0,131	25,0	-76,8	0,133	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:05:44	054:06:04	0,132	25,0	-88,8	0,133	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 173 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	21:06:41	054:07:01	0,132	25,0	-80,6	0,134	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:07:43	054:08:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:07:43	054:08:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:07:43	054:08:03	0,133	25,0	-89,2	0,135	22,0	22,8	24,5	24,8	24,9	22,3	24,7	25,2	
	21:08:44	054:09:04	0,133	25,0	-83,9	0,136	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:09:41	054:10:01	0,135	25,0	-76,7	0,136	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:10:43	054:11:03	0,135	25,0	-88,8	0,137	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:11:40	054:12:00	0,135	25,0	-81,9	0,138	22,0	22,8	24,5	24,8	24,9	22,4	24,6	25,2	
	21:12:41	054:13:01	0,128	25,0	-89,3	0,130	22,0	22,8	24,5	24,8	24,9	22,4	24,6	25,2	
	21:13:44	054:14:04	0,129	25,0	-86,0	0,131	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:14:44	054:15:04	0,130	25,0	-85,6	0,132	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:15:40	054:16:00	0,130	25,0	-89,2	0,132	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:16:40	054:17:00	0,131	25,0	-81,3	0,133	22,0	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:17:40	054:18:00	0,132	25,0	-89,8	0,134	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:18:40	054:19:00	0,132	25,0	-85,4	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:19:42	054:20:02	0,133	25,0	-75,9	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:20:43	054:21:03	0,134	25,0	-88,2	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:21:43	054:22:03	0,134	25,0	-82,1	0,137	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:22:43	054:23:03	0,136	25,0	-91,9	0,138	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:23:44	054:24:04	0,128	25,0	-88,3	0,130	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:24:44	054:25:04	0,130	25,0	-77,8	0,131	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 174 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:25:41	054:26:01	0,129	25,0	-87,8	0,132	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:26:41	054:27:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:26:41	054:27:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:26:41	054:27:01	0,130	25,0	-83,4	0,132	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:27:42	054:28:02	0,131	25,0	-81,6	0,133	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:28:43	054:29:03	0,131	25,0	-88,6	0,134	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:29:44	054:30:04	0,133	25,0	-78,3	0,134	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:30:41	054:31:01	0,133	25,0	-88,5	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:31:41	054:32:01	0,133	25,0	-84,0	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:32:43	054:33:03	0,135	25,0	-86,3	0,137	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:33:43	054:34:03	0,135	25,0	-88,9	0,137	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:34:43	054:35:03	0,130	25,0	-80,3	0,131	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:35:40	054:36:00	0,129	25,0	-89,5	0,131	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:36:41	054:37:01	0,129	25,0	-84,9	0,131	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:37:42	054:38:02	0,130	25,0	-77,3	0,132	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:38:44	054:39:04	0,130	25,0	-88,8	0,133	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:39:40	054:40:00	0,132	25,0	-80,8	0,133	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:40:43	054:41:03	0,133	25,0	-89,2	0,134	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:41:43	054:42:03	0,132	25,0	-84,4	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:42:44	054:43:04	0,133	25,0	-76,5	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:43:40	054:44:00	0,134	25,0	-88,6	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 175 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	21:44:42	054:45:02	0,135	25,0	-81,5	0,137	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:45:42	054:46:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	21:45:42	054:46:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	21:45:42	054:46:02	0,135	25,0	-88,9	0,138	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:46:44	054:47:04	0,129	25,0	-86,9	0,130	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:47:41	054:48:01	0,129	25,0	-78,6	0,131	22,1	22,8	24,5	24,8	24,9	22,4	24,7	25,2	
	21:48:42	054:49:02	0,130	25,0	-89,1	0,132	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:49:42	054:50:02	0,130	25,0	-82,2	0,132	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:50:42	054:51:02	0,131	25,0	-89,0	0,133	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:51:42	054:52:02	0,131	25,0	-85,1	0,134	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	21:52:44	054:53:04	0,132	25,0	-78,3	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	21:53:40	054:54:00	0,133	25,0	-89,0	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:54:40	054:55:00	0,134	25,0	-79,8	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:55:42	054:56:02	0,134	25,0	-88,8	0,137	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:56:43	054:57:03	0,135	25,0	-83,8	0,137	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:57:44	054:58:04	0,129	25,0	-75,7	0,130	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	21:58:41	054:59:01	0,129	25,0	-88,9	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	21:59:41	055:00:01	0,129	25,0	-80,7	0,131	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:00:41	055:01:01	0,130	25,0	-89,9	0,132	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:01:43	055:02:03	0,131	25,0	-84,9	0,133	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:02:44	055:03:04	0,131	25,0	-76,1	0,134	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature:
Document issued by: Page: 176 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	22:03:44	055:04:04	0,132	25,0	-88,3	0,134	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:04:44	055:05:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:04:44	055:05:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:04:44	055:05:04	0,133	25,0	-82,3	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:05:45	055:06:05	0,133	25,0	-91,0	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:06:42	055:07:02	0,134	25,0	-87,8	0,136	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:07:44	055:08:04	0,135	25,0	-78,7	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:08:41	055:09:01	0,128	25,0	-88,9	0,130	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:09:44	055:10:04	0,128	25,0	-82,5	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:10:40	055:11:00	0,129	25,0	-86,1	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:11:40	055:12:00	0,130	25,0	-88,0	0,132	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:12:43	055:13:03	0,130	25,0	-77,3	0,133	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:13:44	055:14:04	0,131	25,0	-89,0	0,133	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:14:40	055:15:00	0,132	25,0	-83,3	0,134	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:15:41	055:16:01	0,132	25,0	-89,6	0,135	22,1	22,8	24,5	24,9	24,9	22,4	24,7	25,2	
	22:16:42	055:17:02	0,133	25,0	-84,5	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:17:40	055:18:00	0,134	25,0	-77,3	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:18:44	055:19:04	0,135	25,0	-89,1	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:19:40	055:20:00	0,135	25,0	-82,3	0,138	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:20:40	055:21:00	0,128	25,0	-86,9	0,130	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:21:41	055:22:01	0,130	25,0	-88,2	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	177 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	22:22:40	055:23:00	0,129	25,0	-77,9	0,132	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:23:41	055:24:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:23:41	055:24:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:23:41	055:24:01	0,130	25,0	-88,1	0,132	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:24:43	055:25:03	0,131	25,0	-82,6	0,133	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:25:43	055:26:03	0,132	25,0	-74,6	0,134	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:26:40	055:27:00	0,132	25,0	-88,6	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:27:40	055:28:00	0,133	25,0	-80,4	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:28:40	055:29:00	0,134	25,0	-89,7	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:29:42	055:30:02	0,134	25,0	-84,4	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:30:44	055:31:04	0,135	25,0	-75,4	0,138	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:31:40	055:32:00	0,128	25,0	-87,8	0,130	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:32:41	055:33:01	0,128	25,0	-82,6	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:33:42	055:34:02	0,129	25,0	-90,1	0,132	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:34:43	055:35:03	0,130	25,0	-87,3	0,132	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:35:43	055:36:03	0,131	25,0	-77,8	0,133	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:36:40	055:37:00	0,131	25,0	-87,7	0,134	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:37:42	055:38:02	0,132	25,0	-82,8	0,134	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:38:43	055:39:03	0,133	25,0	-83,7	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:39:40	055:40:00	0,134	25,0	-89,1	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:40:40	055:41:00	0,135	25,0	-79,5	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 178 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	22:41:42	055:42:02	0,135	25,0	-89,2	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:42:44	055:43:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	22:42:44	055:43:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	22:42:44	055:43:04	0,128	25,0	-82,9	0,130	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:43:40	055:44:00	0,128	25,0	-87,2	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:44:41	055:45:01	0,129	25,0	-87,9	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:45:43	055:46:03	0,131	25,0	-77,3	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,2	
	22:46:44	055:47:04	0,130	25,0	-87,8	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:47:42	055:48:02	0,131	25,0	-83,4	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:48:43	055:49:03	0,132	25,0	-88,5	0,134	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:49:40	055:50:00	0,132	25,0	-88,8	0,135	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:50:42	055:51:02	0,133	25,0	-79,2	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:51:45	055:52:05	0,134	25,0	-88,8	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,2	
	22:52:41	055:53:01	0,134	25,0	-84,5	0,137	22,1	22,8	24,6	24,9	24,9	22,4	24,8	25,3	
	22:53:42	055:54:02	0,132	25,0	-79,8	0,133	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:54:43	055:55:03	0,129	25,0	-88,8	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	22:55:41	055:56:01	0,129	25,0	-80,9	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:56:43	055:57:03	0,131	25,0	-89,3	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:57:44	055:58:04	0,130	25,0	-84,3	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:58:40	055:59:00	0,131	25,0	-92,1	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	22:59:43	056:00:03	0,131	25,0	-88,8	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 179 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:00:43	056:01:03	0,133	25,0	-79,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,7	25,3	
	23:01:44	056:02:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:01:44	056:02:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:01:44	056:02:04	0,133	25,0	-88,3	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,7	25,3	
	23:02:40	056:03:00	0,133	25,0	-85,2	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	23:03:44	056:04:04	0,134	25,0	-88,8	0,137	22,2	22,8	24,6	24,9	24,9	22,4	24,7	25,3	
	23:04:44	056:05:04	0,135	25,0	-88,3	0,138	22,2	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:05:41	056:06:01	0,128	25,0	-79,0	0,130	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:06:44	056:07:04	0,129	25,0	-88,4	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:07:42	056:08:02	0,129	25,0	-84,3	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:08:44	056:09:04	0,130	25,0	-88,3	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:09:41	056:10:01	0,131	25,0	-89,1	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:10:42	056:11:02	0,131	25,0	-80,8	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:11:42	056:12:02	0,132	25,0	-89,1	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:12:43	056:13:03	0,133	25,0	-83,6	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:13:44	056:14:04	0,133	25,0	-76,2	0,136	22,1	22,8	24,6	24,9	24,9	22,4	24,8	25,3	
	23:14:40	056:15:00	0,134	25,0	-88,7	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:15:42	056:16:02	0,135	25,0	-81,9	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:16:43	056:17:03	0,128	25,0	-90,1	0,130	22,2	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:17:42	056:18:02	0,128	25,0	-86,0	0,130	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:18:40	056:19:00	0,129	25,0	-76,4	0,131	22,1	22,8	24,6	24,9	24,9	22,4	24,7	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 180 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	23:19:41	056:20:01	0,130	25,0	-88,1	0,132	22,1	22,8	24,6	24,9	25,0	22,4	24,7	25,3	
	23:20:42	056:21:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:20:42	056:21:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:20:42	056:21:02	0,130	25,0	-82,5	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:21:44	056:22:04	0,131	25,0	-83,5	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:22:41	056:23:01	0,132	25,0	-88,4	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:23:42	056:24:02	0,132	25,0	-78,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,7	25,3	
	23:24:41	056:25:01	0,133	25,0	-89,2	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:25:42	056:26:02	0,134	25,0	-83,7	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:26:42	056:27:02	0,134	25,0	-75,7	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:27:43	056:28:03	0,135	25,0	-88,3	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:28:40	056:29:00	0,129	25,0	-81,9	0,130	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:29:42	056:30:02	0,129	25,0	-91,6	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:30:43	056:31:03	0,129	25,0	-87,9	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:31:43	056:32:03	0,130	25,0	-77,7	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:32:44	056:33:04	0,131	25,0	-87,7	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:33:40	056:34:00	0,131	25,0	-83,8	0,134	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:34:44	056:35:04	0,132	25,0	-85,9	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:35:40	056:36:00	0,133	25,0	-88,8	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:36:42	056:37:02	0,133	25,0	-79,7	0,136	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:37:41	056:38:01	0,134	25,0	-88,6	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 181 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:38:42	056:39:02	0,135	25,0	-85,1	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:39:44	056:40:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:39:44	056:40:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:39:44	056:40:04	0,128	25,0	-90,5	0,130	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:40:44	056:41:04	0,128	25,0	-87,6	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:41:40	056:42:00	0,129	25,0	-78,3	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
	23:42:40	056:43:00	0,130	25,0	-88,2	0,132	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:43:40	056:44:00	0,130	25,0	-83,3	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:44:43	056:45:03	0,131	25,0	-89,7	0,134	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:45:44	056:46:04	0,132	25,0	-87,7	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:46:40	056:47:00	0,132	25,0	-78,8	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:47:41	056:48:01	0,133	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:48:41	056:49:01	0,134	25,0	-83,5	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:49:40	056:50:00	0,134	25,0	-75,6	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:50:41	056:51:01	0,128	25,0	-88,7	0,130	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:51:41	056:52:01	0,128	25,0	-81,4	0,131	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:52:42	056:53:02	0,129	25,0	-89,9	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:53:44	056:54:04	0,129	25,0	-84,2	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:54:42	056:55:02	0,130	25,0	-86,8	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:55:44	056:56:04	0,131	25,0	-88,9	0,133	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	23:56:46	056:57:06	0,132	25,0	-80,4	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus Signature:
Document issued by: Page: 182 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	23:57:43	056:58:03	0,132	25,0	-89,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:58:40	056:59:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	23:58:40	056:59:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	23:58:40	056:59:00	0,133	25,0	-86,8	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	23:59:41	057:00:01	0,133	25,0	-77,4	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,7	25,3	
2021-11-26	00:00:41	057:01:01	0,134	25,0	-88,6	0,137	22,2	22,9	24,6	24,9	24,9	22,5	24,7	25,3	
	00:01:41	057:02:01	0,135	25,0	-82,4	0,138	22,2	22,9	24,6	24,9	24,9	22,5	24,7	25,3	
	00:02:42	057:03:02	0,128	25,0	-90,3	0,130	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:03:40	057:04:00	0,129	25,0	-87,5	0,131	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:04:41	057:05:01	0,129	25,0	-79,9	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:05:43	057:06:03	0,130	25,0	-89,4	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:06:42	057:07:02	0,130	25,0	-81,7	0,133	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:07:42	057:08:02	0,131	25,0	-89,2	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:08:43	057:09:03	0,132	25,0	-84,3	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:09:43	057:10:03	0,132	25,0	-74,5	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:10:43	057:11:03	0,133	25,0	-89,2	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:11:44	057:12:04	0,134	25,0	-80,3	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:12:40	057:13:00	0,134	25,0	-88,9	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:13:42	057:14:02	0,134	25,0	-84,3	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:14:40	057:15:00	0,129	25,0	-76,8	0,131	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:15:44	057:16:04	0,129	25,0	-88,7	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 183 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	00:16:40	057:17:00	0,130	25,0	-80,5	0,132	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:17:43	057:18:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:17:43	057:18:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:17:43	057:18:03	0,130	25,0	-89,4	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:18:41	057:19:01	0,131	25,0	-84,9	0,133	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:19:40	057:20:00	0,132	25,0	-77,4	0,134	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:20:41	057:21:01	0,132	25,0	-88,9	0,135	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:21:41	057:22:01	0,133	25,0	-81,8	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:22:43	057:23:03	0,134	25,0	-89,3	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:23:44	057:24:04	0,135	25,0	-85,2	0,137	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:24:40	057:25:00	0,135	25,0	-88,3	0,138	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:25:41	057:26:01	0,128	25,0	-89,1	0,130	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:26:43	057:27:03	0,129	25,0	-80,2	0,131	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:27:44	057:28:04	0,130	25,0	-88,7	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:28:41	057:29:01	0,130	25,0	-84,6	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:29:44	057:30:04	0,130	25,0	-81,8	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:30:44	057:31:04	0,131	25,0	-88,8	0,134	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:31:44	057:32:04	0,132	25,0	-81,0	0,135	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:32:40	057:33:00	0,132	25,0	-90,2	0,135	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:33:41	057:34:01	0,133	25,0	-86,2	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:34:41	057:35:01	0,134	25,0	-76,9	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 184 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	00:35:43	057:36:03	0,135	25,0	-88,1	0,137	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:36:44	057:37:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:36:44	057:37:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:36:44	057:37:04	0,128	25,0	-82,2	0,130	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:37:41	057:38:01	0,128	25,0	-90,6	0,131	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	00:38:42	057:39:02	0,129	25,0	-87,8	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:39:44	057:40:04	0,131	25,0	-77,3	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:40:40	057:41:00	0,130	25,0	-88,3	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:41:41	057:42:01	0,131	25,0	-84,4	0,133	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:42:42	057:43:02	0,132	25,0	-76,6	0,134	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:43:42	057:44:02	0,132	25,0	-88,7	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:44:42	057:45:02	0,134	25,0	-79,0	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:45:43	057:46:03	0,134	25,0	-88,5	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:46:44	057:47:04	0,134	25,0	-83,7	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	00:47:44	057:48:04	0,128	25,0	-75,3	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:48:40	057:49:00	0,128	25,0	-88,8	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:49:41	057:50:01	0,129	25,0	-80,8	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:50:42	057:51:02	0,130	25,0	-89,2	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:51:43	057:52:03	0,130	25,0	-84,1	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:52:44	057:53:04	0,131	25,0	-77,6	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:53:40	057:54:00	0,131	25,0	-89,0	0,134	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 185 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	00:54:40	057:55:00	0,133	25,0	-80,5	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:55:42	057:56:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	00:55:42	057:56:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	00:55:42	057:56:02	0,133	25,0	-88,7	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	00:56:43	057:57:03	0,133	25,0	-84,8	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:57:43	057:58:03	0,134	25,0	-84,3	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:58:43	057:59:03	0,135	25,0	-88,6	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	00:59:44	058:00:04	0,127	25,0	-78,7	0,130	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:00:40	058:01:00	0,128	25,0	-89,0	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:01:41	058:02:01	0,130	25,0	-83,8	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:02:43	058:03:03	0,130	25,0	-79,0	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:03:44	058:04:04	0,130	25,0	-88,8	0,133	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:04:40	058:05:00	0,132	25,0	-81,2	0,134	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:05:42	058:06:02	0,132	25,0	-89,1	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:06:42	058:07:02	0,132	25,0	-84,9	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:07:43	058:08:03	0,133	25,0	-78,9	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:08:44	058:09:04	0,134	25,0	-88,9	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:09:44	058:10:04	0,134	25,0	-78,8	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:10:44	058:11:04	0,127	25,0	-87,8	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:11:44	058:12:04	0,129	25,0	-83,5	0,131	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:12:44	058:13:04	0,130	25,0	-92,2	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 186 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	01:13:44	058:14:04	0,130	25,0	-88,7	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:14:41	058:15:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:14:41	058:15:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:14:41	058:15:01	0,130	25,0	-80,4	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:15:41	058:16:01	0,132	25,0	-89,6	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:16:42	058:17:02	0,131	25,0	-84,8	0,134	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:17:43	058:18:03	0,132	25,0	-76,5	0,135	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:18:45	058:19:05	0,133	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:19:41	058:20:01	0,133	25,0	-82,0	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:20:42	058:21:02	0,134	25,0	-90,5	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:21:40	058:22:00	0,135	25,0	-87,3	0,138	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:22:42	058:23:02	0,127	25,0	-77,3	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:23:44	058:24:04	0,128	25,0	-88,5	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:24:44	058:25:04	0,129	25,0	-82,4	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:25:44	058:26:04	0,130	25,0	-87,3	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:26:44	058:27:04	0,131	25,0	-88,0	0,133	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:27:41	058:28:01	0,131	25,0	-78,7	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:28:41	058:29:01	0,132	25,0	-89,5	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:29:41	058:30:01	0,132	25,0	-84,7	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:30:41	058:31:01	0,134	25,0	-77,4	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:31:41	058:32:01	0,135	25,0	-88,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 187 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:32:42	058:33:02	0,135	25,0	-81,7	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:33:40	058:34:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:33:40	058:34:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:33:40	058:34:00	0,127	25,0	-88,3	0,130	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:34:40	058:35:00	0,128	25,0	-87,9	0,131	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:35:40	058:36:00	0,129	25,0	-77,7	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:36:41	058:37:01	0,130	25,0	-87,7	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:37:42	058:38:02	0,130	25,0	-83,2	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:38:43	058:39:03	0,131	25,0	-90,1	0,134	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:39:44	058:40:04	0,132	25,0	-85,2	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:40:41	058:41:01	0,132	25,0	-78,2	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:41:43	058:42:03	0,133	25,0	-89,3	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:42:43	058:43:03	0,133	25,0	-81,8	0,136	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:43:43	058:44:03	0,134	25,0	-89,7	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:44:44	058:45:04	0,135	25,0	-84,7	0,138	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:45:41	058:46:01	0,128	25,0	-76,8	0,130	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:46:41	058:47:01	0,128	25,0	-87,4	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:47:41	058:48:01	0,129	25,0	-83,0	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	01:48:43	058:49:03	0,131	25,0	-91,8	0,133	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:49:41	058:50:01	0,130	25,0	-88,6	0,133	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:50:41	058:51:01	0,132	25,0	-78,7	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 188 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	01:51:42	058:52:02	0,132	25,0	-88,3	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	01:52:42	058:53:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	01:52:42	058:53:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	01:52:42	058:53:02	0,132	25,0	-84,6	0,135	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	01:53:42	058:54:02	0,133	25,0	-91,6	0,136	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:54:43	058:55:03	0,134	25,0	-87,7	0,137	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:55:43	058:56:03	0,134	25,0	-77,8	0,137	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:56:40	058:57:00	0,127	25,0	-87,7	0,130	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,4	
	01:57:44	058:58:04	0,129	25,0	-83,2	0,131	22,2	22,9	24,6	25,0	24,9	22,5	24,8	25,3	
	01:58:40	058:59:00	0,129	25,0	-90,8	0,131	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	01:59:41	059:00:01	0,129	25,0	-86,5	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:00:41	059:01:01	0,130	25,0	-76,3	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:01:42	059:02:02	0,131	25,0	-88,0	0,134	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:02:44	059:03:04	0,131	25,0	-81,9	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:03:40	059:04:00	0,132	25,0	-91,2	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:04:42	059:05:02	0,133	25,0	-88,5	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:05:43	059:06:03	0,133	25,0	-78,3	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:06:43	059:07:03	0,134	25,0	-88,1	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:07:40	059:08:00	0,135	25,0	-84,7	0,138	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:08:44	059:09:04	0,127	25,0	-77,0	0,130	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	02:09:44	059:10:04	0,128	25,0	-88,6	0,131	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 189 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:10:40	059:11:00	0,129	25,0	-80,4	0,132	22,2	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	02:11:41	059:12:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:11:41	059:12:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:11:41	059:12:01	0,129	25,0	-89,9	0,132	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	02:12:43	059:13:03	0,130	25,0	-84,9	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:13:42	059:14:02	0,132	25,0	-77,1	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:14:43	059:15:03	0,131	25,0	-88,9	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:15:41	059:16:01	0,132	25,0	-80,5	0,135	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	02:16:44	059:17:04	0,133	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	02:17:40	059:18:00	0,133	25,0	-84,6	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	02:18:41	059:19:01	0,134	25,0	-75,8	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:19:42	059:20:02	0,127	25,0	-89,0	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:20:41	059:21:01	0,129	25,0	-80,7	0,131	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	02:21:42	059:22:02	0,128	25,0	-88,1	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:22:43	059:23:03	0,129	25,0	-83,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:23:43	059:24:03	0,130	25,0	-90,9	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:24:43	059:25:03	0,131	25,0	-88,9	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:25:44	059:26:04	0,132	25,0	-78,7	0,134	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:26:40	059:27:00	0,132	25,0	-88,9	0,135	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:27:40	059:28:00	0,132	25,0	-84,0	0,135	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:28:40	059:29:00	0,134	25,0	-75,8	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 190 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	02:29:42	059:30:02	0,134	25,0	-88,9	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:30:42	059:31:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:30:42	059:31:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:30:42	059:31:02	0,134	25,0	-81,2	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:31:43	059:32:03	0,128	25,0	-90,2	0,130	22,1	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	02:32:44	059:33:04	0,128	25,0	-85,9	0,131	22,1	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	02:33:40	059:34:00	0,129	25,0	-76,9	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:34:41	059:35:01	0,130	25,0	-87,4	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:35:43	059:36:03	0,130	25,0	-83,0	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:36:40	059:37:00	0,132	25,0	-91,8	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:37:41	059:38:01	0,131	25,0	-88,7	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:38:42	059:39:02	0,132	25,0	-78,4	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:39:43	059:40:03	0,133	25,0	-87,9	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:40:40	059:41:00	0,134	25,0	-84,4	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:41:40	059:42:00	0,135	25,0	-79,4	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:42:40	059:43:00	0,134	25,0	-89,0	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:43:41	059:44:01	0,128	25,0	-78,8	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:44:42	059:45:02	0,128	25,0	-88,4	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:45:43	059:46:03	0,129	25,0	-83,7	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	02:46:41	059:47:01	0,130	25,0	-83,9	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:47:42	059:48:02	0,131	25,0	-88,5	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58
Product no.: divers Product name:
Batch no.: Batch name:
Batch size: 20R Batch volume per unit: 5-10
Note: PDA Formulierungen aggressives Rezept
Document issued at: 2021-11-26 10:54:17 LPCplus
Document issued by: Signature: Page: 191 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	02:48:44	059:49:04	0,131	25,0	-78,3	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:49:40	059:50:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	02:49:40	059:50:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	02:49:40	059:50:00	0,132	25,0	-87,9	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:50:40	059:51:00	0,133	25,0	-83,7	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:51:44	059:52:04	0,133	25,0	-88,7	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:52:44	059:53:04	0,134	25,0	-88,5	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:53:41	059:54:01	0,134	25,0	-79,1	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:54:42	059:55:02	0,128	25,0	-88,7	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:55:44	059:56:04	0,128	25,0	-83,9	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:56:44	059:57:04	0,129	25,0	-77,4	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:57:40	059:58:00	0,129	25,0	-88,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:58:41	059:59:01	0,130	25,0	-81,2	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	02:59:43	060:00:03	0,131	25,0	-90,3	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:00:44	060:01:04	0,131	25,0	-86,5	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:01:45	060:02:05	0,132	25,0	-76,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:02:41	060:03:01	0,133	25,0	-87,8	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:03:41	060:04:01	0,133	25,0	-83,0	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:04:42	060:05:02	0,134	25,0	-90,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:05:43	060:06:03	0,135	25,0	-88,3	0,138	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:06:40	060:07:00	0,128	25,0	-79,2	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 192 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	03:07:40	060:08:00	0,129	25,0	-88,3	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:08:40	060:09:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:08:40	060:09:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:08:40	060:09:00	0,129	25,0	-82,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:09:43	060:10:03	0,130	25,0	-90,5	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:10:41	060:11:01	0,130	25,0	-88,1	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:11:41	060:12:01	0,131	25,0	-78,7	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:12:43	060:13:03	0,132	25,0	-88,1	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:13:40	060:14:00	0,133	25,0	-84,1	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:14:43	060:15:03	0,133	25,0	-80,6	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:15:42	060:16:02	0,134	25,0	-88,8	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:16:42	060:17:02	0,134	25,0	-80,4	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:17:43	060:18:03	0,128	25,0	-89,6	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:18:40	060:19:00	0,128	25,0	-85,6	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:19:42	060:20:02	0,129	25,0	-77,8	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:20:40	060:21:00	0,129	25,0	-89,0	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:21:40	060:22:00	0,130	25,0	-81,5	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:22:42	060:23:02	0,131	25,0	-90,1	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:23:43	060:24:03	0,131	25,0	-85,8	0,134	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:24:43	060:25:03	0,132	25,0	-91,4	0,135	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:25:44	060:26:04	0,133	25,0	-89,2	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	193 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:26:44	060:27:04	0,133	25,0	-81,3	0,136	22,1	22,9	24,6	24,9	24,9	22,5	24,8	25,3	
	03:27:41	060:28:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:27:41	060:28:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:27:41	060:28:01	0,134	25,0	-90,0	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:28:42	060:29:02	0,134	25,0	-85,6	0,138	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:29:44	060:30:04	0,127	25,0	-77,1	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:30:41	060:31:01	0,128	25,0	-89,1	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:31:44	060:32:04	0,130	25,0	-82,6	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:32:41	060:33:01	0,129	25,0	-88,7	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:33:43	060:34:03	0,130	25,0	-87,2	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	03:34:46	060:35:06	0,131	25,0	-77,1	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:35:42	060:36:02	0,131	25,0	-87,7	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:36:44	060:37:04	0,132	25,0	-83,3	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:37:40	060:38:00	0,133	25,0	-91,0	0,136	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:38:40	060:39:00	0,133	25,0	-87,0	0,136	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:39:42	060:40:02	0,134	25,0	-77,4	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:40:44	060:41:04	0,130	25,0	-88,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:41:40	060:42:00	0,128	25,0	-82,8	0,130	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:42:41	060:43:01	0,128	25,0	-91,0	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:43:43	060:44:03	0,129	25,0	-87,0	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:44:44	060:45:04	0,130	25,0	-78,4	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 194 / 209

Martin Christ GmbH
USB flash drive
Event 21112303.00



	03:45:40	060:46:00	0,130	25,0	-88,5	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:46:41	060:47:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	03:46:41	060:47:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	03:46:41	060:47:01	0,131	25,0	-84,6	0,134	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:47:43	060:48:03	0,132	25,0	-91,1	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:48:43	060:49:03	0,132	25,0	-87,2	0,135	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:49:45	060:50:05	0,133	25,0	-76,8	0,136	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:50:40	060:51:00	0,134	25,0	-88,3	0,137	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:51:42	060:52:02	0,134	25,0	-82,3	0,137	22,2	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:52:43	060:53:03	0,127	25,0	-91,1	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:53:44	060:54:04	0,128	25,0	-87,2	0,131	22,1	22,9	24,6	24,9	24,9	22,4	24,8	25,3	
	03:54:40	060:55:00	0,128	25,0	-78,4	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:55:41	060:56:01	0,129	25,0	-88,2	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:56:42	060:57:02	0,130	25,0	-84,1	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:57:43	060:58:03	0,131	25,0	-88,9	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:58:45	060:59:05	0,131	25,0	-88,7	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	03:59:41	061:00:01	0,133	25,0	-79,0	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:00:41	061:01:01	0,132	25,0	-88,5	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:01:41	061:02:01	0,133	25,0	-85,0	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:02:41	061:03:01	0,134	25,0	-87,3	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:03:41	061:04:01	0,134	25,0	-89,0	0,138	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 195 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	04:04:42	061:05:02	0,127	25,0	-80,0	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:05:44	061:06:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:05:44	061:06:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:05:44	061:06:04	0,128	25,0	-88,8	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:06:40	061:07:00	0,129	25,0	-85,0	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:07:41	061:08:01	0,129	25,0	-87,6	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:08:44	061:09:04	0,130	25,0	-89,1	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:09:41	061:10:01	0,131	25,0	-79,6	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:10:41	061:11:01	0,131	25,0	-88,1	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:11:43	061:12:03	0,133	25,0	-83,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:12:43	061:13:03	0,133	25,0	-90,3	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:13:40	061:14:00	0,133	25,0	-86,4	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:14:40	061:15:00	0,135	25,0	-78,0	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:15:40	061:16:00	0,128	25,0	-89,3	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:16:40	061:17:00	0,127	25,0	-81,5	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:17:42	061:18:02	0,128	25,0	-89,3	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:18:44	061:19:04	0,129	25,0	-84,0	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:19:41	061:20:01	0,129	25,0	-76,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:20:43	061:21:03	0,130	25,0	-88,9	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:21:41	061:22:01	0,132	25,0	-81,9	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:22:44	061:23:04	0,133	25,0	-90,5	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 196 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	04:23:41	061:24:01	0,132	25,0	-87,7	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:24:42	061:25:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:24:42	061:25:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:24:42	061:25:02	0,133	25,0	-77,7	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:25:43	061:26:03	0,134	25,0	-87,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:26:44	061:27:04	0,134	25,0	-83,5	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:27:44	061:28:04	0,128	25,0	-91,7	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:28:45	061:29:05	0,128	25,0	-88,4	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:29:42	061:30:02	0,129	25,0	-78,4	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:30:43	061:31:03	0,129	25,0	-87,7	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:31:40	061:32:00	0,130	25,0	-82,7	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:32:41	061:33:01	0,130	25,0	-91,8	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:33:42	061:34:02	0,131	25,0	-88,8	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:34:42	061:35:02	0,132	25,0	-79,9	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:35:42	061:36:02	0,133	25,0	-89,0	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:36:40	061:37:00	0,133	25,0	-84,4	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:37:41	061:38:01	0,135	25,0	-76,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:38:43	061:39:03	0,134	25,0	-88,9	0,138	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:39:40	061:40:00	0,127	25,0	-81,8	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:40:42	061:41:02	0,129	25,0	-89,3	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:41:40	061:42:00	0,129	25,0	-87,5	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 197 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	04:42:40	061:43:00	0,129	25,0	-77,9	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:43:41	061:44:01	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	04:43:41	061:44:01	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	04:43:41	061:44:01	0,130	25,0	-89,2	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:44:43	061:45:03	0,131	25,0	-81,8	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:45:43	061:46:03	0,132	25,0	-89,7	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:46:43	061:47:03	0,132	25,0	-85,1	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:47:40	061:48:00	0,133	25,0	-76,3	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:48:43	061:49:03	0,133	25,0	-88,4	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:49:40	061:50:00	0,134	25,0	-82,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:50:41	061:51:01	0,127	25,0	-91,4	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:51:42	061:52:02	0,128	25,0	-87,9	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:52:44	061:53:04	0,128	25,0	-77,1	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:53:40	061:54:00	0,129	25,0	-88,1	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:54:41	061:55:01	0,130	25,0	-84,0	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:55:43	061:56:03	0,131	25,0	-90,7	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:56:43	061:57:03	0,131	25,0	-88,8	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:57:44	061:58:04	0,131	25,0	-78,4	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:58:40	061:59:00	0,132	25,0	-88,5	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	04:59:41	062:00:01	0,133	25,0	-84,5	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	05:00:41	062:01:01	0,133	25,0	-85,6	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 198 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:01:43	062:02:03	0,134	25,0	-88,3	0,138	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	05:02:44	062:03:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:02:44	062:03:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:02:44	062:03:04	0,127	25,0	-77,8	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:03:40	062:04:00	0,128	25,0	-87,7	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:04:40	062:05:00	0,129	25,0	-83,3	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:05:44	062:06:04	0,129	25,0	-78,2	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:06:44	062:07:04	0,130	25,0	-88,7	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:07:41	062:08:01	0,130	25,0	-79,3	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:08:43	062:09:03	0,131	25,0	-88,6	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:09:44	062:10:04	0,133	25,0	-84,5	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:10:44	062:11:04	0,133	25,0	-76,8	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:11:41	062:12:01	0,134	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:12:42	062:13:02	0,135	25,0	-79,2	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:13:44	062:14:04	0,127	25,0	-88,0	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:14:40	062:15:00	0,127	25,0	-83,5	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:15:41	062:16:01	0,128	25,0	-75,7	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:16:41	062:17:01	0,129	25,0	-88,8	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:17:42	062:18:02	0,129	25,0	-79,0	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:18:43	062:19:03	0,130	25,0	-88,3	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:19:45	062:20:05	0,130	25,0	-84,0	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	199 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	05:20:40	062:21:00	0,131	25,0	-77,4	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:21:40	062:22:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:21:40	062:22:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:21:40	062:22:00	0,132	25,0	-88,9	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:22:41	062:23:01	0,133	25,0	-79,2	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:23:43	062:24:03	0,133	25,0	-88,1	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:24:43	062:25:03	0,134	25,0	-84,0	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:25:43	062:26:03	0,128	25,0	-77,5	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:26:44	062:27:04	0,128	25,0	-88,7	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:27:41	062:28:01	0,129	25,0	-79,3	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:28:41	062:29:01	0,130	25,0	-88,9	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:29:43	062:30:03	0,130	25,0	-83,5	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:30:43	062:31:03	0,131	25,0	-77,1	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:31:40	062:32:00	0,131	25,0	-88,9	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:32:41	062:33:01	0,131	25,0	-79,5	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:33:43	062:34:03	0,132	25,0	-88,3	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:34:43	062:35:03	0,133	25,0	-84,1	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:35:43	062:36:03	0,134	25,0	-76,6	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:36:44	062:37:04	0,134	25,0	-88,8	0,138	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:37:41	062:38:01	0,128	25,0	-79,8	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:38:41	062:39:01	0,128	25,0	-89,1	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 200 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	05:39:43	062:40:03	0,128	25,0	-84,0	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:40:43	062:41:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:40:43	062:41:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:40:43	062:41:03	0,129	25,0	-77,2	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:41:44	062:42:04	0,130	25,0	-89,2	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:42:44	062:43:04	0,131	25,0	-79,6	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:43:44	062:44:04	0,131	25,0	-88,1	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:44:44	062:45:04	0,132	25,0	-84,0	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:45:44	062:46:04	0,133	25,0	-77,2	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:46:41	062:47:01	0,133	25,0	-88,8	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:47:42	062:48:02	0,134	25,0	-78,9	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:48:42	062:49:02	0,127	25,0	-88,2	0,130	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:49:42	062:50:02	0,128	25,0	-84,3	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:50:44	062:51:04	0,128	25,0	-90,6	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:51:41	062:52:01	0,129	25,0	-88,5	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:52:41	062:53:01	0,130	25,0	-77,8	0,133	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:53:44	062:54:04	0,130	25,0	-87,8	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	05:54:40	062:55:00	0,131	25,0	-83,3	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:55:40	062:56:00	0,131	25,0	-77,1	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:56:43	062:57:03	0,132	25,0	-88,6	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:57:46	062:58:06	0,133	25,0	-77,4	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 201 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	05:58:42	062:59:02	0,133	25,0	-88,3	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	05:59:43	063:00:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	05:59:43	063:00:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	05:59:43	063:00:03	0,134	25,0	-84,4	0,138	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:00:40	063:01:00	0,128	25,0	-80,1	0,130	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:01:43	063:02:03	0,128	25,0	-88,5	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:02:40	063:03:00	0,128	25,0	-78,8	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:03:46	063:04:06	0,129	25,0	-87,9	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:04:41	063:05:01	0,129	25,0	-84,4	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:05:43	063:06:03	0,130	25,0	-86,8	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:06:44	063:07:04	0,131	25,0	-88,8	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:07:41	063:08:01	0,132	25,0	-78,8	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:08:42	063:09:02	0,133	25,0	-87,8	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:09:43	063:10:03	0,133	25,0	-83,0	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:10:44	063:11:04	0,134	25,0	-84,6	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:11:41	063:12:01	0,134	25,0	-88,9	0,138	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:12:41	063:13:01	0,127	25,0	-78,6	0,130	22,2	23,0	24,6	24,9	25,0	22,5	24,8	25,3	
	06:13:41	063:14:01	0,129	25,0	-88,4	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:14:42	063:15:02	0,128	25,0	-84,1	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	06:15:41	063:16:01	0,130	25,0	-90,4	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:16:43	063:17:03	0,131	25,0	-88,9	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	202 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:17:45	063:18:05	0,131	25,0	-79,8	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:18:40	063:19:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:18:40	063:19:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:18:40	063:19:00	0,131	25,0	-89,3	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:19:40	063:20:00	0,132	25,0	-84,7	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:20:40	063:21:00	0,132	25,0	-76,8	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:21:40	063:22:00	0,133	25,0	-88,7	0,137	22,1	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	06:22:40	063:23:00	0,134	25,0	-81,5	0,137	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	06:23:42	063:24:02	0,127	25,0	-90,3	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:24:40	063:25:00	0,127	25,0	-87,4	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:25:40	063:26:00	0,128	25,0	-77,6	0,131	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:26:43	063:27:03	0,130	25,0	-87,5	0,132	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:27:40	063:28:00	0,129	25,0	-83,4	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:28:44	063:29:04	0,130	25,0	-91,7	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:29:43	063:30:03	0,131	25,0	-88,4	0,134	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:30:44	063:31:04	0,132	25,0	-78,0	0,135	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:31:44	063:32:04	0,132	25,0	-87,8	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:32:41	063:33:01	0,133	25,0	-83,2	0,136	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	06:33:43	063:34:03	0,134	25,0	-82,1	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	06:34:40	063:35:00	0,134	25,0	-89,0	0,138	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:35:43	063:36:03	0,127	25,0	-79,8	0,130	22,1	22,9	24,6	25,0	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10
 Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Signature: Page: 203 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:36:40	063:37:00	0,129	25,0	-88,2	0,131	22,1	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	06:37:42	063:38:02	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:37:42	063:38:02	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:37:42	063:38:02	0,128	25,0	-83,6	0,132	22,1	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	06:38:42	063:39:02	0,129	25,0	-76,9	0,132	22,1	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	06:39:44	063:40:04	0,130	25,0	-88,7	0,133	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	06:40:41	063:41:01	0,130	25,0	-78,7	0,134	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	06:41:42	063:42:02	0,131	25,0	-88,4	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:42:43	063:43:03	0,133	25,0	-82,9	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:43:43	063:44:03	0,132	25,0	-78,2	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:44:41	063:45:01	0,133	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:45:41	063:46:01	0,135	25,0	-79,4	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:46:43	063:47:03	0,127	25,0	-89,0	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:47:41	063:48:01	0,127	25,0	-83,5	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:48:40	063:49:00	0,128	25,0	-88,9	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:49:41	063:50:01	0,129	25,0	-88,7	0,132	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	06:50:45	063:51:05	0,129	25,0	-77,8	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:51:41	063:52:01	0,130	25,0	-88,5	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:52:42	063:53:02	0,131	25,0	-83,8	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:53:43	063:54:03	0,131	25,0	-82,4	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:54:44	063:55:04	0,132	25,0	-88,8	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	

Begin recording: 2021-11-23 14:59:40 End recording: 2021-11-26 8:22:40 Recording time (hhh:mm:ss): 65:22:58

Product no.: divers Product name:
 Batch no.: Batch name:
 Batch size: 20R Batch volume per unit: 5-10

Note: PDA Formulierungen aggressives Rezept

Document issued at: 2021-11-26 10:54:17 LPCplus

Document issued by: Signature: Page: 204 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	06:55:42	063:56:02	0,133	25,0	-80,1	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:56:43	063:57:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	06:56:43	063:57:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	06:56:43	063:57:03	0,133	25,0	-88,7	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	06:57:43	063:58:03	0,134	25,0	-84,2	0,137	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	06:58:40	063:59:00	0,127	25,0	-90,2	0,130	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	06:59:41	064:00:01	0,128	25,0	-88,2	0,131	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:00:42	064:01:02	0,128	25,0	-77,8	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:01:42	064:02:02	0,129	25,0	-87,5	0,132	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	07:02:44	064:03:04	0,129	25,0	-83,0	0,133	22,1	23,0	24,6	24,9	25,0	22,5	24,8	25,3	
	07:03:40	064:04:00	0,130	25,0	-85,4	0,134	22,2	23,0	24,6	24,9	25,0	22,5	24,8	25,3	
	07:04:41	064:05:01	0,131	25,0	-88,8	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:05:40	064:06:00	0,131	25,0	-79,6	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:06:43	064:07:03	0,132	25,0	-88,0	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:07:44	064:08:04	0,133	25,0	-82,9	0,136	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:08:44	064:09:04	0,134	25,0	-75,8	0,137	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	07:09:45	064:10:05	0,127	25,0	-88,7	0,130	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:10:40	064:11:00	0,128	25,0	-80,6	0,131	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:11:40	064:12:00	0,128	25,0	-90,2	0,132	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,4	
	07:12:40	064:13:00	0,129	25,0	-85,8	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:13:40	064:14:00	0,129	25,0	-76,8	0,133	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		
Document issued by:				Page:	205 / 209

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	07:14:43	064:15:03	0,130	25,0	-87,4	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:15:44	064:16:04	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:15:44	064:16:04	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:15:44	064:16:04	0,131	25,0	-82,7	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:16:41	064:17:01	0,131	25,0	-88,8	0,135	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:17:43	064:18:03	0,133	25,0	-88,6	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:18:40	064:19:00	0,133	25,0	-79,8	0,136	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:19:40	064:20:00	0,133	25,0	-88,9	0,137	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:20:44	064:21:04	0,134	25,0	-84,4	0,138	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:21:40	064:22:00	0,127	25,0	-77,1	0,130	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:22:43	064:23:03	0,128	25,0	-88,8	0,131	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:23:40	064:24:00	0,128	25,0	-80,4	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:24:40	064:25:00	0,129	25,0	-89,6	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:25:41	064:26:01	0,130	25,0	-84,6	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:26:41	064:27:01	0,130	25,0	-87,7	0,134	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:27:42	064:28:02	0,131	25,0	-89,2	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:28:42	064:29:02	0,132	25,0	-80,4	0,135	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	07:29:43	064:30:03	0,132	25,0	-88,5	0,136	22,1	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:30:43	064:31:03	0,133	25,0	-84,0	0,136	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:31:40	064:32:00	0,134	25,0	-78,4	0,137	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:32:41	064:33:01	0,135	25,0	-89,0	0,138	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 206 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



	07:33:42	064:34:02	0,127	25,0	-78,6	0,130	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:34:43	064:35:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:34:43	064:35:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:34:43	064:35:03	0,129	25,0	-87,9	0,131	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:35:40	064:36:00	0,128	25,0	-83,1	0,132	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,3	
	07:36:40	064:37:00	0,129	25,0	-77,5	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:37:41	064:38:01	0,130	25,0	-88,9	0,133	22,2	22,9	24,6	24,9	25,0	22,4	24,8	25,4	
	07:38:41	064:39:01	0,130	25,0	-80,3	0,134	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,3	
	07:39:44	064:40:04	0,132	25,0	-88,8	0,135	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	07:40:41	064:41:01	0,132	25,0	-84,0	0,135	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,4	
	07:41:43	064:42:03	0,132	25,0	-76,9	0,136	22,2	23,0	24,6	25,0	25,0	22,5	24,8	25,4	
	07:42:44	064:43:04	0,134	25,0	-89,0	0,137	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,4	
	07:43:44	064:44:04	0,134	25,0	-80,7	0,137	22,2	22,9	24,6	25,0	25,0	22,4	24,8	25,3	
	07:44:42	064:45:02	0,127	25,0	-89,7	0,130	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	07:45:44	064:46:04	0,127	25,0	-84,7	0,131	22,2	23,0	24,6	25,0	25,0	22,5	24,8	25,4	
	07:46:40	064:47:00	0,129	25,0	-77,5	0,131	22,2	23,0	24,6	25,0	25,0	22,5	24,8	25,4	
	07:47:41	064:48:01	0,130	25,0	-88,9	0,132	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	07:48:43	064:49:03	0,129	25,0	-80,8	0,133	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:49:40	064:50:00	0,130	25,0	-90,6	0,133	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:50:42	064:51:02	0,131	25,0	-86,8	0,134	22,2	22,9	24,6	25,0	25,0	22,5	24,9	25,3	
	07:51:43	064:52:03	0,131	25,0	-77,0	0,135	22,2	22,9	24,6	25,0	25,0	22,5	24,9	25,3	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 207 / 209
Document issued by:					

Martin Christ GmbH
USB flash drive
Event 21112303.00



	07:52:44	064:53:04	0,132	25,0	-87,7	0,136	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:53:40	064:54:00	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
	07:53:40	064:54:00	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
	07:53:40	064:54:00	0,133	25,0	-83,5	0,136	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,3	
	07:54:41	064:55:01	0,133	25,0	-92,0	0,137	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:55:42	064:56:02	0,134	25,0	-88,9	0,137	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:56:43	064:57:03	0,127	25,0	-79,5	0,130	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:57:40	064:58:00	0,128	25,0	-89,7	0,131	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:58:42	064:59:02	0,128	25,0	-84,8	0,131	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	07:59:43	065:00:03	0,129	25,0	-76,8	0,132	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	08:00:40	065:01:00	0,130	25,0	-88,9	0,133	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	08:01:41	065:02:01	0,131	25,0	-81,7	0,133	22,2	23,0	24,6	25,0	25,0	22,5	24,8	25,4	
	08:02:41	065:03:01	0,131	25,0	-89,2	0,134	22,2	23,0	24,6	24,9	25,0	22,5	24,8	25,4	
	08:03:40	065:04:00	0,131	25,0	-85,5	0,135	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	08:04:41	065:05:01	0,132	25,0	-78,1	0,135	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	08:05:42	065:06:02	0,133	25,0	-88,8	0,136	22,2	23,0	24,6	25,0	25,0	22,5	24,8	25,4	
	08:06:44	065:07:04	0,133	25,0	-79,2	0,137	22,2	23,0	24,6	24,9	25,0	22,5	24,8	25,4	
	08:07:40	065:08:00	0,134	25,0	-89,8	0,138	22,2	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
	08:08:42	065:09:02	0,127	25,0	-84,6	0,130	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	08:09:42	065:10:02	0,128	25,0	-78,5	0,131	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
	08:10:42	065:11:02	0,129	20,0	-83,8	0,132	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,4	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 208 / 209
Document issued by:					

**Martin Christ GmbH
USB flash drive
Event 21112303.00**



08:11:43	065:12:03	0,129	20,0	-74,1	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
08:12:43	065:13:03	Vac.Pir	SetShel	Ice col	Vac.cap	WTMplu1	WTMplu2	WTMplu3	WTMplu4	WTMplu5	WTMplu6	WTMplu7	WTMplu8	
08:12:43	065:13:03	mBar	Â°C	Â°C	mBar	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	Â°C	
08:12:43	065:13:03	0,130	20,0	-66,5	0,133	22,1	22,9	24,6	24,9	25,0	22,5	24,8	25,4	
08:13:44	065:14:04	0,132	---	-59,8	0,135	22,2	22,9	24,6	25,0	25,0	22,5	24,8	25,4	
08:14:42	065:15:02	1000	20,0	-54,1	1,18	22,2	23,0	24,7	25,0	25,0	22,6	24,9	25,4	
08:15:43	065:16:03	1000	20,0	-48,4	1,18	22,4	23,2	24,8	25,1	25,1	22,7	25,0	25,4	
08:16:44	065:17:04	1000	20,0	-43,5	1,18	22,6	23,4	24,8	25,2	25,1	23,0	25,1	25,4	
08:17:40	065:18:00	1000	20,0	-39,6	1,18	22,8	23,5	24,8	25,1	25,1	23,1	25,2	25,3	
08:18:41	065:19:01	1000	20,0	-35,8	1,18	---	23,6	24,8	25,1	25,1	23,1	25,3	25,3	
08:19:43	065:20:03	1000	20,0	-32,2	1,18	22,9	23,7	24,8	25,1	25,1	23,2	25,3	25,3	
08:20:41	065:21:01	1000	---	-29,2	1,18	22,9	23,7	24,8	25,0	25,2	23,2	25,4	25,2	
08:21:42	065:22:02	1000	---	-26,2	1,18	22,9	23,7	24,9	25,1	25,2	23,3	25,5	25,2	

Begin recording:	2021-11-23 14:59:40	End recording:	2021-11-26 8:22:40	Recording time (hhh:mm:ss):	65:22:58
Product no.:	divers	Product name:			
Batch no.:		Batch name:			
Batch size:	20R	Batch volume per unit:	5-10		
Note:	PDA Formulierungen aggressives Rezept				
Document issued at:	2021-11-26 10:54:17	LPCplus	Signature:		Page: 209 / 209
Document issued by:					