

Life Cycle Assessment (LCA) and Comparison of Five Sterile Injectable Platforms



Poster Presenter:

Laura Weidner, Head of Corporate Development & Quality, ApiJect

Introduction

Third Party, Peer-Reviewed Life Cycle Assessment (LCA) quantifies the environmental impact of five sterile drug delivery solutions, compares results to determine if there is a statistically significant difference in the potential environmental impact of options, and identifies key drivers of sustainability and future opportunities for designs, configurations, and locations.

Five Product Alternatives


- Single-dose glass vial with plastic syringe and two hypodermic needles for drug drawing and drug administration
- Multi-dose glass vials with (per dose) plastic syringe and two hypodermic needles for drug drawing and drug administration
- Glass prefilled syringe with luer needle assembly
- Glass prefilled syringe with staked needle assembly
- Blow-Fill-Seal (BFS) container and needle assembly (Prefilled ApiJect Injector)

International Organization for Standardization, Geneva

- ISO, 2006a. ISO 14040:2006 Environmental management - Life cycle assessment - Principles and framework.
- ISO, 2006b. ISO 14044:2006 Environmental management - Life cycle assessment - Requirements and guidelines.

Methodology and Analysis


- Utilized ISO standard for LCA (14040), defining standard phases to any analysis: (1) Goal and Scope (2); Life Cycle Inventory Analysis; (3) Life Cycle Impact Assessment; and (4) Interpretation.
- Evaluation was based on the delivery of a single-dose of 1 millimeter (mL) of drug, set as the functional unit of comparison for the study.
- Shelf-life, durability during transit, ease of use, and clinical efficacy were assumed to be equivalent based on design specifications and expert judgement.
- Use stage was assumed to be equivalent among all options.
- Drug wastage was chosen as a representative data set from Parvatker et al. (2019).



- Modeling scenarios tested process, location, and transportation assumptions.
- Data assessment quality included: reliability, completeness, temporal correlation, geographic correlation, and technological correlation.
- Analysis included all operational material and energy flows. No universal cut-off thresholds were applied.


Life Cycle Impact Assessment Results

Carbon Footprint Results




- Prefilled ApiJect Injector results of 38 g CO₂e/dose are lowest among options
- Biggest drivers are lightweighting and avoided sterilization

Direct Water Use Results



- Compared to Prefilled ApiJect Injector, other options use 12x-116x more water
- Biggest drivers are avoided washing and sterilization

Overall Environmental Impact Results



- Prefilled ApiJect Injectors has lowest environmental impacts across all EPA LCA categories

Conclusion

This peer-reviewed, independent study finds a substantial difference in the environmental impact comparing the Prefilled ApiJect Injector and other injection formats.

- Blow-Fill-Seal has lower environmental footprint than glass options in all stages of production.
- Compared to the Prefilled ApiJect Injector result of 38 g CO₂-eq per dose:
 - The estimated single-dose glass vial with plastic syringe impacts are ~125% higher per dose.
 - The glass prefilled syringe with needle assembly impacts are ~100% higher per dose.
 - The multi-dose glass vials with plastic syringe impacts are 65-75% higher per dose.
- Single dose glass vials use more than 100x the water than the Prefilled ApiJect Injector - 2 gallons versus 2 ounces per dose.
- The Prefilled ApiJect Injector uses less total plastic than the glass-based injection alternatives.
- The Prefilled ApiJect Injector is designed to improve product safety by preventing microbial contamination risks associated with multiple withdrawals from a vial.
- A typical single-dose glass vial requires 750 mL of direct water for cleaning and sterilization for a single 1 mL dose. In contrast, the Prefilled ApiJect Injector requires approximately 6.4 mL of water for the manufacturing process.
- The Prefilled ApiJect Injector has the lowest impact in all 10 Life-Cycle Assessment categories set by the EPA.

Acknowledgments

Authors

Matthew Eckelman, Associate Professor, Civil and Environmental Engineering, Northeastern University and Adjunct Associate Professor, Yale School of Public Health
Robert Litan, PhD, Economist and Attorney in Law, Economic Research, and Policy

Independent Review Panel

Terrie Boguski, Harmony Environmental (panel chair)
Cassandra Thiel, New York University
C. Jason Pierce, LCA Consultant and Certified Practitioner

For Further Information:
apiject.com



apiject™